A200

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 Instructions for use

Quality for life

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Instructions for use

Table of content

1	General information	77	4	Delivery / Preparing the wheelchair for use	93
1.1	Preface	77	4.1	Delivery	93
1.2	Intended use	78	4.2	Initial operation	94
1.3	Field of application	78	5	Transportation and storage	95
1.4	Service	79	5.1	Transport in a wheelchair accessible motor vehicle	95
2	Safety	79	5.1.1	Required accessories	96
2.1	Explanation of symbols	79	5.1.2	Using the product in the vehicle	96
2.2	Standards and directives	80	5.1.3	Restrictions for use	98
2.3	General safety instructions	80	5.2	Disassembling the power wheelchair	100
2.4	Safety requirements for transportation, storage and assembly	81	5.3	Storage	104
2.5	Safety requirements for operation	83	6	Operation	105
2.6	Safety requirements for care, maintenance and disposal	88	6.1	Adjustment possibilities	105
2.7	Requirements for the user	89	6.1.1	Backrest	105
2.8	Safety functions	89	6.1.2	Seat angle	106
2.9	Warning symbols and type plates	90	6.1.3	Side panel with armrest	106
3	Product description	93	6.1.4	Control panel	107
			6.1.5	Footrest	108
			6.2	Getting into and out of the power wheelchair	109
			6.2.1	From the side	109

6.2.2 From the front.....	110	8 Malfunctions / Troubleshooting	127
6.3 Control unit	110	8.1 Warning	128
6.3.1 Control panel	111	8.2 Errors	128
6.3.2 Switching on and off.....	112	9 Maintenance, cleaning and care	131
6.3.3 Driving function.....	112	9.1 Maintenance intervals.....	131
6.3.4 "Charge level" LED display	114	9.2 Changing the fuse	134
6.3.5 Range	115	9.3 Changing tyres/inner tube	134
6.3.6 Drive-away lock.....	115	9.4 Cleaning and care	136
6.4 Releasing and locking the brake.....	116	9.4.1 Disinfection	136
6.5 Batteries	118	10 Disposal	137
6.5.1 Charging.....	120	10.1 Disposal information.....	137
6.5.2 Charger	121	10.2 Information on re-use.....	137
7 Accessories	122	11 Legal Information	138
7.1 Mechanical back angle adjustment	123	11.1 Service life	138
7.1.1 Specialties for preparing power wheelchair for transport	124	11.2 Liability.....	139
7.2 Mechanically elevating footrest.....	125	11.3 CE conformity	139
7.3 Swing-away control panel holder.....	125	11.4 Warranty.....	139
7.4 Lap belt.....	126	11.5 Trademarks.....	139
7.5 Other options	127	12 Technical data	140

1 General information

INFORMATION

Last update: 2015-05-18

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

INFORMATION

- New information regarding product safety and product recalls can be obtained from the Customer Care Center (CCC) at oa@ottobock.com or from the manufacturer's service department (see inside back cover or back page for addresses).
- You can request this document as a PDF file from the Customer Care Center (CCC) at oa@ottobock.com or from the manufacturer's service department (see inside back cover or back page for addresses). It is possible to increase the display size of the PDF document.

- For further questions about the instructions for use, please contact the authorised personnel who issued the product to you.

1.1 Preface

These instructions for use provide the user as well as his or her attendants with all the required knowledge on the design, function, operation and maintenance of the A200 power wheelchair from Otto Bock Mobility Solutions GmbH. The instructions contain all the information that is needed to use the power wheelchair safely, to determine the possible causes of a malfunction and to help eliminate it.

Knowledge of these instructions for use is absolutely necessary for ensuring that the power wheelchair is used safely. Therefore, the user and his or her attendants must read the instructions for use thoroughly, especially the section on "Safety", before using the power wheelchair. This will ensure the complete utilisation of all the features of the power wheelchair.

1.2 Intended use

The power wheelchair is designed solely for people who are unable to walk or who have mobility problems in order to provide them with individual mobility in indoor areas.

The power wheelchair may only be combined with the options mentioned in these instructions for use. Ottobock assumes no liability for combinations with medical products and/or accessories from other manufacturers outside of the modular system. An option for the power wheelchair is to have it controlled by an attendant with the use of an attendant control.

Any other use is considered improper use. The manufacturer is not liable for any personal injury or damage to property resulting from improper use; in such cases, the user has sole liability.

The power wheelchair may only be used by persons trained to use it. Training users and attendants to use the power wheelchair is required for protecting persons from danger and for ensuring that the power wheelchair is operated safely and correctly.

The operational safety of the power wheelchair can only be ensured if it is used properly in accordance with the informa-

tion contained in these instructions for use. The user is ultimately responsible for accident-free operation.

1.3 Field of application

The power wheelchair is suitable for patients who have walking impediments or walking disabilities due to:

- Paralysis
- Loss of limbs (lower limb amputation)
- Defective or deformed limbs
- Joint contractures / defects
- Other diseases

The power wheelchair was specially designed for users who are able to move independently in such a power wheelchair.

Fitting considerations:

- Body size and body weight (max. load capacity 100 kg; standard seat Junior: 75 kg)
- Physical and psychological limitations
- Age of the patient
- Home conditions
- Environment

1.4 Service

INFORMATION

Service and repairs on the power wheelchair are always carried out by experts at specialist dealers who have been authorised and trained by Ottobock. Should any problems arise, please contact your power wheelchair supplier.

Should you have any questions or a problem that cannot be resolved despite using the instructions for use, please contact Ottobock Service (see back side of cover for address).

Ottobock endeavours to support their customers in all respects in order to keep them satisfied with their product for a long time.

2 Safety

2.1 Explanation of symbols

WARNING

Warnings regarding possible risks of severe accident or injury.

CAUTION

Warnings regarding possible risks of accident or injury.

NOTICE

Warnings regarding possible technical damage.

INFORMATION

Information regarding operation. Information for service personnel.

2.2 Standards and directives

All information on safety contained in these instructions for use refers to the currently valid national laws and regulations of the European Union. In other countries, compliance with the applicable laws and national regulations is required.

In addition to the safety instructions contained in these instructions for use, the user is required to observe and comply with the BGV (Employer's Liability Insurance Association regulations), UVV (accident prevention regulations), and environmental protection regulations. All information contained in these instructions for use must be complied with at all times without restrictions.

The power wheelchair has been constructed in accordance with the currently valid technical rules and is safe to operate. The safety of the power wheelchair is confirmed by the CE symbol and the declaration of conformity.

2.3 General safety instructions

WARNING

Risk of suffocation. Packaging materials must be kept out of reach of children.

CAUTION

Risk of accident and injury in case of failure to follow the safety instructions.

All safety instructions contained in these instructions for use and in all other applicable documents must be observed and complied with. The instructions for use must be available to the user at all times.

CAUTION

Risk of accident and injury due to improper use. The may only be used properly. The may only be used by persons trained to use it.

CAUTION

Risk of burns when near to fire. The back upholstery and seat cushion of the power wheelchair are not highly flammable, but can however catch fire. Therefore utmost caution near any sources of open flame or sparks, especially lit cigarettes, is required.

⚠ CAUTION**Risk of hypothermia or burns on wheelchair parts.**

Parts can become extremely hot or cold due to extreme temperatures.

Do not expose the product to any extreme temperatures (e.g. direct sunlight, sauna, extreme cold) in order to prevent injuries by coming into contact with the parts.

INFORMATION

Only use original manufacturer's options. The optional components may be mounted only as described here. Failure to comply will void the warranty.

2.4 Safety requirements for transportation, storage and assembly

Transportation and storage

⚠ WARNING**Risk of accident and injury due to improper use for transportation in wheelchair accessible motor vehicles.**

The power wheelchair is permissible for transporting passengers in motor vehicles when using safety components (such as lap belt) as well as appropriate restraint systems offered by Ottobock.

Observe the limitations regarding installed options (see section 5.1.3).

The power wheelchair may be used for transporting only one person. We recommend that wherever and whenever possible, power wheelchair users transfer in the seats installed in the motor vehicle and use the corresponding vehicle restraint system, because this is the only way to ensure optimal protection of the passengers in case of an accident.

For more information, please refer to our brochure on using Ottobock products in wheelchair accessible vehicles, order number 646D158.

The power wheelchair may generally be used in wheelchair accessible motor vehicles without seat tilt.

⚠ WARNING

Risk of accident and injury due to incorrect transportation in aeroplanes. In aeroplanes, the power wheelchair has to be transported in compliance with IATA regulations (International Air Transport Association). To this end, the fuse must always be removed and the battery connectors must be insulated so that they are short-circuit-proof prior to check-in of the power wheelchair as luggage.

Note that those batteries in particular which may leak or will not be transported upright must be removed and packaged so they cannot leak or short circuit.

For more information please visit the www.iata.org website. Ottobock recommends contacting the airline directly before every flight to obtain information regarding special transport regulations.

⚠ WARNING

Risk of damage as a result of improper transportation.

For transportation of the power wheelchair, only use hoisting devices with a sufficient capacity. The A200 must be secured in accordance with the regulations for the transportation device used. Only attach the tensioning straps to the corresponding attachment eyes.

During transport on lifting platforms or in lifts, the wheelchair controls must be switched off. Lock the brake. Ensure that the power wheelchair is centred on the lifting platform.

None of the power wheelchair's components, such as its anti-tippers or other components, may be in the danger area.

NOTICE

Battery damage due to total discharge. Remove the fuse for shipping or when the power wheelchair is not being used for an extended period of time.

INFORMATION

The tyres of the power wheelchair contain chemical substances that can react to other chemical substances (such as cleaning agents, acids).

For this reason, place a suitable mat underneath the wheelchair, if it is not going to be used over an extended period of time.

Assembly

CAUTION

Risk of injury due to unsecured screw connections. After loosening screw connections with thread lock, replace these screw connections with new ones or secure them with medium-strength thread lock substance (e.g. Euro-Lock A24.20).

After the power wheelchair has been set and adjusted, the attachment screws and/or nuts must be firmly retightened. During tightening observe torques when specified.

CAUTION

Risk of tipping as a result of incorrectly mounted anti-tipper. To ensure safe driving operation, the anti-tipper must have been mounted correctly and must be in proper condition.

INFORMATION

Prior to using the power wheelchair, all the necessary mechanical adaptations (e. g. mounting special controls) and software settings (e. g. programming the control) must be made to comply with the individual requirements and abilities of the user. The settings may only be made by trained specialists who have been authorised by Ottobock.

2.5 Safety requirements for operation

WARNING

Risk of accident and injury as a result of incorrect safety functions. The user is obliged to make sure that the power wheelchair and its safety functions are in safe and proper condition every time before using it.

The power wheelchair may only be operated if all safety functions, e.g. the automatic brakes, are functional. Brake failure can result in serious accidents with fatal injuries.

WARNING

Risk of accidents and injury due to incorrect configuration settings. Modified parameter settings in the configuration can lead to changes in driving characteristics. In particular, changes to the speed, acceleration, braking or joystick settings can lead to unexpected and therefore uncontrollable operating performance with a risk of accidents. Always test the driving characteristics of the power wheelchair after configuration / programming is complete. Programming must only be completed by authorised personnel. Neither Ottobock nor the control unit manufacturer are liable for damages caused by programming that was not properly / professionally adapted to the abilities of the wheelchair user.

⚠ WARNING

Risk of accident and injury if brake function is unavailable or unlocked. Take into account the missing brake function during unlocked brake especially when driving on slopes or inclines.

⚠ WARNING

Risk of injury if the power wheelchair tips over during driving. The power wheelchair has been approved for driving on inclines and slopes of no more than **7° (12 %)**. Navigating inclines above this percentage value is not permitted.

The critical obstacle height of the A200 is **15 mm**. It is not permitted to cross obstacles higher than **15 mm**.

Operating the power wheelchair on stairs is not permitted.

⚠ WARNING

Risk of accident if the power wheelchair tips over during driving. Reduce the driving speed when driving downhill (e.g. set it to speed level 1). When driving on inclines or slopes, it is not permissible to negotiate obstacles. Avoid getting into and out of the power wheelchair on inclines and slopes.

Before driving up slopes or over obstacles, a tilted seat must be lowered to its basic position and the back must be in the upright position.

Obstacles such as steps may only be crossed at a reduced speed (max. of 3 km/h). Always approach obstacles at a right angle and cross over them without stopping.

⚠ WARNING

Risk of tipping over when driving on inadequate surface. Using the power wheelchair on very slippery ground (such as icy surfaces) or on very coarse-grained surfaces (gravel or pebbles) is not permitted.

⚠ WARNING

Risk of tipping when using lifting platforms. During transport on lifting platforms or in lifts, the wheelchair controls must be switched off. Lock the brake.

⚠ WARNING

Risk of tipping due to centre of gravity shifting. Getting to know how the power wheelchair reacts when the centre of gravity shifts, for example, on slopes or inclines or when clearing obstacles like steps and curbs should be done only with assistance from another person before using the wheelchair for the first time.

⚠ WARNING

Risk of accident and injury as a result of lifting the wheelchair incorrectly. Attendants must lift the wheelchair only by parts that are firmly attached and not by the footrests or armrests.

⚠ CAUTION

Risk of accident when driving without experience. Driving without experience can lead to falls and other dangerous situations. Familiarise yourself with the new wheelchair by practising on even, straightforward terrain first.

⚠ CAUTION

Risk of accident due to uncontrolled driving behaviour. Uncontrolled movements can occur during the operation of the power wheelchair as a result of malfunction. In this case, please contact your authorised dealer immediately. If any faults, defects or other dangers that can lead to personal injury are detected, the power wheelchair must be put out of operation immediately.

⚠ CAUTION

Risk of accident and injury as a result of incorrect getting into and out of the power wheelchair. Switch off the wheelchair control before getting into and out of the power wheelchair. The footrests and armrests are not capable of bearing full body weight, and therefore must not be used for getting into or out of the wheelchair.

⚠ CAUTION

Risk of accident and injury if the wheelchair starts rolling. Releasing the brake may result in uncontrolled rolling of the power wheelchair. Make sure that the brake is engaged after parking the power wheelchair.

⚠ CAUTION

Risk of accident as a result of wearing improper clothing. When driving in the dark, Ottobock strongly recommends that the user wears bright clothes or clothing with reflectors.

⚠ CAUTION

Risk of accident as a result of bad tyres. Check visually before every use that the tread depth is sufficient and that the tyres are inflated to the correct tyre pressure.

Whenever possible, avoid parking the wheelchair outdoors and on direct sunlight (UV light) as it causes the tyres to age prematurely. As a result, the tread surface hardens and corner pieces break out of the tread. Ottobock recommends to replace the tyres every 2 years regardless of wear and tear.

When the wheelchair is not used for extended periods of time, or if the tyres are heated strongly (e. g. near radiators or by sunlight shining through a window), permanent deformation of the tyres will result. Therefore always make sure that the wheelchair has sufficient distance from sources of heat, move your wheelchair from time to time, or jack up the wheelchair when storing it.

NOTICE

Risk of tyre damage. Too high an air pressure in your tyres may result in tyre defect. Please observe the data in Section "Technical data". The maximum air pressure must not be exceeded.

NOTICE

Damage due to exposure to extreme temperatures. The power wheelchair may be used within a temperature range of -15 °C to +40 °C. It must not be operated outside this temperature range.

NOTICE

Damage caused by overload. The maximum load for the power wheelchair is 100 kg / 220 lb (standard seat Junior: 75 kg). This load must not be exceeded.

NOTICE

Electric interference due to electromagnetic fields. The power wheelchair has been tested according to EMC regulations. Take care of the following particularities during the operation:

The driving characteristics of the power wheelchair can be affected by electromagnetic fields (mobile phones or other radiating devices). Please switch off all mobile devices when driving.

The power wheelchair can generate electromagnetic fields that can cause interference with other devices. Therefore, switch off the controls whenever you do not need them. Despite of compliance with all applicable EMC directives and standards, the power wheelchair can be affected by interference from other electric devices (e.g. department store EAS systems) or cause interference to such devices. If you notice such a behaviour, move your power wheelchair outside the interference range.

INFORMATION

When the wheelchair is used on the street, its user is required to obey the public traffic regulations.

INFORMATION

If the driving mode is changed while driving, the power wheelchair will accelerate or decelerate.

INFORMATION

After each emergency stop, the controls of the power wheelchair must be turned on again. In the event of communication problems in the bus system of the controls, the system triggers an emergency stop and thus prevents any uncontrolled functions. If the driving function is still not available after switching the controls on again, unlock the brakes to activate the push mode. In this case, contact a specialist dealer as soon as possible.

INFORMATION

If the controls do not emit an error signal as soon as the brake lever is activated, this indicates a malfunction. The settings must be checked by a specialist dealer.

INFORMATION

Manoeuvring is only allowed at reduced speeds.

2.6 Safety requirements for care, maintenance and disposal

CAUTION

Risk of accident and injury as a result of incorrect maintenance, repair or adjustment. Only trained staff who have been authorised by Ottobock may do maintenance work on the power wheelchair. This also applies to all repairs and settings on the brakes. Incorrect settings can lead to brake failure.

WARNING

Risk of injury due to explosive gases. Explosive gases can develop while the batteries are charging. For this reason, the following safety precautions are to be taken when charging the battery:

Switch the control off.

Ensure sufficient ventilation when charging the batteries in closed rooms.

Smoking and open flames are not permitted. Sparks must be avoided.

The air vents in the trim must not be covered.

CAUTION

Risk of injury due to uncontrolled movements. The fuse must always be removed for any maintenance work where the battery cover is open.

NOTICE

Unauthorised battery replacement. The battery may only be replaced by a specialist dealer. The characteristic curve of the battery charger set at the factory corresponds to the battery provided and must not be changed. Setting the characteristic curve incorrectly can result in permanent damage to the battery.

NOTICE

Damage of electronics as a result of penetrating moisture. Never use a water jet or high-pressure cleaning apparatus for cleaning the power wheelchair. Avoid direct contact of water with electronics, motors, and batteries under all circumstances.

INFORMATION

The power wheelchair must be inspected for functional reliability and driving safety at least once a year by an authorised specialist.

INFORMATION

Defective batteries must be disposed of properly in accordance with country-specific regulations.

2.7 Requirements for the user

⚠ CAUTION

Risk of accident and injury as a result of improper use.
 The power wheelchair may only be used by informed users. For this purpose, the user and maybe any attendant must receive instructions on how to use the power wheelchair from trained specialists who have been authorised by Ottobock. The operator must have read and understood all the information in the instructions for use.
 The wheelchair may not be operated in cases of exhaustion or under the influence of alcohol or medications.
 The operator must not have any mental limitations which can temporarily or permanently restrict attentiveness and judgment.

2.8 Safety functions



INFORMATION

In dangerous situations, the power wheelchair can be turned off at any time using the on/off button. When the button is pressed, the power wheelchair brakes immediately and the electric functions will be stopped.
 If any malfunctions occur such as an insufficient supply of energy to the brake, the software will recognise them and trigger the emergency stop brake or reduce the speed of the power wheelchair. At the same time, a warning signal will be emitted.

2.9 Warning symbols and type plates

Fig. 1 Signage on the A200

Label/Type plate	Explanation
 The type plate is located on the side of the frame below the seat.	A Type designation B European article number (EAN) C Maximum load capacity (see section “Technical data”) D Maximum climbing ability (see section “Technical data”) E Maximum speed (see section “Technical Data”) F Symbol for separate collection of electric and electronic devices. Components of the power wheelchair and the batteries must not be disposed of like regular domestic waste. G CE marking – product safety in conformity with EC Directives H Serial number I Read the instructions for use prior to using the product. Observe the safety instructions in the instructions for use. J Manufacturer / address K Product reference number
	Close the lock bar prior to use. Observe the safety instructions in the instructions for use manual.

Label/Type plate	Explanation
	Components of the power wheelchair and the batteries must not be disposed of like regular domestic waste.

Table 1 Signage on the A200

3 Product description

The power wheelchair can only be used inside. It is designed to be compact and manoeuvrable indoors.

The power wheelchair is controlled by a VR2 wheelchair control unit. The latter includes a control panel to enter driving commands and display the current status as well as a controller that operates the drive motors based on the inputs. Data transmission is realised via a bus system.

The VR2 control unit is programmable so that the controls can be adjusted to individual user requirements. The speed, acceleration and deceleration values, for example, can be adjusted.

The special features of the power wheelchair include:

- Compact design.
- Easy to disassemble for transport.
- Easy to service due to modular component group concept.
- Easy access to all component groups.

4 Delivery / Preparing the wheelchair for use

4.1 Delivery

INFORMATION

The options included in the delivery depend on the product configuration that has been selected for the power wheelchair.

Scope of delivery includes:

- Adapted power wheelchair with main components (see fig. 2)
- Charger
- Instructions for use
- Options (see section 7)

Upon delivery by the specialist dealer the power wheelchair is ready for use. All the settings correspond to the indications on the order form or are adjusted directly on site by the dealer. The power wheelchair is adapted to the personal requirements of the individual.

The functions of the individual components can be tested by following the instructions in section 6. Possible malfunctions are described in section 8.

The tool kit, an optional accessory, contains a set of Allen wrenches (3, 4, 5 and 6 mm) and one open-end 13 mm wrench.

Fig. 2 Main components

- | | |
|-----------------------------|---|
| 1 Backrest | 5 Footrest |
| 2 Armrest (side panel) | 6 Joystick |
| 3 Joystick and control unit | 7 Anti-tipper (fig. 3) |
| 4 Seat cushion | 8 Connector with battery cable (fig. 3) |

Fig. 3 Anti-tipper / Connector with battery cable

4.2 Initial operation

WARNING

Risk of suffocation. Packaging materials must be kept out of reach of children.

Prior to putting the power wheelchair into operation, the completeness (see fig. 2) and function of all its components must be checked. Insert and charge batteries if need be (see section 6.5 / 6.5.1).

Check the firm seat of the battery cable in connector (fig. 3, item 8) or insert the cable into the connector.

5 Transportation and storage

5.1 Transport in a wheelchair accessible motor vehicle

WARNING

Risk of accident and injury due to improper use for transportation in wheelchair accessible vehicles. The power wheelchair is permissible for transporting passengers in motor vehicles when using safety components (such as lap belt) as well as appropriate restraint systems offered by Ottobock.

Observe the limitations regarding installed options (see section 5.1.3).

The A200 may be used for transporting only one person. We recommend that wherever and whenever possible, A200 users transfer in the seats installed in the motor vehicle and use the corresponding vehicle restraint system, because this is the only way to ensure optimal protection of the passengers in case of an accident.

For more information, please refer to our brochure on using Ottobock products in wheelchair accessible vehicles, order number 646D158.

The power wheelchair may generally be used in wheelchair accessible vehicles without seat tilt.

WARNING

Risk of injury when using the belt system as a restraint system in a wheelchair accessible vehicle, which is forbidden.

- Under no circumstances may the belt system (e.g. lap belt) be used as part of a restraint system for transportation in a wheelchair accessible vehicle.
- Note that the belt system is only intended to help support the user sitting in the power wheelchair.

WARNING

Prohibited transportation of the passenger with activated back angle adjustment.

Loss of safe restraint in the product.

- Move the backrest to a nearly upright position prior to travel. Observe the related information in the product instructions for use.

Ensure that the wheelchair is fully secured with attachment straps for transport in a wheelchair accessible vehicle. To attach straps, 4 eyebolts are provided on the frame of the A200 (fig. 4a/4b). These are installed at the factory.

⚠ CAUTION

Risk of injury as a result of not engaging the locks sufficiently.

Secure the power wheelchair during transport in another vehicle by sufficiently tightening the tensioning straps.

Fig. 4a **Front eyebolts**

Fig. 4b **Rear eyebolts**

5.1.1 Required accessories

To use the power wheelchair as a seat in a wheelchair accessible vehicle, additional accessories have to be mounted (491S00=SK028 anchor point kit).

The qualified personnel which fit the wheelchair can provide more information.

5.1.2 Using the product in the vehicle

The power wheelchair has been tested according to ANSI/RESNA and ISO 7176-19.

Securing the product in the wheelchair accessible vehicle

1. Position the power wheelchair in the wheelchair accessible vehicle. For further information see section 5 in the brochure "Using your product for transportation in a wheelchair accessible vehicle", order number 646D158.
2. Turn the control unit off (see section 6.3.2).
3. Engage the brake (see section 6.4).
4. Engage the attachment straps at the front and rear before tightening them (fig. 4c).

Fig. 4c Engaged attachment straps at the front and rear

Placement of the restraint lap belt

1. Pull each end of the restraint lap belt from the inner side of the seat through to the outside.
2. Hook the end of the restraint lap belt onto the pin (fig. 4d).


Fig. 4d Restraint lap belt hooked onto pin

5.1.3 Restrictions for use

 WARNING

Risk of accidents and injury due to using the product with certain settings and/or installed options.

- Before using the product as a seat in a wheelchair accessible vehicle, remove options that need to be taken off for safe transportation. Please observe the following table.
- Stow all dismantled components securely in the wheelchair accessible vehicle.
- Please note that certain settings on the product exclude the use of the product in a wheelchair accessible vehicle.

Option	Transportation in a wheelchair accessible vehicle NOT possible	Option must be removed	Secure option on product
Height-adjustable control panel holder, swing-away			X
Headrest with mounting kit		X	
Lap belt with buckle			X*
Tray, swing-away to the side		X	

* The lap belt can be used to position the passenger during transportation. Using the personal restraint system is nevertheless required. For more information please refer to section 5 of the brochure on using Ottobock products in wheelchair accessible vehicles, order number 646D158.

5.2 Disassembling the power wheelchair

NOTICE

Risk of cable damage. Make sure that no cables are pinched during disassembly / reassembly.

You can reduce the folding size of the power wheelchair for transportation by easy dismantling. Carefully lay down the removed components.

Please observe the following order for disassembly:

1. Remove the control panel from the armrest (see section 6.1.4).
2. Remove the side panels (see section 6.1.3).
3. Remove the footrests (see section 6.1.5).
4. Pull the backrest release strap and fold the backrest down to the front (fig. 5).

INFORMATION

To proceed when using the option of “Mechanical back angle adjustment” see section 7.1.1.

5. Release the seat lock by simultaneously pulling the release strap and folding the locking bar to the rear (fig. 6). The battery packs are free.

Fig. 5

Backrest folded down

Fig. 6

Seat lock release

6. Remove the battery packs (fig. 23).
7. Separate the frame from the drive unit sustainer. Press down with your foot on the anti-tipper above the wheels until the wheels touch the ground (fig. 7).
8. Now you can remove the frame from the drive unit sustainer by slightly lifting the frame and pulling it backwards (fig. 8).

Fig. 7 Lowered anti-tipper

Fig. 8 Lifting the frame

9. Fold the seat bottom onto the frame: To do so you must loosen the two quick clamps on the front seat brackets. You can now fold down the seat surface by slight lifting the seat brackets and simultaneously pushing them backwards (fig. 9).

To reduce the space needed, the battery packs can be replaced into the drive unit sustainer. Take care that the red arrows on the battery packs point in driving direction. To prevent the battery packs from falling out of the drive unit sustainer, the locking lever must be relocated (fig. 10).

Fig. 9 Loosening the quick clamp

Fig. 10 Securing the battery packs

The control panel can either be laid onto the battery packs or put into the holder provided for this purpose. Make sure that the control panel remains switched off and that no cables are pinched.

The power wheelchair is now disassembled (fig. 11) and can be packed into a car.

NOTICE

Risk of damage as a result of not engaging the locks sufficiently. Sufficiently secure the disassembled power wheelchair components against slipping during transport.

Fig. 11 A200 disassembled

Perform the steps listed above in reverse order to reassemble the power wheelchair.

Ensure that the locking bar is pushed back into its proper position (see fig. 11a, item 1) and both locking pins engage correctly (see fig. 11a, item 2).

⚠ CAUTION

Risk of injury if the locking mechanism is not properly engaged.

Check that the locking pins are properly engaged at both sides on the outside of the base support. The bolts must be clearly visible so that the frame cannot unexpectedly disengage from the drive unit sustainer. (see fig. 11a, item 3).


Fig. 11a Locking bar and locking pins, engaged

Fig. 11b Locking pin, engaged

5.3 Storage

The power wheelchair must be stored in a dry place.

For transport and storage, an ambient temperature range from -15 °C to +40 °C (14 °F to +104 °F) must be observed.

INFORMATION

Tyres contain chemical substances that can react to other chemical substances (such as cleaning agents, acids, etc.).

INFORMATION

If your power wheelchair is not moved for several days, permanent colour changes might result on the ground where the wheelchair comes into contact with it. For this reason, place a suitable mat underneath the wheelchair, if it is not going to be used over an extended period of time.

INFORMATION

Black tyres can leave black traces on the ground under some circumstances. For this reason, we recommend choosing grey tyres when using the power wheelchair mainly indoors.

INFORMATION

Direct sunlight/UV light causes the tyres to age prematurely. As a result, the tread surface hardens and corner pieces break out of the tread. Avoid parking the wheelchair outdoors whenever possible.

INFORMATION

The tyres should be replaced every 2 years regardless of wear and tear.

INFORMATION

Remove the fuse for shipping or when the power wheelchair is not being used for an extended period of time.

6 Operation

6.1 Adjustment possibilities

⚠ CAUTION

Risk of injury due to unsecured screw connections. Risk of injury due to unsecured screw connections. After loosening screw connections with thread lock, replace these screw connections with new ones or secure them with medium-strength thread lock substance (e.g. EuroLock A24.20). After the power wheelchair has been set and adjusted, the attachment screws and/or nuts must be firmly retightened. During tightening observe torques when specified.

Various adjustments can be made to the power wheelchair. The seat height and seat position have been set in accordance with the indications on the customer order form and may only be changed by the specialist dealer.

The following items can be adapted by the user:

- Back angle
- Seat angle
- Armrest height
- Armrest position
- Lower leg length

If need be, the footrests and side panels can be removed.

6.1.1 Backrest

The backrest can be set to four different angles by pulling the release strap at the base of the backrest (fig. 12).

Once the backrest is in the desired position, release the strap to engage the locks.

Fig. 12 Release strap for the backrest

6.1.2 Seat angle

⚠ CAUTION

Risk of injury as a result of improper adjustment. During seat angle adjustments the user must not sit in the wheelchair. Ensure the quick release levers are snugly tightened after all seat angle adjustments!

The power wheelchair offers seat angle settings in three different positions (0°, 3° and 6°). For setting the position, two quick release levers are used which are located at the front below the seat.

To change the seat angle, loosen the two quick release levers (see fig. 9). Adjust the desired seat angle by simultaneously lifting the seat slightly and pushing the brackets to the rear. Ensure that the brackets completely engage in the new position.

6.1.3 Side panel with armrest

The side panels can be removed by loosening the thumb screws and lifting the side panels (fig. 13). The thumb screws are located at the lower end of the armrest holders (fig. 13, item A).

Both side panels can be removed, however, only the side without the control panel can be removed directly. Before removing the side panel on the side with control panel, the control panel must be removed (see section 6.1.4).

The height of the armrest can be adjusted by loosening the set screw at the upper end of the armrest holder (fig. 13, item B). Ensure the set screw is firmly retightened after making any adjustments.

Fig. 13 Removing the side panel

6.1.4 Control panel

NOTICE

Cable damage. Positioning the cable incorrectly can lead to pinching and thus damage to the cable. The cable must not be attached too tightly or too loosely. Avoid bending or squeezing the cable.

Remove the control panel

The control panel can be removed by pulling it to the front.

In case of a swing-away control panel (see section 7.3), pull off the control panel from the ball head of the control panel holder (fig. 14).

Fig. 14 Pulling off the swing-away control panel

Adapting the control panel to arm length

To adapt the control panel to the arm length of the user, you must loosen the three screws on the underside of the arm-rest (size 3) (see fig. 15).

You can now slide the rail with the control panel backwards and forwards. Firmly re-tighten the screws after making adjustments.

INFORMATION

If the control panel rail is too long, the protruding part can simply be cut off with a saw.

Fig. 15 Adapting the control panel

Changing the control panel side

INFORMATION

The change of the control panel side may only be carried out by authorised personnel.

The standard configuration has the control panel mounted on the right side. If desired by the user, it can be also mounted on the other side of the wheelchair.

6.1.5 Footrest

⚠ CAUTION

Risk of pinching. Make sure your fingers are not in the danger area when flipping the footplates up or down. Complete removal of the footrest is permitted only for transportation of the wheelchair.

Removal

1. Remove the calf band.
2. Flip the footplate up, release the footrest lock (see fig. 16) and swing the footrest either to the inside or outside.
3. Pull the footrest up to remove it.

Reassembly

1. Reinsert the footrest from above by placing it in the holding device and swinging it forward until the footrest lock engages.
2. Reattach the calf band to the holding device.


Fig. 16 Footrest lock

Adjusting the lower leg length (see fig. 17)

1. Loosen the screws on the footrest bar.
2. Move the footplate up or down to adapt the height to the individual lower leg length and seat cushion thickness.
3. Retighten the screws.

INFORMATION

The footrest bar must not be pulled out of the holding device by more than **160 mm**.


Fig. 17 Adjusting the lower leg length

6.2 Getting into and out of the power wheelchair

CAUTION

Risk of injury if the power wheelchair starts rolling. When getting into or out of the power wheelchair, the wheelchair controls must always be switched off. This will automatically engage the motor brake.

CAUTION

Risk of breakage due to overloading. The user must not apply his or her entire body weight to the footrests and armrests for support when getting into or out of the wheelchair.

The modular design of the power wheelchair and the ease with which you can remove the side panels and footrests make it easy to get into and out of the wheelchair from the side or from the front.

Getting into and out of the wheelchair can be done by the user individually in a way that suits him or her best.

6.2.1 From the side

Getting into the wheelchair from the side requires the user to remove the right or left side panel depending on the side where the user will get into the wheelchair. Bring the power wheelchair as close as possible to the place where its user is sitting.

If the control panel is on the side the user wants to use to get into/out of the wheelchair, then undo the hook and loop closures that fix the cable for the control panel to the wheelchair and pull off the control panel. Remove the side panel as described in section 6.1.3. If necessary, remove the footrest (see section 6.1.5).

The user can then slide onto the wheelchair. Use of a transfer board will make this even easier.

6.2.2 From the front

Flipping up the two footplates makes it possible for the user to get into and out of the wheelchair from the front. The space for getting into/out of the wheelchair can be increased by swinging the footrests to the side (fig. 18/19).

The assistance of an attendant or a transfer lifter make it easy for the user to get into or out of the power wheelchair. Use of a rotation plate is also possible.

Fig. 18 Footplates flipped up

Fig. 19 Footrests swung to the side

6.3 Control unit

NOTICE

Risk of impairing the driving performance of the power wheelchair. The driving characteristics of the power wheelchair can be affected by electromagnetic fields (cell phones or other radiation emitting devices). For this reason, all mobile devices are to be switched off when driving.

NOTICE

Damage of other devices. The power wheelchair can generate electromagnetic fields that can cause interference with other devices. Therefore, switch off the controls whenever you do not need them.

6.3.1 Control panel

The control panel is used to control the power wheelchair.

The control panel consists of a keypad, LED display and joystick. The charging/programming receptacle is on the underside. The control panel is used to switch the power wheelchair on and off, to enter driving commands and to display the current state of certain functions and components.


Fig. 20 Control panel

- | | |
|---------------------------------------|--------------------------------------|
| 1 Joystick | 5 Horn |
| 2 Speed button slow | 6 On/Off button |
| 3 Speed button fast | 7 Battery charge level (LED display) |
| 4 Speed level indicator (LED display) | 8 Charging receptacle underside |

On/Off key

The on/off key is used for switching the power wheelchair on, for activating the drive-away lock and for switching the power wheelchair off.

Speed buttons

Pressing the speed buttons shortly increases or reduces the speed level. After reaching the maximum speed level, the acoustic signal is changed.

Joystick

The joystick controls the direction and speed of travel.

Horn

The horn will sound as long as the horn key is being pressed.

“Charge level” LED display

The LEDs show the battery charge level (see section 6.3.4)

Speed Level indicator

The LEDs show the currently selected speed level.

6.3.2 Switching on and off

WARNING

Danger to life if the brake fails to function. The brake release lever must be locked when using the power wheelchair. The automatic brakes must be operational and functional.

CAUTION

Risk of accident as a result of a wrong tyre pressure. The power wheelchair must be visually checked every time before it is used to make sure that the tread depth is sufficient and that the tyres are inflated to the correct tyre pressure. Incorrect tyre pressure reduces the tyres' service life and makes the wheelchair's driving behaviour worse.

The control panel of the power wheelchair is turned on or off by pressing the on/off key (see fig. 20, item 6). If the control panel is not used for an extended period of time, the wheelchair turns off automatically. It is also possible to switch off the power wheelchair with the on/off key during driving. In this case, the wheelchair brakes immediately until it stops.

INFORMATION

Every time you switch on the control unit, it will return to the previously selected speed level.

6.3.3 Driving function

CAUTION

Risk of injury if the power wheelchair gets out of control. Switch off the controls on the power wheelchair whenever you do not need them. This will prevent you from unintentionally activating the joystick.

CAUTION

Risk of injury if the power wheelchair tips over during driving. Observe the following safety instructions when driving the power wheelchair:

- Do not drive on inclines and slopes of more than **7° (12 %)**.
- When driving downhill, reduce speed in accordance with the degree of the incline.
- Do not cross obstacles with differences in height of more than **15 mm (0.6 in)**.
- Do not cross steps without reducing the speed.

The power wheelchair has been approved for driving on inclines and slopes of no more than **7° (12 %)**. You are not allowed to drive on inclines or slopes steeper than this.

To ensure safe downhill driving, the driving speed must be reduced in accordance with the degree of incline (e.g. set it to speed level 1).

The critical obstacle height of the power wheelchair is **15 mm**. Obstacles exceeding **15 mm** in height may not be crossed. Obstacles such as steps must be crossed at a reduced speed.

If there are obstacles in the travel path, it is important to drive around them with a large margin of safety. On uneven ground, the driving behaviour of the wheelchair may get out of control. Therefore the speed must always be adapted to the ground conditions.

The joystick (see fig. 20, item 1) is used for driving. The further the joystick is moved away from the mid-position, the faster the power wheelchair will drive in this direction.

The maximum speed with full deflection of the joystick depends on the selected speed level.

Releasing the joystick automatically activates the brake function, which brings the wheelchair to a halt. When stand-

ing still, the mechanical brakes are automatically active so the power wheelchair cannot move.

The power wheelchair has five speed levels. The speed buttons (see fig. 20, items 2/3) are used to increase or decrease the speed level. The current speed level is indicated by lighting the “Speed Level” LED. The signal tone will change when the highest/lowest speed level was achieved.

Adapting the driving characteristics

WARNING

Risk of accidents and injury due to incorrect configuration settings. Modified parameter settings in the configuration can lead to changes in driving characteristics. In particular, changes to the speed, acceleration, braking or joystick settings can lead to unexpected and therefore uncontrollable operating performance with a risk of accidents. Always test the driving characteristics of the power wheelchair after configuration / programming is complete. Programming must only be completed by authorised personnel. Neither Ottobock nor the control unit manufacturer are liable for damages caused by programming that was not properly / professionally adapted to the abilities of the wheelchair user.

A hand programming device, which is to be connected to the control panel, can only be used by service personnel to adapt the speed, acceleration and deceleration values to the individual wishes of the user.

6.3.4 “Charge level” LED display

INFORMATION

Immediately after switching the power wheelchair on, the battery indicator shows the battery capacity saved before the wheelchair was switched off the last time. The exact battery status is displayed after driving a short distance.

INFORMATION

At temperatures of **<0 °C** the battery capacity drops by up to **35 %** in relation to the capacity for an outside temperature of **20 °C**. This shortens the range of the power wheelchair accordingly. Moreover at low temperatures the charge level displayed on the control panel can differ more significantly from the actual battery capacity.

The battery indicator in the LED display is divided into 10 segments and shows the current charge level.

All 10 LED segments are illuminated when the battery is fully charged. As the remaining battery charge decreases, the LED segments turn off one by one (see Table 2).

If the last three segments flash, the battery must be charged immediately.

The charging process is indicated by LED running light. When the battery is charging, the driving function is blocked.

Display	Information
	Battery is charged
	Charge the battery, if possible
	Low battery, recharge urgently
 Running light	Battery charging
 Flashing light	Battery undervoltage
 Flashing light	Battery overvoltage

Table 2 Battery indicator on the control panel

6.3.5 Range

INFORMATION

Information about the range of the power wheelchair can be found in the chapter 'Technical data'.

Note that the batteries reach their full capacity only after **20–30** charge cycles, and the power wheelchair can achieve the stated range only then.

Please note that the specified range was determined under defined conditions in accordance with ISO 7176-4. In practice the range can be reduced by up to **50 %**.

The following factors influence the range of a power wheelchair:

- Battery capacity
- Age of the batteries
- Ambient temperature
- Driving conditions (e.g. terrain profile, condition of surface)
- Charging method
- Use of power options
- Body weight of user

6.3.6 Drive-away lock

The control unit of the power wheelchair features an electronic drive-away lock. This function is activated/deactivated via control panel.

Activating the drive-away lock:

1. With the control unit switched on, press and hold the On/Off switch.
2. Once you hear a signal sound (approx. 1 sec.), release the On/Off switch.
3. Move the joystick to the front until a signal sound is heard.
4. Now move the joystick to the rear until a further signal sound is heard.
5. A long beep confirms that the driving function has been locked. The control unit then turns off.

The activation of the drive-away lock is indicated by a running light bar on the LED indicator "Speed levels" after switching the control unit on (see Table 3).

Display	Information
LED running light Speed level indicator	Drive-away lock

Table 3 Display of the drive-away lock on the control panel

Deactivating the drive-away lock:

1. After switching on, the battery capacity indicator is off and the LED indicator "Speed levels" is in running light mode.
2. Move the joystick to the front until a signal sound is heard.
3. Now move the joystick to the rear until a further signal sound is heard.
4. A long signal sounds to confirm that the driving function has been unlocked and the LED indicator "Battery charge level" is on.

Now the drive-away lock is deactivated and the wheelchair can be driven.

INFORMATION

If the joystick is not moved correctly, the drive-away lock remains activated. For a new deactivation of the drive-away lock, the you must switch off the control unit. Afterwards you can switch on the power wheelchair again and deactivate the drive-away lock.

6.4 Releasing and locking the brake

WARNING

Danger to life due to brake failure. Incorrect brake setting can lead to brake failure and therefore to serious bodily injuries or death. Repairs and adjustments to the brake may only be carried out by authorised service staff.

WARNING

Risk of accident and injury if brake function is unavailable.

No brake functionality is available when the brake is unlocked (push mode).

- When moving the power wheelchair on an incline, the person pushing must provide the appropriate brake force.
- The brake function may only be released in the presence of an attendant.
- Should the user be unable to release the brake himself, the brake can be released by the attendant.

NOTICE

Damage due to parking without locked brake. Releasing the brake may result in uncontrolled rolling of the power wheelchair. Therefore, make sure that the brake is engaged after parking the power wheelchair.

It is possible to push the power wheelchair in case of control unit failure or insufficient battery charge level. To do so, the brake is released by mechanically unlocking it. The brake release mechanism is located on the left-hand and right-hand sides on the drive motors.

Unlocking the brake

For activation of the push mode, pull out both red brake release levers located on the drive motors on the left-hand and right-hand sides. In this position, the control unit will recognise that the brake has been released and automatically deactivate the driving function. The unlocked brake will be indicated by the flashing "Battery charge level" LED.

INFORMATION

Once the brake unlocking lever has been released, all braking systems are deactivated.

Display	Information
 Flashing light	Unlocked brake for push mode

Table 4 Display of the brake release on the control panel


Fig. 21 Brake release

Locking the brake

To lock the brake, push the red brake release lever inwards. To reactivate the driving function, you must switch off the wheelchair controls and then switch them on again.

6.5 Batteries

NOTICE

Risk of damage to the battery packs. The batteries are factory-encapsulated in battery packs. Battery packs should not be opened because this could cause irreversible damage to the encapsulated battery or battery pack.

INFORMATION

Regularly checking of the remaining battery capacity and charging of the battery early enough are indispensable preconditions for the operational safety of your power wheelchair.

INFORMATION

Before doing any maintenance work on the batteries, please read the battery manufacturer's warnings thoroughly!

Please observe the disposal instructions (see section 11).

The power wheelchair comes factory-equipped with two maintenance-free AGM batteries with a total capacity of **28 Ah (C5)**. The batteries are located in two battery packs under the seat of the power wheelchair (fig. 22).

See section 6.3.4 for battery charge level display on the control panel.

Fig. 22 Battery pack

The applied AGM technology fulfils the latest provisions regarding safety and environmental tolerance through close construction.

For increased safety and to facilitate handling, the batteries are integrated in separate battery packs. Thanks to special lead plates, these batteries are particular well-suited for electric drive systems and have been designed with a high cycle solidity.

For exchange or disposal of the batteries (batteries within battery pack), please contact your authorised dealer. The batteries are to be replaced as complete units (battery packs).

To remove / replace the battery packs proceed as follows:

1. Release the seat lock by simultaneously pulling the release strap and folding the locking bar to the rear (fig. 6).
2. The battery packs are free. Remove the battery packs (fig. 23).
3. Place the new battery packs into the drive unit sustainer. Take care that the red arrows on the battery packs point in driving direction and that the plug connections (fig. 24) engage into the battery contacts.
4. To prevent the battery packs from falling out of the drive unit sustainer, the locking lever must be relocated (fig. 10).

Fig. 23 Removing battery pack

Fig. 24 Bottom side of battery pack

6.5.1 Charging

WARNING

Risk of injury due to explosive gases. Explosive gases can develop while the batteries are charging. The following safety instructions must be followed under all circumstances:

- Ensure sufficient ventilation when charging the batteries in closed rooms.
- Smoking and fire are not permitted.
- Sparks must be avoided.
- Do not cover the air vents in the trim.

NOTICE

Risk of battery damage. Driving for longer periods with low battery charge level will result in deep discharge and damage to the battery. There is a risk that the power wheelchair may stop due to zero battery charge level and bring the user into a dangerous situation.

The remaining battery capacity determines the distance range of the power wheelchair. The battery capacity is influenced by many factors. In addition to the temperature, battery age, and amount of use, the charging cycle has a pronounced effect on the capacity and therefore on the range.

The batteries reach their full capacity only after **20–30** charge cycles. Only if the full capacity of the batteries has been reached can the power wheelchair achieve the stated range.

At temperatures of **<0 °C** the battery capacity drops by up to **35 %** in relation to the capacity for an outside temperature of **20 °C**. This shortens the range of the power wheelchair accordingly. For an optimal charging frequency the following applies:

- The batteries can be charged at any time, regardless of the remaining charge.
- It takes about 12 hours until a discharged battery (only one, flashing segment) is completely charged. When the charging process is complete, the battery charger can remain connected with no risk of overcharging or damaging the battery. The battery charger features a programmed recharging phase that will maintain the battery capacity at the previously reached level.
- For daily use, it is recommended to connect the charger overnight so that the full capacity is available every day.
- If the wheelchair is not used for extended periods of time, the batteries gradually discharge. If the power wheelchair is not used for a longer period of time, a

charging cycle must be carried out at least once a week to maintain the battery capacity. When the power wheelchair is not used for an extended period of time, we recommend removing the fuse.

- The batteries should never be fully discharged (deep discharge).
- When the batteries are charging, the controls on the power wheelchair must be switched off to allow the charging current to be fed completely into the battery.

The following must be observed when charging the batteries:

- Only the battery charger provided by Ottobock may be used for charging. Failure to comply will void the warranty.
- The voltage settings on the battery charger must correspond to the voltage used in your country.

6.5.2 Charger

WARNING

Explosion hazard as a result of sparks! Always switch off the battery charger and remove the mains plug before disconnecting the battery.

NOTICE

Unauthorised battery replacement. The battery may only be replaced by a specialist dealer. The characteristic curve of the battery charger set at the factory corresponds to the battery provided and must not be changed. Setting the characteristic curve incorrectly can result in permanent damage to the battery.

The battery charger is designed for gel batteries (maintenance-free batteries). The charger is powered according to the appropriate characteristic curve allowing for optimal charging of the battery type. If the battery charger is to be used with another power wheelchair or if new batteries are installed, the characteristic curve has to be checked.

NOTICE

Risk of battery damage. Setting the characteristic curve incorrectly can result in permanent damage to the battery.

When handling the battery charger, the following safety instructions must be observed:

- D Always place the rubber feet of the battery charger on a level surface.

- Protect the battery charger from direct sunlight to prevent the charger from getting any warmer than it already is.
- The location where the battery charger is installed must be dry and well ventilated. Avoid letting dust and dirt enter the battery charger. It can impair the functionality of the battery charger.
- To clean the charger, use a dry piece of cloth.

To charge the batteries, proceed as follows:

1. Turn off the control unit on the power wheelchair.
2. Connect the plug of the battery charger to the charging receptacle on the control panel of the power wheelchair (see fig. 20).
3. Connect the battery charger to a wall outlet. The batteries will start charging automatically, and the current battery capacity is indicated on the LEDs on the battery charger.
4. When the charging process is complete: Disconnect the battery charger from the mains supply.
5. Turn on the control unit; the power wheelchair is now ready for operation.

Please see the instructions for use supplied with the battery charger for further details on use and on the LED indicators.

7 Accessories

WARNING

Risk of pinching due to exposed pinch points in the seat adjustment area

- Note that when manual seat options are used, inherent pinch and shear points are located between the seat frame and the power wheelchair frame.
- Ensure that no body parts, such as hands or feet, are in the danger area while manual seat options are used.
- Ensure that no interfering objects, such as clothing or other obstacles, are in the danger area while manual seat options are used.

CAUTION

Risk of injury due to unsecured screw connections. After loosening screw connections with thread lock, replace these screw connections with new ones or secure them with medium-strength thread lock substance (e.g. EuroLock A24.20).

After the power wheelchair has been set and adjusted, the attachment screws and/or nuts must be firmly retightened. During tightening observe torques when specified.

INFORMATION

Only original options provided by the manufacturer may be used. The optional components may be mounted only as described here. Failure to comply will void the warranty.

INFORMATION

All available optional components are contained on the order form and in the wheelchair accessories catalogue.

INFORMATION

In case of power wheelchairs which feature gas compression springs, attention must be paid to oil leaking from the piston rod. Leakages reduce the functional capabilities of the gas compression springs and result in their failure. Defective gas compression springs (e. g. due to scratches or damage to the piston rod caused by impact; deformation of the cylinder tube) must be replaced immediately. Piston rods do not require lubrication. They are maintenance-free.

The power wheelchair has been designed as a modular system. Certain component groups can be exchanged and other accessory components can be added to the power wheelchair. The complete range of accessory components

is listed on the order form and in the wheelchair accessories catalogue. Some important options are described in the following.

7.1 Mechanical back angle adjustment

The backrest of the power wheelchair with standard or contoured seat can be equipped with a mechanical back angle adjustment.

In the sitting position, the backrest can be tilted back continuously by up to **30°**.

To adjust the backrest proceed as follows (see fig. 25):

1. Use the release lever on the armrest (to unblock the gas compression spring).
2. Move the backrest to the desired position.
3. Let go of the release lever (to again block the gas compression spring).

Fig. 25 Mechanical back angle adjustment

7.1.1 Specialties for preparing power wheelchair for transport

⚠ CAUTION

Risk of accident and injury due to incorrectly installed backrest. Make sure that the locking on cross bolt has engaged (fig. 26). Check if the cross bolt firmly locks in place.

To fold the backrest down proceed as follows:

1. Release the cross bolt at the lower end of the gas compression spring (by pushing up the lever, see fig. 26).
2. Remove the bolt out of the holding device.

3. Remove the side panels and fold down the back part onto the seat bottom.

Fig. 26 Cross bolt locked in the holding device

- | | |
|--------------|-------------------|
| 1 Cross bolt | 2 Lock with lever |
|--------------|-------------------|

To put a wheelchair with mechanical back angle adjustment into operation:

1. Unfold the backrest.
2. Insert the cross bolt at the end of the gas compression spring into the holding device (fig. 26, item 2).
3. Lock the cross bolt (push down the lever, see fig. 26, arrow). The locking will engage on the cross bolt.

7.2 Mechanically elevating footrest

The power wheelchair can be equipped with mechanically elevating footrests (see fig. 27).

To elevate the footrest:

1. Use the release lever on the footrest (to unblock the gas compression spring, see arrow in fig. 27).
2. Move the footrest to the desired position.
3. Let go of the release lever (to again block the gas compression spring).

Fig. 27 Mechanically elevating footrest

7.3 Swing-away control panel holder

To allow the user to drive the power wheelchair closer to an object or under the edge of a table, the control panel can be swung to the side with a special control panel holder (see fig. 28).

1. Apply a little bit of pressure to press the control panel holder to the side and release the pivot element.
2. Swing the control panel holder to the side.
3. When bringing the control panel holder back to its original position, the pivot element automatically locks into place.

Fig. 28 Swing-away control panel holder

7.4 Lap belt

⚠ WARNING

Risk of injury due to improper use of lap belt. The lap belt helps to additionally stabilise the person sitting in the power wheelchair. The lap belt must never be used as part of a restraint system for transportation of the person in wheelchair accessible motor vehicles.

⚠ CAUTION

Risk of injury due to insufficient stabilisation of the seated person. Always wear the safety belt / lap belt when driving in public.

⚠ CAUTION

Pressure points/constriction due to putting the lap belt on incorrectly.

- Ensure that the buckle lies in the middle of the body.
- Ensure that the lap belt is not too tight against the body.
- Remove objects or clothing which get caught.

A lap belt can be added to the power wheelchair. The lap belt provides additional support and keeps the user from sliding out of the seat.

To apply the lap belt, insert the two buckle halves into each other until they lock in place (fig. 29). The belt buckle must snap in place audibly. Then, verify that the belt has locked by trying to pulling it apart.

To open the lap belt, press the red release button.

The belt length can be adjusted on both sides. To set the length, position both buckle halves in the middle of the body. After positioning the buckle halves at a right angle in relation to the belt, their position can be varied. Any excessive belt length will be held in place by the plastic slides.

Fig. 29 Applying the lap belt

7.5 Other options

- Curb climbing assist: The curb climbing assist is used to clear curbs and steps that have a maximum height of **80 mm**.
- Puncture-proof tyres: Solid rubber tyres
- Contoured seats, small/large: for better sitting comfort and more adaptation possibilities
- Standard seat, junior
- Armrest accessories: Special adapter for the armrests of our Accessories Catalogue
- Adapter for headrest mounting kit:
For attachment to the back tube
- Additional lights

These and other optional add-on components are included in the order form and in the accessories catalogue.

8 Malfunctions / Troubleshooting

INFORMATION

Should you encounter problems while troubleshooting or if you do not manage to completely eliminate a problem by following the measures described here, please contact your specialist dealer.

Malfunctions are indicated by the flashing "Battery charge level" LED on the control panel. Table 5 describes the individual indications with the corresponding sources of the problem as well as the possible causes and measures to be taken.

If you do not manage to completely eliminate the problems with the described measures, your specialist dealer may be able to read the exact error code with the hand programming device and undertake a more systematic analysis of the system.

All problems that have ever occurred are saved in a list and can be retrieved, for example, in case of a general overhaul of the power wheelchair. The saved data can be used to determine future service and maintenance intervals, for example.

8.1 Warning

A warning indicates a status or malfunction of one or several components of the power wheelchair. The function of components that have no errors is not restricted.

8.2 Errors

CAUTION

Risk of injury if the power wheelchair stops abruptly. In the event of communication problems in the bus system of the controls, the system triggers an emergency stop to prevent any uncontrolled functions. Depending on the respective malfunction, it is possible to drive the wheelchair out of a danger zone by turning the control unit back on.

If the driving function is still not available after switching on the controls again, the brakes must be released (see section 6.4) to switch over to the pushing mode. After the power wheelchair has triggered an emergency stop, contact a specialist dealer as soon as possible!

CAUTION

Risk of accident due to uncontrolled driving behaviour.

Uncontrolled movements can occur during the operation of the power wheelchair as a result of malfunction. In this case, please contact your authorised dealer immediately.

NOTICE

Malfunction on the control unit. When the joystick is activated while the brakes are unlocked, the control system emits an error signal on the control panel. If this is not the case, there is a malfunction that must be corrected immediately by a specialist dealer.

An error affects one or several functions of the power wheelchair. The system is not fully operational until the error has been corrected.

Flashing LED(s)	Error/Warning	Cause	Possible corrective action
	Battery undervoltage	Battery deep discharge Battery cable malfunctioning / faulty connection to the battery	Charge the battery; check the connection to the battery (in case of proper con- nection, charge the battery)
	Incorrect/defective cabling of the left motor Defective motor	E.g. faulty plug connection	Check the connections to the left motor Check the motor
	Short circuit on the battery connection to the left motor	E.g. broken cable	Check the battery connection to the left motor
	Incorrect/defective cabling of the right motor Defective motor	E.g. faulty plug connection	Check the connections to the right motor Check the motor
	Short circuit on the battery connection to the right motor	E.g. broken cable	Check the battery connection to the right motor
	Driving function is blocked due to external influences	Maybe the battery charger is connected	Disconnect the battery charger
	Joystick error	Joystick not in zero position when switching on	Bring joystick to zero position prior to switching on


Flashing LED(s)	Error/Warning	Cause	Possible corrective action
	Controller error	Defective controller	Check all connections
	Brake release	Brake release is open	Check the motor brakes/ Check the connections to the controller
	Battery overvoltage	Voltage too high Battery contacts loose	Continue driving slowly Check the plug contacts
	Communication error between control panel (joystick) and controller	Poor cable, loose plug connection	Check the connections

Table 5 Status and error messages

9 Maintenance, cleaning and care

INFORMATION

A spare parts catalogue from Ottobock is available upon request for ordering spare parts. Only spare parts supplied by Ottobock may be used. Failure to comply with this requirement will render the warranty null and void.

INFORMATION

Should you encounter problems during maintenance, contact your authorised dealer. The power wheelchair is to be checked and serviced once a year by an authorised dealer.

9.1 Maintenance intervals

The correct function of the power wheelchair should be checked every time before using it. The items listed in table 7 must be checked by the user at the indicated intervals.

Component	Activity	Prior to every use	Monthly	Weekly
Armrest and side panel	Fastening screws tightened Armrest and control panel secured Check armrest for damage	X	X	X
Drive wheels	Wheels must rotate freely and without axial run out Central nut on the driving shaft tightened Check if wheel mounts are seated securely Directional stability of the entire wheelchair		X	X X X
Tyres	Air pressure (printed on the sidewall of the tyre) Sufficient tread depth, at least 1 mm Check for damage			X X X
Batteries	Check battery charge level	X		
Electronics	Control system free of errors Battery charger does not show any error messages on the LEDs Check plug connections	X	X	X
Brake	Activate brake lever while control system is switched on Brake function active with engaged brake	X		X

Component	Activity	Prior to every use	Monthly	Weekly
Mechanically elevating footrest	Check ratchet mechanism for functionality and secure seating			X
	Check footplates for damage			X
	Visual check for scratches on the piston rod and oil leakage			X
Swivelling wheels / casters	Fork seated in the receiver without play			X
	Wheels must rotate freely and without axial run out			X
	Fastening nut tightened			X
Padding and belts	Proper condition of padding			X
	No wear on the seat belts			X
	Check belt buckle for functionality		X	
Seat attachment	Check if attachment screws are seated securely			X
Gas compression spring (option) or actuator	Visual check for scratches on the piston rod and oil leakage			X

Table 6 Maintenance measures and intervals

9.2 Changing the fuse

The fuses (for values see section “Technical Data” > table “Electric equipment”) are located respectively in one of the fuse holders on the bottom of the battery packs (see fig. 30).

1. Remove the batteries (see section 6.5).
2. Simply pull out the fuse on the underside.
3. Insert the new fuse between the contacts into the fuse holder on the battery underside. Make sure that the fuse is pressed into the centre of the spring contacts and that it is not at an angle.
4. Reinstall the batteries (see section 6.5).

Fig. 30 Battery pack with fuse

9.3 Changing tyres/inner tube

INFORMATION

Direct sunlight (UV light) causes the tyres to age prematurely. As a result, the tread surface hardens and corner pieces break out of the tread.

INFORMATION

Avoid parking the wheelchair outdoors whenever possible. Regardless of wear and tear, the tyres should be replaced every 2 years.

If the wheelchair is parked for an extended period of time or the tyres overheat (e.g. in the vicinity of radiators or in case of exposure to strong sunlight behind glass), the tyres may become permanently deformed. Therefore, sufficient clearance should always be maintained from sources of heat, and the wheelchair should be moved frequently or placed on blocks for storage.

To change the tyre/inner tube of a drive wheel, proceed as follows:

1. Secure the wheelchair to prevent it from tilting to the side by placing a suitable base under the drive unit sustainer.

2. To dismount a drive wheel, loosen the four Allen head screws using the Allen wrench size 8 in the middle of the wheel (fig. 31) and remove the wheel.
3. If the tyre/inner tube of a drive wheel needs to be replaced, loosen screws on the inner side of the rim using the Allen head wrench size 8 and pull apart the two-piece rim.
4. The defective inner tube is now freely accessible and can be replaced.

Fig. 31 **Disassembling drive wheel**

To change the front wheel, proceed as follows:

1. To dismount the front wheel, loosen the screw of the axle using a 6 mm Allen wrench (see fig. 32, item 1) and pull out the axle.
2. The defective wheel is now freely accessible and can be replaced.

Fig. 32 **Disassembling front wheel**

- 1 Axle screw

9.4 Cleaning and care

NOTICE

Damage of electronics as a result of penetrating moisture. When cleaning the power wheelchair, the electronic components, the motor and the batteries must not come into contact with water to avoid malfunctions.

NOTICE

Damage to power wheelchair components. Do not use any aggressive cleansing agents, solvents, or hard brushes for cleaning the power wheelchair. Never use a water jet or high-pressure cleaning apparatus for cleaning the power wheelchair.

INFORMATION

Prior to disinfection, clean the seat and back upholstery as well as the seat cushion, the control panel and the armrest.

The power wheelchair must be cleaned regularly, depending on the amount of use and the degree of soiling.

The control panel, battery charger, armrests, and trim can be cleaned with a damp cloth and mild cleaning solution.

Use a dry brush to clean the seat and back upholstery as well as the seat cushion.

Use a damp plastic brush to clean the wheels and frame.

9.4.1 Disinfection

Wipe all parts of the wheelchair down with disinfectant.

Important information about disinfection

- Water based disinfectants should be used. Observe the usage information provided by the manufacturer.
- Prior to disinfection, clean the seat and back upholstery, seat cushion, control panel, and armrest.

10 Disposal

10.1 Disposal information

NOTICE

Environmental pollution due to incorrect disposal of the battery

- Observe the information printed on the batteries by the manufacturer.
- Note that the batteries may not be disposed of as household waste.

INFORMATION

If a wheelchair is to be disposed of, all components and materials of the power wheelchair must be recycled or disposed of properly.

If the power wheelchair is no longer in use, it must be disposed of properly in accordance with national regulations.

Please return defective batteries (in battery pack) to your dealer when buying new ones

10.2 Information on re-use

The power wheelchair is suitable for re-use.

Similar to second-hand machines or cars, products that are being re-used are subject to increased strain. Features and functions must not change in a way that could endanger patients or other persons within the product's life cycle.

Based on market observations and the current state of technology, the manufacturer has calculated that the power wheelchair can be used for a period of **5 years**, provided that it is used properly and that the service and maintenance instructions are observed. Periods during which the wheelchair is stored at the dealer or with the cost bearer are not included in this period. It should be clearly pointed out, however, that the wheelchair is a reliable product far beyond this defined period of time, provided that it is cared for and maintained appropriately.

In cases of re-use, the corresponding product must first be thoroughly cleaned and disinfected.

Afterwards, the condition of the product must be examined by an authorised dealer to check for wear and tear as well as any damage.

Any worn and damaged components as well as components which do not fit or are not suitable for the new user must be replaced.

The service manual includes a service schedule for each model, detailed information, and a list of the required tools.

11 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

11.1 Service life

Based on market observations and the current state of technology, the manufacturer has calculated that the product can be used for a period of **5 years**, provided that it is used properly and that the service and maintenance instructions are followed. Storage times at the dealer or with paying parties are not included in this period.

However, we would like to emphasise that the product is reliable far beyond this defined period of time, provided that it is cared for and maintained properly.

If the service life is reached, the user or a responsible attendant should contact the qualified personnel who fitted the product or the manufacturer's service (see inside or outside of rear cover for address).

11.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

11.3 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

The product meets the requirements under the RoHS Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

11.4 Warranty

Further information on the warranty terms and conditions can be obtained from the qualified personnel that has fitted this product or the manufacturer's service (see inside back cover for addresses).

11.5 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

12 Technical data

General information	
Application class (according to DIN EN 12184)	Class A
Dimensions and weights	
Seat width	Standard seat: 340 – 480 mm Contour seat: 380 – 480 mm
Seat depth	Standard seat: 340 – 500 mm Contour seat: 360 – 480 mm
Seat height	Standard seat: 430/480 mm Contour seat: 510/560 mm
Armrest height	225 – 350 mm Standard seat, junior: 205 – 275 mm
Armrest length	260 mm
Lower leg length	Standard seat: 250 – 470 mm Contoured seat: 280 – 540 mm
Back height	550/500/450 mm
Back angle	Manually, in 10° increments: -9°/1°/11°/21° or 0°/10°/20°/30°
Back angle adjustment	Mechanical: adjustable by up to 30°

Dimensions and weights	
Overall width	570 mm
Total height	1030 mm
Total length	1000 mm
Turning area	1290 mm
Turning radius	870 mm
Tyre size	8" 12.1/2x2.1/4"
Caster	
Drive wheel	
Air pressure	printed on the sidewall of the tyre
Weight empty**	66 kg (dependent on options)
Shipping weight**	see weight when empty**, of which: Side panel: < 1 kg (2.2 lbs) Footrest (standard equipment): approx. 1 kg (2.2 lbs) Manual elevating footrest: 1.8 kg (4 lbs) Battery, removable: approx. 11 kg (24 lbs) (per unit)
Max. load capacity	100 kg (patient weight); standard seat Junior: 75 kg

Electric equipment***	
IP protection rating (according to DIN EN 60529)	IPX4
Operating voltage	24 V
Batteries (AGM batteries)	2 x 12 V, 28 Ah (C5) / 33 Ah (C20)
Control: Model Operating voltage Max. output current per motor	VR2 with controller and control panel 24 V DC 60 A
Fuse	60 A in each battery pack
Battery charger	For more information, see the included battery charger instructions for use

Driving data	
Speed	6 km/h [3,7 mph]
Maximum safe inclination	Climbing ability****: 7° (12 %)
Obstacle height that can be cleared	15 mm
Range (on level surfaces)*****	approx. 15 km
Braking distance (according to DIN EN 12184:2009)*****	for 6 km/h [3,7 mph]: 1000 mm (horizontal); 2000 mm (inclination)
Operating temperature range	-15 °C to +40 °C

Driving data	
Transport and storage temperature range	-15 °C to +40 °C

Corrosion protection	
Corrosion protection	Powder-coated frame

* 3-point turn by 180°

** The specified weights vary according to the selected options and model.

*** The product meets all requirements under ISO 7176-14.

**** Depending on the individual, the continuous climbing ability can be significantly lower than the value specified.

***** The specified range was determined under defined conditions according to ISO 7176-4. In practice the range can be reduced by up to 50 %. For information on this, see the section "Range" in the instructions for use (user).

***** The braking distance can be correspondingly longer due to user weight, installed options and condition of the tyres, and due to weather and surface conditions.

Technical changes reserved.

Table 7 Technical data

This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and extend across the width of the page. There are no margins, text, or other markings on the paper.

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Otto Bock HealthCare Canada
Burlington, Ontario, L7L 5N5
Canada
T +1 289 288-4848
F +1 289 288-4837
infocanada@ottobock.com

Otto Bock HealthCare
Andina Ltda.
Bogotá · Colombia
T +57 1 8619988
F +57 1 8619977
info@ottobock.com.co

Otto Bock de Mexico S.A. de C.V.
C.P. 01180 México, D.F. · Mexico
T +52 55 5575 0290
F +52 55 5575 0234
info@ottobock.com.mx

Otto Bock HealthCare
Minneapolis, MN 55447 · USA
T +1 763 553 9464
F +1 763 519 6153
usa.customerservice
@ottobockus.com

Asia/Pacific

Otto Bock Australia Pty. Ltd.
Baulkham Hills NSW 2153
Australia
T +61 2 8818 2800
F +61 2 8814 4500
healthcare@ottobock.com.au

Beijing Otto Bock
Orthopaedic Industries Co., Ltd.
Beijing, 100015, P.R. China
T +8610 8598 6880
F +8610 8598 0040
news-service@ottobock.com.cn

Otto Bock Asia Pacific Ltd.
Kowloon, Hong Kong · China
T +852 2598 9772
F +852 2598 7886
info@ottobock.com.hk

Otto Bock HealthCare
India Pvt. Ltd.
Mumbai 400 021 · India
T +91 22 2274 5500 / 5501 / 5502
information@indiaottobock.com

Otto Bock Japan K. K.
Tokyo, 108-0023 · Japan
T +81 3 3798-2111
F +81 3 3798-2112
ottobock@ottobock.co.jp

Otto Bock Korea HealthCare Inc.
137-070 Seoul · Korea
T +82 2 577-3831
F +82 2 577-3828
info@ottobockkorea.com

Otto Bock
South East Asia Co., Ltd.
Bangkok 10900 · Thailand
T +66 2 930 3030
F +66 2 930 3311
obsea@ottobock.co.th

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Ottobock has a certified Quality Management System in accordance with ISO 13485.

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Otto Bock HealthCare GmbH