BRAVO Racer

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EN**Instructions for use – BRAVO Racer rigid-frame active wheelchair**

Table of contents	Page
1 General information	41
1.1 Preface	41
1.2 Indications for use	41
1.3 Field of application.....	41
1.4 Legal information	42
1.5 Technical service.....	42
2 Safety instructions	43
2.1 Explanation of symbols	43
2.2 General Safety instructions	43
2.3 Further notes.....	47
2.4 Warning symbols and nameplates.....	48
3 Delivery/preparing the wheelchair for use	48
3.1 Storage.....	48
3.2 Preparation for use	48
4 Transportation	50
4.1 Transfer	50
4.2 Use in a wheelchair accessible vehicle.....	52
5 Options.....	53
5.1 Anti-tipper.....	53
5.2 Crutch holder	54
5.3 Transport wheels.....	54
5.4 Tip-assist.....	54
5.5 Backpack	55
5.6 Footrest, independently angle adjustable	55
5.7 Lateral heel block.....	56
5.8 Seat upholstery, adaptable	56
5.9 Back upholstery, breathable.....	56
5.10 Fold-down push handles	56
5.11 Push bar	57
5.12 Designer clothing protectors.....	57
5.13 Side panel with height adjustable armrest	57
5.14 Tray	57
5.15 Designer casters	58
5.16 Luminous wheels.....	58

5.17 Wheelbase extension with second axle	58
5.18 Infinity Ultralight rim	58
5.19 Rear wheel with drum brake	59
5.20 Lap belt	59
6 Adjustment/assembly instructions	60
6.1 Adjusting the height of the footrest	60
6.2 Longitudinal adjustment of the footrest	61
6.3 Adjusting the seat height.....	61
6.4 Adjusting the seat angle	62
6.5 Fine adjustment of seat/footplate angle	63
6.6 Longitudinal adjustment of the seat unit.....	64
6.7 Adjusting the back unit.....	64
6.8 Changing the back angle.....	65
6.9 Readjusting seat upholstery tension	65
6.10 Adjusting the wheel lock	66
6.11 Installing the frame accessories	66
7 Maintenance, cleaning, and disinfection	67
7.1 Maintenance	67
7.2 Cleaning and disinfection.....	68
8 Changing tyres	68
9 Re-use / disposal	70
9.1 Information on re-use.....	70
9.2 Information on disposal	70
10 Technical data	71

1 General information

INFORMATION

Last update: 2015-04-08

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

INFORMATION

- You can request this document as a PDF file from the Customer Care Center (CCC) at oa@ottobock.com or from the manufacturer's service department (see inside back cover or back page for addresses).
- It is possible to increase the display size of the PDF document.
- For further questions about the instructions for use, please contact the authorised personnel who issued the product to you.

1.1 Preface

With the purchase of the BRAVO Racer active wheelchair, you have selected a quality product for versatile use during everyday life, both indoors and outdoors.

The chapter "Options" presents various additional components of the BRAVO Racer, which can expand its field of application and improve comfort for the wheelchair user. The chapter "Adjustment/assembly instructions" provides an overview of the ways the wheelchair can be customized to your child's requirements. The model as described in these Instructions for Use is subject to technical changes without notice.

1.2 Indications for use

The BRAVO Racer active wheelchair is designed solely for individual indoor and outdoor use by persons who are unable to walk or who have a walking impediment, and can be operated by the patient or by another person.

The BRAVO Racer active wheelchair may only be combined with the options specified in these instructions. Ottobock assumes no liability for combinations with medical products and/or accessories from other manufacturers outside of the modular system.

1.3 Field of application

The versatility of this power wheelchair and its modular design make it suitable for patients who have walking impediments or walking disabilities, for example due to:

- Paralysis (paraplegia / tetraplegia or tetraparesis)
- Loss of limbs (dysmelia/lower limb amputation)
- Infantile/spastic cerebral palsy
- Spina bifida
- Muscle and nerve disorders

- Imperfect osteogenesis
- Poliomyelitis

The BRAVO Racer was specially designed for individuals who are able to independently use a wheelchair or those of certain age groups who require an attendant.

Fitting considerations:

- Body height and weight (maximum load is 60 kg)
- Physical and mental condition as well as degree of development
- Age of user
- Home and living conditions
- Environment

1.4 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

1.4.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

1.4.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

1.4.3 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

1.5 Technical service

Service and repairs on the Ottobock product may only be carried out by an authorised dealer. Should any problems arise, please contact the supplier of your product. Any necessary repairs will be made exclusively with authentic Ottobock spare parts there. Your product requires regular maintenance / regular service (see 647G502=GB Service Instructions).

Your authorised Ottobock dealer:

2 Safety instructions

2.1 Explanation of symbols

 WARNING	Warnings regarding possible risks of severe accident or injury.
 CAUTION	Warnings regarding possible risks of accident or injury.
 NOTICE	Warnings regarding possible technical damage.
 INFORMATION	Information regarding operation. Information for service personnel.

2.2 General Safety instructions

INFORMATION

Please read the Instructions for Use first! Before using the product, you should become familiar with the handling, function and use of the product. You may compromise your safety if you do not observe these instructions.

However, it is impossible to protect against and avoid all unforeseeable circumstances and situations.

Danger when overcoming stairs/obstacles

WARNING

Risk of falling in case of lacking assistance. Do not ascend or descend stairs without the assistance of attendants. If devices such as ramps or elevators are available, they should be used. If wheelchair-friendly access is missing, two attendants must carry the wheelchair over the obstacle.

WARNING

Risk of falling as a result of incorrectly set anti-tipper. If only one attendant is available when ascending or descending stairs, an incorrectly set anti-tipper (if mounted at all) can lead to severe falls. Adjust the anti-tipper so that it does not come in the way of the steps during transport. Afterwards, swing the anti-tipper back to its operational position.

WARNING

Risk of falling as a result of lifting the wheelchair incorrectly. Attendants must lift the wheelchair only by parts that are firmly attached (e.g. not at the footrest or rear wheels). If your wheelchair is equipped with height-adjustable push handles, ensure that the clamping levers are always firmly tightened.

WARNING

Risk of tipping when using lifting platforms. Ensure that the anti-tipper (if mounted) is outside the danger-area.

Danger when getting into the wheelchair

CAUTION

Risk of injury due to incorrect boarding.

- Activate the wheel lock each time before you get into, out of, or transfer to or from the wheelchair.
- Always get into the wheelchair from the side, where possible.

⚠ CAUTION

Risk of injury as a result of using the footplate for getting into the wheelchair. Depending on footplate settings and wheelchair geometry, the wheelchair may tip over if the user boards the chair using the footplate. First practise boarding the chair with the child and an attendant who can secure the wheelchair, and modify footplate and seat height settings if the chair has a tendency to tip over. In addition, turn the caster fork to the front prior to using the footplate for getting into the wheelchair; this increases the wheelbase and thereby the wheelchair's stability against tipping.

Danger while driving**⚠ CAUTION**

Risk of accident when driving without experience. This can lead to falls and other dangerous situations. Therefore, practise with the new wheelchair on even, straightforward terrain first, together with the child. Together with the child, learn how the wheelchair reacts when the centre of gravity shifts; for example on slopes or inclines or when clearing obstacles such as steps and curbs. This should be done only with assistance from another person. Using an anti-tipper is strongly recommended for inexperienced wheelchair users.

⚠ CAUTION

Risk of tipping when ascending too steep slopes. When ascending slopes or ramps and when crossing obstacles on upward slopes, always lean the wheelchair user's upper body far forward.

⚠ CAUTION

Risk of tilting or tipping over. A displacement of the centre of gravity will result in a loss of stability. The wheelchair could be prone to tilting and tipping.

- Do not lean out of the wheelchair too far.
- Prior to activities that require you to bend forward in the wheelchair (e.g. tying your shoes), maximise the stability of the wheelchair. In order to do so, push the wheelchair backwards until the steering casters turn forward.
- Always ensure two anti-tippers are deployed when using a backpack

⚠ CAUTION

Risk of tipping when driving against obstacles (steps, curbs) without braking. Use your wheelchair properly. Do not "jump" the chair down from higher surfaces.

⚠ CAUTION

Risk of accident as a result of incorrect use of the wheel lock. All brakes acting on the tyres do not serve as service brake but are only designed as parking brake (wheel lock). The wheel locks must not be used as driving brakes for slowing down the wheelchair, as in extreme cases, the abrupt stopping of the wheelchair can lead to falls.

⚠ CAUTION

Risk of accident if the wheelchair starts rolling. Engage the wheel locks to prevent the wheelchair from moving on uneven ground or during transfers (e.g. into a car).

⚠ CAUTION

Risk of accident as a result of driving in the dark. In the dark, the user should wear light clothing or clothing with reflectors in order to improve visibility. Ensure that the reflectors installed on the sides and rear of the wheelchair are easily visible. We also recommend installing active illumination.

⚠ CAUTION

Risk of tipping when descending slopes. When descending slopes, do not drive without braking and reduce your speed. Reduced load on the casters due to centre of gravity shifting can cause the casters to flutter.

⚠ CAUTION

Risk of falling due to lack of tipping resistance on public transport

- Always observe currently applicable legal requirements when using public transport.
- Always ensure that you are held in place securely when travelling on public transport. To do so, use the wheelchair areas, wheelchair bays and restraint systems provided. Firmly tighten the wheel locks.
- Please note that the anti-tipper, if activated, may be exposed to heavy loads when the public transport vehicle starts to move. In order to avoid damage, the manufacturer recommends that the wheelchair is positioned at a right angle to the direction of travel if no restraint system for passengers with reduced mobility is available.

⚠ CAUTION

Risk of pinching between certain components. Avoid reaching into the pinching area

- between the side panel and backrest when folding the backrest;
- between the wheel lock lever and side panel; and
- between the brake block and rear wheel.

Danger as a result of incorrect assembly or settings**⚠ CAUTION**

Risk of tipping as a result of incorrectly set wheels. Extreme settings (e.g. rear wheels mounted in the foremost position) combined with an unfavourable body posture can cause the wheelchair to tip even on level ground.

⚠ CAUTION

Risk of tipping as a result of missing or incorrectly mounted anti-tipper. Under no conditions should the anti-tipper assume the function of transport wheels, for example to transport a person in the wheelchair with the rear wheels removed. The anti-tipper must audibly lock in place, before it is able to bear loads. Firm seating must be verified by the user or by an attendant.

⚠ CAUTION

Risk of tipping as a result of changed diameters/mounting positions of the wheels. If the size and position of the casters or the rear wheel size are modified, the casters may flutter when driving at higher speeds. This can lead to blocking of the casters and cause the wheelchair to tip over. For this reason, you should always maintain the factory settings or, if modifications are required, ensure that the wheelchair frame is aligned horizontally (see "Adjustment / assembly instructions").

⚠ CAUTION

Risk of falling as a result of unintentional loosening of the rear wheels. Always make sure that the quick-release axles are correctly set on the rear wheel. It must not be possible to remove the rear wheel unless the button on the quick-release axle is pushed.

Danger as a result of improper use of the tyres

⚠ CAUTION

Risk of accident as a result of insufficient or excessive tyre pressure. The effectiveness of the wheel lock and the overall driving quality are dependent on adequate air pressure. The wheel locks are only effective with sufficient air pressure and proper adjustment (approx. **3 mm** distance, subject to technical changes).

Before starting to use your wheelchair, check that the tyres are inflated correctly. The required air pressure is printed on the side of the tyre. For rear wheels, it should be **at least 7.5 bar**. With properly inflated rear wheels and even tyre pressure on both sides, your wheelchair is much easier to operate and manoeuvre.

⚠ CAUTION

Risk of accident as a result of bad tyres. Insufficient tread depth of the tyres reduces adhesion. Please note that when you operate your wheelchair in public road traffic, you are subject to road traffic regulations.

Risk of skin damage

⚠ CAUTION

Be careful in case of unhealthy skin. Using the wheelchair with unhealthy skin can cause medically indicated complications such as reddening of the skin and pressure sores. Over several hours of use, the skin on the buttocks, back and rear thighs is subjected to considerable strain. Therefore, always check these areas carefully before using the wheelchair.

No liability will be assumed for health injuries caused by use of this wheelchair with damaged skin.

Danger caused by fire/heat and cold

⚠ CAUTION

Risk of burns when near to fire. The seat and back upholstery might catch fire. Keep the wheelchair away from all sources of ignition, especially lit cigarettes.

⚠ CAUTION

Be careful in case of extreme temperatures. The wheelchair can heat up significantly when in the sun or in the sauna. In extreme cold, there is a risk of hypothermia.

Risk of hand injuries

⚠ CAUTION

Risk of hand injuries. To avoid hand injuries, the user must not grasp between the rear wheel and the wheel lock when driving the wheelchair.

⚠ CAUTION

Risk of pinching between wheelchair components. The wheel lock lever and the clothing protector side panel can represent a risk of crushing with various adjustment options. Teach the child not to place fingers between the clothing protector lip and the rear wheel or in the recess for the wheel lock operating lever.

⚠ CAUTION

Risk of friction burns. Slowing down from high speeds or when descending longer slopes tends to heat up the hands and fingers, especially if using aluminium push rings. When using the wheelchair outdoors, leather gloves should be worn. Gloves provide the wheelchair user with a better grip and protect his or her fingers from dirt and hot metal.

Warnings regarding wheelchair damage**NOTICE**

Damage caused by overload. The maximum load of the BRAVO Racer is 60 kg/132 lbs.

NOTICE

Damage caused by objects lying upon the appliance. If appliances are transported in the folded state, the weight of objects lying upon the appliance can cause deformations and result in problems when unfolding the appliance. Never put heavy objects on a folded rehabilitation device.

NOTICE

Damage caused by wear and tear. If your seat and back upholstery become worn out or damaged, replace it immediately.

NOTICE

Damage caused by incorrect packaging. Only use original packaging for shipment of the product.

2.3 Further notes**INFORMATION**

Even in the event of compliance with all applicable guidelines and standards, it is possible that alarm systems (e.g., in department stores) may respond to your product. Should this happen, remove your product from the area where the alarm was triggered.

2.4 Warning symbols and nameplates

Label / Nameplate	Explanation
	<p>A Type designation B Manufacturer's article number C Maximum load capacity (see Section „Technical data“) D Manufacturer information / address / country of manufacture E Serial number F International Article Number G Read the Instructions for Use before using the product. H CE marking – product safety in conformity with EU Directives</p>
	<p>Anchoring point/securing point or securing the wheelchair in motor vehicles for the transportation of disabled persons.</p>

3 Delivery/preparing the wheelchair for use

3.1 Storage

The wheelchair must be stored in a dry place.

The ambient temperature must be maintained between -10 °C and +40 °C during transportation and storage.

During extended storage, the knee lever wheel lock on wheelchairs with PU tyres must be released since tyre deformation may otherwise result.

INFORMATION

Tyres contain chemical substances that can react with other chemical substances (such as cleaning agents, acids, etc.).

3.2 Preparation for use

⚠ CAUTION

Risk of pinching. Only grip the components shown when folding or collapsing the wheelchair.

Your wheelchair is usually delivered fully assembled, with the backrest folded and the rear wheels removed.

Three simple steps prepare it for use:

- 1) Unlock the backrest by pulling the cable control (Figure 1).
- 2) Fold the backrest to the rear (Figure 2) until it contacts the stop. The angle between the seat and backrest is pre-set to 90° at the factory. The side panels are automatically moved to the

correct position via a guide (Figure 3).

3) Attach the quick-release rear wheels to the quick-release axle housings. To do so, reach into the spokes close to the hub of the rear wheel and press the button on the quick-release axle with your thumb. Then you can easily slide the rear wheel into place (Figure 4). Ensure that the quick-release axle is securely locked in the axle housing!

The overhang of the back upholstery can be easily pushed under the seat upholstery. We recommend keeping the overhang stored like that at all times.

⚠ CAUTION

Risk of tipping. The rear wheels have been pre-positioned according to your indications. Before using the wheelchair, this pre-adjusted position must be checked by the wheelchair user for function and stability against tipping with the assistance of qualified personnel.

4 Transportation

In a few simple steps, the BRAVO Racer can be disassembled into several units which are easier to store. Remove the rear wheels by pressing the button on the quick-release axle as described in the previous chapter (Figure 4). Unlock the backrest using the cable control (Figure 5). Remove the seat padding, if necessary. Now fold the backrest forward until it locks (Figure 6).

CAUTION Risk of accident. During assembly, it is essential to ensure that the backrest is folded back to the stop and the quick-release axles for the rear wheels are securely locked in the axle housings. It must not be possible to remove the wheel unless the button is pressed. Check the function of the wheel locks.

4.1 Transfer (getting into and out of the wheelchair)


There are basically two ways of getting into and out of the wheelchair. Children should use the footplate as a step for getting into and out of the wheelchair.

CAUTION Risk of injury when using the footplate for getting into the wheelchair. Depending on footplate settings and wheelchair geometry, the wheelchair may tip over if the user boards the chair using the footplate. First practise boarding the chair with the child and an attendant who can secure the wheelchair, and modify footplate and seat height settings if the chair has a tendency to tip over. In addition, turn the caster fork to the front prior to using the footplate for getting into the wheelchair; this increases the wheelbase and thereby the wheelchair's stability against tipping.

CAUTION Risk of pinching between the wheel lock lever and side panel. Reaching into the recess of the clothing protector side panel must be avoided.

First activate the wheel locks (Figure 7). Have the child practise getting into and out of the chair until the child can do so confidently. Whether the child requires assistance from an attendant should be decided on an individual basis. The front frame can be used as a grip or support for stepping on the footplate (Figure 8).

CAUTION Risk of pinching between the brake block and rear wheel. When getting into the wheelchair using the front frame, the user must grasp the front frame in front of the brake block.



If the child sits or kneels on the footplate, the child can use the side panels and wheel lock levers as a grip or support when getting into the seat (Figure 9). If entering from the front, the child must turn around in the seat. In this case, the child can use the push handles for support (Figure 10). The wheel locks may not be released until the child is fully seated in the wheelchair!

Reverse steps should be followed for getting out of the chair.

For adolescents of a certain age, weight and physical constitution, it may be more appropriate to get in and out of the wheelchair via the side panels.

To do this, manoeuvre the wheelchair at a 45° angle to the furniture or wheelchair from which the user will be transferred (Figure 11). Activate the wheel locks. The seat/backrest, push rings and side panels can be used for support while getting into and out of the wheelchair (Figure 12).

The user should practise getting into and out of the wheelchair with an attendant until he or she can do so confidently. The user should be able to get into and out of the chair in one step without stopping. Next, the user's feet must be positioned on the footplate (Figure 13) and the seating position can be adjusted, if necessary. Finally, the wheel locks (Figure 14) can be released and the wheelchair can be used. Reverse steps should be followed for getting out of the chair.

4.2 Use in a wheelchair accessible vehicle

⚠ WARNING

Use as a seat for transportation in a wheelchair accessible vehicle

Risk of injury as a result of unauthorised use in a wheelchair accessible vehicle. Whenever possible while travelling in a wheelchair accessible vehicle, passengers should use the seats and vehicle restraint systems already installed in the vehicle. This is the only way to provide optimal protection for the passengers in case of accident. You can use the product as a seat for transporting passengers in wheelchair accessible vehicles when using the relevant Ottobock safety components as well as the appropriate restraint systems.

For further information please refer to our instructions for use „Using your product for transportation in a wheelchair accessible vehicle“, order number 646D158.

The product may be used as a seat in a wheelchair accessible vehicle.

⚠ WARNING**Risk of accident and injury due to improper use as seat during transportation in a wheelchair accessible vehicle.**

Use in a wheelchair accessible vehicle is not permissible if the following options are used:

- 30° back angle adjustment
- folding back.

A list of further options not permitted in a wheelchair accessible vehicle can be found in our instructions for use „Using your product for transportation in a wheelchair accessible vehicle“, order number 646D158.

5 Options

The wheelchair is designed as a modular component system. This means that certain options and accessories can be adapted to your wheelchair. We recommend using a seat cushion at all times. The seat upholstery of the BRAVO Racer has a loop texture so that seat cushions can be attached with hook and loop connection. All seat cushions listed on the order form have the required hook texture.

Furthermore, we would like to introduce you to a selection of our versions and accessories that can make the operation of your wheelchair easier.

5.1 Anti-tipper (Fig. 15, 16)

Prevents the wheelchair from tipping too far backward. The anti-tipper is highly recommended for all children and untrained or unconfident adolescents.

If the anti-tipper is temporarily not required, it can be folded upwards by pressing the button (Figure 16).

⚠ CAUTION Risk of tipping. The anti-tipper must audibly lock in place before it is able to bear loads. The fixed position is to be verified by the user/attendant.

15

16

5.2 Crutch holder (Fig. 17)

The crutch holder makes it possible to transport crutches on the wheelchair.

5.3 Transport wheels (Fig. 18)

With the rear wheels removed, the wheelchair can be used as a rolling base. This allows for manoeuvring through very narrow passage ways (such as bathroom doorways or airplane aisles).

⚠ CAUTION Risk of accident as a result of unavailable wheel lock function. Please note that once the rear wheels are removed, the wheel locks can no longer be activated!

5.4 Tip-assist (Fig. 19)

The tip-assist makes it easier for the attendant to tip the wheelchair (e.g., for pushing the chair over a step).

5.5 Backpack (Fig. 20)

The backpack allows the user to carry objects. It can be easily hung on the push handles.

⚠ CAUTION

Risk of tilting or tipping over. A backpack and its contents attached to the backrest may cause the wheelchair to become unstable and tilt/tip. Therefore Ottobock recommends that two anti-tippers are always deployed when using a backpack.



INFORMATION

Regarding the maximum weight of the backpack and its contents that can be used; this depends on several factors such as adjustments made to the wheelchair and/or weight and sitting position of the user. However the total combined weight of both the user and all accessories including the backpack and its contents should not exceed the maximum load capacity (60Kgs) of the chair. Before use the user/attendant must ensure the stability of the chair is checked.


5.6 Footrest, independently angle adjustable (Fig. 21)

The standard footrest of the BRAVO Racer matches the angle of the seating unit. If a solution independent of the seating unit is required, the angle of the footrest can be individually adjusted by loosening the clamping screws (Fig. 22).

5.7 Lateral heel block (Fig. 23, 24)

This option prevents the feet from uncontrolled slipping off the footplate. The position of the lateral heel block is adjustable.

5.8 Seat upholstery, adaptable (Fig. 25)

This option makes it possible to readjust the tension of the seat upholstery at any time without the use of tools. You can also set the desired seat contour.

5.9 Back upholstery, breathable (Fig. 26)

The standard back upholstery of the BRAVO racer can be adjusted. The option shown here features padding made of spacer fabric, which allows for air circulation. This prevents heat and perspiration buildup during warm weather and keeps body heat in the fabric during cold weather.

5.10 Fold-down push handles (Fig. 27)

A push button is located underneath the push handles; when it is activated, the push handles can be folded down. When they are folded up into the active position, the handles lock automatically.

5.11 Push bar (Fig. 28)

The removable push bar offers maximum comfort for the attendant while controlling the wheelchair. The bar can be extended so that a more ergonomic position for pushing and tipping can be set without the need for tools.

CAUTION Risk of falling as a result of lifting the wheelchair incorrectly. Do not use the push bar for lifting the wheelchair.

5.12 Designer clothing protectors (Fig. 29)

This option enhances the wheelchair without hindering the functionality of the standard clothing protector.

5.13 Side panel with height adjustable armrest (Fig. 30)

The armrests on the side panel offer the user additional forearm support. Their height can be adjusted without tools.


5.14 Tray (Fig. 31)

This option enables activities, games and therapy on a removable tray mounted directly on the wheelchair.

5.15 Designer casters (Fig. 32)

These casters improve the appearance of the wheelchair.


5.16 Luminous wheels (Fig. 33)

These wheels improve the appearance of the wheelchair and also provide a passive safety feature. The self-charging LEDs make the wheelchair more visible in low light or darkness.

5.17 Wheelbase extension with second axle (Fig. 34)

This option is designed especially for the adaptation of hand bikes. When adapting hand bikes, the rear wheels must be moved to the rear axle.

5.18 Infinity Ultralight rim (Fig. 35)

Weight reduction always improves the manoeuvrability of a wheelchair. It also makes loading easier. When every last bit of weight matters, these ultralight rims can reduce weight significantly.

5.19 Rear wheel with drum brake (Fig. 36)

Choose this option if a brake for the attendant is required instead of or in addition to the standard user wheel lock. The drum brake is activated from the push handle and safely locks the rear wheels.

5.20 Lap belt (Fig. 37)

Use of this option effectively protects users with limited upper body control from sliding forward. Alternatively, the lap belt provides useful support during sporting activities.

6 Adjustment/assembly instructions

All wheelchair users have their own expectations regarding wheelchair type, accessories, seating and handling. For this reason, users can adjust the BRAVO Racer active wheelchair in a variety of different ways to suit their needs. We recommend that you and the user consult with your dealer or therapist to determine which wheelchair settings will work best.

INFORMATION

The following tools are required for repair and maintenance:


- Allen wrenches, sizes 3 mm, 4 mm and 5 mm
- Open-end wrenches, sizes: 10 mm, 17 mm, 19 mm and 24 mm
- Set of screwdrivers
- Torque wrench
- Tyre levers

6.1 Adjusting the height of the footrest (Fig. 38, 39)

To adapt to the lower leg length, the footrest bar is continuously adjustable. Loosen the inner screws on the back of the clamps that guide the footrest bar and cross brace (Figure 39, marking A) until you are able to move the footrest bar in the clamp. Set the desired lower leg length by moving the footrest forward or back. Ensure that the footrest bar is in the same position on both sides.

If you are unable to reach the desired lower leg length, remove the footrest bar completely and turn it around. Then proceed as described in 6.2. After adjustment, snugly retighten the screws.

INFORMATION Ensure that all screws are properly re-tightened before using the footrest. The correct torque for the screw (Figure 39, marking A) is **8 Nm**.

INFORMATION Make sure that the footrest bars are guided completely in the clamp profiles to achieve a sufficient clamping effect.


6.2 Longitudinal adjustment of the footrest (Fig. 40, 41)

To adjust the footplate longitudinally, loosen the four screws on the top of the footplate until you can move the footplate on the footrest bars (Figure 40). Tighten the screws once you have reached the desired position. If this adjustment is not sufficient, completely unscrew the four screws and move the footplate all the way to the front or rear (Figure 41). Then reinsert the screws, position the footplate as desired, and tighten the screws.

If you turned the footrest bar around to reach the correct lower leg length, completely remove the four screws. Pull the footplate forward away from the footrest bars, turn it around and put it back into place. Then make settings as described above.

If you cannot reach the desired lower leg length because neither assembly method allows the adjustment range of footrest bar to be positioned correctly, use the spacers provided. The spacers can be installed between the footplate and the footrest bar.

INFORMATION The correct torque for the screw (Figures 40 and 41) is **10 Nm**.



6.3 Adjusting the seat height (Fig. 42, 43)

Completely open and remove the two nuts on the front of the adjustment tube (Figure 42, marking A). Pull the U-shaped part back (Figure 43) until you can freely move the adjustment tube up and down. Proceed in the same way on both sides.

Loosen the two outer screws on the back of the clamps that guide the footrest bar and cross brace until the cross brace moves easily in the clamp (Figure 39, marking B).

Reposition the adjustment tube over the adjustment hole channel (Figure 42, marking B) in the desired location and reinstall it by repeating the steps described above in the reverse sequence. Firmly retighten the screws previously loosened.

INFORMATION Ensure that all screw connections are properly re-tightened before using the wheelchair. The correct torque for the screw (Figure 42, marking A) is **10 Nm**.



6.4 Adjusting the seat angle (Fig. 44, 45, 46)

Slightly loosen the two nuts on the front of the adjustment tube (Figure 44, marking A). Proceed in the same way on both sides.

Loosen the two outer screws on the back of the clamps that guide the footrest bar and cross brace as well as the central recessed screws until the cross brace moves easily in the clamp (Figure 45, marking A/B).

Now change the seat angle by moving the cross brace up and down in the clamps. The driving axle is the centre of rotation (Figure 46). The angle of the footrest moves automatically with the seat angle. Once you have reached the desired angle, retighten all screws.

INFORMATION The rear wheel adapter's hexagon exclusively serves the form-fit attachment of frame accessories. Never try to loosen this connection using a wrench. The correct torque for the screw (Figure 44, marking A) is **10 Nm**.

INFORMATION Ensure that the vertical position of the adjustment tube on the axle retainer/base frame is always the same on both sides.

INFORMATION Ensure that all screw connections are properly re-tightened before using the wheelchair.

6.5 Fine adjustment of seat/footplate angle (Fig. 47, 48)



A profile rail is located on the underside of the seat tube (Figure 48). It is attached on both sides by four nut sliders.

On both sides, loosen the screws of the second nut slider from the front (Figure 47, Marking A/ Figure 48, marking B). These sliders connect the cross brace to the seat tubes.

Loosen the two central recessed screws on the back of the clamps that guide the footrest bar and cross brace until the cross brace moves easily in the clamp (Figure 47, marking B).

Now you can make fine adjustments to the seat and footplate angle by moving the cross brace on the seat tube. Once you have found the desired position, retighten all screws.

INFORMATION Ensure that all screw connections are properly re-tightened before using the wheelchair. The correct torque for the screw (Fig. 47, marking A/B) is **10 Nm**.

6.6 Longitudinal adjustment of the seat unit (centre of gravity) (Fig. 48)

A profile rail is located on the underside of the seat tube. It is attached on both sides by four nut sliders.

First loosen the counter nut at the upper end of the adjustment tube. Then loosen the screws on both sides of the nut sliders (Fig. 48, marking A). You can move the seat tubes lengthwise and thus, change the centre of gravity. Different adjustment ranges, marked as a grid at the bottom of the tube, help to move both tubes evenly.

Once you have reached the desired position, retighten all screws.

INFORMATION Ensure that the adjustments at both tubes are parallel and that all screw connections are properly re-tightened before using the wheelchair.

6.7 Adjusting the back unit (Fig. 49)

A profile rail is located on the underside of the seat tube. It is attached on both sides by four nut sliders.

Loosen each of the rear screws on both sides (Figure 49, marking A). Next, loosen the screws under the clothing protector lip in the rear wheel housing (Figure 49, marking B).

Now you can freely move the back unit on the seat tube. This allows you to adjust seat depth or position the tube for a firm back shell.

Use the increment guide (Figure 49, marking C) on the bottom of the seat tube to ensure parallel adjustment.

INFORMATION Ensure that all screw connections are properly tightened before using the wheelchair. The correct torque for the screw (Figure 49, marking A) is **10 Nm**.


6.8 Changing the back angle (Fig. 50)

It may be useful to customize the preset angle between the seat and the backrest for therapeutic reasons. To do this, loosen the screws under the clothing protector lip in the rear wheel housing (Figure 50, marking A).


Next, remove the rear screw from the perforated plate of the back unit (Figure 50, marking B). Remove the screws on both sides from their current position, move the backrest to the desired angle, and replace the screws in the corresponding holes. For reassembly, proceed in reverse order. Finally, retighten all screws.

INFORMATION The correct torque for the screw (Figure 50, marking B) is **10 Nm**.

6.9 Readjusting seat upholstery tension (Fig. 51, 52)

Should the seat upholstery stretch after extended use, it can be retightened. To do this, remove the seat cushion and open the hook-and-loop strap underneath the seat upholstery on the left side (relative to driving direction) (Figure 51, marking A). Remove the protective cap from the profile tube (Figure 51, marking B) and pull out the upholstery bar (Figure 51, marking C). Now remove the seat upholstery from the profile rail (Figure 52, marking A) and use the hook-and-loop closure to tighten it, as necessary (Figure 52, marking B). Next, thread the seat upholstery back into the profile rail and reinsert the upholstery rail and protective cap. Then replace the seat cushion, ensuring that it is properly aligned.

6.10 Adjusting the wheel lock

The wheel locks are activated via cables. If braking function is insufficient due to tyre or material wear, the wheel lock must be readjusted. To do this, open the nuts on the wheel lock lever component group on the lower frame and set the desired cable tension by loosening or tightening the nuts. There should be approximately 3 mm between the open brake block and the outer wall of the tyre. It is important to use the tension setting to position the brake block so that the rear wheels are fully blocked when the wheel lock is locked.

CAUTION **Risk of accident as a result of insufficient or excessive tyre pressure.** Before starting to use your wheelchair, check that the tyres are inflated correctly. The required air pressure is printed on the side of the tyre. For rear wheels, it should be at least **7.5 bar**. The wheel locks are only effective with sufficient air pressure and proper adjustment (approx. **3 mm** distance, subject to technical changes).

INFORMATION Check the tyre air pressure and the effectiveness of the locking devices regularly. The air pressure is indicated on the side wall of the tyre.


INFORMATION Ensure that all screws and nuts are properly retightened after making any adjustments.

6.11 Installing the frame accessories (Fig. 53, 54)

INFORMATION The rear wheel adapter's hexagon exclusively serves the form-fit attachment of frame accessories. Never try to loosen this connection using a wrench.

For subsequent installation of mounting accessories (anti-tipper, tip-assist, crutch holder, transport wheels), first remove the rear wheels. First, proceed as described in section 6.3 until you have loosened the U-shaped part. Next, completely remove the U-shaped part by pulling it back (Figure 53). Place the holder for the frame accessories with the hexagonal clamp on the hexagonal profile of the rear wheel axle (Figure 54, marking A) and tighten the clamp. Reinstall all parts in reverse order. Now you can install the desired frame accessories on the holder, for example: anti-tipper, crutch holder, tip-assist. The anti-tipper is the only frame accessory that can be swung out of the way when it is not needed (see Section 5.1).

INFORMATION The correct torque for the screw (Figure 53) is **10 Nm**.

7 Maintenance, cleaning, and disinfection

7.1 Maintenance

Your wheelchair has been provided with the CE marking. The manufacturer herewith guarantees that this medical product as a whole conforms to the requirements of the European Directive for Medical Products 93/42/EEC.

The proper function of the wheelchair, especially of the wheel locks, should be checked before every use. Safety nuts should be used only once. If they have been loosened several times, they must be replaced.

In the case of defects, contact your authorized dealer immediately to eliminate them. This applies in particular when the driving behaviour changes or the wheelchair becomes unstable. **We also recommend that you have the wheelchair serviced every 12 months by your Ottobock authorized dealer.** The following table lists functions which the user must check at the indicated intervals.

Check	before each use	monthly	quarterly
Function test of the wheel locks/brakes	X		
Sagging of the seat or back upholstery		X	
Check stability of footplate		X	
Visual inspection of wear parts (e.g. tyres, bearings)		X	
Dirt on bearings		X	
Push rings for damages		X	
Air pressure (see indication on tyre)		X	
Folding mechanism wear		X	
Spoke tension of the rear wheels			X
Screw connections			X

To ensure smooth operation at all times, the following precautions and cleanings should be carried out. Users may carry these out by themselves.

- Hair and dirt particles generally accumulate between the caster wheel and fork. This can restrict the caster wheels from rotating smoothly. Remove the caster and thoroughly clean the fork and caster using a mild detergent.
- The rear wheels feature a quick-release system. To keep this system operational, ensure that no dirt adheres to the quick-release axle or axle housing. The quick-release axle should also be lightly lubricated regularly with resin-free sewing machine oil.
- If your wheelchair gets wet, towel-dry it as soon as possible.
- Do not use the wheelchair in salt water and keep sand or other particles from damaging the wheel bearings.

- Screw connections should be checked frequently, in particular when beginning to use the wheelchair and after any adjustment. If a screw connection becomes loose repeatedly, consult your dealer.

7.2 Cleaning and disinfection

Clean the padding and upholstery with warm water and a washing-up liquid. Remove spots with a sponge or a soft brush. Rinse with clear water and let the components dry.

INFORMATION Do not use any aggressive cleansing agents, solvents, or hard brushes etc.

INFORMATION No wet-cleaning. No machine-washing. For disinfection, water-based agents (e.g., Sagrotan original concentrate) should be used. The manufacturer's instructions for use must be observed.

INFORMATION Before disinfecting, the pads and handles must be cleaned. Plastic parts, frame parts, and the chassis and wheels can be cleaned with a mild detergent and a damp cloth. Dry thoroughly.

8 Changing tyres (Fig. 55–58)

Changing flat tyres requires only the necessary tools and minimal technical skills. Users may change tyres themselves if they wish. We recommend carrying an emergency repair kit and tyre pump at all times (if using pneumatic tyres).

Suitable tyre pumps are listed on the order form and are supplied with the product. Alternatively, users can keep a spray can of hardening foam (available at bicycle shops) to fill a tyre in case of punctures.

- In case of a flat tyre, carefully remove the tyre from the rim using the appropriate tools.
- Take care not to damage the rim or the inner tube.
- Repair the inner tube according to the directions in the repair kit, or, replace the damaged inner tube with a new one.
- Before re-mounting the tyre, inspect the inner rim surface and the inside of the tyre for any object that may have caused the puncture.
- Be sure that the rim bands are not damaged. They protect the tube from being damaged by the ends of the spokes.

Mounting

- Place the protective rim band over the valve. Then push the valve through the rim. Screw on the valve nut. The rim band can then be easily applied.
- Ensure that all spoke ends are covered.

Tyres

- Beginning behind the valve, place the lower side of the tyre over the edge of the wheel rim. Pump the tube with a small amount of air until it assumes its round shape, and then insert it into the tyre.
- Check to ensure the inner tube is properly positioned all around and is free of folds; if this is not the case, release some of the air. Now you can easily install the upper side of the tyre starting across from the valve, using both hands and working towards the valve.

Inflating

- Check around the tyre on both sides to ensure that the tube is not clamped between the edge of the tyre and the rim.
- Push the valve back slightly and then pull it out again to ensure that the tyre is properly seated in the valve area.
- Inflate the tyre until it can still be pressed easily with your thumb. If the control line on both sides of the tyre is equally distanced to the rim, the tyre sits centrally. If the tyre does not sit centrally, let out some air and realign the tyre. Then inflate the tyre to maximum operating pressure (see sidewall) and tighten the dust cap.

9 Re-use / disposal

9.1 Information on re-use

⚠ CAUTION

Danger due to improper use. Seat paddings that come directly into contact with the skin can cause functional and hygienic risks if the product is used by another person. They must be replaced for re-use.

The product is suitable for re-use.

Similar to second-hand machines or cars, products that are being re-used are subject to increased strain. Features and functions must not change in a way that could endanger patients or other persons within the product's life cycle.

Based on market observations and the current state of technology, the manufacturer has calculated that the product can be used for a period of **4 years**, provided that it is used properly and that the service and maintenance instructions are observed. Periods during which the product is stored at the dealer or with the third party payer are not included in this period. It should be clearly pointed out, however, that the product is reliable far beyond this defined period of time, provided that it is cared for and maintained appropriately.

In cases of re-use, the corresponding product must first be thoroughly cleaned and disinfected. Subsequently, the product must be examined by an authorized specialist to check the condition, wear and possible damages. Any worn and damaged components as well as components which do not fit or are not suitable for the user must be replaced.

The service manual includes a service schedule, detailed information and a list of the required tools.

9.2 Information on disposal

If a product is to be disposed of, all components of the product must be disposed of properly in accordance with the applicable environmental protection regulations in the respective country.

10 Technical data

Seat width:	20 – 36	
Seat depth:	24 – 38	
Back height:	20 – 40	
Back angle:	-15° – +15° adjustable, can additionally be engaged in position +15°	
Lower leg length:	14.5 – 44.5	
Seat height:	31.5 – 46	with 20" rear wheel
	34 – 48.5	with 22" rear wheel
	36.5 – 51	with 24" rear wheel
Overall width:	47 – 63	
	48 – 66	with 9° camber
Overall length:	66.5 – 81.5 with 20" rear wheel	
	69 – 84 with 22" rear wheel	
	71.5 – 86.5 with 24" rear wheel	
Overall height:	51.5 – 86	with 20" rear wheel
	54 – 88.5	with 22" rear wheel
	56.5 – 91	with 24" rear wheel
	} without push handles	
Weight:	approx. 8.5 kg	
Load capacity:	60 kg	

Seat height (in cm)¹⁾

Rear wheel size		20"			22"			24"			
		3°	6°	9°	3°	6°	9°	3°	6°	9°	
Camber											
	Position										
Adjustment tube, long	Adjustment tube, short	1	32.2	32.0		34.7	34.5	33.9	37.2	37.0	36.4
		2	33.9	33.7		36.4	36.2	35.6	38.9	38.7	38.1
		3	35.6	35.5		38.1	38.0	37.3	40.6	40.5	39.8
		4	37.3	37.2		39.8	39.7	39.0	42.3	42.2	41.5
	5	39.1	38.9		41.6	41.4	40.8	44.1	43.9	43.3	
	6	40.8	40.6		43.3	43.1	42.5	45.8	45.6	45.0	
	7	42.5	42.4		45.0	44.9	44.2	47.5	47.4	46.7	
	8	44.2	44.1		46.7	46.6	45.9	49.2	49.1	48.4	
	9	46.0	45.8		48.5	48.3	47.7	51.0	50.8	50.2	

1) Specified at 0° seat angle, without seat cushion

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Otto Bock South Africa (Pty) Ltd
Building 3 Thornhill Office Park · 94 Bekker Road
Midrand · Johannesburg · South Africa
T +27 11 312 1255
info-southafrica@ottobock.co.za
www.ottobock.co.za

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Otto Bock HealthCare Canada
5470 Harvester Road
Burlington, Ontario, L7L 5N5, Canada
T +1 289 288-4848 · F +1 289 288-4837
infocanada@ottobock.com · www.ottobock.ca

Otto Bock HealthCare Andina Ltda.
Calle 138 No 53-38 · Bogotá · Colombia
T +57 1 8619988 · F +57 1 8619977
info@ottobock.com.co · www.ottobock.com.co

Otto Bock de Mexico S.A. de C.V.
Prolongación Calle 18 No. 178-A
Col. San Pedro de los Pinos
C.P. 01180 México, D.F. · Mexico
T +52 55 5575 0290 · F +52 55 5575 0234
info@ottobock.com.mx · www.ottobock.com.mx

Otto Bock HealthCare
Two Carlson Parkway North, Suite 100
Minneapolis, MN 55447 · USA
T +1 763 553 9464 · F +1 763 519 6153
usa.customerservice@ottobockus.com
www.ottobockus.com

Asia/Pacific

Otto Bock Australia Pty. Ltd.
Suite 1.01, Century Corporate Centre
62 Norwest Boulevard
Baulkham Hills NSW 2153 · Australia
T +61 2 8818 2800 · F +61 2 8814 4500
healthcare@ottobock.com.au · www.ottobock.com.au

Beijing Otto Bock Orthopaedic Industries Co., Ltd.
B12E, Universal Business Park
10 Jiuxianqiao Road, Chao Yang District
Beijing, 100015, P.R. China
T +8610 8598 6880 · F +8610 8598 0040
news-service@ottobock.com.cn
www.ottobock.com.cn

Otto Bock Asia Pacific Ltd.
Unit 1004, 10/F, Greenfield Tower, Concordia Plaza
1 Science Museum Road, Tsim Sha Tsui
Kowloon, Hong Kong · China
T +852 2598 9772 · F +852 2598 7886
info@ottobock.com.hk · www.ottobock.com

Otto Bock HealthCare India Pvt. Ltd.
20th Floor, Express Towers
Nariman Point, Mumbai 400 021 · India
T +91 22 2274 5500 / 5501 / 5502
information@indiaottobock.com · www.ottobock.in

Otto Bock Japan K. K.
Yokogawa Building 8F, 4-4-44 Shibaura
Minato-ku, Tokyo, 108-0023 · Japan
T +81 3 3798-2111 · F +81 3 3798-2112
ottobock@ottobock.co.jp · www.ottobock.co.jp

Otto Bock Korea HealthCare Inc.
4F Agaworld Building · 1357-74, Seocho-dong
Seocho-ku, 137-070 Seoul · Korea
T +82 2 577-3831 · F +82 2 577-3828
info@ottobockkorea.com · www.ottobockkorea.com

Otto Bock South East Asia Co., Ltd.
1741 Phaholyothin Road
Kwaeng Chatuchark · Khet Chatuchark
Bangkok 10900 · Thailand
T +66 2 930 3030 · F +66 2 930 3311
obsea@ottobock.co.th · www.ottobock.co.th

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Ottobock has a certified Quality Management System in accordance with ISO 13485.

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Otto Bock HealthCare GmbH
Otto Bock Mobility Solutions GmbH