Motus CV/CS

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EN**Instructions for use for the Motus adaptive wheelchair**

Contents	Page
1 General information	43
1.1 <i>Preface</i>	43
1.2 <i>Intended use</i>	43
1.3 <i>Field of application</i>	43
1.4 <i>Legal information</i>	44
1.5 <i>Service</i>	45
2 Safety instructions.....	45
2.1 <i>Explanation of symbols</i>	45
2.2 <i>General safety instructions</i>	45
2.3 <i>Further notes</i>	50
2.4 <i>Warning symbols and type plates</i>	50
3 Delivery/Preparing the wheelchair for use	51
4 Transportation	52
4.1 <i>Transfer</i>	53
5 Options	55
5.1 <i>Anti-tipper</i>	55
5.2 <i>Transport wheels</i>	55
5.3 <i>Crutch holder</i>	56
5.4 <i>Tip-assist</i>	56
5.5 <i>Plug-on side panels</i>	56
5.6 <i>Desk side panels</i>	57
5.7 <i>Plug-on wheel lock lever extension</i>	57
5.8 <i>Knee lever wheel lock for one-handed operation</i>	57
5.9 <i>Rear wheel with drum brake</i>	58
5.10 <i>Tray</i>	58
5.11 <i>Elevating footrest</i>	58
5.12 <i>Motus CS</i>	58
5.13 <i>Push handles, height-adjustable</i>	58
5.14 <i>Wheelbase extension</i>	59
5.15 <i>Angle adjustable back</i>	59
5.16 <i>Spoke protector</i>	59
6 Adjustment/Assembly instructions	60
6.1 <i>Adjusting seat height and seat inclination</i>	60
6.2 <i>Changing the wheelbase</i>	61

6.3	<i>Setting the wheel camber and toe-in</i>	61
6.4	<i>Caster trailing adjustment</i>	63
6.5	<i>Adjusting the wheel lock</i>	64
6.6	<i>Adjusting the brake force - drum brake</i>	64
6.7	<i>Angle adjustment of the footplate</i>	65
6.8	<i>Adjusting the lower leg length</i>	65
6.9	<i>Adjusting the angle of the elevating footrest</i>	65
6.10	<i>Adjusting the height of the armrest for desk side panel</i>	66
6.11	<i>Adjusting the back height</i>	66
6.12	<i>Installing the anti-tipper</i>	67
6.13	<i>Adapting the back upholstery</i>	68
6.14	<i>Changing the back angle</i>	68
6.15	<i>Push handles, height-adjustable</i>	68
6.16	<i>Attaching the safety belt</i>	69
7	<i>Maintenance, cleaning, and disinfection</i>	69
7.1	<i>Maintenance</i>	69
7.2	<i>Cleaning and disinfection</i>	70
8	<i>Changing tyres</i>	71
9	<i>Technical data</i>	73
10	<i>Re-use/Disposal</i>	75
10.1	<i>Information on re-use</i>	75
10.2	<i>Information on disposal</i>	75

1 General information

INFORMATION

Date of the last update: 2015-10-06

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

INFORMATION

- You can request this document as a PDF file from the Customer Care Center at oa@ottobock.com or from the manufacturer's service department.
- It is possible to increase the display size of the PDF document.
- For further questions about the instructions for use, please contact the authorised personnel who issued the product to you.

1.1 Preface

With the purchase of the Motus adaptive wheelchair, you have selected a quality product for versatile use during everyday life, both indoors and outdoors.

The chapter "Options" presents various additional components of the Motus, which can expand its field of application and improve comfort for the wheelchair user. The chapter "Adjustment/Assembly information" provides an overview of the ways the wheelchair can be customised to your requirements. The model as described in these instructions for use is subject to technical changes without notice.

1.2 Intended use

The Motus adaptive wheelchair is designed solely for individual indoor and outdoor use by persons who are unable to walk or who have a walking impediment, and can be operated by the patient or by another person.

The Motus adaptive wheelchair may only be combined with the options mentioned in these instructions for use and vice versa. Ottobock assumes no liability for combinations with medical products and/or accessories from other manufacturers outside of the modular system.

1.3 Field of application

The **Motus CV** features advanced functions. Thanks to its frame geometry, it enables precise leg guidance. The footrests can be swung away and removed.

Thanks to its closed, stiff frame geometry, particular active driving is possible with the **Motus CS**. The footrest remains on the wheelchair.

The versatility of this power wheelchair and its modular design make it suitable for patients who have walking impediments or walking disabilities due to:

- Palsies/paralysis

- Loss of limbs
- Defective or deformed limbs
- Joint contractures or defects
- Muscle and nerve disorders
- Diseases such as cardiac or circulatory insufficiency, balance disorders or cachexia as well as for geriatric patients who still have usable residual strength in the upper limbs.

The Motus was specially designed for individuals who are able to independently move in a wheelchair.

Fitting considerations:

- Body size and body weight
- Physical and psychological limitations
- Age of the user
- Home and living conditions, and
- Environment.

1.4 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

1.4.1 Service Life

Expected service life: **4 years**.

The design, manufacturing and requirements for the intended use of the product are based on the expected service life. These also include the requirements for maintenance, ensuring effectiveness and the safety of the product.

Using the product beyond the specified expected service life leads to increased residual risk and should only take place subject to the due diligence and deliberations of qualified personnel.

If the service life is reached, the user or a responsible attendant should contact the qualified personnel who fitted the product or the manufacturer's servicing department. Here the user can obtain information about known risks and the current options for refurbishing the product.

1.4.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

1.4.3 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

1.4.4 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.


1.5 Service

Service and repairs on the product may only be carried out by authorised dealers. Should any problems arise, please contact the supplier of your product. Any necessary repairs will be made exclusively with authentic Ottobock spare parts there. Your product requires regular maintenance/ regular service .

Your authorised Ottobock dealer:

2 Safety instructions

2.1 Explanation of symbols

 WARNING	Warnings regarding possible risks of severe accident or injury.
 CAUTION	Warnings regarding possible risks of accident or injury.
 NOTICE	Warnings regarding possible technical damage.
 INFORMATION	Information regarding operation. Information for service personnel.

2.2 General safety instructions

INFORMATION

Please read the instructions for use first! Before using the product, you should become familiar with the handling, function and use of the product. You may compromise the safety of your child if you do not observe these instructions.

However, it is impossible to protect against and avoid all unforeseeable circumstances and situations.

⚠ WARNING**Use as a seat for transportation in a wheelchair accessible**

Risk of injury as a result of unauthorised use in a wheelchair accessible. We recommend that, wherever and whenever possible, Motus users transfer to the seats installed in the motor vehicle and use the corresponding vehicle restraint systems, because this is the only way to ensure optimal protection of the passengers in case of an accident. Only this can ensure optimal protection of the passengers in case of an accident. Your Motus adaptive wheelchair is permissible for transporting passengers in motor vehicles when using Ottobock safety components as well as appropriate restraint systems.

For more information, please refer to our instructions for use manual "Using your wheelchair/mobility base with seating shell or buggy as a seat for transportation in motor vehicles", order no. 646D158.

⚠ WARNING**Serious injuries due to exceeding the service life.**

- Using the product beyond the specified expected service life leads to increased residual risk and should only take place subject to the due diligence and deliberations of qualified personnel.
- If the service life is reached, the user or a responsible attendant should contact the qualified personnel who fitted the product or the manufacturer's servicing department. Here the user can obtain information about known risks and the current options for refurbishing the product.

Danger when overcoming stairs/obstacles**⚠ WARNING**

Risk of falling in case of lacking assistance. Do not ascend or descend stairs without the assistance of attendants. Make use of ramps or elevators whenever available. If wheelchair-friendly access is missing, two attendants must carry the wheelchair over the obstacle.

⚠ WARNING

Risk of falling as a result of incorrectly set anti-tipper. If only one attendant is available when ascending or descending stairs, an incorrectly set anti-tipper can lead to severe falls. Adjust the anti-tipper so that it does not come in the way of the steps during transport. Afterwards, swing the anti-tipper back to its operational position.

⚠ WARNING

Risk of falling as a result of lifting the wheelchair incorrectly. Attendants must lift the wheelchair only by parts that are firmly attached. If your wheelchair is equipped with height-adjustable push handles, ensure that the clamping levers are firmly tightened.

⚠ WARNING

Risk of tipping when using lifting platforms. Ensure that the anti-tipper is outside the danger-area.

Danger when getting into the wheelchair**⚠ CAUTION**

Risk of injury as a result of using the footplate for getting into the wheelchair. If possible, the footplates should be flipped up before getting into or out of the wheelchair.

Danger while driving**⚠ CAUTION**

Risk of accident when driving without experience. This can lead to falls and other dangerous situations. Therefore, familiarise yourself with the new wheelchair by practicing on even, straightforward terrain first. Get to know how the wheelchair reacts when the centre of gravity shifts; for example on slopes or inclines or when clearing obstacles such as steps and curbs. This should be done only with assistance from another person. Using the anti-tipper is strongly recommended for untrained wheelchair users.

⚠ CAUTION

Risk of tipping when ascending too steep slopes. When ascending slopes or ramps and when crossing obstacles on upward slopes, always lean your upper body far forward.

⚠ CAUTION

Risk of tilting or tipping over as a result of centre of gravity shifts. When reaching for objects in front, to the side or behind the wheelchair, do not lean out of the wheelchair too far.

⚠ CAUTION

Risk of tipping when driving against obstacles without braking. Use your wheelchair properly. Do not “jump” the chair down from higher surfaces.

⚠ CAUTION

Risk of accident as a result of incorrect use of the wheel lock. All brakes acting on the tyres do not serve as operating brake but are only designed as parking brake. The wheel locks must not be used as driving brakes for slowing down the wheelchair, as in extreme cases, the abrupt stopping of the wheelchair can lead to falls.

⚠ CAUTION

Risk of accident if the wheelchair starts rolling. Engage the wheel locks to prevent the wheelchair from moving on uneven ground or during transfers .

⚠ CAUTION

Risk of accident as a result of driving in the dark. In the dark, wear light clothing or clothing with reflectors in order to improve your visibility. Ensure that the reflectors installed on the sides and rear of the wheelchair are easily visible. We also recommend installing lights.

⚠ CAUTION

Risk of tipping when descending slopes. Especially on slopes, do not drive without braking and reduce your speed.

⚠ CAUTION

Risk of falling due to lack of tipping resistance on public transport

- Always observe currently applicable legal requirements when using public transport.
- Always ensure that you are held in place securely when travelling on public transport. To do so, use the wheelchair areas, wheelchair bays and restraint systems provided. Firmly tighten the wheel locks.
- Please note that the anti-tipper, if activated, may be exposed to heavy loads when the public transport vehicle starts to move. In order to avoid damage, the manufacturer recommends that the wheelchair is positioned at a right angle to the direction of travel if no restraint system for passengers with reduced mobility is available.

Danger as a result of incorrect assembly or settings**⚠ CAUTION**

Risk of tipping as a result of incorrectly set wheels. Extreme settings combined with an unfavourable body posture can cause the wheelchair to tip even on level ground. Avoid extreme settings. For transfemoral amputees, the rear wheels must be shifted backwards. Use of an anti-tipper is required.

⚠ CAUTION

Risk of tipping as a result of missing or incorrectly mounted anti-tipper. Under no condition should the anti-tipper be used as transport wheels, e.g. when wanting to push the wheelchair while the rear wheels are removed. The anti-tipper must audibly lock in place, before it is able to bear loads. Firm seating must be verified by the user or by an attendant.

⚠ CAUTION

Risk of tipping as a result of changed diameters/mounting positions of the wheels. If the size and position of the casters or the rear wheel size are modified, the casters may flutter when driving at higher speeds. This can lead to blocking of the casters and cause the wheelchair to tip over. For this reason, you should always maintain the factory settings or, if modifications are required, ensure that the wheelchair frame is aligned horizontally.

⚠ CAUTION

Risk of falling as a result of unintentional loosening of the rear wheels. Always make sure that the quick-release axles are correctly set on the rear wheel. The rear wheel must not be removable without the button of the quick-release axle being depressed.

Danger as a result of improper use of the tyres**⚠ CAUTION**

Accident risk due to insufficient or excessive tyre pressure! The effectiveness of the wheel lock and the overall driving quality are dependent on adequate air pressure. The wheel locks are only effective with sufficient air pressure and proper adjustment.

Before starting to use your mobility base with seating shell, check that the tyres are inflated correctly. The required air pressure is printed on the side of the tyre. For rear wheels, it should be **at least 7.5 bar**. Properly inflated rear wheel tyres and equal air pressure in both tyres considerably improve the manoeuvrability of your wheelchair.

⚠ CAUTION

Risk of accident as a result of bad tyres. Insufficient tread depth of the tyres reduces adhesion. When using your wheelchair in public traffic, you are subjected to public traffic regulations.

Risk of skin damage**⚠ CAUTION**

Possible complications in case of unhealthy skin. Using the wheelchair with unhealthy skin can cause medically indicated complications such as reddening of the skin and pressure sores. Over several hours of use, the skin on the buttocks, back and rear thighs is subjected to considerable strain. Therefore, always check these areas carefully before using the wheelchair.

No liability will be assumed for health injuries caused by use of this wheelchair with damaged skin.

Danger caused by fire/heat and cold**⚠ CAUTION**

Risk of burns when near to fire. The seat and back upholstery may catch fire. Keep the wheelchair away from all sources of ignition, especially lit cigarettes.

⚠ CAUTION

Be careful in case of extreme temperatures. The wheelchair can heat up significantly due to solar radiation or in the sauna. In extreme cold, there is a risk of hypothermia.

Risk of hand injuries**⚠ CAUTION**

Risk of hand injuries. To avoid hand injuries, do not grasp between the rear wheel and the wheel lock when driving your wheelchair.

⚠ CAUTION

Risk of pinching between wheelchair components. The wheel lock lever and the side or frame sections can represent a risk of crushing with various adjustment options.

⚠ CAUTION

Risk of friction burns. Slowing down from high speeds or when descending longer slopes tends to heat up the hands and fingers, especially if using aluminum push rings. When using your wheelchair outdoors, wear leather gloves. Gloves provide you with a better grip and protect your fingers from dirt and hot metal.

Warnings regarding wheelchair damage**NOTICE**

Damage caused by overload. The maximum load capacity of the Motus adaptive wheelchair is 125 kg or 140 kg in connection with double crossbrace.

NOTICE

Damage caused by objects lying upon the device. If devices are transported in the folded state, the weight of objects lying upon the device can cause deformations and result in problems when unfolding the device. Never put heavy objects on a folded rehabilitation device.

NOTICE

Damage caused by wear and tear. If the upholstery gets worn or damaged, replace it immediately.

NOTICE

Damage caused by incorrect packaging. Only use original packaging for shipment of the product.

2.3 Further notes

INFORMATION

Even in the event of compliance with all applicable guidelines and standards, it is possible that alarm systems may respond to your product. Should this happen, remove your product from the area where the alarm was triggered.

2.4 Warning symbols and type plates

Label/Type plate	Explanation
	A Type designation B Manufacturer's article number C Maximum load capacity D Manufacturer information/address/country of manufacture E Serial number F International article number G Read the instructions for use before using the product. H CE marking – product safety in conformity with EU directives
Achtung! Die Antriebsräder sind Ihren Angaben entsprechend vorpositioniert. Die voreingestellte Position muß vor Inbetriebnahme vom Rollstuhlfahrer mit sichernder Unterstützung durch Fachpersonal auf Kippsicherheit und Funktion überprüft werden. Attention! The rear wheels are positioned according to your specifications. Before putting the wheelchair into operation the pre-adjusted position must be checked for tilting safety and proper functioning by the user supported safely by specialist staff.	Attention! The rear wheels are positioned according to your instructions. Before putting the wheelchair into operation, the pre-adjusted position must be checked for tilting safety and proper functioning by the user supported safely by the specialist staff.
	Anchoring point/securing point for securing the wheelchair in wheelchair accessible vehicles.

3 Delivery/Preparing the wheelchair for use

CAUTION

Risk of pinching. Only grip the components shown when folding or collapsing the wheelchair.

Your wheelchair is usually delivered fully assembled with the backrest folded and with the rear wheels removed.

A few simple steps are sufficient to prepare it for use:

1. Attach rear wheels. To attach the rear wheel, press the button on the quick-release axle, slide four fingers between the spokes and press the button of the quick-release axle with your thumb . Then attach the rear wheels to the fittings and the axles will automatically engage.
2. While standing next to the wheelchair, tip it slightly toward you. Press apart the crossbrace with flat of the hand .
3. Pull the bib of the back upholstery upwards and attach the bib to the wheelchair's seat upholstery .
4. Insert the seat cushion .
5. Finally, fold down the footplates .

⚠ CAUTION Accident risk. When remounting, be sure the quick-release axle securely locks into the quick-release axle housing! The rear wheel must not be removable when the button is not being pressed!

⚠ CAUTION Risk of tipping. The rear wheels have been pre-positioned according to your indications. Before using the wheelchair, this pre-adjusted position must be checked by the wheelchair user for function and stability against tipping with the assistance of skilled staff.

4 Transportation

Pull the footrest up and remove the seat cushion.

Grab the seat upholstery at the front and the back in the middle and pull it upwards: The wheelchair will fold .

For more manageable transportation, e.g. in cars, the footrests can be swung outwards or be entirely removed in just a few easy steps.

To move or remove the footrest, pull the swivel handle to the rear .

INFORMATION During reassembly, ensure that the footrest locking device has engaged.

The standard equipment includes removable quick-release axles on the rear wheels which facilitate transportation of your wheelchair.

To remove the rear wheels, press the button on the quick-release axle as described in chapter 3. The rear wheels can then be removed. .

CAUTION **Accident risk.** During assembly, confirm that all bolts engage in the corresponding bore holes and that the quick-release rear wheel axles and casters securely lock into their axle housings. The rear wheels and casters must not be removable without the button being depressed! Check the function of the brakes.

4.1 Transfer

Transfer into the wheelchair is described below. The Motus is designed for persons who normally move around in the wheelchair independently. Therefore, the independent transfer is described. Should you require an attendant, the steps described here are carried out correspondingly but with support.

Maneuver the wheelchairs so that they are positioned at an angle of approx. 45° degrees to each other .

Then activate the wheel locks of the wheelchair you want to transfer into .

Transfer into the second wheelchair. The method illustrated is one possible way to support yourself . Find the most suitable path for you by practicing with an attendant.

Push the footplates of the Motus down again . Now you can use the wheelchair.

5 Options

The wheelchair is designed as a modular assembly system. This means that certain accessories can be adapted to your wheelchair.

We recommend using a seat cushion at all times. The seat upholstery of the Motus has a loop texture so that seat cushions can be attached with hook and loop connection. All seat cushions listed on the order form have the required hook texture.

Furthermore, we will present you a selection of our variants and accessories which improve the use and comfort of your wheelchair.

5.1 Anti-tipper

Prevents the wheelchair from tipping too far backward. The anti-tipper can be adjusted in length without need of tools depending on seat height and seat angle .

5.2 Transport wheels

With the rear wheels removed, the wheelchair can still be used as a rolling base; this makes it possible to manoeuvre through narrow passages .

⚠ CAUTION Accident risk. Please note that once the rear wheels are removed, the wheel locks can no longer be activated!

5.3 Crutch holder

Allows crutches to be attached to the wheelchair.

5.4 Tip-assist

Makes it easier for an attendant to tip up the wheelchair, e.g. to cross a step. The tip-assist must be mounted so that the distance to the ground is at least **4 cm**.

5.5 Plug-on side panels

These side panels can be removed for getting into or out of the wheelchair. To remove the side panels, pull the out of the side panel attachment. Additionally, the armrest can be adjusted in height.

5.6 Desk side panels

The side panels can be swung away to the rear for transfer. For this purpose, press the half-moon push buttons and swing the side panel away to the rear. Remove the side panel in vertical position by pulling it upward. Moreover, the armrest is adjustable in height.

5.7 Plug-on wheel lock lever extension

The extension facilitates the wheel lock lever control to the users with limited hand function.

5.8 Knee lever wheel lock for one-handed operation

This wheel lock is recommended especially for hemiplegics. It can be activated either on the right or left side and ensures the safe locking of both rear wheels with a pull-wire system.

INFORMATION

Please note that the wheel lock lever extension must always be attached when releasing or engaging the wheel lock.

5.9 Rear wheel with drum brake

The drum brakes allow the attendant to comfortably and safely brake the wheelchair. The rear wheels remain removable via the quick-release axle system.

5.10 Tray

Serves as a supporting surface during meals, when working or when playing. The clear material allows visibility of the legs and correction of the sitting posture.

5.11 Elevating footrest

Allows positioning of the leg at different angles. Upon manipulation of the release lever which is integrated into the swivel segment the footrest angle can be adjusted in small increments.

5.12 Motus CS

Version with closed, stiff frame geometry for particularly active driving. To make getting into the wheelchair easier, the rigid footplate can be flipped up.

5.13 Push handles, height-adjustable

Can be adjusted vertically to a more comfortable height for the attendant.


5.14 Wheelbase extension

This option makes the wheelchair particularly stable against tipping and at the same time increases the turning radius.

INFORMATION For transfemoral amputees the long wheelbase setting is indispensable.

5.15 Angle adjustable back

The seat design allows the backrest angle adjustment by 30°. The unlocking of the adjustable back is activated via the pull cable.

5.16 Spoke protector

Protect the fingers from getting caught in the wheel spokes. Also improve the appearance of the wheelchair.

6 Adjustment/Assembly instructions

INFORMATION

The following tools are required for repair and maintenance:

- Allen wrenches, sizes 4 and 5 mm
- Open-end wrenches, sizes: 10 mm, 11 mm, 13 mm, and 24 mm
- Set of screwdrivers
- Torque wrench
- Tyre levers
- Tyre pump

All wheelchair users have expectations regarding wheelchair type, accessories, seating and handling. Your Motus adaptive wheelchair offers many options that permit the wheelchair to be adapted according to your requirements. We recommend consulting with your dealer or therapist to determine which wheelchair settings will work best for you.

6.1 Adjusting seat height and seat inclination

The higher the rear wheel is attached to the frame, the more the seat inclines downward. This has two effects:

- First, the wheelchair has a greater tendency to tilt backwards.
- Secondly, the user sits deeper in the wheelchair and thus has more stability.


A careful and synchronised height adjustment of both the rear wheels and the casters will allow to adapt the seat height to the user's individual requirements.

The Motus offers the possibility to change vertical position of the rear wheel by shifting the rear wheel adapter along the frame. Loosen all four screw connections and remove the screws. Shift the adapter to the desired position and firmly retighten the screws.

To achieve additional rear wheel positions, the rear wheel adapter can be turned by 180°.

INFORMATION When the rear wheel position is changed, the angle of the caster axle in relation to the ground also changes. Ensure that this angle is always approx. 90° wide .

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments. The proper torque for the screw connection of the rear wheel adapter is **8 Nm**. The proper torque for the screw connection of the fitting is **50 Nm**.

6.2 Changing the wheelbase

Shifting the rear wheel backwards lengthens the wheelbase and thus provides for greater stability of the wheelchair. Shifting the rear wheel forwards relieves the load on the casters. This increases the manoeuvrability of the wheelchair. The caster wheels can then also be lifted more easily to clear obstacles such as curbs or steps.

Loosen the screw connections from the fitting and remove the fitting. Shift the fitting to the desired position and firmly retighten the screws.

The distance between rear wheel and side panel can be continuously adjusted by adjusting the counter nut on the outer side of the rear wheel adapter.


When changing the adjustments originally made in the factory, two other component groups of the wheelchair must be readjusted, refer to chapters 6.3 and 6.4.

CAUTION Risk of tipping. Please note: If the rear wheels are in a more forward position and the user's body is not appropriately positioned, the user may tip backward - even on level ground!

CAUTION Risk of tipping. For untrained wheelchair users and when there are extreme adjustments to the rear wheel, an anti-tipper is strongly recommended!

CAUTION Risk of tipping. For transfemoral amputees it is essential to shift the rear wheels backwards!

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments. The proper torque for the screw connection of the rear wheel adapter is **8 Nm**. The proper torque for the screw connection of the fitting is **50 Nm**.

6.3 Setting the wheel camber and toe-in

The rear wheel adapter of the Motus allows for different camber angle settings of the rear wheels. The greater the negative camber angle of the rear wheels, the more manoeuvrable the wheelchair becomes.

The rear wheel adapter itself has a wheel camber of 0° or 2.5°, depending on the arrangement of the wedge-shaped camber washers set at the factory. The wheel camber can be changed by turning round the mounted camber washers. Proceed as follows:

First disassemble the nuts on the rear wheel fitting.




Setting the rear wheel camber to 2.5°: Position the camber washer on the outside of the rear wheel adapter with the thinner side pointing upwards and into the desired engaging position .

Setting the rear wheel camber to 0°: Remove the camber washer from the outside of the rear wheel adapter and place it against the washer on the inside of the adapter such that the thinner side of the washers is next to the thicker side of the other washer .

Further proceeding: Once you have applied the retaining washer from the inside , slightly tighten the nuts on the rear wheel fitting. Check the space between the rear wheel and the side panel. You can adjust the space by rotating the fitting.

Toe-in correction: For toe-in correction, the nuts on the rear wheel fitting remain loosened. Then also loosen the oval head screw on the set plate . Then firmly retighten all screw connections.

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments. Observe the values indicated in chapter 6.1.

6.4 Caster trailing adjustment

After the rear wheels have been mounted in the appropriate position, the caster wheel journal angle must be readjusted. The threaded axle should be as horizontal as possible in relation to the ground to ensure optimal driving behaviour of the wheelchair. The caster adapter allows for continuous angle adjustment.

Remove the cap from caster head bearing with a screwdriver .

Loosen the rear screw connection as shown in fig. 43.

Open the screw connection of the eccentric. Place the bubble level on the caster wheel journal and adjust the eccentric with a screwdriver . Retighten all screw connection and put on the cover caps.

INFORMATION The proper torque for tightening the screws is **8 Nm**.

6.5 Adjusting the wheel lock

To change the position of the rear wheels, first loosen the attachment screws of the clamping adapters of the wheel lock and push the wheel lock forward. Then reposition the rear wheels. The wheel lock should then be remounted such that the distance between the tyres and wheel lock bolts is max. **5 mm** when the wheel lock is not activated .

CAUTION **Accident risk.** The wheel locks are only effective with sufficient air pressure and proper adjustment. The proper air pressure is indicated on the side wall of the tyre. It should be at least **7 bar** . To ensure sufficient braking effect, use only original rear wheels with an approved radial excursion of maximally $\pm 1 \text{ mm}$.

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments.

INFORMATION Tighten the wheel lock attachment with a torque of **10 Nm**.

6.6 Adjusting the brake force - drum brake

In order to achieve optimal braking effect, the brake force can be adjusted by means of the adjustment screw . The brake force is increased by turning the screw counter-clockwise and decreased by turning the screw clockwise.

Turn the adjustment screw counter-clockwise until frictional noise occurs at the turning wheel. Then tighten the adjustment screw until the frictional noise disappears. The wheel should run freely.

After finishing the adjustment, fasten the adjustment screw by tightening the counter nut


INFORMATION Ensure that the brake force of both rear wheels is adjusted equally.

INFORMATION Make sure that the drum brake engages with sufficient force only when the manual brake lever is set to the second ratchet position.

6.7 Angle adjustment of the footplate

Prior to adjusting the footplate angle, loosen the screw connections as indicated in the figure . Subsequently tilt the footplate until the desired angle is achieved.

INFORMATION Make sure to re-tighten the loosened screw connections. The proper torque is **6 Nm**.

6.8 Adjusting the lower leg length

By loosening the set screw, the footplate can be adapted to your lower leg length and to the thickness of your seat cushion .

Make sure when adjusting the footrest height, that the footrest bar is slid into the swivel segment at least **40 mm**. A marking on the footrest bar indicates how far the bar should be slid in.

INFORMATION After having made adjustments, retighten the set screw to a torque of **6 Nm**.

6.9 Adjusting the angle of the elevating footrest

To adjust the footrest angle, turn the release lever up to the stop. While propping up the footrest bar, set the desired angle. Then carefully let the release lever turn back. At the next free position, the footrest will snap into place.

6.10 Adjusting the height of the armrest for desk side panel

The armrest can be adjusted in height by bringing it into various positions without use of any tools. For this purpose, pull the button in the round opening in the side panel upward and shift the armrest into the desired position. The armrest automatically locks in place after releasing the button.

6.11 Adjusting the back height

The back height is adjustable by **10 cm**. The back tube which is integrated into the posterior frame allows for the adjustment of the back height in increments of 2.5 cm.

For this purpose, remove the screws from the rear wheel adapter and set the desired back height.

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments. The proper torque is **8 Nm**.

6.12 Installing the anti-tipper

1. Mount the receiver tube on the frame in the position as shown in the figure .
2. Slide the anti-tipper from below into the mounted receiver tube. First screw the anti-tipper in the utmost position of the receiver tube as shown in the figure .
3. Adjust the anti-tipper using the longitudinal adjustment must always be facing up.
4. Verify that the setting is correct.

As a minimum, the anti-tipper rollers must project fully beyond the largest diameter of the rear wheel. The distance between the anti-tipper rollers and the ground must be max. **50 mm** .

If required: After releasing the screw connections in the receiver tube, the anti-tipper can be adjusted vertically.

If required: After pressing the push-button, the length of the anti-tipper can be readjusted.

Note: Combine the two adjustments. Find the appropriate position with the assistance of a helper.

6.13 Adapting the back upholstery

You can adjust the back upholstery in segments according to your needs.


To do this, remove the back upholstery padding and undo the hook-and-loop fasteners of the straps .

Reposition the straps as desired, attach them again with the hook-and-loop fasteners and pull on the back upholstery padding.

6.14 Changing the back angle

It may be useful to customise the preset angle between the seat and the backrest for therapeutic reasons. The angle can be adjusted between 90° and 120°.

Pull the release cable of the back release until the angle adjustment is released. Bring the backrest into the desired position and release the release cable. Ensure that the locking bolts are safely locked on both sides.

6.15 Push handles, height-adjustable

To adjust the suitable height of the push handles, use the adjustable clamping lever to loosen the clamping screws.

INFORMATION Be sure to snugly re-tighten the clamping levers after all adjustments!

6.16 Attaching the safety belt

Attach the safety belt with the screws provided to the uppermost bore hole of the posterior frame part.

7 Maintenance, cleaning, and disinfection

7.1 Maintenance

Your wheelchair has been provided with the CE sign. The manufacturer herewith guarantees that this medical product as a whole conforms to the requirements of the European Directive for Medical Products 93/42/EEC.

The proper function of the wheelchair, especially of the brakes, should be checked before every use. Safety nuts should be used only once. If they have been loosened several times, they must be replaced.

Should any problems or wear be encountered, please contact your authorised specialist dealer immediately in order to have them corrected. This is particularly important if you notice that the driving characteristics change or the wheelchair becomes less stable. **We also recommend that you have the wheelchair serviced every 12 months by your Ottobock authorised dealer.** The following table lists functions which the user must check at the indicated intervals.

Check	before each use	monthly	quarterly
Function test of the wheel locks and drum brake	X		
Sagging of the seat or back upholstery		X	
Footplate stability		X	
Sight-check of wearing parts		X	
Dirt on bearings		X	
Push ring damage		X	

Check	before each use	monthly	quarterly
Air pressure		X	
Wear on folding mechanism		X	
Spoke tension of the rear wheels			X
Screw connections			X

To ensure smooth operation at all times, the following precautions and cleanings should be carried out. Users may carry these out by themselves.

- Hair and dirt particles generally accumulate between the caster wheel and fork. This can restrict the caster wheels from rotating smoothly. Remove the caster and thoroughly clean the fork and caster using a mild detergent.
- The rear wheels feature a quick-release system. To keep this system operational, ensure that no dirt adheres to the quick-release axle or axle housing. The quick-release axle should also be lightly lubricated regularly with resin-free sewing machine oil.
- If your wheelchair gets wet, towel dry it as soon as possible.
- Do not use your wheelchair in salt water and keep sand or other particles from damaging the wheel bearings.
- Screw connections should be checked frequently, in particular when beginning to use the wheelchair and after any adjustment. If a screw connection becomes loose repeatedly, consult your dealer.

7.2 Cleaning and disinfection

Clean the padding and upholstery with warm water and a washing-up liquid. Remove spots with a sponge or a soft brush. Rinse with clear water and let the components dry.

INFORMATION Do not use any aggressive cleansing agents, solvents, or hard brushes etc.

INFORMATION No wet-cleaning. No machine-wash.

For disinfection, water-based agents should be used. The manufacturer's instructions for use must be observed.

INFORMATION Before disinfecting, the pads and handles must be cleaned.

Plastic parts, frame parts as well as chassis and wheels can be cleaned with a mild detergent. Dry thoroughly.

8 Changing tyres

Changing flat tyres requires only the necessary tools and minimal technical skills. Users may change tyres themselves if they wish. We recommend carrying an emergency repair kit and tyre pump at all times .

Suitable tyre pumps are listed on the order form and are supplied with the product. Alternatively, users can keep a spray can of hardening foam to fill a tyre in case of punctures.

- In case of a flat tyre, carefully remove the tyre from the rim using the appropriate tools.
- Take care not to damage the rim or the inner tube.
- Repair the inner tube according to the directions in the repair kit, or, replace the damaged inner tube with a new one.
- Before re-mounting the tyre, inspect the inner rim surface and the inside of the tyre for any object that may have caused the puncture.
- Be sure that the rim bands are not damaged. They protect the tube from being damaged by the ends of the spokes.

Mounting

- Place the protective rim band over the valve. Then push the valve through the rim. Screw on the valve nut. The rim band can then be easily applied.
- Ensure that all spoke ends are covered.

Tyres

- Beginning behind the valve, place the lower side of the tyre over the edge of the wheel rim. Pump the tube with a small amount of air until it assumes its round shape, and then insert it into the tyre.
- Check to ensure the inner tube is properly positioned all around and is free of folds; if this is not the case, release some of the air. Now you can easily install the upper side of the tyre starting across from the valve, using both hands and working towards the valve.

Inflating

- Check around the tyre on both sides to ensure that the tube is not clamped between the edge of the tyre and the rim.
- Push the valve back slightly and then pull it out again to ensure that the tyre is properly seated in the valve area.
- Inflate the tyre until it can still be pressed easily with your thumb. If the control line on both sides of the tyre is equally distanced to the rim, the tyre sits centrally. If the tyre does not sit centrally, let out some air and realign the tyre. Then inflate the tyre to maximum operating pressure and tighten the dust cap.

9 Technical data

Seat width:	35.5 – 55.5
Seat depth:	36.0 – 54.0
Seat angle:	adjustable from -5° to 15°
Seat height:	42.5 – 49.0 caster fork short
	46.0 – 53.0 caster fork long
Back angle:	adjustable up to -30°
Lower leg length:	34.0 – 54.0
Armrest height:	21.0 – 30.0
Anterior support area of the armrest:	30.0
Overall width:	52.0 – 72.0 with standard rear wheels*
	55.5 – 75.5 with rear wheels with drum brake*
Overall length:	80.5 – 108.5 with 22" rear wheel
	83.0 – 111.0 with 24" rear wheel
Horizontal position of the rear wheel axle:	5.0
Overall width, folded	32.0 – 35.5
Overall length, folded:	80.5 – 111.0
Overall height, folded:	73.0 – 109.0
Weight:	approx. 13 kg
Static stability in forward direction:	10°
Static stability in sideward direction:	10°
Static stability in rearward direction:	0 – 13°, depending on adjustment. Observe safety instructions!
Weight of the heaviest part:	from 8 kg
Load capacity:	125/140 kg

* Applicable to push ring attachment, narrow and a 0° camber of the rear wheels.

Anterior seat height

Caster fork short		
Position		Range of adjustment
Caster diameter	3"	16,5" – 17"
	4"	16,5" – 17,5"
	5"	17,5" – 18"
	5,5"	18" – 18,5"
	6"	18" – 19"
	7"	19"
Caster fork long		
Position		Range of adjustment
Caster diameter	3"	18" – 18,5"
	4"	18" – 19"
	5"	18" – 19,5"
	5,5"	18,5" – 20"
	6"	18,5" – 20"
	7"	19" – 20,5"

Posterior seat height

Rear wheel size	Range of adjustment	Range of adjustment z50	Range of adjustment e-Support
22"	15" – 19,5"	17,5", 18,5", 19,5", 20"	17,5"
24"	15" – 20,5"	18,5", 19,5", 20"	18,5"

INFORMATION The values indicated are measurements which have been theoretically determined. Not all adjustment possibilities can be used for all wheelchair models. Further, the adjustment combinations are limited due to the sturdy frame geometry. We reserve the right to technical alterations and deviations by 5 %.

10 Re-use/Disposal

10.1 Information on re-use

CAUTION

Risk due to improper use. Seat padding which come into direct contact with the skin may cause functional or hygienic risks for another patient using one and the same product. They should be replaced in case of re-use.

The product is suitable for re-use.


Similar to second-hand machines or cars, products that are being re-used are subject to increased strain. Features and functions must not change in a way that could endanger patients or other persons within the product's life cycle.

In cases of re-use, the corresponding product must first be thoroughly cleaned and disinfected. Subsequently, the product must be examined by an authorised specialist to check the condition, wear, and possible damages. Any worn and damaged components as well as components which do not fit or are not suitable for the user must be replaced.

The service manual includes a service schedule, detailed information, and a list of the required tools.

10.2 Information on disposal

If a product is to be disposed of, all components of the product must be disposed of properly in accordance with the applicable environmental protection regulations in the respective country.

[illegible]

Lined area for writing or editing content.

Kundenservice/Customer Service

Europe

Otto Bock HealthCare Deutschland GmbH
Max-Näder-Str. 15 · 37115 Duderstadt · Germany
T +49 5527 848-3455 · F +49 5527 848-1510
healthcare@ottobock.de · www.ottobock.de

Otto Bock Healthcare Products GmbH
Brehmstraße 16 · 1110 Wien · Austria
F +43 1 5267985
service-admin.vienna@ottobock.com · www.ottobock.at

Otto Bock Adria Sarajevo D.O.O.
Omladinskih radnih brigada 5
71000 Sarajevo · Bosnia-Herzegovina
T +387 33 766200 · F +387 33 766201
obadria@bih.net.ba · www.ottobockadria.com.ba

Otto Bock Bulgaria Ltd.
41 Tzar Boris III' Blvd. · 1612 Sofia · Bulgaria
T +359 2 80 57 980 · F +359 2 80 57 982
info@ottobock.bg · www.ottobock.bg

Otto Bock Suisse AG
Luzerner Kantonsspital 10 · 6000 Luzern 16 · Suisse
T +41 41 455 61 71 · F +41 41 455 61 70
suisse@ottobock.com · www.ottobock.ch

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Otto Bock South Africa (Pty) Ltd
Building 3 Thornhill Office Park · 94 Bekker Road
Midrand · Johannesburg · South Africa
T +27 11 312 1255
info-southafrica@ottobock.co.za
www.ottobock.co.za

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Otto Bock HealthCare Canada
5470 Harvester Road
Burlington, Ontario, L7L 5N5, Canada
T +1 289 288-4848 · F +1 289 288-4837
infocanada@ottobock.com · www.ottobock.ca

Otto Bock HealthCare Andina Ltda.
Calle 138 No 53-38 · Bogotá · Colombia
T +57 1 8619988 · F +57 1 8619977
info@ottobock.com.co · www.ottobock.com.co

Otto Bock de Mexico S.A. de C.V.
Prolongación Calle 18 No. 178-A
Col. San Pedro de los Pinos
C.P. 01180 México, D.F. · Mexico
T +52 55 5575 0290 · F +52 55 5575 0234
info@ottobock.com.mx · www.ottobock.com.mx

Otto Bock HealthCare LP
11501 Alterra Parkway Suite 600
Austin, TX 78758 · USA
T +1 800 328 4058 · F +1 800 655 4963
usa.customerservice@ottobockus.com
www.ottobockus.com

Asia/Pacific

Otto Bock Australia Pty. Ltd.
Suite 1.01, Century Corporate Centre
62 Norwest Boulevard
Baulkham Hills NSW 2153 · Australia
T +61 2 8818 2800 · F +61 2 8814 4500
healthcare@ottobock.com.au · www.ottobock.com.au

Beijing Otto Bock Orthopaedic Industries Co., Ltd.
B12E, Universal Business Park
10 Jiuxianqiao Road, Chao Yang District
Beijing, 100015, P.R. China
T +8610 8598 6880 · F +8610 8598 0040
news-service@ottobock.com.cn
www.ottobock.com.cn

Otto Bock Asia Pacific Ltd.
Unit 1004, 10/F, Greenfield Tower, Concordia Plaza
1 Science Museum Road, Tsim Sha Tsui
Kowloon, Hong Kong · China
T +852 2598 9772 · F +852 2598 7886
info@ottobock.com.hk · www.ottobock.com

Otto Bock HealthCare India Pvt. Ltd.
20th Floor, Express Towers
Nariman Point, Mumbai 400 021 · India
T +91 22 2274 5500 / 5501 / 5502
information@indiaottobock.com · www.ottobock.in

Otto Bock Japan K. K.
Yokogawa Building 8F, 4-4-44 Shibaura
Minato-ku, Tokyo, 108-0023 · Japan
T +81 3 3798-2111 · F +81 3 3798-2112
ottobock@ottobock.co.jp · www.ottobock.co.jp

Otto Bock Korea HealthCare Inc.
4F Agaworld Building · 1357-74, Seocho-dong
Seocho-ku, 137-070 Seoul · Korea
T +82 2 577-3831 · F +82 2 577-3828
info@ottobockkorea.com · www.ottobockkorea.com

Otto Bock South East Asia Co., Ltd.
1741 Phaholyothin Road
Kwaeng Chatuchark · Khet Chatuchark
Bangkok 10900 · Thailand
T +66 2 930 3030 · F +66 2 930 3311
obsea@ottobock.co.th · www.ottobock.co.th

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Ottobock has a certified Quality Management System in accordance with ISO 13485.

© Ottobock · 647G232=D/GB-17-1510

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Otto Bock Mobility Solutions GmbH