**NXT™ BACK SUPPORTS**

*DOSSIERS NXT<sup>MD</sup>*

**INSTRUCTION MANUAL**  
*MANUEL D'INSTRUCTION*

## IMPORTANT:

The installation and set up instructions should be carried out by a qualified technician.

It is important to refer to the wheelchair owner's manual for additional safety information related to the safe operation of the wheelchair.

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Complies with:

- ISO 16840-4 Wheelchair Seating - part 4: Seating Systems for use in motor vehicles.
- RESNA WC-4:2012 Section 20: Seating Systems for use in motor vehicles.

It is recommended that a headrest be used when the back rest is used for transport in a motor vehicle.

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Advena Limited | Tower Business Centre, 2nd Flr., Tower Street Swatar, BKR 4013 Malta

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

## IMPORTANT INFORMATION: READ BEFORE INSTALLATION

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**The Quickfit mounting hardware can accommodate different size back canes.**

3/4"/19mm, 7/8"/22mm and 1"/25mm back canes can be accommodated without inserts.

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## IMPORTANT INFORMATION: FOLLOW BEFORE INSTALLATION

**WEIGHT CAPACITY:** The weight capacity of the **NXT** Back is 300lbs/136kg – 400lbs/182kg on 22"/55cm. 23"/58cm – 29"/73cm is 550lbs/250kg with 4 point hardware.


Please note that the Quickfit mounting hardware on the **NXT** Back are pre-assembled with the cane clamp in the rear most mounting position to maximize seat depth of the wheelchair. See *Fig. 1* below.

The configuration can be changed by orienting the Quickfit cane clamp forward, providing a forward mounting position. See *Fig. 2* below.

**Note:** The Quickfit mounting hardware provides 3"/7.6cm of depth adjustment, in front of and behind the back canes.

To set the Quickfit mounting hardware in the forward most position:

- Remove the cane clamp from the mounting plate
- Rotate the cane clamp 180° and re-attach it to the mounting plate
- Install clamp as shown in *Fig. 2*

**Please Note:** Setting the desired position of the **NXT** back may compromise the stability of the wheelchair.

Be aware that the forward and aft position for the backrest will change the user's centre of gravity within the wheelchair. A back set in the most rearward position can significantly reduce the rearward stability of the wheelchair. In contrast, a more forward back position may reduce the wheelchair's forward stability.

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

# NXT BACK INSTALLATION INSTRUCTIONS

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**WARNING: DO NOT** substitute hardware provided. Use only high strength hardware provided.

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**Fig. 4**

**Cane clamps must be installed at the same height and parallel to each other**

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
**Fig. 5**

**Loosen 4 cap screws to adjust back height**

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7. To install the backrest on the mounting plates:

**Step 1:** Insert the top of the backrest slot onto the top adjustment pins

**Step 2:** Rotate the lower section of the backrest rearward to lock it into the locking latch to secure the backrest in place.

**Step 3:** Ensure that the latches lock completely by pushing down on the release handle or latches if necessary.

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

ADJUST MOUNTING WIDTH  
AJUSTER LA LARGEUR  
DE MONTAGE

MAXIMUM WIDTH  
LARGEUR MAXIMALE

8. To ensure the backrest is properly fitted you may need to adjust the mounting hardware to match the back cane width.

**Step 1:** Loosen the screws to adjust the mounting width and the one-hand release.

**Step 2:** Make sure the backrest is centered horizontally and vertically between the back canes and the locking latch on the bottom is on the inside of the post

**Step 3:** When satisfied, tighten screws.

**Step 4:** The width adjustable mounting hardware allows the backrest to fit 3 sizes of wheelchair, ie a 16"/40cm wide back rest can fit a 16"/40cm, a 17"/43cm or 18"/45cm wheelchair.

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
## CONTINUED: NXT BACK INSTALLATION INSTRUCTIONS

9. Make sure that the overall fit and alignment of the back is correct, if not make adjustments as required.
10. Secure the cane clamp screws onto each back cane.
11. Setting the depth and recline settings it is best done with the client seated in the wheelchair. While the client is seated and supported in the wheelchair, loosen the depth screws on the mounting plate and set to the desired position. Tighten all screws to lock the backrest in position. The backrest can be mounted asymmetrically, so be sure that it is symmetrical if desired.


### QUICK RELEASE MOUNT

TUBE CLAMPING  
SCREWS.



TO LOCK LATCH MOVE THIS  
SCREW TO THE ADJACENT  
HOLE, THIS WILL PREVENT  
THE LATCH FROM LIFTING.  
REPLACE SCREW FOR  
NORMAL OPERATION.

NOTE: EXTRA HOLE IS NOT  
THREADED BUT THE SCREW  
WILL MAKE ITS OWN THREAD.

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
## REMOVING THE BACKREST *RETRAIT DU DOSSIER*

To remove the backrest from the wheelchair

**Step 1:** Lift on the one hand release lever to release the latches, push the bottom of the back forward.

**Step 2:** Lift the backrest off the upper pins to remove completely.

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1. Remove the upholstery from the backrest.
2. For light cleaning of the backrest fabric may be wiped down with a damp cloth.
3. Follow laundering instructions provided with the cover.

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**Undo velcro closure & remove foam pad before washing outer cover**

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1. The Pelvic-lumbar support pads can be uniquely positioned to provide optimal pelvic positioning and lumbar support. The supports come uninstalled, allowing them to be specifically fit to the person at the time of the back installation. Fit as shown between fabric and inner shell of backrest. Hook Velcro on the front of the pelvic pad attached to the Velcro compatible fabric of the cover.
2. Where additional lateral support is required, optional lateral supports are. Fit pads between the fabric and the shell of the backrest.
3. After positioning support pads, seal velcro closure, reposition cover on the backrest shell, attach the upholstery.

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## INSTALLING HEAD/NECK SUPPORTS

### INSTALLATION DE SUPPORTS DE TÊTE/COU

Attachment slots are stamped into the back shell for attaching the NXT Multifit Head Support or after market head/neck supports.

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# Armadillo™ ARRIÈRE

## 1. Setting Back Height and Angle

The Armadillo back installs, like all NXT backs, using Universal Clamps attached to the back canes. The back height and angle of the middle shell segment are adjusted using these clamps as shown in the NXT INSTALLATION/USERS MANUAL.

NOTE: When installing a back, always ensure the latches are closed (down) on the lower pins.

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## 2. Adjusting Back to match Chair Width


The Armadillo back is available in sizes to fit chair widths 14"-24". To adjust the mounting hook width, loosen screws A and B and slide the mounting hooks to the correct width and tighten all screws. Screw A fits into a fixed threaded hole in the middle shell segment and screws B fit into sliding nuts which also hold the upper laterals in place if installed. The hooks are infinitely adjustable between the minimum and maximum width and are normally placed symmetrically on the back but may be asymmetrical if required.

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NOTE: Screws A & B on the left and right side of the back must be used to make width adjustments.

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
### 3. Head Rest Mounting

The Upper and Middle segments of the Armadillo shell have slots to accept mounting hardware for the NXT MultiFit Head Rest. Choose the location which provided the best positioning for the head rest.

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### 4. Lateral Installation and Adjustment

Upper laterals are held in place by screws B and the lower laterals by screws E. These screws thread into sliding nuts in slots on the laterals allowing the laterals to be infinitely adjustable between their most inner and outer positions only limited by the postural requirements of the occupant of the wheelchair and any structural components of the wheelchair. Loosening screws B or E allowing the laterals to slide will not affect the positioning of the mounting hooks.

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
## 5. Angle Adjustments

The Upper and Lower segments of the Armadillo shell are individually angle adjustable with respect to the Middle segment. Loosening screws G (on both sides of hinge joint) allows the Upper segment to be tilted forward up to 40 degrees or backward up to 20 degrees. Loosening screws F (on both sides of hinge joint) allows the Lower segment to be tilted forward or backward up to 30 degrees. Tightening the screws locks the segment preventing a change in angle. Ensure all 8 screws are tightened.

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
# üFit™ CUSTOMIZABLE CONTOUR SUPPORT BACK

## 1. Customization

üCube™ sets can be removed completely, or cut for greater customization. To remove/cut, first unzip both zipper pulls on the cover and fold the front of the cover down. Remove the center foam overlay, to expose the üCube™. To remove the üCube™, grasp the cube set by the darker coloured closed cell foam base and gently pull up – the cube will release from the back shell.

If cutting, cut the üCube™ once it has been removed, and re-install on back shell by aligning the studs (white pins on the underside of the üCube™) with the holes and pushing back into place until an audible click is heard. Once all desired customization is complete, re-install the center foam overlay and zip-up the cover.

## 2. Installing Head Support

To install a head support on the chair, you must remove the cover fabric overlays the slots which the head support is installed in. Conveniently, there is a punched template to follow. Simply connect the ends of the existing cuts with scissors or a knife to remove the rectangle of fabric.


Next remove the back cover and then the üCube™ that are in the area around the mounting slots. Install the mounting hardware and adjust the angle as desired by rotating the bracket slightly either clockwise or counter clockwise.

Reassemble the üCube™ and cover ensuring the mounting hardware fits through the hole cut into the cover.

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The **NXT™** Back Support has a Limited Warranty for 24 months from the date of delivery to the original purchaser. If any defect in material or workmanship is found, the manufacturer will repair or replace the component, at our discretion. This limited warranty does not cover daily wear and tear or damage that is a result of mishandling, misuse, neglect, or not following the care and maintenance instructions. Your back support cover is under warranty for 90 days. The warranty includes damage due to workmanship, material defects or errors caused by the manufacturer. This warranty does not apply to cigarette burns, damage to the cover by sharp objects that may cause tears, or damage as a result of not following the washing instructions. Claims and repairs should be processed through your nearest authorized dealer. A copy of the invoice issued to the purchaser of the product may be required prior to processing a warranty claim.

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