



CE

Kimba Neo

EN Instructions for Use (User)	3
---	---

Table of contents

EN

1	Foreword	5
2	Product Description	5
2.1	Function	5
2.2	Product Overview	6
3	Safety	6
3.1	Explanation of Warning Symbols	6
3.2	General Safety Instructions	7
3.3	Safety Instructions for Use	7
3.3.1	Further information	8
3.4	Nameplate and warning labels	8
3.4.1	Nameplate	8
4	Delivery	8
4.1	Scope of Delivery	8
4.2	Storage	8
5	Preparation for Use	8
5.1	Assembly	8
6	Use	9
6.1	Information on use	9
6.2	Basic functions of the Kimba Neo seating unit	10
6.2.1	Removing and mounting the seating unit	10
6.2.2	Rotating the seating unit	11
6.2.3	Adjusting the tilt-limiting mechanism after rotating the seat	11
6.2.4	Adjusting the back angle	11
6.2.5	Adjusting the seat angle	12
6.2.6	Adjusting the footrest	13
6.2.6.1	Adjusting the knee angle	13
6.2.6.2	Folding up the footrest	13
6.3	Additional options	14
6.3.1	Headrests	14
6.3.1.1	Headrest supports	14
6.3.1.2	Integrated headrests	14
6.3.2	Safety belts	14
6.3.2.1	Lap belt	15
6.3.2.2	Four-point belt	15
6.3.2.3	Five-point belt	16
6.3.3	Positioning aids	17
6.3.3.1	Thoracic supports	17
6.3.3.2	Abductor	18
6.3.3.3	Pelvic band	18
6.3.3.4	Chest and shoulder support	18
6.3.3.5	Fixation vest	19
6.3.3.6	Groin strap	20
6.3.3.7	Neck support	21
6.3.3.8	Lumbar padding	21
6.3.3.9	Hook-and-loop foot straps	22
6.3.4	Tray	22
6.3.5	Canopy	23
6.3.6	Rain canopy	24
6.3.7	Rain cover	25
6.3.8	Summer/winter slip sack	25
6.3.9	Seat depth reducer	26
6.3.10	Incontinence cover	27
6.3.11	Armrests	28

6.3.12	Grab rail	29
6.4	Disassembly/transport	29
6.5	Use in a wheelchair accessible vehicle	30
6.5.1	General Safety Instructions	30
6.5.2	Permitted use	30
6.5.3	Important information.....	30
6.6	Care.....	31
6.6.1	Detaching/attaching the padding for cleaning	31
6.6.1.1	Detaching/attaching the seat padding.....	31
6.6.2	Detaching/attaching the back padding.....	31
6.6.3	Detaching/attaching the covers.....	32
6.6.4	Detaching/attaching the grab rail cover (option)	33
6.7	Cleaning	33
7	Maintenance/Repair	34
7.1	Maintenance.....	34
7.1.1	Maintenance intervals	34
8	Disposal	34
8.1	Disposal Information.....	34
9	Legal Information	34
9.1	Service Life	34
9.2	Liability	35
9.3	CE Conformity.....	35
9.4	Warranty	35
9.5	Trademarks	35
10	Technical data	35

1 Foreword

INFORMATION

Last update: 2014-05-09

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

INFORMATION

- ▶ You can request this document as a PDF file from the Customer Care Center (CCC) at oa@ottobock.com or from the manufacturer's service department (see inside back cover or back page for addresses).
- ▶ It is possible to increase the display size of the PDF document.
- ▶ For further questions about the instructions for use, please contact the authorised personnel who issued the product to you.

You have purchased a high-quality product which can be put to versatile, daily use at home and outdoors.

Please read these instructions for use before you start using the product. Familiarise yourself with the handling, functions and use of the product before using it. Observe all safety information and instructions in order to prevent the possibility of injury. Store these instructions for use in a safe place for future reference.

Please note the following in particular:

- Attendants should also be instructed in the use of the product by qualified personnel with the help of these instructions for use.
- The product has been adapted to the needs of the user. Further changes may be made only by qualified personnel. We recommend checking the product settings regularly in order to assure an optimum fit over the long term. For growing children and youths in particular, fitting should be performed every six months.
- In case of questions or problems, please consult the qualified personnel that adapted the product or the manufacturer's service (see inside back cover or back page for addresses).
- The product may be combined only with the options listed here. The manufacturer assumes no liability for combinations with third-party medical devices and/or accessories not included in the modular system. Please also observe the information in the section "Liability".
- Service and repairs to the product may only be carried out by qualified personnel. If you have any problems, please contact your specialist dealer. This ensures that any necessary repairs will be made exclusively with Ottobock spare parts.
- Your product may differ from the models shown. In particular, not all the options described in these instructions for use will be installed on your product.
- The manufacturer reserves the right to make technical changes to the model described in these instructions for use.

2 Product Description

2.1 Function

This seating unit for a rehab buggy is intended solely for transporting children with walking disabilities or who are unable to walk.

This seating unit for a rehab buggy has only been approved for transporting one child aged approx. 1 to 10 years (depending on the size and weight of the child). It is operated by an attendant (conveyance by another person).

This seating unit for a rehab buggy can be used indoors and outdoors. It is fully adjustable to provide an optimised, individual fitting.

2.2 Product Overview



Kimba Neo seating unit

1	Angle adjustable backrest
2	Removable seat with lateral supports
3	Adjustable footrest
4	Back angle adjustment
5	Seat angle adjustment

3 Safety

3.1 Explanation of Warning Symbols

⚠ WARNING	Warning regarding possible serious risks of accident or injury.
⚠ CAUTION	Warning regarding possible risks of accident or injury.
ℹ NOTICE	Warnings regarding possible technical damage.

3.2 General Safety Instructions

CAUTION

Lack of familiarisation with instructions for use

Injuries due to user error

- ▶ Read the instructions for use before using the product.
- ▶ Familiarise yourself with the handling, functions and use of the product before using it. The safety of the user may be compromised if you do not follow these instructions.

3.3 Safety Instructions for Use

INFORMATION

The seating unit is intended for use with the Kimba Neo outdoor mobility base. Observe safety information in the 647G891 Instructions for Use.

Hazards during preparation for use

WARNING

Independent modification of settings

Serious injuries to the user due to unallowable changes to the product

- ▶ Do not modify the settings established by the qualified personnel. Only the settings described in the section "Usage" in these instructions for use may be adjusted independently.
- ▶ In case of problems with the settings (unsatisfactory seating position, front wheels wobble, etc.) please contact the qualified personnel who adjusted your product.
- ▶ Discuss all changes to the settings with the qualified personnel/therapist in advance in order to minimise health hazards and avoid putting the therapy results at risk.

WARNING

Improper handling of packaging materials

Risk of suffocation due to neglect of the duty to supervise

- ▶ Packaging materials must be kept out of the reach of children.

Hazard in case of broken skin

CAUTION

Prolonged contact with broken skin

Reddening of the skin or pressure points due to contamination with germs or incorrectly adapted seat cushion

- ▶ Check before using the product that the skin is not broken in areas subject to prolonged pressure (e.g., buttocks, back and rear side of the thighs).
- ▶ In case of problems, please contact the qualified personnel who adjusted the product.
- ▶ No liability is assumed for injuries caused by the use of the product with broken skin.

Hazards due to fire, heat and cold

CAUTION

Seat cushion or back upholstery may ignite

Burns due to user error

- ▶ The seat cushion and back upholstery are flame retardant but nevertheless flammable. Therefore the utmost caution must be exercised in the vicinity of an open flame.
- ▶ Keep away from all ignition sources, especially lit cigarettes.

CAUTION

Extreme temperatures

Hypothermia or burns through contact with components, failure of components

- ▶ Do not expose the product to any extreme temperatures (e.g. direct sunlight, sauna, extreme cold).

NOTICE**Use under incorrect environmental conditions**


Damage to the product

- Do not use the product in salt water.
- Avoid damage to the product due to sand or other particles of dirt.

3.3.1 Further information**INFORMATION**

Even in the event of compliance with all applicable guidelines and standards, alarm systems (e.g. in department stores) may respond to your product. Should this happen, remove your product from the area where the alarm was triggered.

3.4 Nameplate and warning labels**3.4.1 Nameplate**

Label	Meaning	
	A	Type designation
	B	Manufacturer article number
	C	Maximum load capacity (see section "Technical Data")
	D	Manufacturer information/address/country of origin
	E	Serial number
	F	European Article Number/Global Trade Item Number
	G	Read the instructions for use before using the product.
	H	CE marking – product safety in accordance with EU directives
<p>The nameplate "Seating unit" is located on the backrest of the seat.</p>		

4 Delivery**4.1 Scope of Delivery**

The scope of delivery includes:

- Kimba Neo seating unit
- Options according to the order
- Adjustable footrest
- Instructions for use (user)

4.2 Storage

The seating unit must be stored in a dry place. An ambient temperature between **-10° C** and **+40° C** must be maintained. Direct exposure to sunlight or UV radiation causes the product to age more quickly. This causes the material surface to harden.

5 Preparation for Use**5.1 Assembly****⚠ CAUTION****Exposed pinch points**

Crushing, pinching due to incorrect handling

- When folding out the backrest, only grip by the specified components.


⚠ CAUTION**Improper assembly**

Tipping, falling over of the user due to components coming off

- When attaching the seating unit, make sure you hear it engage into the mobility base.
- Secure the seating unit using the seat lock.

Mounting the seat

- 1) Place the seat on the rear tube of the seat adapter at an angle of approximately 45° (see Fig. 3).
- 2) Press the front edge of the seat onto the mobility base until the gripper jaw audibly snaps on to the front tube of the seat adapter (see Fig. 4).
- 3) Slide the seat lock to the left to engage it (see Fig. 5).
- 4) Check that the seat and the mobility base are securely connected. You can do this, for example, by pulling back the backrest.



6 Use

6.1 Information on use

⚠ WARNING

Independent modification of settings

Serious injuries to the user due to unallowable changes to the product

- ▶ Do not modify the settings established by the qualified personnel. Only the settings described in the section "Usage" in these instructions for use may be adjusted independently.
- ▶ In case of problems with the settings (unsatisfactory seating position, front wheels wobble, etc.) please contact the qualified personnel who adjusted your product.
- ▶ Discuss all changes to the settings with the qualified personnel/therapist in advance in order to minimise health hazards and avoid putting the therapy results at risk.

⚠ WARNING

Impermissible use

Tipping over, rolling over of the user due to failure to observe specifications

- ▶ Using a rehab buggy beyond normal conditions may be dangerous.
- ▶ Do not exceed the maximum load capacity (see the section "Technical Data").

⚠ CAUTION**Failure to check for proper functioning**

Falling, tipping over of the user due to lack of maintenance

- ▶ In the interest of the user's safety, check the product for correct function before each use.
- ▶ Please also read the instructions in the section "Maintenance" > "Maintenance intervals".

INFORMATION

Please also observe the safety instructions in the section "Safety" > "Safety instructions for use".

In the interest of user safety, the following inspection must be conducted prior to use:

- 1) Verify firm locking of the seat and back angle adjustment mechanism.
- 2) Verify the functionality of all belt buckles on the safety belts and/or positioning aids.

6.2 Basic functions of the Kimba Neo seating unit

INFORMATION

Adjusting the seat and back angle is not possible in combination with the Kimba Cross and Kimba Inline mobility bases.

6.2.1 Removing and mounting the seating unit

⚠ CAUTION**Incorrect mounting of the seating unit**


Falling out of the user if the seating unit is not secured

- ▶ Secure the seating unit using the seat lock.
- ▶ When attaching the seating unit, make sure you hear it engage into the mobility base.

To facilitate transportation of the rehab buggy in a passenger vehicle or to divide the weight into two units for lifting, the seating unit can be easily detached from and reattached to the outdoor mobility base.

Detaching the seating unit

- 1) Slide the seat lock to the right to open it (see Fig. 6).
- 2) Stand next to the seating unit and hold the seat by the top of the backrest with one hand.
- 3) Pull the release lever located under the seat up to the edge of the seat (see Fig. 7).
 - The clutch will now release the seat adapter.
- 4) Tip the seat back at an angle of approx. 45° and then lift it up (see Fig. 3).





Mounting the seating unit

See the section "Assembly: Mounting the Seating Unit".

Rotating and mounting the seating unit

If the seat is rotated before mounting, the tilt-limiting mechanism then has to be adjusted for safety reasons. See the section "Assembly: Mounting the Seating Unit" and the two sections that follow.

6.2.2 Rotating the seating unit

CAUTION

Non-activated tilt-limiting mechanism

Falling out of the user due to user error

- Before rotating the seating unit, always adjust the tilt-limiting mechanism.

The seating unit can be mounted onto the mobility base in two directions (facing forward or backward) (see Fig. 9, see Fig. 10).

To remount the seating unit after rotating it 180°, the tilt limit must be adjusted (see the following section).



6.2.3 Adjusting the tilt-limiting mechanism after rotating the seat

Adjusting the tilt-limiting mechanism after rotating the seat is described in the mobility base instructions for use.

6.2.4 Adjusting the back angle

CAUTION

Change in centre of gravity following back angle adjustment

Falling over, falling out of the user due to adjustment errors

- After each modification of the back angle, check the tilt stability. To do so, hold the product with both hands by the push handle/push bar.

The backrest can be adjusted as needed in a range from 90° to 180°:

- 1) Twist the lower knob counterclockwise to release the back angle adjustment (see Fig. 11).

- 2) Move the backrest forward/backward to the desired position. The maximum position is lying flat (see Fig. 12)
INFORMATION: Please note that in the maximum position (lying flat), the footrest must be adjusted so that the user's feet can comfortably reach it.
- 3) Turn the knob clockwise to fixate the back angle.



6.2.5 Adjusting the seat angle

CAUTION

Incorrect use of the seat tilt feature

Tipping over, falling from the seating shell/seating system due to user error

- ▶ Practise the operation of the seat angle adjustment mechanism without the seated user at first.
- ▶ Only operate the seat angle adjustment mechanism on a level, firm surface.
- ▶ Only negotiate slopes and obstacles in the horizontal neutral position of the seat. Tilting the seat angle back slightly when driving downhill is recommended.
- ▶ Only negotiate slopes where the holding forces remain manageable.
- ▶ When adjusting the seat tilt feature, secure the user from falling out to the front or rear. In order to do so, firmly hold on to the push handle/push bar while operating the seat angle adjustment mechanism.
- ▶ Do not reach into the adjustment mechanism when changing the seat tilt settings.

INFORMATION

The function of the upper knob for the adjustment mechanism is not available on all outdoor mobility bases. The procedure described below refers to the functions of the Kimba Neo seating unit when it is ordered and delivered together with the Kimba Neo mobility base. Deviating from this, the seat angle adjustment release can be mounted on the outdoor mobility base. Please observe the information given in the outdoor mobility base instructions for use.

The seat can be adjusted as needed within a range of **-10° to +35°**:

- 1) Twist the upper knob clockwise to release the adjustment mechanism (see Fig. 13).
NOTICE! Damage to the seat adjustment. Turn the knob with normal manual force to the stop, until you feel a clear resistance. Excessive turning of the knob will damage the adjustment mechanism.
- 2) Set the seat angle to the desired position (see Fig. 14).
CAUTION! With one hand, secure the occupant against falling out (e.g. by holding the upper body).
- 3) Release the knob.
 → The adjustment mechanism will lock automatically.
- 4) Press the backrest gently to the front/back to make sure that the seat angle is securely set. If this is not the case, the position of the upper knob must be checked.



6.2.6 Adjusting the footrest

⚠ CAUTION

Getting into buggy when the footrest is folded down is dangerous

Falling, tipping over of the user due to user error

- Fold up the footrest before getting in.

⚠ CAUTION

Improper lifting by attendants

Tipping over, falling of the user due to lifting on removable components

- Attendants may lift the buggy only by welded or permanently attached components. It is particularly dangerous to lift the buggy by the footrest, push bar, lap bar, or armrests.

6.2.6.1 Adjusting the knee angle

The knee angle of the footrest can be adjusted within a range of 0° to +90° as needed:

- 1) Open the adjustment lever (see Fig. 15).
 - 2) Adjust the knee angle.
 - 3) Firmly close the adjustment lever.
 - 4) Prior to use, check the secure position of the footrest. To do so, press the footrest down.
- The knee angle must not change.



6.2.6.2 Folding up the footrest

The footrest can be folded up for transport.


Standard footrest

- 1) Hold the edge of the footrest (see Fig. 16, item 1) and fold it up (see Fig. 16, item 2).
- 2) The footrest can be unfolded again by pressing down on it.

Angle-adjustable footrest (optional)

- 1) Pull the release strap (see Fig. 17).

- 2) Fold up the footrest.
- 3) The footrest can be unfolded again by pressing down on it. It locks into place automatically.



6.3 Additional options

6.3.1 Headrests

The seat of the rehab buggy can be equipped with headrests that provide additional support for certain indications.

6.3.1.1 Headrest supports



In this option (see Fig. 18), the free space in the head region is adjusted by shifting the supports in the vertical slots (large adjustment range).

The option has been preset by qualified personnel.

6.3.1.2 Integrated headrests

With this option (see Fig. 19), the free space in the head region is adjusted through the insertion of foam pads into the lateral supports (adjustment range: **30 mm** in **1 mm** increments).

The option has been preset by qualified personnel.



6.3.2 Safety belts

⚠ WARNING

Use of the belt system as a restraint system in a wheelchair accessible vehicle is forbidden

Falling out of the product, serious injuries due to improper handling of the product

- Under no circumstances may the belt system be used as part of a restraint system when transporting the user in a vehicle for transporting the disabled.
- Note that the belt system is only intended as additional stabilisation for the user sitting in the rehab buggy.

⚠ WARNING**Incorrect application of the belt system/positioning aid**

Strangulation, suffocation due to incorrectly positioned safety belts/positioning aids

- ▶ Follow the instructions for proper application so that the user cannot slip into a position that would endanger him or her.
- ▶ Note that the user's back must be up against the back upholstery when securing the belt system/positioning aid and during use of the product.
- ▶ Do not use safety belts or positioning aids if the user is able to sit without them.

The product can be equipped with different types of safety belts.

The option has been installed and preset by qualified personnel. If the option has been adjusted too tightly or cannot be engaged, the safety belt can be subsequently adjusted.

6.3.2.1 Lap belt

The seating unit can be equipped with a lap belt. It provides additional support for the pelvis and keeps the user from sliding out of the seat (see Fig. 20).

Adjusting the belt

- 1) Pull the belt out of the belt restraint (see Fig. 21, item 1).
- 2) Close the belt buckle firmly (see Fig. 21, item 2) and pull the belt in the desired direction.
Make sure that the belt extends at least **20 mm** beyond the buckle so that the belt restraint can be pushed over the loose end.
- 3) Thread the belt back into the belt restraint (see Fig. 21, item 1).

Using the belt

- 1) To open the belt buckle, push the "press" button (see Fig. 20, item 1).
- 2) Sit the user into the seat.
- 3) To close the belt, slide the metal tab into the buckle until you hear it click.

**6.3.2.2 Four-point belt**

The seating unit can be equipped with a four-point belt (see Fig. 23). It additionally stabilises the pelvis and shoulders and prevents the user from sliding out of the seat. The buckle makes it easy to get in and out.



Adjusting the belt

- 1) The setting of the lap belt part corresponds to the setting of the lap belt (see Page 15)
INFORMATION: If the lap belt was run under the padding, the closures of the backrest pad must be opened.
 - 2) Adjust the shoulder straps using the clip buckles at the base of the backrest (see Fig. 27).
- The four-point belt is adjusted correctly when the closed belt buckle lies against the child's abdomen (see Fig. 24).

Using the belt

- 1) Press the release button to open the belt (see Fig. 23, item 1).
- 2) Move the shoulder/lap belt straps to the side.
- 3) Sit the user into the seat.

- 4) Place the belt over his/her shoulders (see section "Five-point belt").
- 5) Make sure the two belt tabs click into the buckle (see Fig. 23, item 2).




6.3.2.3 Five-point belt

The seating unit can be equipped with a five-point belt (see Fig. 24). It additionally stabilises the pelvis and shoulders and prevents the user from sliding out of the seat. The buckle makes it easy to get in and out.


Adjusting the belt

- 1) Lengthen/shorten the crotch strap at the belt buckle (see Fig. 25).
- 2) Adjust the lap-belt part behind the seat using the adjustment buckle (see Fig. 26).
INFORMATION: If the lap belt was run under the padding, the closures of the backrest pad must be opened.
- 3) Adjust the shoulder straps using the clip buckles at the base of the backrest (see Fig. 27, item 1 and item 2).
 → The five-point belt is adjusted correctly when the closed belt buckle lies against the child's abdomen (see Fig. 24).



Using the belt

- 1) Press the red button to unfasten the belt (see Fig. 24, item 1).
- 2) Move the shoulder/lap belt straps to the side.
- 3) Sit the user into the seat.
- 4) Place the belt over his/her shoulders.
- 5) Make sure the two belt tabs click into the buckle (see Fig. 29).



6.3.3 Positioning aids

⚠ WARNING

Use of the belt system as a restraint system in a wheelchair accessible vehicle is forbidden

Falling out of the product, serious injuries due to improper handling of the product

- Under no circumstances may the belt system be used as part of a restraint system when transporting the user in a vehicle for transporting the disabled.
- Note that the belt system is only intended as additional stabilisation for the user sitting in the rehab buggy.

⚠ WARNING

Incorrect application of the belt system/positioning aid

Strangulation, suffocation due to incorrectly positioned safety belts/positioning aids

- Follow the instructions for proper application so that the user cannot slip into a position that would endanger him or her.
- Note that the user's back must be up against the back upholstery when securing the belt system/positioning aid and during use of the product.
- Do not use safety belts or positioning aids if the user is able to sit without them.

The option has been installed and preset by qualified personnel. If the option has been adjusted too tightly or cannot be engaged, the safety belt can be subsequently adjusted.

The seating unit can be equipped with positioning systems.

6.3.3.1 Thoracic supports

The seating unit can be equipped with thoracic supports that provide additional support to the upper body for certain indications (see Fig. 30).

The option has been preset by qualified personnel.



6.3.3.2 Abductor

The seating unit can be equipped with an abductor (see Fig. 31) that ensures that the legs are abducted (spread to the sides) from the centre of the body for certain indications.

The option has been preset by qualified personnel.

6.3.3.3 Pelvic band

INFORMATION

For further information, please see the instructions for use included with the product.

The seating unit can be equipped with a pelvic band (see Fig. 32) that secures the pelvis in the seat for certain indications.


Opening the pelvic band

- 1) Press on the back end of a buckle with the thumb/forefinger (see Fig. 33).
- 2) Pull out the plastic belt while holding the buckle down.

Fastening the pelvic band

- 1) Sit the user into the seat. Make sure that the user's buttocks are up against the backrest padding.
- 2) Press on the back end of a buckle with your forefinger (see Fig. 33).
- 3) Insert the plastic belt while holding the buckle down.

CAUTION! If the pelvic band is too loose, the pelvis can slide forward and lead to incorrect positioning or cause the user to slide out of the buggy. Make sure the user is correctly positioned in the seat.



6.3.3.4 Chest and shoulder support



INFORMATION

For further information, please see the instructions for use included with the product.

The seating unit can be equipped with a chest and shoulder support that secures the trunk in the seat for certain indications (see Fig. 34). It is used to achieve upright symmetrical position of the pelvis and trunk.

Adjusting the chest and shoulder support



The adjustment of the final belt tension takes place when putting it on by adjusting the belt length at the two pairs of clip buckles on the back of the seat (see Fig. 35, item 1).



Putting on the chest and shoulder support

- 1) If needed: Before putting the user into the seat, unfasten the clip buckles at the back of the seat (see Fig. 35, item 1) and pull the belt straps of the chest and shoulder support forward somewhat.
- 2) Unfasten the bottom belt straps of the chest and shoulder support at the clip buckles (see Fig. 36/see Fig. 37, item 1), pull the belt straps through and move the chest and shoulder support back.
- 3) Sit the user into the seat. Make sure that the user's back is up against the backrest padding.
- 4) Pull the chest and shoulder support over the user's head onto his or her chest.
- 5) Put the lower belt straps of the chest and shoulder support into the clip buckles (see Fig. 36/see Fig. 37, item 1).
- 6) If needed: Tighten the belt straps in the clip buckles at the back of the seat and fasten them (see Fig. 35, item 1). Make sure that the chest and shoulder support lies symmetrically on the pelvis.

CAUTION! If the chest and shoulder support is too loose, the user can slide forward. This may lead to incorrect positioning or cause the user to slide out of the buggy. Please ensure that the child is positioned correctly in the seat. The lower support belts must lie firmly against the spinae.



6.3.3.5 Fixation vest



INFORMATION

For further information, please see the instructions for use included with the product.

The seating unit can be equipped with a fixation vest that secures the trunk in the seat for certain indications (see Fig. 38). The fixation vest does not provide any support for the pelvis.


Adjusting the fixation vest

The adjustment of the final belt tension takes place when putting it on by adjusting the belt length at the six clip buckles on the back of the seat (see Fig. 39).



Putting on the fixation vest

- 1) Unfasten the hook-and-loop tabs over the zipper (see Fig. 38, item 1).
- 2) Unzip the zipper (see Fig. 38, item 2).
- 3) If needed: Before putting the user into the seat, unfasten the clip buckles at the back of the seat and pull the belt straps of the fixation vest forward somewhat (see Fig. 39).
- 4) Sit the user into the seat. Make sure that the user's back is up against the backrest padding.
- 5) Close the zipper (see Fig. 38, item 2).
- 6) If needed: Tighten the belt straps in the clip buckles at the back of the seat (see Fig. 39, item 3) and fasten them.
- 7) Close the hook-and-loop tabs over the zipper (see Fig. 38, item 1). **In daily use, the fixation vest can also be put on/removed by opening/closing the clips on the front (see Fig. 40, item 1/2).**



6.3.3.6 Groin strap

INFORMATION

For further information, please see the instructions for use included with the product.

The seating unit can be equipped with a groin strap that secures the pelvis in the seat for certain indications (see Fig. 41).

Adjusting the groin strap



The adjustment of the final belt tension takes place when putting it on as follows:

- Adjusting the belt length at the two clip buckles on the back of the seat (see Fig. 42, item 1).
- Adjusting the belts in the strap adjustment at the leg straps (see Fig. 43, item 1).

Putting on the groin strap

- 1) If needed: Open the clip buckles and pull the belt straps of the leg straps to the front somewhat (see Fig. 42, item 1).
- 2) Open the two clips at the leg straps by pressing on the fastener tabs (see Fig. 43, item 1).
- 3) Place the user on the groin strap and thread the leg straps between the legs and over the thighs to the outside (see Fig. 41).

- 4) Close the two clips at the leg straps (see Fig. 43, item 1).
- 5) If needed: Tighten the belt straps in the clip buckles at the back of the seat and fasten them (see Fig. 42, item 1).
- 6) By pulling on the free belt ends in the belt adjustment, the length of the leg straps can be adjusted, if needed (see Fig. 43, item 2).



6.3.3.7 Neck support

The shaped padding helps to immobilise the user's head in the seating unit.
The option has already been installed by qualified personnel.



6.3.3.8 Lumbar padding

The shaped padding offers additional support for the lumbar spine.
The option has been positioned and fastened below the seat padding by qualified personnel.



6.3.3.9 Hook-and-loop foot straps

The seating unit can be equipped with hook-and-loop foot straps that secure the child's feet firmly to the footrest for certain indications (see Fig. 47).

Using the hook-and-loop foot strap

- 1) Open the hook-and-loop fastener.
- 2) Place the foot on the footrest.
- 3) Pull the strap firmly around the shoe and close the hook-and-loop fastener fully.



6.3.4 Tray

⚠ CAUTION

Improper lifting by attendants

Falling, tipping over of the user due to lifting by removable components

- The product must not be lifted by the tray.

NOTICE

Overloading

Damage to the product through user error

- Trays must not be loaded with heavy objects.

INFORMATION

- Please note that the tray cannot be combined with the options "armrests" or "grab rail", as the same mounts are used.
- Please note that the tray can be used only when the seat is facing forward.

The seating unit can be equipped with a tray on which objects can be placed directly next to the seat.

The tray is mounted on the seating unit and secured against rotation.

Mounting the tray



- 1) Plug the tray into the adapters.
- 2) Push in the tripod spring (see Fig. 48, item 1).
- 3) Press down the tray until the tripod springs lock into place.

Removing the tray

- 1) Push in the tripod spring (see Fig. 48, item 1).
- 2) Pull the tray out and up.

Adjusting the tray position

- 1) **Adjust height:** loosen the Allen screws on the frame with a size 4 Allen wrench (see Fig. 48, item 2) and adjust the frame in parallel on both sides. Retighten the Allen head screws.
- 2) **Adjust depth:** release the adjustment lever under the tray and set the depth by moving the tray (see Fig. 49, item 1). Then firmly close the adjustment lever.
- 3) **Adjust the angle:** release the adjustment lever under the tray and set the tray angle by turning the tray at the joint (see Fig. 49, item 2). Then firmly close the adjustment lever.



6.3.5 Canopy

⚠ CAUTION

Exposed pinch points

Crushing, pinching due to incorrect handling

- When folding out the backrest, only grip by the specified components.

NOTICE

Improper use of the canopy

Damage to the fabric through user error

- Remove the canopy before folding the rehab buggy.
- Remove the canopy from the rehab buggy if it is stored for a long period to protect it from damage.

The rehab buggy can be equipped with a canopy to protect the user from sun or rain (see Fig. 50).

Mounting the canopy

- 1) Move the backrest to an upright position.
- 2) Move the locking tongues for the canopy into parallel mounting position (see Fig. 55).
- 3) Press the locking tongues for the canopy into the buckles (see Fig. 55, item 1) until they click into place audibly (see Fig. 55, item 2).


Adjusting the canopy angle

- 1) Fold up the canopy.
- 2) Hold the canopy at the front support and change the angle by pressing or pulling on the support. The setting joint locks into place automatically (see Fig. 55, item 3).

Removing the canopy

- 1) Move the backrest to an upright position.

- 2) Release the locking tongues for the canopy from the buckle by pressing lightly (see Fig. 55, item 2).
- 3) Pull the locking tongues out and remove the canopy.



Special features of the canopy

- Opening the window allows for permanent visual contact with the child (see Fig. 52).
- Removal of the back part allows for improved ventilation (e.g. in case of accumulating heat in summer; see Fig. 53).
- The rain cover bag contains the rain cover (see Fig. 57).



6.3.6 Rain canopy

The seating unit can be equipped with a rain canopy to protect the user from sun or rain (see Fig. 54).

Mounting the rain canopy


- 1) Move the backrest to an upright position.
- 2) Move the locking tongues for the rain canopy into parallel mounting position (see Fig. 55).
- 3) Press the locking tongues for the rain canopy into the buckles (see Fig. 55, item 1) until they click into place audibly (see Fig. 55, item 2).

Adjusting the rain canopy angle

- 1) Fold up the rain canopy.
- 2) Hold the rain canopy at the front support and change the angle by pressing or pulling on the support. The setting joint locks into place automatically (see Fig. 55, item 3).

Removing the rain canopy

- 1) Move the backrest to an upright position.
- 2) Release the locking tongues for the rain canopy from the buckle by pressing lightly (see Fig. 55, item 2).
- 3) Pull the locking tongues out and remove the rain canopy.



Special features of the rain canopy

- The transparent rain canopy offers protection against rain from the side and permits permanent visual contact with the child (see Fig. 54).

6.3.7 Rain cover

NOTICE

Incorrect attachment

Damage to product due to incorrect installation

- During the attachment, make sure that the product cannot grind against the wheels.

The rehab buggy can be equipped with a rain cover that prevents the user from getting wet (see Fig. 56). The rain cover can be stored in the side bag of the canopy (see Fig. 57).

Using the rain cover

- 1) Take the rain cover from the side bag of the canopy (see Fig. 57).
- 2) Put the rain cover over the canopy and the footrest. Make sure that the air intake is placed open over the user's nose (see Fig. 56, item 1).



6.3.8 Summer/winter slip sack

NOTICE

Incorrect attachment

Damage to product due to incorrect installation

- During the attachment, make sure that the product cannot grind against the wheels.

The rehab buggy can be equipped with a summer/winter slip sack (see Fig. 58). The warming winter part can be inserted/removed using snap fasteners.

Fastening the slip sack

- 1) Place the slip sack in the rehab buggy.

- 2) If needed, place the lap belt (see Fig. 59, item 1), the abductor/five-point belt (see Fig. 59, item 2) or the groin strap (see Fig. 59, item 3) through the openings underneath the slip sack.
- 3) Fasten the attachment cords (see Fig. 59, item 4) depending on the attachment situation:
 - Back: e.g., at the seat plate (see Fig. 60, item 1)
 - Front: e.g., at the footrest tube (see Fig. 60, item 2)

Using the slip sack

- 1) If needed: insert/remove winter part using snap fasteners.
- 2) Open the zipper and place the user inside.
- 3) Close the zipper and the closures.



6.3.9 Seat depth reducer

⚠ CAUTION

Exposed pinch points

Crushing, pinching due to incorrect handling

- When folding out the backrest, only grip by the specified components.

The seat depth of the Kimba Neo seating unit in size 1 can be reduced by **50 mm** to **140 mm** with the seat depth reducer.

If the seat depth needs to be increased, the seat depth reducer can be removed. Please contact the qualified personnel who adjusted the product.




6.3.10 Incontinence cover

The incontinence cover (see Fig. 63) protects the padding of the Kimba Neo seating unit. The incontinence cover can be removed for washing. Close the hook-and-loop straps before washing.



Attaching the incontinence cover

- 1) Remove the seat padding and the back padding from the seating unit (see Page 31 et seq.).
- 2) Lay out the incontinence cover ready for the seat and back padding. To do this, open all the hook-and-loop straps (see Fig. 64).
- 3) Place the seat padding in the incontinence cover.
- 4) Secure the incontinence cover to the back padding:
 - Place the back padding on the incontinence cover (see Fig. 65).
 - Guide the lower 3 hook-and-loop straps around the back padding and strap them together (see Fig. 66, item 1).
 - Thread the vertical hook-and-loop straps through the back padding (see Fig. 66, item 2).
- 5) Secure the back padding to the seating unit:
 - Pull the back padding with its cap over the back insert.
 - Guide the vertical hook-and-loop straps towards the rear over or through the back plate (see Fig. 67, item 1).
 - Guide the upper horizontal hook-and-loop strap around the back plate and strap it together (see Fig. 67, item 2). In doing so, the vertical hook-and-loop straps must be secured with the horizontal hook-and-loop strap (see Fig. 67, item 3).
- 6) Push the back padding towards the rear and secure it to the back of the seating unit with the twist lock.



Removing and cleaning the incontinence cover

- 1) Remove the seat padding and the back padding from the seating unit (see Page 31 et seq.).
- 2) Remove the incontinence cover from the seat padding and lay it down.
- 3) Remove the incontinence cover from the back padding. To do this, open all the hook-and-loop straps and remove the incontinence cover from the back padding.
- 4) Wash the incontinence cover according to the washing symbols on the label with care instructions.

INFORMATION: Fully close the hook-and-loop straps before washing.

6.3.11 Armrests

⚠ CAUTION

Improper lifting by attendants

User may fall or tip over

- The armrests must not be used to lift the seating unit or the entire rehab buggy.


INFORMATION

Please note that the armrests cannot be combined with the options "grab rail" or "tray", as the same mounts are used.

The rehab buggy can be equipped with armrests that support the arms for certain indications. The removable armrests have been attached to the seating unit and secured against rotation.

Removing/mounting the armrests

- 1) Press in the tripod spring (see Fig. 69).
- 2) Pull the armrest up and out or insert the armrest into the adapter.



6.3.12 Grab rail

⚠ CAUTION

Improper lifting by attendants

User may fall or tip over

- The grab rail must not be used to lift the seating unit or the entire rehab buggy.

INFORMATION

Please note that the grab rail cannot be combined with the options "armrests" or "tray", as the same mounts are used.

The rehab buggy can be equipped with a grab rail that prevents the sitting child falling forward for certain indications.

The grab rail can be removed to facilitate getting into the buggy.

Removing the grab rail

- 1) Press in the tripod springs (see Fig. 68, item 1).
- 2) Pull out and up the grab rail.

Mounting the grab rail

- 1) Insert the proper sides of the grab rail into the adapter (see Fig. 68).
- 2) Press in the tripod springs (see Fig. 68, item 1).
- 3) Press down the grab rail until the tripod springs lock into place.

6.4 Disassembly/transport

⚠ CAUTION

Exposed pinch points

Crushing, pinching due to incorrect handling

- When folding out the backrest, only grip by the specified components.

NOTICE

Deformation when folded

Damage to the product, problems unfolding due to unallowable loads

- Never place heavy objects on the folded product.

INFORMATION

The product can be equipped with the Kimba Neo seating unit or alternative seating systems. In this case, please observe the information in the separate instructions for use for the seating system.

The seating unit has to be removed from the mobility base (see Page 10) for transport.

6.5 Use in a wheelchair accessible vehicle

6.5.1 General Safety Instructions

WARNING

Use as a seat in a wheelchair accessible vehicle

Serious injuries in case of accidents due to user error

- ▶ Always use the seats and restraint systems in the wheelchair accessible vehicle first. This is the only way to provide optimal protection for the occupants in case of an accident.
- ▶ The product was tested for use in a wheelchair accessible vehicle in combination with a Kimba outdoor mobility base. For more information please refer to our brochure on using Ottobock products in wheelchair accessible vehicles, order number 646D158.
- ▶ The product may only be used as a seat in a wheelchair accessible vehicle with the securing elements offered by the manufacturer and with suitable restraint systems.
- ▶ Please note that you must always remove certain options before using the product in a wheelchair accessible vehicle.
- ▶ As a matter of principle, only use the product in wheelchair accessible vehicles with the seat facing in the driving direction.
- ▶ As a matter of principle, only use the product in wheelchair accessible vehicles with the seating unit locked horizontally and the backrest locked in the upright position.

WARNING

Incorrect transport weight when using in a wheelchair accessible vehicle

Serious injuries in case of accidents due to exceeding the allowable load capacity

- ▶ During use in a wheelchair accessible vehicle, the maximum user weight in the Kimba Neo seating unit is limited to **40 kg**.

WARNING

Use of the belt system as a restraint system in a wheelchair accessible vehicle is forbidden

Serious injuries due to improper handling of the product

- ▶ Under no circumstances may the safety belts and positioning aids offered with the product be used as part of a restraint system for transportation in a wheelchair accessible vehicle.
- ▶ Note that the safety belts and positioning aids offered with the product are only intended to help support the user sitting in the product.

INFORMATION

For transporting children up to **36 kg**, Ottobock urgently recommends using a restraint system for children (known as a child car seat) tested according to ECE 44.03 or higher, or a restraint system according to ECE 44.04 developed especially for disabled children (e.g. "Lars" child car seat from Ottobock). Observe the legal regulations for your country/state in all cases.

6.5.2 Permitted use

The product may be used for seating in a wheelchair accessible vehicle if a suitable restraint system is used. The safety belts and positioning aids offered by the manufacturer only serve to further stabilise the person using the product.

The following options must be detached from the rehab buggy prior to transport in the wheelchair accessible vehicle:

- Sun canopy
- Grab rail
- Tray
- Diaper bag
- Backpack

6.5.3 Important information

The following important information is contained in the 647G891 Instructions for Use of the Kimba mobility bases:

- Adjusting and securing the seat angle
- Adjusting the seating unit

- Adjusting and securing the back angle
- Securing the product in the vehicle
- Securing the user in the product

6.6 Care


- Do not allow sand or other particles to get into the product. They could degrade the wheel bearings and the locking mechanism.
- Use in saltwater is not permitted.

6.6.1 Detaching/attaching the padding for cleaning

The padding can be removed for care and maintenance of the product.

6.6.1.1 Detaching/attaching the seat padding

- 1) Lift the seat padding and remove it from the hook-and-loop fastener on the seat plate (see Fig. 70).
- 2) Clean the seat padding.
- 3) Fit the seat padding onto the hook-and-loop fastener.
- 4) Push the seat padding under the hip supports (see Fig. 71).



6.6.2 Detaching/attaching the back padding

INFORMATION

The manufacturer offers two versions of the back padding:

- Back padding for standard back insert (see Fig. 72)
- Back padding for back insert with headrests (see Fig. 73)

- 1) Release the twist lock of the back padding behind the seat bottom (see Fig. 74).
- 2) Unzip the back padding zipper (see Fig. 75) and remove/clean the back padding.
- 3) Put the back padding in place and pull it through between the back and seat bottom.
- 4) Close the twist lock (see Fig. 74).
- 5) Slide the back padding under the lateral supports (no illustration).







6.6.3 Detaching/attaching the covers


INFORMATION

The manufacturer offers the covers in the following versions:

- ▶ Covers for headrest supports (see Fig. 76, item 1)
- ▶ Covers for lateral supports (see Fig. 76, item 2)
- ▶ Covers for thoracic supports (no illustration)
- ▶ Cover for the grab rail (see Fig. 76, item 3; see the following section for replacement)
- ▶ Covers for hip supports (see Fig. 76, item 4)
- ▶ Cover for the abductor (see Fig. 76, item 5)


- 1) Remove the seat padding or back padding if required.
- 2) Open the hook-and-loop closure on the cover.
- 3) Remove the cover from the support or abductor and clean it.
- 4) Slide the cover over the support with proper alignment. Note the following:
 - If present: position the side reinforced with foam material on the **outside** (over the screw heads) (see Fig. 77, item 1).
 - Pull the cover over the support starting with the side that does not have a lug (see Fig. 78, item 1).
 - Pull the cover over the support with the help of the lug (see Fig. 78, item 2).
 - Before fastening the hook-and-loop closure, push and fold down the lug (see Fig. 79).
- 5) Put the seat padding or back padding back into place.





6.6.4 Detaching/attaching the grab rail cover (option)

- 1) Push in the tripod spring and remove the grab rail (see Fig. 80, item 1).
- 2) Pull the grab rail pad (see Fig. 76, item 3) off the grab rail and clean it.
- 3) Slide the grab rail pad onto the grab rail (see Fig. 80).
- 4) Reattach the grab rail.



6.7 Cleaning

Cleaning by hand

- 1) Clean all frame components and plastic parts using mild detergents only.
- 2) Wash covers made of spacer fabric and microfibre at a maximum temperature of **40 °C**. To do this, remove the covers from the padding and wash the covers only.

INFORMATION: You can wash covers in the washing machine – but only in a laundry bag or pillow case. However, wiping with a damp cloth is sufficient in most cases.

Important information on cleaning

- Do not use any aggressive cleaners, solvents or hard brushes, etc.

Cleaning in a washing chamber

- 1) Remove the padding and covers, and follow the washing instructions for each article.
- 2) Wash the seating unit at max. **60 °C** for max. 10 minutes.
- 3) **If necessary:** Remove all components where water has accumulated.
- 4) Allow the seating unit to dry. In order to do so, set the folded rehab buggy vertically onto the back wheels and push handles so that water can run out (see Fig. 81).
INFORMATION: We recommend drying with compressed air.
- 5) **If necessary:** Reinstall all components that were removed.
- 6) Before storing or using the seating unit again, verify that it is dry. Wipe wet components with a cloth as needed.
- 7) Ensure that all labels (nameplate and warnings) are still legible. Otherwise replace these labels.



7 Maintenance/Repair

7.1 Maintenance

7.1.1 Maintenance intervals

The functions described below can be checked by the user at home at the specified intervals. During these tasks, the user must not be sitting in the seating unit:

Component	Activity	Before each use	Weekly	Monthly
Seating unit	Check attachment on mobility base			X
	Check locking of seat tilt feature	X		
	Verify the function of knobs			X
Footrest	Check for damage			X
	Check footrest and angle adjustment mechanism for stability			X
Adapters for optional features	Check adapters for damage			X
Padding/belts	Check padding for proper condition			X
	Check belts for wear		X	
	Verify the belt buckle functionality	X		

8 Disposal

8.1 Disposal Information

Return the product to the specialist dealer for disposal.

All components of the product must be disposed of properly in accordance with the respective national environmental regulations.

9 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Service Life

Based on market observations and the current state of technology, the manufacturer has calculated that the product can be used for a period of **4 years**, provided that it is used properly and that the service and maintenance instructions are observed. Periods of storage at the specialist dealer or with the paying party are not included in this period.

However, we would like to emphasise that the product is reliable far beyond this defined period of time, provided that it is cared for and maintained properly.

If the service life is reached, the user or a responsible attendant should contact the qualified personnel who fitted the product or the manufacturer's service (see inside rear cover or back page for address).

9.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.3 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

9.4 Warranty

Further information on the warranty terms and conditions can be obtained from the qualified personnel that has fitted this product or the manufacturer's service (see inside back cover for addresses).

9.5 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

10 Technical data

Kimba Neo seating unit	Size 1	Size 2
Seat width	200 – 300	260 – 350
Seat depth	190 ²⁾ – 310	240 – 400
Backrest height	410 – 610	560 – 700
Backrest angle	80° – 180°	80° – 180°
Height of shoulder strap slots (min./max.)	300/450	300/450
Width of shoulder strap slots	50	50
Lower leg length	190 – 310	200 – 370
Max. load capacity ¹⁾	40	40
Max. load capacity for use in a wheelchair accessible vehicle ³⁾	40	40
Weight	7	8.5
Folded size with lateral supports (LxWxH)	580x370x420	590x430x580

¹⁾ The load capacity differs between the Kimba outdoor mobility base and the seating unit:

- Seating unit (max. 40 kg) + own weight of the seating unit, size 1 (7.5 kg) = 47.5 kg
- Outdoor mobility base (max. 55 kg) – fully loaded seating unit (47.5 kg) = 7.5 kg for other accessories + weight in storage bag

²⁾ Possible with seat depth reducer option from **140 mm**.

³⁾ Observe the information in the section "Use in a wheelchair accessible vehicle".



A series of horizontal lines for writing, spanning the width of the page.

[illegible]



A series of horizontal lines for writing, spanning the width of the page.

Kundenservice/Customer Service

Europe

Otto Bock HealthCare Deutschland GmbH
Max-Näder-Str. 15 · 37115 Duderstadt · Germany
T +49 5527 848-3411 · F +49 5527 848-1414
healthcare@ottobock.de · www.ottobock.de

Otto Bock HealthCare Deutschland GmbH
Max-Näder-Str. 15 · 37115 Duderstadt · Germany
T +49 5527 848-3411 · F +49 5527 848-1414
healthcare@ottobock.de · www.ottobock.de

Otto Bock HealthCare Deutschland GmbH
Max-Näder-Str. 15 · 37115 Duderstadt · Germany
T +49 5527 848-3455 · F +49 5527 848-1510
healthcare@ottobock.de · www.ottobock.de

Otto Bock HealthCare Deutschland GmbH
Max-Näder-Str. 15 · 37115 Duderstadt · Germany
T +49 5527 848-3433 · F +49 5527 848-1460
healthcare@ottobock.de · www.ottobock.de

Otto Bock HealthCare
Deutschland GmbH
37115 Duderstadt · Germany
T +49 5527 848-3433
F +49 5527 848-1460
healthcare@ottobock.de

Otto Bock Healthcare Products GmbH
Kaiserstraße 39 · 1070 Wien · Austria
T +43 1 5269548 · F +43 1 5267985
vertrieb.austria@ottobock.com · www.ottobock.at

Otto Bock Healthcare
Products GmbH
1070 Wien · Austria
T +43 1 5269548
F +43 1 5267985
vertrieb.austria@ottobock.com

Otto Bock Adria Sarajevo D.O.O.
Omladinskih radnih brigada 5
71000 Sarajevo · Bosnia-Herzegovina
T +387 33 766200 · F +387 33 766201
obadria@bih.net.ba · www.ottobockadria.com.ba

Otto Bock Adria Sarajevo D.O.O.
71000 Sarajevo
Bosnia-Herzegovina
T +387 33 766200
F +387 33 766201
obadria@bih.net.ba

Otto Bock Bulgaria Ltd.
41 Tzar Boris III' Blvd. · 1612 Sofia · Bulgaria
T +359 2 80 57 980 · F +359 2 80 57 982
info@ottobock.bg · www.ottobock.bg

Otto Bock Bulgaria Ltd.
1612 Sofia · Bulgaria
T +359 2 80 57 980
F +359 2 80 57 982
info@ottobock.bg

Otto Bock Suisse AG
Pilatusstrasse 2 · CH-6036 Dierikon
T +41 41 455 61 71 · F +41 41 455 61 70
suisse@ottobock.com · www.ottobock.ch

Otto Bock Suisse AG
CH-6036 Dierikon
T +41 41 455 61 71
F +41 41 455 61 70
suisse@ottobock.com

Otto Bock ČR s.r.o.
Protetická 460 · 33008 Zruč-Senec · Czech Republic
T +420 377825044 · F +420 377825036
email@ottobock.cz · www.ottobock.cz

Otto Bock ČR s.r.o.
33008 Zruč-Senec
Czech Republic
T +420 377825044
F +420 377825036
email@ottobock.cz

Otto Bock Iberica S.A.
C/Majada, 1 · 28760 Tres Cantos (Madrid) · Spain
T +34 91 8063000 · F +34 91 8060415
info@ottobock.es · www.ottobock.es

Otto Bock Iberica S.A.
28760 Tres Cantos (Madrid)
Spain
T +34 91 8063000
F +34 91 8060415
info@ottobock.es

Otto Bock France SNC
4 rue de la Réunion · CS 90011
91978 Courtaboeuf Cedex · France
T +33 1 69188830 · F +33 1 69071802
information@ottobock.fr · www.ottobock.fr

Otto Bock France SNC
91978 Courtaboeuf Cedex
France
T +33 1 69188830
F +33 1 69071802
information@ottobock.fr

Otto Bock Healthcare plc
32, Parsonage Road · Englefield Green
Egham, Surrey TW20 0LD · United Kingdom
T +44 1784 744900 · F +44 1784 744901
bocuk@ottobock.com · www.ottobock.co.uk

Otto Bock Healthcare plc
Egham, Surrey TW20 0LD
United Kingdom
T +44 1784 744900
F +44 1784 744901
bocuk@ottobock.com

Otto Bock Hungária Kft.
Tatai út 74. · 1135 Budapest · Hungary
T +36 1 4511020 · F +36 1 4511021
info@ottobock.hu · www.ottobock.hu

Otto Bock Hungária Kft.
1135 Budapest · Hungary
T +36 1 4511020
F +36 1 4511021
info@ottobock.hu

Otto Bock Adria d.o.o.
Dr. Franje Tuđmana 14 · 10431 Sveta Nedelja · Croatia
T +385 1 3361 544 · F +385 1 3365 986
ottobockadria@ottobock.hr · www.ottobock.hr

Otto Bock Adria d.o.o.
10431 Sveta Nedelja · Croatia
T +385 1 3361 544
F +385 1 3365 986
ottobockadria@ottobock.hr

Otto Bock Italia Srl Us
Via Filippo Turati 5/7 · 40054 Budrio (BO) · Italy
T +39 051 692-4711 · F +39 051 692-4710
info.italia@ottobock.com · www.ottobock.it

Otto Bock Italia Srl Us
Via Filippo Turati 5/7 · 40054 Budrio (BO) · Italy
T +39 051 692-4711 · F +39 051 692-4710
info.italia@ottobock.com · www.ottobock.it

Otto Bock Italia Srl Us
Via Filippo Turati 5/7 · 40054 Budrio (BO) · Italy
T +39 051 692-4711 · F +39 051 692-4720
info.italia@ottobock.com · www.ottobock.it

Otto Bock Italia Srl Us
Via Filippo Turati 5/7 · 40054 Budrio (BO) · Italy
T +39 051 692-4711 · F +39 051 692-4720
info.italia@ottobock.com · www.ottobock.it

Otto Bock Italia Srl Us
40054 Budrio (BO) · Italy
T +39 051 692-4711
F +39 051 692-4720
info.italia@ottobock.com

Otto Bock Benelux B.V.
Ekkersrijt 1412 · 5692 AK
Son en Breugel · The Netherlands
T +31 499 474585 · F +31 499 476250
info.benelux@ottobock.com · www.ottobock.nl

Otto Bock Benelux B.V.
5692 AK Son en Breugel
The Netherlands
T +31 499 474585
F +31 499 476250
info.benelux@ottobock.com

Industria Ortopédica Otto Bock Unip. Lda.
Av. Miguel Bombarda, 21 - 2º Esq.
1050-161 Lisboa · Portugal
T +351 21 3535587 · F +351 21 3535590
ottobockportugal@mail.telepac.pt

Industria Ortopédica
Otto Bock Unip. Lda.
1050-161 Lisboa · Portugal
T +351 21 3535587
F +351 21 3535590
ottobockportugal@mail.telepac.pt

Otto Bock Polska Sp. z o. o.
Ulica Korolowa 3 · 61-029 Poznań · Poland
T +48 61 6538250 · F +48 61 6538031
ottobock@ottobock.pl · www.ottobock.pl

Otto Bock Polska Sp. z o. o.
61-029 Poznań · Poland
T +48 61 6538250
F +48 61 6538031
ottobock@ottobock.pl

Otto Bock Romania srl
Șos de Centura Chitila · Mogoșoia Nr. 3
077405 Chitila, Jud. Ilfov · Romania
T +40 21 4363110 · F +40 21 4363023
info@ottobock.ro · www.ottobock.ro

Otto Bock Romania srl
077405 Chitila, Jud. Ilfov
Romania
T +40 21 4363110
F +40 21 4363023
info@ottobock.ro

OOO Otto Bock Service
p/o Pultikovo, Business Park „Greenwood”,
Building 7, 69 km MKAD
143441 Moscow Region/Krasnogorskiy Rayon
Russian Federation
T +7 495 564 8360 · F +7 495 564 8363
info@ottobock.ru · www.ottobock.ru

OOO Otto Bock Service
143441 Moscow
Region/Krasnogorskiy Rayon
Russian Federation
T +7 495 564 8360
F +7 495 564 8363
info@ottobock.ru

Otto Bock Scandinavia AB
Koppargatan 3 · Box 623 · 60114 Norrköping · Sweden
T +46 11 280600 · F +46 11 312005
info@ottobock.se · www.ottobock.se

Otto Bock Scandinavia AB
60114 Norrköping · Sweden
T +46 11 280600
F +46 11 312005
info@ottobock.se

Otto Bock Slovakia s.r.o.
Röntgenova 26 · 851 01 Bratislava 5 · Slovak Republic
T +421 2 32 78 20 70 · F +421 2 32 78 20 89
info@ottobock.sk · www.ottobock.sk

Ihr Fachhändler | Your specialist dealer



Otto Bock Mobility Solutions GmbH
Lindenstraße 13 · 07426 Königsee-Rottenbach/Germany
www.ottobock.com

Ottobock has a certified Quality Management System in accordance with ISO 13485.