



Bariatric Shower Commode Wheeler

We have aimed at designing the best shower commode wheeler on the market. It is designed for easy access, independence, safety and on the caregiver terms. The 24" rear wheels provide more directional stability and smooth pushing.

The inside width of the chair is 61 or 71 cm, but using a five cm high PU foam seat and seven cm free space to the armrests allows the body mass to fall down, instead of fading out to the sides. The impossible "wide inside but small outside" is obtained.

The footrest is always the weak point when it comes to bariatric products. Normally the producer makes a short cut here, by using a standard solution. We did not only design the footrest to carry 100 kg each, we did more!

The footrest is also side and width adjustable, so access to intimate hygiene is possible and people with extra tissue placed on the inside of their leg can sit with legs apart (up to 110 cm) and still have foot support. The footrests are covered with PU foam.

One of the unique new features is a split toilet seat, a long open hole: 44 cm and 22 cm at the widest point. This to enable a correct use of the chair, it eases access for caregivers and avoids incidents where material ends up in places it wasn't meant to end up!

The ability for bariatric users to reach far back to reach the armrest is usually not present. Therefore we have made the armrests with a length of 64 cm and let it stick out seven cm from the seat base, ready to welcome and support the hands of mobile users. It enables the user to get into the chair by him/ herself. Another option is the butterfly armrests. Butterfly arm-rests add 14 cm to the seating width.

The space between backrest and seat has a clearance of 32 cm. This gives space for a possible bulbous gluteal region, but also enables the caregivers to reach for the sling and position the user better.



- Shower Commode Wheeler
with 24" rear wheels is easy to drive
for the attendant

MAX 325 kg

Bariatric Shower Commode Wheeler Specifications

Bariatric Shower Commode Wheeler	
Maximum weight capacity is	325 kg
Seat width with standard armrests	61 / 71 cm
Total outer width with standard armrests	68 / 78 cm
Seat depth	56 cm
Seat height	55 cm
Backrest	32-52 cm
Seat to footrest	40-42-44-46-48 cm
Width of footrest	58-110 cm
Pushbar from floor	100 cm
Width between armrests	54 cm
Height of backrest	48 cm
Height of footrests	40-48 cm
Distance from floor to seat underside	44.5 cm
Total length	129 cm
Total weight	27 / 30 kg
Front wheels	125 mm
Rear wheels	24"
Seat	Pressure moulded PUR foam
Frame	Powder-lacquered chromium molybdenum steel
ISO - classification	ISO 09 12 03

Ideal use for.

- Bariatric users in institutions, private homes and hospitals.



Wheel lock



Calf support or amputation support

Model	Seat width cm	Seat depth cm	Order No
XXL-Rehab Shower commode wheeler incl. standard armrests, standard footrests and commode bucket	61	56	0152-061-000
XXL-Rehab Shower commode wheeler incl. standard armrests, standard footrests and commode bucket	71	56	0152-071-000
XXL-Rehab Shower commode wheeler incl. butterfly armrests, standard footrests and commode bucket	61 + 14	56	0153-061-000
XXL-Rehab Shower commode wheeler incl. butterfly armrests, standard footrests and commode bucket	71 + 14	56	0153-061-000

Accessories	
Butterfly armrests left (add 14 cm to seat width)	0150-999-012L
Butterfly armrest right (add 14 cm to seat width)	0150-999-012R
Commode bucket	0140-999-001
Standard footrest right	0156-999-090H
Standard footrest left	0156-999-090V
Full size seat for shower commode with lid for bucket	0150-999-023
Calf support / amputation support right	0156-999-030
Calf support / amputation support left	0156-999-031



Seat with lid



Split toilet seat