

# 4FRONT™



The 4Front® is a quiet, more responsive front-wheel drive power chair that features automotive-grade suspension resulting in unprecedented comfort and ride quality. The 4Front® features standard speeds of up to 6 mph in the lowered position, and speeds up to 3.2 mph while elevated at 10" with our optional safe seat elevation.



ELECTRIC BLUE

  
QUANTUM®



QUANTUMREHAB.CO.UK

## Specifications

<b>WEIGHT CAPACITY</b> <sup>4</sup>	21 stone 6 lbs. (136 kg)
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 6 mph Up to 3.2 mph w/10" lift
<b>GROUND CLEARANCE</b> <sup>2</sup>	48 mm (1.9") at anti-tip wheels 67 mm (2.6") at motor pivot 95.25 mm (3.75") at battery tray
<b>TURNING RADIUS</b> <sup>2</sup>	622.5 mm (24.5")
<b>OVERALL SIZE</b> <sup>2, 6</sup>	Length: 1163 mm (45.7") with front riggings 1014 mm (39.9") without front riggings Width: 616 mm (24.2")
<b>DRIVE WHEELS</b>	355.6 mm (14") flat free
<b>CASTER WHEELS</b>	228.6 mm (9") flat free
<b>ANTI-TIP WHEELS</b>	101.6 mm (4") flat free
<b>SUSPENSION</b>	SRS (Smooth Ride Suspension)
<b>COMPONENT WEIGHTS</b> <sup>5, 7</sup>	Base: 140 lbs (63.5 kg) Batteries: 51.5 lbs (23.36 kg) each (standard) 38.5 lbs (17.46 kg) each (optional)
<b>DRIVETRAIN</b>	Two motor, front-wheel drive
<b>BRAKES</b>	"Intelligent Braking" electronic regenerative, disc park brake
<b>BATTERIES</b>	Two 12-volt, deep-cycle 73 Ah (standard)
<b>BATTERY CHARGER</b>	8-amp, Off-board
<b>MAXIMUM OBSTACLE CLIMBING ABILITY</b>	85 mm (3.34")
<b>SEATING SIZES</b>	
<b>TRU-Balance® 3 Static Seat</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>TRU-Balance® 3 Power Tilt</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>TRU-Balance® 3 Lift &amp; Tilt</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>TRU-Balance® 3 Static Seat</b> <sup>5</sup>	15.25-17.25"
<b>TRU-Balance® 3 Power Tilt</b> <sup>5</sup>	17.5-19.5"
<b>TRU-Balance® 3 Lift &amp; Tilt</b> <sup>5</sup>	17.75-19.75"

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.  
 2) Without front riggings  
 3) Battery capacity may change when adding power positioning systems.  
 4) Weight capacity is dictated by the selection of seating.  
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 21 stone 6 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details.  
 Listed measurements are in the lowest position; does not include elevating system unless specified.  
 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.  
 7) Battery weight may vary +/- 3 lbs. based on manufacturer.  
 8) All Quantum Rehab® units can be ordered with specialty controls.  
 9) Call for availability.

## Features

- SRS (Smooth Ride Suspension)
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions
- Q-Logic 3 EX Advanced Drive Controls
- LED front and rear fender lights
- Quiet ride technology

## Options

- 10" of power adjustable safe seat elevation at speeds of up to 3.2 mph
- TRU-Balance® 3 Seating
- TRU-Balance® 3 Power Tilt
- TRU-Balance® 3 Power Recline
- Powered Elevating leg rests

## Colours\*



RED  
CORVETTE



ELECTRIC  
BLUE



LIME  
LIGHT



PRETTY  
IN PINK



NIGHT  
TRAIN

\* Please contact your representative for the full range of options and colours.