

ORDER FORM

Quickie® All Court Ti



Failure to use an order form for a Quickie chair order may result in delay or wrong shipment

Please send all orders to:

Fax: +44 (0)845 605 66 89

Email: orders@sunmed.co.uk

Sunrise Medical Ltd


Thorns Road, Brierley Hill, West Midlands, DY5 2LD

www.SunriseMedical.co.uk

Standard features:


- ◆ Base frame single integral anti-tip
- ◆ Choice of 28 frame colours
- ◆ Choice of 4 frame lengths
- ◆ Choice of seat widths ranging from 12" to 20"
- ◆ Fully adjustable frame
- ◆ Patented centre of mass adjustment
- ◆ Choice of seat depths ranging from 10" to 20"
- ◆ 24" spoked rear wheels
- ◆ Choice of folding and non-folding back
- ◆ 3" Micro castors
- ◆ Quick release axles
- ◆ Choice of 16° or 20° camber
- ◆ Maximum user weight limit 112kg (18st)

How to measure:




Seat Depth

From front of seat sling to back of seat sling




Seat Width

From outside seat rail tube to outside seat rail tube



Frame Length

From outside rear frame to outside front frame



Back Height

From seat sling to top of backrest sling

Prices:

All prices are in Pounds Sterling/Euros, exclude VAT. No responsibility held for technical changes or printing failures. Please consult Sunrise Medical with any queries you may have.

Quickie® All Court Ti

Purchase Order (Repeat on all pages):

		£			£	TOTAL
FRONT SEAT HEIGHT (Mark one box)						
154SH14	<input type="radio"/> Seat height 16" (40.5cm)	Std	154SH21	<input type="radio"/> Seat height 19" (48cm)	Std	
154SH16	<input type="radio"/> Seat height 17" (43cm)	Std	154SH24	<input type="radio"/> Seat height 20" (51cm)	Std	
154SH18	<input type="radio"/> Seat height 18" (45.5cm)	Std	154SH27	<input type="radio"/> Seat height 21" (53.5cm)	Std	
BACKREST (Mark one type and choices underneath)						
Type:			Back Rest Upholstery:			
154BT1	<input type="radio"/> Folding backrest	N/C	154BU1	<input type="radio"/> Standard upholstery	N/C	
154BT3	<input type="radio"/> Non-folding backrest	N/C	154BU13	<input type="checkbox"/> Adjustable sports upholstery	150	
Back Height:			154BU14	<input type="radio"/> Omit back upholstery	N/C	
154B128	<input type="radio"/> 9" (23cm)	N/C				
154B129	<input type="radio"/> 10" (25.5cm)	N/C				
154B130	<input type="radio"/> 11" (28cm)	N/C				
154B131	<input type="radio"/> 12" (30.5cm)	N/C				
154B33	<input type="radio"/> 13" (33cm)	N/C				
154B132	<input type="radio"/> 14" (35.5cm)	N/C				
154B35	<input type="radio"/> 15" (38cm)	N/C				
154B36	<input type="radio"/> 16" (40.5cm)	N/C				
154B37	<input type="radio"/> 17" (43cm)	N/C				
154B38	<input type="radio"/> 18" (45.5cm)	N/C				
FOOTRESTS (Mark one box)						
154F234	<input type="radio"/> Platform	N/C	154F39	<input type="radio"/> High mount platform	N/C	
154IG20	<input type="radio"/> Impact Guard Neoprene	N/C				
ADAPTIVE PLATFORM SPACER (Mark one box)						
154AD1	<input type="checkbox"/> 2" Spacer (5)	108	154AD3	<input type="checkbox"/> 6" Spacer (5)	100	
154AD2	<input type="checkbox"/> 4" Spacer (5)	108				
CASTORS (Mark one box)						
154CW1	<input type="radio"/> 3" Micro	N/C				
REAR WHEEL SIZE (Mark one box)						
154RW53	<input type="radio"/> 24" X/S Spoke wheel (540)	N/C	154RW59	<input type="checkbox"/> 24" Spinergy Sport Wheel (559) Black	623	
154RW69	<input type="checkbox"/> 24" X/S QP Heat treated tennis wheel (559)	355	154RW61	<input type="checkbox"/> 26" Spinergy Sport Wheel (590) Black	623	
154RW70	<input type="checkbox"/> 24" X/S QP Heat treated basketball wheel (559)	355				
154RW66	<input type="checkbox"/> 26" X/S QP Heat treated basketball wheel (590)	355				
154RW72	<input type="checkbox"/> 26" X/S QP Heat treated basketball wheel (700C)	355	154RW63	<input type="radio"/> Omit rear wheel and axle (540)	N/C	
154RW57	<input type="checkbox"/> 24" Spinergy Sport Wheel (540) Black	672	154RW64	<input type="radio"/> Omit rear wheel and axle (559)	N/C	
			154RW65	<input type="radio"/> Omit rear wheel and axle (590)	N/C	
			154RW71	<input type="radio"/> Omit rear wheel and axle (700C)	N/C	
REAR WHEEL TYRES (Mark one box)						
154RT31	<input type="radio"/> High pressure clincher	N/C				
154RT5	<input type="radio"/> Turbo tread	N/C				
HANDRIMS (Mark one box)						
154HR1	<input type="radio"/> Aluminium anodised	Std	154HR7	<input type="checkbox"/> Omit handrims	-24	
154HR4	<input type="checkbox"/> Plastic coated	92	Handrim Tabs:			
			154TB1	<input type="radio"/> Regular tab	N/C	
AXLES (Mark one box)						
154AX4	<input type="radio"/> Quick-release axles	N/C				
154AX11	<input type="checkbox"/> Quad quick-release axles	92				
CAMBER (Mark one box)						
154CI10	<input type="radio"/> 16°	N/C	154CI11	<input type="radio"/> 20°	N/C	
REAR ANTI-TIP (Mark one box)						
154AT39	<input type="radio"/> Rear anti-tip integral	N/C				
154AT40	<input type="checkbox"/> Rear anti-tip quick release	80				
ACCESSORIES (Optional)						
Back Pack:			Positioning Belt:			
154A1	<input type="checkbox"/> Back pack adult	63	154A17	<input type="checkbox"/> Hook and loop closure	43	
154A2	<input type="checkbox"/> Back pack kids	63	154A18	<input type="checkbox"/> Nylon buckle closure	43	
Seat Pouch:			154A190	<input type="checkbox"/> Ratchet positioning belt (Large)	81	
154A8	<input type="checkbox"/> Seat pouch	63	154A191	<input type="checkbox"/> Ratchet positioning belt (Extra large)	81	
Spoke Guard:			Miscellaneous:			
154A24	<input type="checkbox"/> Spoke guard (Clear)	118	154A103	<input type="checkbox"/> Chest strap 2" (length 48")	22	
154A23	<input type="checkbox"/> Spoke guard (Black)	118	154A106	<input type="checkbox"/> Chest strap 4" (length 48")	22	
154A25	<input type="checkbox"/> Spoke guard sport (High profile) (11)	138	154A90	<input type="checkbox"/> Pads for backrest crossbar (10)	22	
154A26	<input type="checkbox"/> Spoke guard (Low profile) (11)	138	154A98	<input type="checkbox"/> Tyre pump	85	
Side Guard:			154A180	<input type="checkbox"/> Toe clips	36	
154SG3	<input type="checkbox"/> Plastic (Regular)	174	154A181	<input type="checkbox"/> 5th Wheel support (8)	81	
154SG4	<input type="checkbox"/> Plastic (Kids)	174				
154SG5	<input type="checkbox"/> Aluminium (Regular) (7)	240				

Base Price:	£ 3,737
Options:	£
Total Price:	£

(Not Including 5% Surcharge)

All prices are in Pounds Sterling/Euros, exclude VAT. No responsibility held for technical changes or printing failures. Please consult Sunrise Medical with any queries you may have.

Quickie All Court Ti

FOOTNOTES:

- 1 1/2" seat width is built by selecting the Aluminium side guards
- 2 Not available with 10" to 11" seat depths
- 3 Not available with 17" to 20" seat depths
- 4 Not available with 19" to 20" seat depths
- 5 Not available with high mount platform footrest
- 6 Not available with Spinergy wheels
- 7 Aluminium side guards required for 1/2" seat width option
- 8 Only available with base frame single integral anti-tip
- 9 Not available on Spinergy, 700C rear wheels
- 10 Required for tournament play
- 11 Black only

Seat Matrix:

	DEPTH	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"
WIDTH												
12"		X	X	X	X	X	X	X				
13"			X	X	X	X	X	X	X			
14"				X	X	X	X	X	X			
15"					X	X	X	X	X	X		
16"					X	X	X	X	X	X		
17"						X	X	X	X	X	X	X
18"						X	X	X	X	X	X	X
19"						X	X	X	X	X	X	X
20"						X	X	X	X	X	X	X

All prices are in Pounds Sterling/Euros, exclude VAT. No responsibility held for technical changes or printing failures. Please consult Sunrise Medical with any queries you may have.