



Bariatric Wheelchair Cushion Anatomical

Strongly engineered and functionally designed exclusively for bariatric users.

Anatomical Cushion is designed exclusively for bariatric users at risk of developing pressure sores.

Anatomical Cushion is made of a molded base 11 cm firm foam and with a three cm viscoelastic foam on the top.

The molded base provides increased stability and seat support and provides unique pressure relief on the Tuber Ischii.

Anatomical Cushion is 14 cm high and has a two-way stretch vapour permeable incontinence and flame resistant cover.

- Anatomically shaped
- Two-way stretch incontinence cover

MAX 325 kg



Anatomical Cushion

Maximum weight capacity is	325 kg
Weight Approx.	3 kg
Basis cushion	HD PUR foam / Viscoelastic PUR foam
Cover	PU
Incontinence cover	Yes
Washable cover	Yes
Flame retardant cover	Yes

Ideal use for...

- Bariatric users at risk of developing pressure ulcers.

Model	Seat width cm	Seat depth cm	Height cm	Maximum weight capacity	Order No
XXL-Rehab Anatomical cushion	56	50	14	325	0233-056-050
XXL-Rehab Anatomical cushion	61	50	14	325	0233-061-050
XXL-Rehab Anatomical cushion	66	50	14	325	0233-066-050
XXL-Rehab Anatomical cushion	71	50	14	325	0233-071-050
XXL-Rehab Anatomical cushion	81	50	14	325	0230-081-050
XXL-Rehab Anatomical cushion	91	50	14	325	0233-091-050
XXL-Rehab Anatomical cushion	96	50	14	325	0233-096-050
XXL-Rehab Anatomical cushion	100	50	14	325	0233-100-050

Accessories

Cover for XXL-Rehab Anatomical cushion

0233-999-099

