

Revolutionary Suspension Comfortable Ride



There's smooth, and then there's the Apex Rapid with CTS! This scooter has revolutionary suspension for the smoothest, most comfortable ride. From the park to your backyard, this feature-rich scooter offers best-in-class ride quality.



Specifications


Weight Capacity	21 st 6 lbs (136 kg)
Maximum Speed ¹	Up to 4 mph
Ground Clearance ³	38 mm (1.5") at motor
Turning Radius ³	1162 mm (44.75)
Overall Length ³	1041 mm (41")
Overall Width ³	521 mm (20.5")
Max. Obstacle Climb	50 mm (2")
Drive Tyres	203 mm (8") mm front & rear - solid
Suspension	Comfort-Trac Suspension - front & rear
Total Weight Without Batteries ^{3, 5}	39.6 kg (87.4 lbs)
Heaviest Piece When Disassembled ³	15 kg (33 lbs)
Drivetrain	24 volt DC motor
Dual Breaking System	Regenerative & electromechanical
Battery Requirements ³	2 x 12 or 17 Ah deep cycle
Battery Weight ⁴	12 Ah 8.9 kg (19.6 lbs) 17 Ah 13.8 kg (30.4 lbs)
Battery Charger	2-amp off-board
Colour	Blue & red

Features

- Compact, foldable seat with armrests
- Easy to use wrap around delta tiller
- Feather-touch disassembly into 5 pieces
- Comfort-Trac Suspension
- High level tiller charger
- LED lighting
- Multiple colour choices including optional colour shrouds

Also available in:

- Red
- Blue



¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

² Includes standard seat weight.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.