



COBI REHAB®
The Positioning Company

Bariatric Shower Commode Rise N Tilt

Rise N Tilt is a combination chair with built-in lift and tilt function.

It is the world's first shower commode for bariatric users to lift and tilt the user to a 30° degree. Rise N Tilt Shower Commode has a weight capacity of 325 kg.

Making tasks easier and safer

Shower Commode Rise N Tilt provides better access for caregivers in relation to hygiene, skin and tissue care. Nursing and dressing of feet may, for caregivers, take place in a standing position, thus achieving an ergonomically sound working position. The shower commode also allows the caregiver and the user to have eye contact, which provides far better communication.

Improved and more efficient care

This shower commode raises the efficiency of care because it enables a single caregiver to perform the whole showering procedure, as well as a wide range of everyday hygiene tasks.

The high level of support provided by the shower commode results in streamlined routines that offer considerable time savings compared with traditional methods.

Shower Commode Rise N Tilt is designed to give bariatric users a comfortable and safe bathing experience and can also be used where the user's condition is poor with reduced postural control.

Enhancing quality of life

Shower Commode Rise N Tilt can enhance the quality of life for a bariatric user by making assisted showering and other hygiene routines not only safer and more comfortable, but also a more dignified experience.

- Lift and tilt function up to 30°
- Possibility of extra wide butterfly armrests
- Brakes on all 4 wheels
- Unique lift and tilt function
- Hand control gives the caretaker full control
- High backrest and sturdy armrests

MAX 325 kg



In lowest position with standard footrests

Bariatric Shower Commode Rise N Tilt Specifications

Bariatric Shower Commode Rise N Tilt	
Maximum weight capacity	325 kg
Seat depth adjustable	42 / 48 / 54 cm
Total length (lowest position)	140 cm
Total width	74 cm
Seat width	66 cm
Seat height	51-69 cm
Electric tilt	0-30°
Tilt of seat independent from height adjustment	Yes
Width between standard armrests	61 cm
Width between butterfly armrests	76 cm
Width toilet opening	Diameter 24 cm, 35 cm long
Seat height electric adjustable	55-100 cm
Distance from floor to seat underside	44,5 cm
Pushbar from floor	100 cm
Seat to footrest	39-42-45-48 cm
Height seat to armrest	22 cm
Weight	46 kg
Swing-away armrests	Yes
Swing-away footrests	Yes
Removable footrests	Yes
Height adjustable footrests	Yes
Headrest	Yes
4 wheels with individual locks	Yes
Control	Via the handset
Capacity up/down fully charged battery	Approx. 15 times
Turning radius	156 cm
Frame	Powder lacquered chromium molybdenum steel
Wheels	Diameter 150 mm
Seat / back / armrests	Polyurethane
IP - classification	IPX6
IPX6	Washable, also in washing tunnel, 20°C IEC 60601-2-52
ISO - classification	ISO 09 12 03 - 23

Ideal use for.

- Bariatric users.
- Users with reduced muscle tone, pressure sores, spasm or amputations.



Commode bucket with lid



Seat with lid



Split toilet seat

Model	Seat width cm	Seat depth cm	Order No
XXL-Rehab Shower Commode Rise N Tilt incl. standard footrests, charger, headrest and commode bucket	66	42 / 48 / 54	0157-066-000
Accessories			
Butterfly armrest right (add 14 cm to seat width)			0150-999-012R
Butterfly armrest left (add 14 cm to seat width)			0150-999-012L
Headrest adjustable			0156-999-020
Full size seat for shower commode with lid for bucket			0150-999-023
Commode bucket			0140-999-001
Standard footrest right			0156-999-090H
Standard footrest left			0156-999-090V
Calf support / amputation support right			0156-999-030
Calf support / amputation support left			0156-999-031
Battery for XXL-Rehab Rise N Tilt			0157-999-050

The handset gives the caregiver control over battery-powered functions

