

Boogie Fix & Boogie Swing



Standing device for children and adolescents (User height from 80 up to 150cm)

Boogie Fix (Center column fixed)

Boogie Swing (Center column tiltable forward by gas pressure spring)

Included in the base price:

- Standard knee truss pads with kneecap cut out
- Standard therapy table with framed edge
- Posterior and chest truss pads
- Standard foot plates

Delivery address / Company stamp:		Billing address: <input type="checkbox"/> same as delivery address	
Name:		Name:	
Address:		Address:	

Ordering information:	
Reference/ internal description:	Date:
Contact person:	Signature:

**Quotation
Order**

Weight capacity: 60 kg


Prices excl. VAT, valid from 01.03.2021

Base price from £2,375

**Please choose
Model:**



<input type="checkbox"/> Boogie Fix German aid index N° without casters 28.29.01.0xxx with casters 28.29.01.1xxx Base price from £2,375	<input type="checkbox"/> Boogie Swing German aid index N° without casters 28.29.01.0xxx with casters 28.29.01.1xxx Base price from £2,599
--	--



We mount the mobile standing device according to your information			Please do not use measuring tape for your measurements, but rather a meter. Please measure:
A	body height:	_____cm	from the shoe sole to the parting
B	knee height:	_____cm	from the shoe sole to the middle of the kneecap
C	knee width:	_____cm	the widest part of the knee
D	knee distance:	_____cm	the distance between the middle of both kneecaps
E	crotch height:	_____cm	from the shoe sole to the crotch
F	pelvis width:	_____cm	the widest part of the pelvis (trochanter major)
G	sternum height:	_____cm	from the shoe sole to the middle of the sternum
H	chest width:	_____cm	the width of the ribcage at the height of the armpits

Boogie Fix



Boogie Swing







Measurements base models Boogie Fix / Swing

	Size 1	Size 2	Size 3
Overall length (without therapy table)	67 cm	73 cm	80 cm
Overall width	54 cm	58 cm	64 cm
Table height adjustment from foot plate	61 - 81 cm	75 - 95 cm	88 - 108 cm
Maximum user size	80 - 110 cm	100 - 130 cm	120 - 150 cm
Maximum payload	60 kg		

Frame sizes					
Size 1		Size 2		Size 3	
<input type="checkbox"/> for body height 80-110 cm, Measures (LxW): 67x54cm		<input type="checkbox"/> for body height 100-130 cm, Measures (LxW): 73x58 cm		<input type="checkbox"/> for body height 120-150 cm, Measures (LxW): 80x64cm	
Tilt adjustment if Boogie Swing is ticked: 30° to the front with gas pressure spring		Tilt adjustment if Boogie Swing is ticked: 30° to the front with gas pressure spring		Tilt adjustment if Boogie Swing is ticked: 30° to the front with gas pressure spring	
Boogie Fix					
N° 5911		£2,375	N° 5912		£2,500
N° 5913		£2,600			
Boogie Swing					
N° 5914		£2,599	N° 5915		£2,700
N° 5916		£2,800			

Frame colours					
Coating options		Slightly metallic (except <i>Black & King's Blue</i>)		Special colours (without additional charge, except for Nightblue)	
<input type="checkbox"/> shiny		<input type="checkbox"/> Creme White (17)		strong metallic effect:	
		<input type="checkbox"/> Zinc Yellow(02)		<input type="checkbox"/> Silver Metallic (36)	
		<input type="checkbox"/> Chrome Yellow (06)		<input type="checkbox"/> Anthracite Metallic (22)	
		<input type="checkbox"/> Pure Orange (21)		coarse structured surface:	
<input type="checkbox"/> dull (not possible with special colours!)		<input type="checkbox"/> Crimson (04)		<input type="checkbox"/> Coarse Black (27)	
		<input type="checkbox"/> Girly Pink (19)		Dormant colours:	
		<input type="checkbox"/> Pastel Blue (20)		<input type="checkbox"/> Apple Green (18)	
		not metallic:		<input type="checkbox"/> Firered (45)	
		<input type="checkbox"/> Black (31)		<input type="checkbox"/> Violet (50)	
		<input type="checkbox"/> King's Blue (32)		Dormant colour with AC:	
N° 1561				<input type="checkbox"/> Nightblue (47) £65	
£ 0					
Truss pad covers (synthetic leather)					
Colours (N° 15..)		<input type="checkbox"/> Black (90) <input type="checkbox"/> Firered (92) <input type="checkbox"/> Violet (93) <input type="checkbox"/> Turquoise (94) <input type="checkbox"/> Powder Blue (95) <input type="checkbox"/> Lime (96) <input type="checkbox"/> Green (97)			

Casters	
	<input type="checkbox"/> Standing feet with transport casters
	N° 286101 £ 0
	<input type="checkbox"/> 4 Steerable and braked casters Ø 75 mm
	N° 286102 £145

Chest truss pad			
	Chest truss pad stepless adjustment in height and in tilting angle. Incl. safety belt		
	<input type="checkbox"/> Size 1 N° 683107		
	Padsize: 20x13 cm (WxH) recommended for breast width 16-24 cm		
	<input type="checkbox"/> Size 2 N° 683108		
	Padsize: 26x17 cm (WxH) recommended for breast width 22-32 cm		
	<input type="checkbox"/> Size 3 N° 683109		
	Padsize: 32x21cm (WxH) recommended for breast width 28-42 cm		
	£ 0		
	<input type="checkbox"/> Side chest truss pads stepless adjustment in width and depth. Incl. safety belt		
	Depending on selected chest truss pad		
		for chest truss pad size 1 N° 684307, Padsize 9x9cm(WxH)	
		for chest truss pad size 2 N° 684308, Padsize 12x12cm(WxH)	
		for chest truss pad size 3 N° 684309, Padsize 12x12cm(WxH)	
		£215	

Posterior truss pad standard		
	Posterior truss pad stepless height adjustment, incl. crotch truss pad	Side posterior truss pads stepless adjustment in width and depth. incl. belt
	<input type="checkbox"/> Size 1 pelvis depth 13-22cm N° 684201, padsize 16x16cm (WxH) <input type="checkbox"/> Size 2 pelvis depth 18-31,5cm N° 684202, padsize 22x18cm (WxH) <input type="checkbox"/> Size 3 pelvis depth 18-31,5cm N° 684203, padsize 30x19,5cm (WxH)	<input type="checkbox"/> Size 1 pelvis width 16-24 cm N° 684211, padsize 9x9cm (WxH) <input type="checkbox"/> Size 2 pelvis width 22-32 cm N° 684212, padsize 12x12cm (WxH) <input type="checkbox"/> Size 3 pelvis width 28-42 cm N° 684213, padsize 12x12cm (WxH)
£ 0		£215

Posterior truss pad with swivel bracket		
	Posterior truss pad with swivel bracket, for stepless raising in the pelvic area, crank operable from the front and removable, incl. side posterior truss pads.	<input type="checkbox"/> Size 1 N° 684261 pelvis width 18-26cm pelvis depth 13-24cm posterior truss pad: 16x16cm (WxH) side posterior truss pads: 9x9cm (WxH) <input type="checkbox"/> Size 2 N° 684262 pelvis width 23-36 cm pelvis depth 18-32cm posterior truss pad: 22x18cm (WxH) side posterior truss pads: 12x12cm (WxH) <input type="checkbox"/> Size 3 N° 684263 pelvis width 29-42 cm pelvis depth 18-32cm posterior truss pad: 30x19,5cm (WxH) side posterior truss pads: 12x12cm (WxH)
	<p><i>The side pelvic truss pads are not possible.</i></p> <p><i>Pelvis width = Distance between pad and traverse</i></p> <p><i>Depth adjustment: max. 8 cm</i></p>	£360

Posterior truss pad with swivel bracket and crank		
	Posterior truss pad with crank and swivel bracket, for stepless raising in the pelvic area, crank operable from the front and removable, incl. side posterior truss pads.	<input type="checkbox"/> Size 1 N° 684221 pelvis width 18-26cm pelvis depth 13-24cm posterior truss pad: 16x16cm (WxH) side posterior truss pads: 9x9cm (WxH) <input type="checkbox"/> Size 2 N° 684232 pelvis width 23-36 cm pelvis depth 18-32cm posterior truss pad: 22x18cm (WxH) side posterior truss pads: 12x12cm (WxH) <input type="checkbox"/> Size 3 N° 684242 pelvis width 29-42 cm pelvis depth 18-32cm posterior truss pad: 30x19,5cm (WxH) side posterior truss pads: 12x12cm (WxH)
	<p><i>The side pelvic truss pads are not possible.</i></p> <p><i>Pelvis width = Distance between pad and traverse</i></p> <p><i>Depth adjustment with crank: max. 10 cm</i></p>	£420

Side pelvis truss pad and back truss pad		
	Side pelvis truss pad stepless adjustment in width, height and depth	Back truss pad with swivel bracket
	<input type="checkbox"/> Size 1 N° 684401 pelvis width 16-24cm Padsize 9x9cm (WxH) <input type="checkbox"/> Size 2 N° 684402 pelvis width 22-32cm Padsize 12x12cm (WxH) <input type="checkbox"/> Size 3 N° 684403 pelvis width 28-42cm Padsize 12x12cm (WxH)	<input type="checkbox"/> Size 1 N° ..05 Padsize 12x16 cm (WxH) <input type="checkbox"/> Size 2 N° ..06 Padsize 16x19 cm (WxH)
<p><i>Not possible in combination with posterior truss pad with crank OR side posterior truss pads</i></p>		<p>Only in combination with side chest truss pads</p>
£225		N° 6848.. £235

Spine truss pad			
	Spine truss pad two-part, for mounting on pelvis truss pad or posterior truss pad with crank. <i>Selection of pelvis truss pad or posterior truss pad is necessary.</i>		
	<input type="checkbox"/> Spine truss pad size 1 (...01) Pad size 9x10cm only with selection pelvis-/posterior truss pad size 1 <input type="checkbox"/> Spine truss pad size 2 (...02) Pad size 13x13cm only with selection pelvis-/posterior truss pad size 2 or 3 Only in combination with compatible pelvis truss pad or posterior truss pad with crank		
	N° 6847..		£175
Knee truss pads			
	Standard knee truss pad with kneecap cut-out, stepless, 3-dimensional setting, individual width setting.		
	<input type="checkbox"/> Size 1 (...01) knee width 4,5-7,5 cm <input type="checkbox"/> Size 2 (...02) knee width 6-11 cm		
	N° 683601/02	£ 0	
	Knee truss pad with resting brackets individual, stepless, 3-dimensional adjustment, also recommended for Foot plate 3-dimensional angle (N° 408500)		
	<input type="checkbox"/> Size 1 (...11) knee width 6,5 cm <input type="checkbox"/> Size 2 (...12) knee width 8,5 cm <input type="checkbox"/> Size 3 (...13) knee width 10,5 cm		
	N° 6836..		£125
Foot plate			
	Standard foot plates <input type="checkbox"/> Divided standard (...400) individually abductable, adjustable in height and width and can be turned horizontally. <input type="checkbox"/> Undivided Standard (...600) adjustable in height		
	N° 408...	£ 0	
	<input type="checkbox"/> Foot plate 3-dimensional angle adjustment can be set individually in all directions, as well as in the tilting angle, the abduction, height and width and can be turned horizontally. The maximum possible body size is reduced by 3 cm.		
	N° 408500		£225
Therapy table (absolutely necessary for the Boogie Swing)			
	<input type="checkbox"/> Therapy table standard with framed edge Beechwood Multiplex Size depends on selected chest truss pad for chest truss pad size 1 N° 921021 for chest truss pad size 2 N° 921022 for chest truss pad size 3 N° 921023		
		£ 0	
	<input type="checkbox"/> Therapy table with therapy box and lid, with framed edge Beechwood Multiplex with framed edge Size depends on selected chest truss pad for chest truss pad size 1 N° 921031 for chest truss pad size 2 N° 921032 for chest truss pad size 3 N° 921033		
			£175
Therapy table accessories			
	<input type="checkbox"/> Upholstery for therapy table made of polyethylene foam for therapy table Standard (...01) for therapy table with toy box (...02)		
	N° 9210..	£110	
	<input type="checkbox"/> Arm rest adjustable arm rest for chosen therapy table		
	N° 921010	£135	
	<input type="checkbox"/> Upholstery for arm rest made of polyethylene foam		
	N° 921011	£60	
Communications helper Topper Orderable via www.nickylinedesign.de (not in combination with toy box or upholstery for therapy table)			
example from nickylinedesign.de			
Options			
	<input type="checkbox"/> clamp lever Instead of screws, for easier adjustment of the truss pads, especially useful for frequent adjusting of the mobile standing device.		
	N° 570210	£75	
	Foot shell Shell for max. shoe width <input type="checkbox"/> size 1 (...11) 6,5 cm <input type="checkbox"/> size 2 (...12) 8,5 cm <input type="checkbox"/> size 3 (...13) 9,5 cm		
	N° 4020..	£85	
	<input type="checkbox"/> Belts for foot shells Matching belts for your chosen foot shells		
	N° 402015	£95	
<input type="checkbox"/> Leg length compensation for divided standard foot plates (N° 408400) Foot plate increased on one side: Side: <input type="checkbox"/> left <input type="checkbox"/> right Difference: <input type="checkbox"/> 1 cm <input type="checkbox"/> 1,5 cm <input type="checkbox"/> 2 cm <input type="checkbox"/> 2,5 cm <input type="checkbox"/> 3 cm			N° 408410 £28/Set

Terms and conditions

1. Preface and Scope

(1) Our terms and conditions apply exclusively. Any deviating, opposing, or complementing terms and conditions on the part of the business partner do not form integral part of contract unless explicitly approved of by us. There are no exceptions to the consent requirement, which also applies in cases where we conduct deliveries without restrictions while aware of the business partner's terms and conditions.

(2) Agreements made in the individual case with the purchaser (including side-agreements, supplements, and alterations) always take priority over these terms and conditions. A contract or our confirmation - both in written form - are authoritative for the content of such arrangements.

(3) Any legally relevant declarations and notifications the purchaser submits to us after conclusion of the contract (e.g. setting of deadlines, notice of defects, declarations of rescission or abatement) require the written form to be valid.

(4) Any information on the application of legal regulations serves the mere purpose of clarification. All regulations are legally binding without any such clarifications unless explicitly altered or excluded in these terms and conditions.

2. Prices

All prices are ex works, excluding insurance and packaging. The prices are net prices; the value-added tax is stated and charged separately according to the legally binding rate at the time of invoicing. The prices valid on the day of delivery are binding unless explicitly arranged otherwise.

3. Conclusion of Contract

SORG's offers are invitatio ad offerendum (invitations to treat). Ordering goods counts as an obligatory offer of contract on the part of the purchaser. Unless the order suggests otherwise, we have the right to accept this contract offer within two weeks time after reception. The acceptance of an offer of contract is effected by the confirmation of the order.

4. Terms of Payment

(1) Provided that no other agreements exist, payment is to be made by paying in advance - 50% at placing of order, 50% pre-delivery of goods. The purchaser bears any costs arising from transactions or discounting.

(2) If a 30-day payment deadline is exceeded, a delay rate of 9% on top of the respective base rate of the European Central Bank falls due without the need for prior notices. Beyond that, the purchaser is charged 3,- € for any reminder sent. We reserve the right to make use of the commercial due rate according to § 353 HGB. A higher damage for delay may be enforced via a suitable verification.

5. Offers and Validity

All current price lists are generally valid until December 31 of the respective year, or until new price lists are published. The latest version of the price list will always be provided on our homepage. Offers older than 3 months at that point expire and must be replaced with new offers and prices.

6. Time of Delivery and Delay

We state the terms of delivery when accepting the order. The delivery date is non-binding and has mere guideline value - especially with regard to customised products. The terms for delay in delivery are defined by legal regulations. In any case, a writing notice issued by the purchaser is necessary for them to take effect.

7. Force Majeure (Act of God)

(1) SORG is not to be held responsible for the following delays of delivery. This also applies if these causes occur with respect to legal representatives, auxiliary persons, or deliverers of SORG: Circumstances of force majeure as well as other extraordinary events that are unpredictable and unavoidable for, and not caused by, SORG, and that only occur after ordering or reserving, and that SORG could not have known about. Moreover, this applies in case of strikes, fire, floods, labour disputes, business disruptions, altered magisterial approbation and legal situation, and magisterial orders which are not due to operational risks.

(2) In the cases of paragraph 1, SORG is eligible to delay delivery/performance for the time of hindrance as well as an appropriate start-up time. If SORG thus causes unacceptable complications for the purchaser, both SORG and the purchaser are eligible to withdraw from the contract. Circumstances that only lead to a temporal and hence acceptable delay of delivery and performance are excluded from this right. The legal rights of the client remain unaffected.

8. Dispatch

We generally dispatch uncovered by insurance and at the recipient's risk. We transfer the risk of coincidental loss or deterioration of the goods as well as delay on to the client at the moment of passing the goods on to the forwarder, bearer, or any other person/institute responsible for delivery. SORG will make any effort to respect the client's wishes in terms of path and kind of delivery. Any costs caused by this (including agreed carriage free delivery) are to be covered by the client.

In the case of carriage free delivery conducted by SORG, the recipient is required to put forward any carriage costs for a refund. If the client wishes the dispatch to be delayed, any risk of loss or deterioration of the goods is carried by the recipient from the first day of readiness for dispatch. The goods dispatched are only insured upon special demand and at the cost of the recipient.

9. Damage through Delivery Angelieferte Gegenstände sind, auch wenn sie unwesentliche Mängel. Upon delivery, the goods need to be checked for obvious damage in the presence of the bearer. Any damage is to be documented in a written form and co-signed by the bearer. This protocol is to be forwarded to us immediately. Delayed reclamations, or reclamations not signed by the bearer, cannot be taken into account.

10. Warranted Features, Receipt, and Technical Alterations

Delivered objects, even with minor damage, are to be accepted by the client, even if § 433 paragraph 1 BGB is not fulfilled. Partial deliveries are permissible. Features only count as warranted if explicitly designated as such. A warranty of condition can only be asserted explicitly and in a written form. SORG reserves the right for technical alterations, especially for solutions of equal or higher value.

11. Warranty for Defects

(1) SORG warrants for any defects of delivered new goods for 12 months as long as the purchaser is not the user (§§ 474, 13 BGB). Should the purchaser be the user then the warranty period is two years for new goods, for used goods 12 months. Should the purchaser also be an entrepreneur, then the liability for used items is void, as long as there is no claim for compensation according to N° 12.

(2) The purchaser's claim for defects presupposes that he has attended his obligation to give notice of defects on the legal grounds for his reasons of inquiry (§§ 377, 381 HGB). If defects are found during the inspection or at a later point, we are to be notified in writing without delay. The notice counts as not delayed if issued within two weeks. Posting the notice on time suffices for the keeping of the term. Independently from this obligation to inspect and give notice, the purchaser is to report obvious defects (including wrong and short delivery) within two weeks upon delivery (again, posting the notice on time suffices). If the purchaser misses the proper inspection and/or notice of defects,

we do not warrant for the reported defects. If the delivered good is defected, we may choose to conduct supplementary performance either by repairing the defect (subsequent improvement) or by delivering flawless goods (compensation delivery). Our right to refuse the chosen method of supplementary performance on the respective legal grounds remains untouched.

(3) Our liability for defects expires when alteration or maintenance works are done by the client or a third-person-party without our prior written permission. Consequences of excessive strain, negligent or incorrect handling, and violent damage are excluded from the warranty for defects.

12. Further Liability

(1) With respect to compensation - regardless of the legal cause - SORG is only liable in case of malice and act of gross negligence. In case of slight negligence, SORG is only liable for damage of life, body, health, and for damage caused by the negligence of an essential contractual obligation (an obligation whose fulfilment make the proper realisation of the contract possible in the first place and in whose fulfilment the contract partner regularly trusts and can trust). In this case, however, our liability does not exceed the covering of predictable, typically occurring damage.

(2) A free right to termination on the part of the purchaser (acc. to § 651, 649 BGB in particular) is precluded. Otherwise, the usual prerequisites and legal consequences apply.

13. Costs of Supplementary Performance

Costs of supplementary performance cannot be charged on SORG's account.

14. Voluntary Withdrawal of Goods

SORG is not obliged to take back any delivered end products or parts thereof. If SORG voluntarily decides otherwise in individual cases (especially for therapeutically altered requirements or deaths), SORG refunds a maximum of 75% of the delivery costs for unused parts whose delivery does not date back more than three months. There is no legal obligation concerning withdrawal in this respect. Deliveries which date back more than three months as well as sewed parts and customised products are not returnable. A copy of the delivery note or the bill serves as proof of age of delivery.

15. Reservation of Ownership

Goods delivered by SORG remain the property of SORG until any obligations from the parties' business contract on the part of the purchaser have been redeemed. This also applies if SORG has adapted individual or all claims into a running account which has been balanced and acknowledged. The purchaser has the right to dispose of the goods in the proper course of business. The purchaser hereby yields all claims to SORG which arise from the resale against the purchaser or third parties. The purchaser is authorised to withdraw these claims even upon conveyance. SORG's right to withdraw the claims remains untouched. However, SORG is obliged not to withdraw the claims as long as the purchaser complies with his duty of payment and any other duties. SORG has the right to demand that the purchaser makes known the yielded claims and their debtors, that he provides all information necessary for conveyance, and that he hands in all associated documents. As long and insofar the reservation of ownership is intact, the purchaser must not transfer any goods for securing or bailment without SORG's permission. The purchaser may process or alter the goods for the benefit of SORG exclusively. If processed with other goods which do not belong to SORG, SORG obtains co-ownership of these new goods in relation to the invoice value to the cost price of the other processed goods. If a third party interferes with SORG's ownership, the purchaser is to notify SORG immediately. The purchaser must not make any agreements with other purchasers which may exclude or limit SORG's rights. In particular, the purchaser must not make any agreements which destroy or limit the advance assignments to SORG. If the estimated value of the existing guarantees exceeds the claims to be ensured by more than 20%, SORG is insofar obliged to release by request of the purchaser.

16. Copyright Protection

If SORG conveys any images, drawings, indications of weight or measurements, technical documentations, or other product descriptions and documents (including electronic ones), the respective indications only apply if explicitly stated that they do. SORG reserves ownership and copyright of all these documents. They may only be made available to third parties upon SORG's explicit permission.

17. Place of Fulfilment, Place of Jurisdiction, and Severability Clause

(1) The law of the Federal Republic of Germany applies for these Terms and Conditions and all legal relations between us and our business partners. All international and supranational (contractual) legal systems, the UN Convention on Contracts for the International Sale of Goods in particular, are hereby excluded. Prerequisites and effects of the reservation of ownership acc. to no. 15 are subject to the right of the respective storage location of the good if the choice of law made in favour of the German law is illegitimate or ineffective. The only place of jurisdiction (also on an international basis) for any disputes arising directly or indirectly from the contractual relationship is, according to German Commercial Code, our office in Oberhausen-Rheinhausen, Germany. However, we also have the right to file a suit concerning the general jurisdiction of the purchaser.

(2) If any clause in these Terms of Conditions should turn out to be ineffective, this does not affect the effectiveness of all other clauses.

Technical alterations and misprints reserved. As of: January 2017

SORG Rollstuhltechnik GmbH + Co. KG
Benzstr. 3-5, 68794 Oberhausen-Rheinhausen / Germany
Phone +49 7254 / 9279-0
Fax +49 7254 / 9279-10
info@sorgrollstuhltechnik.de, www.sorgrollstuhltechnik.de