

BRAVOracer

The Bravo Racer is a rigid frame wheelchair with various positioning aids for children with severe disabilities. The highly stable frame with the wide extension arms provides maximum space for hip abduction.

What is special about me?

- Stable, torsion-resistant design for outstanding driving manoeuvrability
- Frame design allows for easy standing transfer into and out of the wheelchair while standing on the footrest and offers a wide hip abduction
- Grows with the child: setting options for seat height, seat width, back position, back angle and centre of gravity.

Standard Model

- Upholstery colours: red, black, blue
- Aluminium frame with rigid chassis
- Front frame with wide abduction
- Angle-adjustable and foldable backrest, adaptable back upholstery
- 20", 22" and 24" rear wheels for all frame lengths
- Wheel lock lever integrated into the side panel

Options

- Height-adjustable back upholstery, short push handle
- Height-adjustable back upholstery, long push handle
- Push bar
- Frame bar
- Adjustable abduction on side panel
- Wheelbase extension
- Luminous caster wheel
- Anti-tipper
- Spoke protectors, various designs
- Schwalbe RightRun light
- Tire-a-light
- Designer clothing protector
- Angle-adjustable footrest
- Non-skid footrest insert

Quality for life



Technical Data

Max. user weight	60kg
Seat width	20 – 36cm
Seat depth	24 – 38cm
Seat height	32 – 51cm
Total width	49 – 65cm
Total length	61.5 – 81.5cm
Back height	20 – 40cm
Lower leg length	15 – 45cm
Armrests	Optional
Back angle	-15° / +15°
Weight	From 9kg