

Date

User name

Company stamp

Responsible

Rigid frame wheelchair

£5,900

Carbon fibre frame

Wheelchair measures

Seat width (LS) ☐ 33 ☐ 36 ☐ 39 ☐ 42 ☐ 45 ☐ 48

The seat width is measured from the outer sides of the frame, while the distance between the two side guards is 1 cm wider.

Seat depth (PS) ☐ 35 ☐ 37,5 ☐ 40 ☐ 42,5 ☐ 45 ☐ 47,5 ☐ 50

The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.

Footplate distance (DP) (Real distance between seat and footplate)

IMPORTANT NOTICE! From the footplate to the floor minimum 4cm.!! (Minimum Front height = DP + 4cm.)

Setting ☒ Standard ☐ Prudential ☐ Active ☐ Extreme

Backrest Height (TS) Std.Value 33,5(min 24.5-max 47)

camber ☒ 0° ☐ 2° ☐ 4°


Front height STD 49

Rear height STD 42

24" Rear wheel (min 38-max 46)

25" Rear wheel (min 39-max 47)

26" Rear wheel (min 40-max 48)



☒ Version mudguard (Available only with rear wheels 24")

☐ Version straight sideguard

Seat

☒ Seat cloth with pocket (under seat)

☐ Seat with tension straps 08883 **£ 196**

☐ Carbon rigid seat 08336 **£ 415**

☐ Pocket (under seat) [25x11x4,5] 08281 **£ 65**

STYLE

☒ CARBON Essence


☐ ALA ROSSA 07947 **£ 109**

☐ ALA BIANCA 07948 **£ 109**

Footplate

☒ Carbon fibre footplate - adjustable distance +/- 1,5cm. SPECIAL FLIP-UP!

Minimum footplate distance 35cm.



English
Date Order User name Address

Backrest with rake adjustable from -4° to +16°

- ☆ Backrest with double locking (upright and folded)
- ☐ Backrest with upright locking system ONLY (no locking when folded)
- ☐ Fixed back(non-folding back)

☆ Aluminium backrest

☐ Titanium backrest

☐ Carbon fibre backrest



06811 £ 164

(Lightest configuration)

04389 £ 655

- ☐ **A I R - Integrated posture backrest (not removable) - Lateral retaining 6,5cm.** 06815 £ 1015
not available for seat width 48 - Std. Backrest upholstery not provided - Handle tube aluminium with short tubes

Lateral retaining wing B=low

Carbon shell H25



Backrest Height (TS) min 29 max 41

Carbon shell H33



Backrest Height (TS) (37 + 46)

Lateral retaining wing A=high

Carbon shell H30



Backrest Height (TS) (32 + 42,5)

AIR Cover

- ☆ Cover Standard - fabric insert: GREY, stitchings GREY

☐ Cover SPORT


09195 £ 87

- ☐ fabric insert: BLACK - Stitchings GREY
- ☐ fabric insert: BLACK - Stitchings LIME
- ☐ fabric insert: BLACK - Stitchings ORANGE
- ☐ fabric insert: BLACK - Stitchings RED
- ☐ fabric insert: BLACK - Stitchings BLUE
- ☐ fabric insert: BLACK - Stitchings BLACK
- ☐ fabric insert: BLACK - Stitchings FUCHSIA
- ☐ fabric insert: BLACK - Stitchings ANTHRACITE

☐ Cover ACTIVE

09195 £ 87

- ☐ fabric insert: RED - Stitchings BLACK
- ☐ fabric insert: LIME - Stitchings BLACK
- ☐ fabric insert: BLUE - Stitchings BLACK



Transpiring fabric stitchings

Handle tube

☆ Handle tube STD

- ☐ Handle tube with short tubes
- ☐ Handle tube SPORT Without push handles

☐ Carbon fibre handle (with short tubes)

(Lightest configuration)

04390 £ 131

Not available with: Posture backrest PHYSIO - Posture backrest A I R

☐ Folding handles aluminium



02958 £ 218

Casters

☆ 4" Soft

☐ 3" Sport

(Lightest configuration)

☐ 5" Soft

☐ 3" Active Soft (wide)

(Aluminium alloy rim)

06824 £ 109

☐ 4" Active Sport (narrow)

(Aluminium alloy rim)

01398 £ 109

☐ 4" Active Soft (wide)

(Aluminium alloy rim)

03943 £ 109

☐ 5" Active Soft (wide)

(Aluminium alloy rim)

07518 £ 109

English

Date Order User name Address

Rear wheels

Available with : VERSION MUDGUARD - VERSION STRAIGHT SIDEGUARD

☆ 24" Ultra-light Progeo - Standard!

☐ 24" ULTRA-LIGHT PROGEO KEY SPORT	(Lightest configuration)	04383	£ 764
☐ 24" Ultra-light "Spinergy LX"		01685	£ 764
☐ 24" Ultra-light "Spinergy LXL"	☐ Black hub and white blade spokes ☐ Hub and blade spokes black	07846	£ 764
☐ 24" Off-road (aluminium handrim black)		01392	£ 196

Available only with: VERSION STRAIGHT SIDEGUARD

☐ 22" Ultra-light Progeo - Standard!			
☐ 25" Ultra-light Progeo - Standard!	(colour: black)		
☐ 26" Progeo - Standard!	(colour: black)		
☐ 25" Ultra-light "Spinergy LX"		02101	£ 764
☐ 25" Ultra-light "Spinergy LXL"	☐ Black hub and white blade spokes	07847	£ 764
☐ 26" Ultra-light "Spinergy LX"		02102	£ 764

Rear wheel tyre

☆ High pressure tyre	(22"-24"-25"-26")	(Lightest configuration)		
☐ High pressure profile tyre	(24")			
☐ Marathon Plus	(24"-25")	with anti-puncture protection	01669	£ 55
☐ Anti-puncture - Colour Grey	(22")		04181	£ 87
☐ Solid tyre Kik Mako - Colour Black	(24")		01041	£ 109
☐ Solid tyre Kik Mako - Colour Grey	(24")		01751	£ 109

Handrim - Wheel distance

☆ 3 cm

☐ 4 cm

☆ Aluminium	(22"-24"-25"-26")	(Lightest configuration)		
☐ Aluminium black SOFT-TOUCH	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	09170	£ 109
☐ Chromed Steel	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	01684	£ 131
☐ Titanium	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00016	£ 240
☐ Ultra Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00017	£ 273
☐ Max Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00834	£ 305
☐ Ergo Para	(24")		00632	£ 305
☐ SURGE LT	(24"-25")	Oval countour with a rubber insert	06427	£ 316
☐ Quad (24")	(24")		00693	£ 175
☐ Silicon cover for standard aluminium handrim		Available for measures (22"- 24" - 26")	04050	£ 142

Wheel locks

☆ Without brake	(Lightest configuration)			
☐ Bent lever push to lock				
☐ Straight lever Push to lock				
☐ Straight lever Pull to lock				
☐ Sport brake			06428	£ 109
☐ Brake with extension lever push to lock			01949	£ 55
☐ Brake with extension lever pull to lock			03372	£ 55

Armrests

☐ Swing away armrests	(Swing away, removable and height adjustable armrest)		00694	£ 218
	☆ Pad STANDARD			
	☐ Ergonomic pad (30x6 cm.)		04104	£ 55

Accessories

☐ Pad for front frame protectors		07058	£ 44
☐ Backrest-seat connection cloth		10241	£ 16
☐ Removable storage bag (with quick coupling)		09811	£ 109

Date Order English User name Address

Optional 1

- ☐ Transit wheels  03290 £ 131
- ☐ Anti-tipper with rapid hookup (inside rotation to reduce the encumbrance)
- ☐ Right  06821 £ 109
- ☐ Left 06822 £ 109

Optional 2

- ☐ Lumbar support (to insert in the backrest upholstery) 00013 £ 33
- ☐ Tetraclip (Quick release axle for quadriplegic) 03947 £ 153

Additional Kit

- ☐ 24" Off-road rear wheels Off-Road tyre - aluminium handrim (black colour) 03378 £ 436
- ☐ 24" FAT rear wheels FAT tyre 3" - aluminium handrim (black colour) 09357 £ 764

PHYSIO - Posture backrest

Std. Backrest upholstery not provided

- ☐ Aluminum shell - Lateral retaining STANDARD 7,5cm
- Height ☐ 25 ☐ 33 ☐ 41 05308 £ 884
- Height ☐ 52 05309 £ 938
-  Lateral retaining wing ☐ A=high ☐ B=low
- ☐ Aluminum shell - Lateral retaining DEEP 15cm
- Height ☐ 25 ☐ 33 ☐ 41 05310 £ 938
- H25 not available for seat width 39-42-45-48
- Height ☐ 52 05311 £ 993
-  Lateral retaining wing ☐ A=high ☐ B=low
- ☐ Carbon shell - Lateral retaining STANDARD 7,5cm
- Height ☐ 25 ☐ 33 ☐ 41 05312 £ 1015
-  Lateral retaining wing ☐ A=high ☐ B=low


COVER

☆ Cover with fabric-cover of backrest's rear side and high refraction fabric hems




- ☐ SWING AWAY - Lateral support
- Not available with Physio backrest with lateral retaining DEEP and in carbon backrest

Internal thoracic support to the backrest width, adjustable (0 ÷ -1,5cm)




- Large (17,5x11,5)
- ☐ Left 05646 ☐ Right 05647
- Small (14x9,5)
- ☐ Left 05642 ☐ Right 05643
- ☐ Additional Kit for internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)

Internal thoracic support to the backrest width, adjustable (-4 ÷ -5,5cm)



- Large (17,5x11,5)
- ☐ Left 05646 ☐ Right 05647
- Small (14x9,5)
- ☐ Left 05642 ☐ Right 05643

Internal thoracic support to the backrest width, adjustable (-8,5 ÷ -10cm)



- Large (17,5x11,5)
- ☐ Left 05646 ☐ Right 05647 £ 175
- Small (14x9,5)
- ☐ Left 05642 ☐ Right 05643 £ 175
- (with additional kit)

- ☐ Headrest adjustable - Double ball mounting system
- Not available for Physio backrest's height 25 and with carbon shell

- ☐ Large (34x14)  05648 £ 393
- ☐ Small (28x14)

English

Date Order User name Address

Rehateam s.r.l. Vicolo Negrelli, 5 - 31040 Castagnole di Paese (TV) - Tel. +39.0422.484657 / +39.0422.480387 - Fax +39.0422.484684
P.IVA IT 03182520266 - C.C.I.A.A. TV n.225867 - Iscr.Trib. TV n.41690 Capitale Sociale I.V. Euro 100.000,00

Email: info@rehateamprogeo.com - web: http://www.rehateamprogeo.com

REFERENCE TABLES

1	Rear wheels	
	Model	Features
	24" Ultra-light Progeo - Standard!	(available for all coulors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white blade spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)

Italia Certified company normative UNI EN ISO 9001

Important notice! s page only if any optional have been chosen

Pag. 5 / 5

Wheelchair : CARBOMAX - Rev. 2020.3

Prices in EURO - VAT excluded
to be added to the price list

☒ **standard** (supplied if no requests are indicated)
☐ **on request** (without extra charge)
☐ **optional** (cost to be added at the basic price)