

Embrace your potential! The Colt Deluxe 2.0 is fully loaded with style and advanced features, including frosted LED headlights and an integrated storage space with cup holder. Drive at speeds up to 6.25 mph and see what life has in store for you!





## **Specifications**

Weight Capacity	28 st 8 lbs (181 kg)
Maximum Speed¹	Up to 6.25 mph
Ground Clearance <sup>3</sup>	64 mm (2.5") at motor
Turning Radius³	1334 mm (52.5")
Overall Length <sup>3</sup>	1200 mm (47.25")
Overall Width <sup>3</sup>	563 mm (22.25")
Max. Obstacle Climb	80 mm (3.25")
Drive Tyres	269 mm front & rear - solid* or pneumatic tyres available
Suspension	Front & rear suspension
Total Weight Without Batteries <sup>3, 4</sup>	77.25 kg (170.25 lbs)
Heaviest Piece When Disassembled <sup>3</sup>	30.25 kg (66.75 lbs)
Drivetrain	24 volt DC motor
Dual Breaking System	Regenerative & electromechanical
Battery Requirements <sup>3</sup>	2 x 40 Ah deep cycle
Battery Weight <sup>4</sup>	14.75 kg (32.75 lbs)
Battery Charger	5-amp off-board
Colour	Titanium Grey, Ocean Blue, Candy Apple Red



## **Features**

- Comfortable high back seat
- Integrated USB charging point for devices
- Feather-touch disassembly
- Comfort-Trac Suspension
- High level tiller charger
- Frosted LED headlight package
- Convienient storage/cup holder in tiller

## Available in:

- Candy Apple Red
- Ocean Blue
- Titanium Grey







\* Solid tyres – available in grey only

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or – 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.



<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%

be subject to a variance of (+ or -) 10%. 2 Includes standard seat weight.

<sup>3</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

this specification can be subject to a variance of (+ or -) 3%.

4 Battery weight may vary based on battery manufacturer.