

PRIDE[®]
MOBILITY PRODUCTS LTD.

LIVE YOUR BEST[®]

Large Scooter Range 2020



“Having tried various other brands, I opted for Pride! The colt range is the most comfortable option on the market and excellent quality! Highly recommend!”





PRIDE[®]

MOBILITY PRODUCTS LTD.

LIVE YOUR BEST[®]

Pride Mobility Products is a world leader in the design, development and manufacture of mobility products - powerchairs, scooters and lift chairs - for people with disabilities and mobility impairments. With flagship product ranges such as the Jazzy Range, the British manufactured Hudson Range and our innovative Zero Turn technology – Pride Mobility are committed to helping improving people's lives through quality, innovative products.

Available on
 Motability



WEIGHT CAPACITY

25 stone (158.75 kg)

MAX SPEED

4 mph (6.4 kmh)

TURNING RADIUS

1370 mm (54")

OVERALL LENGTH

1190 mm (47")

OVERALL WIDTH

560 mm (22")

FRONT/REAR TYRES

Front: 10"
Rear: 10"

**TOTAL WEIGHT
WITHOUT BATTERIES**

9 stone 8 lbs (62.5 kg)

BATTERIES REQUIRED

(2) 32 AH (deep cycle)

25stone
Weight Capacity

Up To
4mph
Maximum Speed

1370mm
Turning Radius



Available on
 Motability

WEIGHT CAPACITY
25 stone (158.75 kg)
MAX SPEED
4 mph (6.4 kmh)
TURNING RADIUS
1067 mm (42")
OVERALL LENGTH
1160 mm (46")
OVERALL WIDTH
560 mm (22")
FRONT/REAR TYRES
Front: 9" Rear: 10"
TOTAL WEIGHT WITHOUT BATTERIES
9 stone 8 lbs (62.5 kg)
BATTERIES REQUIRED
(2) 35 AH (deep cycle)



25stone
 Weight Capacity

Up To
4mph
 Maximum Speed

1067mm
 Turning Radius

Available on
 Motability



- Red
- Matte Blue
- Grey



Up To
28stone 8lbs | **6.25mph** | **1334mm**
Weight Capacity | Maximum Speed | Turning Radius

WEIGHT CAPACITY
28 stone 8 lbs (181.4 kg)
MAX SPEED
6.25 mph (10 kmh)
TURNING RADIUS
1334 mm (52.5")
OVERALL LENGTH
1200 mm (47.25")
OVERALL WIDTH
563 mm (22.2")
FRONT/REAR TYRES
Front: 10.6" Rear: 10.6"
TOTAL WEIGHT WITHOUT BATTERIES
12 stone 1 lbs (77.2 kg)
BATTERIES REQUIRED
(2) 40 AH (deep cycle)



Available on
 Motability

WEIGHT CAPACITY
25 stone (158.75 kg)
MAX SPEED
8 mph (12.8 kmh)
TURNING RADIUS
1333 mm (52.5")
OVERALL LENGTH
1194 mm (47")
OVERALL WIDTH
657 mm (25.9")
FRONT/REAR TYRES
Front: 10" Rear: 10"
TOTAL WEIGHT WITHOUT BATTERIES
10 stone 7 lbs (66.6 kg)
BATTERIES REQUIRED
(2) 40 AH (deep cycle)



25stone
Weight Capacity

Up To
8mph
Maximum Speed

1333mm
Turning Radius



Our patent-pending, intelligent-turning *iTurn Technology™* ensures effortless navigation in small spaces and around corners.

Available on  Motability



- Red
- Blue
- White



Up To
28stone 8lbs | **6.25mph** | **1092mm**
 Weight Capacity | Maximum Speed | Turning Radius

WEIGHT CAPACITY
28 stone 8 lbs (181.4 kg)
MAX SPEED
6.25 mph (10 kmh)
TURNING RADIUS
1092 mm (43")
OVERALL LENGTH
1556 mm (61.25")
OVERALL WIDTH
613 mm (24.1")
FRONT/REAR TYRES
Front: 10.2" Rear: 10.2"
TOTAL WEIGHT WITHOUT BATTERIES
12 stone 9 lbs (82.5 kg)
BATTERIES REQUIRED
(2) 30 AH (deep cycle)

FINESSE SPORT

Available on


WEIGHT CAPACITY
28 stone 8 lbs (181.4 kg)
MAX SPEED
8 mph (12.8 kmh)
TURNING RADIUS
1566 mm (61.25")
OVERALL LENGTH
1190 mm (47")
OVERALL WIDTH
648 mm (25.5")
FRONT/REAR TYRES
Front: 10" Rear: 10"
TOTAL WEIGHT WITHOUT BATTERIES
11 stone 5 lbs (72 kg)
BATTERIES REQUIRED
(2) 50 AH (deep cycle)



28stone 8lbs
 Weight Capacity

Up To
8mph
 Maximum Speed

1566mm
 Turning Radius

Available on
 Motability



Black

Silver

25stone
Weight Capacity

Up To
8mph
Maximum Speed

1727mm
Turning Radius

WEIGHT CAPACITY

25 stone (158.75 kg)

MAX SPEED

8 mph (12.8 kmh)

TURNING RADIUS

1727 mm (68")

OVERALL LENGTH

1359 mm (53.5")

OVERALL WIDTH

635 mm (25.9")

FRONT/REAR TYRES

Front: 13"
Rear: 13"

**TOTAL WEIGHT
WITHOUT BATTERIES**

11 stone 5 lbs (73 kg)

BATTERIES REQUIRED

(2) 50 AH (deep cycle)



Available on
 Motability

WEIGHT CAPACITY
28 stone 8 lbs (181 kg)
MAX SPEED
8 mph (12.8 kmh)
TURNING RADIUS
1675 mm (66")
OVERALL LENGTH
1350 mm (53")
OVERALL WIDTH
635 mm (25.9")
FRONT/REAR TYRES
Front: 13" Rear: 13"
TOTAL WEIGHT WITHOUT BATTERIES
15 stone 3 lbs (97.5 kg)
BATTERIES REQUIRED
(2) 50 AH (deep cycle)



28stone 8lbs
 Weight Capacity

Up To
8mph
 Maximum Speed

1675mm
 Turning Radius

Available on
 Motability



28stone 8lbs
Weight Capacity

Up To
8mph
Maximum Speed

1973mm
Turning Radius

WEIGHT CAPACITY

28 stone 8 lbs (181 kg)

MAX SPEED

8 mph (12.8 kmh)

TURNING RADIUS

1973 mm (76.25")

OVERALL LENGTH

1428 mm (56.25")

OVERALL WIDTH

692 mm (27.25")

FRONT/REAR TYRES

Front: 13"
Rear: 13"

TOTAL WEIGHT WITHOUT BATTERIES

18 stone 4 lbs (117 kg)

BATTERIES REQUIRED

(2) 100 AH (deep cycle)



Available on
 Motability

WEIGHT CAPACITY
25 stone (158.75 kg)
MAX SPEED
8 mph (12.8 kmh)
TURNING RADIUS
2407 mm (94.75")
OVERALL LENGTH
1556 mm (61.25")
OVERALL WIDTH
829 mm (32.5")
FRONT/REAR TYRES
Front: 14.5" Rear: 14.5"
TOTAL WEIGHT WITHOUT BATTERIES
19 stone 5 lbs (124.25 kg)
BATTERIES REQUIRED
(2) 75 AH (deep cycle)



Stone Grey

Desert Sand

25stone
Weight Capacity

Up To
8mph
Maximum Speed

2407mm
Turning Radius

Available on
 Motability

MAXIMA



WEIGHT CAPACITY
36 stone (227 kg)
MAX SPEED
4 mph (6.4 kmh)
TURNING RADIUS
3 wheel: 1321 mm 4 wheel: 1372 mm
OVERALL LENGTH
3 wheel: 1320 mm 4 wheel: 1370 mm
OVERALL WIDTH
648 mm (25.5")
FRONT/REAR TYRES
Front: 10" Rear: 10"
TOTAL WEIGHT WITHOUT BATTERIES
3 wheel: 77.5 kg 4 wheel: 80 kg
BATTERIES REQUIRED
(2) 55 AH (deep cycle)

36stone
Weight Capacity

Up To
4mph
Maximum Speed

1321mm
Turning Radius

ACCESSORIES



XLR USB Charger

Charge your phone on the go!
This can be used while your powerchair or scooter is in motion.
Compatible with most mobile phones.



Baskets & Bags

Our handy and versatile range of additional large baskets and storage bags make transporting your essentials hassle free.



Rear Oxygen Tank Holder

Available with most models, featuring 2 thumb screws for safe transportation.



Rear View Mirror

Optional rear view mirror for enhanced visibility and rider safety.



Crutch & Walker Holders

Quick and easy access for your cane, crutch or walking stick with our selection of holders.



Phone Mounts

Keep your phone safe and secure while on the move with our phone mount.
Compatible with most mobile phones.



Cup Holders

Satisfy your thirst on the move with our handy drinks holders, compatible with most of our range of products.



Replacement Keys

We offer a replacement key service.

*Please contact your local Pride dealer for our full range of accessories and further information

PRIDE[®]

MOBILITY PRODUCTS LTD.

LIVE YOUR BEST[®]

Also from Pride:



Travel Scooter Range



Powerchair Range



Lift Chair Range

Pride Mobility LTD
 32 Wedgwood Road
 Bicester, Oxon
 OX26 4UL, England

Tel: 01869 324600
 Email: sales@pride-mobility.co.uk
 Web: www.pride-mobility.co.uk



Pride Mobility LTD reserve the right to alter, at any time, without prior notice, the design, specification, packaging or price without the incurring any obligation.