

Embrace your potential! The Colt Sport is fully loaded with sporty style and advanced features, including frosted LED headlights and an advanced suspension system. Drive at speeds up to 8 mph and see what life has in store for you!





Specifications

Weight Capacity	25 st (159 kg)
Maximum Speed ¹	Up to 8 mph
Ground Clearance ³	57 mm (2.25")
Turning Radius ³	1333.5 mm (52.5")
Overall Length ³	1194 mm (47")
Overall Width ³	657 mm (25.9")
Max. Obstacle Climb	50mm (2")
Drive Tyres	254 mm (10") front & rear - solid
Suspension	Front & rear suspension
Total Weight Without Batteries ^{3, 4}	68 kg (150 lbs)
Heaviest Piece When Disassembled ³	Front section: 29.48 kg (65 lbs)
Drivetrain	24 volt DC motor
Dual Breaking System	Regenerative & electromechanical
Battery Requirements ³	2 x 40 Ah deep cycle
Battery Weight ⁴	11 kg (24.5 lbs) each
Battery Charger	Off-board
Colour	Red & blue



Features

- Comfortable high back seat with armrests
- High visability front and rear LED lighting
- Wrap around delta tiller
- Automatic brake lights & LED kerb light
- Rear view mirror
- Exclusive & stylish low profile tyres
- Convienient storage in tiller

Also available in:

- Red
- Blue





NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or -3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.



¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.
2 Includes standard seat weight.

³ Due to manufacturing tolerances and continual product improvement,

this specification can be subject to a variance of (+ or -) 3%.

Battery weight may vary based on battery manufacturer.