

Dash Carbon Max HD

The lightest HD powerchair on the market



USER
WEIGHT
UPTO
200kg

ONLY
17kg



Frog leg front castors

Frog Legs shock-absorbing caster forks reduce shock for a smoother ride.



Carbon fibre for strength

Made from Japanese T1000 carbon for extra strength in the frame and also keeps the overall carry weight down, making it the lightest HD powerchair on the market.



Folding frame

Compact form factor when folded for easy transport and storage.



Weight info

Chair weight 17kg
(excluding battery)

Max user weight 200kg

DIMENSIONS

Seat width	535 mm
Seat depth	460 mm
Seat height	520 mm
Backrest height	600 mm
Overall length	1055 mm
Overall width	720 mm
Overall height	1050 mm
Overall length (folded)	850 mm
Overall width (folded)	440 mm
Overall height (folded)	720 mm
Front wheel size	8 inch
Rear wheel size	12 inch
Total weight (exc. battery)	17 kg

OTHER TECHNICAL DATA

Max. speed	6km/h
Braking distance	1.2 m
Turning radius	1.2 m
Climbing angle	9°
Battery weight	1.8 kg

KEY FEATURES

- ▶ Two year warranty
- ▶ Shock absorbing front and rear wheels with suspension
- ▶ Lightweight and strong carbon fibre frame
- ▶ 21" seat width
- ▶ Flip up armrest for side transfer
- ▶ Flip up footplate with depth adjustment
- ▶ 10 amp lithium battery as standard, extra 20 amp battery available

For Enquiries: