

E-Drive Blackline, Electric drive unit, specifications:

E-Drive, onboard programmable or programmable with Bluetooth Windows laptop/tablet:

- Speed 0-6km/h, braking distance, turning speed, acceleration, engine power
- 3 programs pre-programmed, Soft, Standard, Sport
- Program and joystick individually programmable
- Motors maximum 2 x 120W
- Wheel self-shutdown time 10, 60 minutes or not shutdown
- Weight per wheel 24" only 6.4 kg
- Max. user weight 125 kg
- Normal charge about 2-3 hours, Refresh charge 7-8 hours
- Maximum camber 3 degrees
- 3-in-1: drive yourself with pushrims, drive yourself with joystick or drive with attendant control



- Standard USB-Anschluss im Joystick zum Laden von z.B. Handy, 5V, 500 mA maximal

[Reset form](#)

[Save form](#)

Company name:

P.O No.

Date:

E-Drive 16" (No heavy duty option available)

MED316 ✓ **£5,695.00 16" E-Drive Blackline** wheels with QR axle, nickel-hydride battery, charger, battery bag

16" E-Drive Available wheelchair compatibility - please choose from the following:

Date of manufacture:

- | | | |
|---|-----------------------|--|
| ☐ Braunability | Carony | MED3200 |
| ☐ Etac | Transit 1 & 3 | MED3039 |
| ☐ Hoggi | Sherpa | MED3324 |
| ☐ HD Rehab | HD Balance 16" | MED3072 |
| ☐ HD Rehab | HD Motion 16" | MED3366 |
| ☐ Icon | A1 | MED3174 |
| ☐ InterCo | Economy | MED3387 |
| ☐ InterCo | Simply Light | MED3137 |
| ☐ Invacare | Rea Dahlia Assist 16" | MED3123 |
| ☐ Invacare | Rea Azalea Assist | MED3053 |
| ☐ Invacare | Rea Spirea4 NG 16" | MED3097 |
| ☐ Ki Mobility | Catalyst 5 & 5 Ti | MED3131 |
| ☐ Ki Mobility | Focus CR | MED3354 |
| ☐ Life & Mobility | Roxx & Roxx 2 | MED3143 |
| ☐ LIW Care Technology | Buffin Trio size 1 | MED3184 |
| ☐ LIW Care Technology | Buffin Trio size 2 | MED3177 |
| ☐ Otto Bock | Discovery | MED3371 (Includes 200mm x 50mm black front castor) |
| ☐ Progeo | Tekna Tilt 2.0 | MED3316 |
| ☐ Rehatec | Sunny | MED3373 |
| ☐ R82 | Kudu 16" | MED3181 |
| ☐ Sorg | Dynamis TSD 16" | MED3365 |
| ☐ Sunrise | IBIS XC | MED3141 |

Joystick control selection -

- JOY-N ☐ Without standard joystick
- JOY-R ☐ Mount joystick on right
- JOY-L ☐ Mount joystick on left

Joystick programme

Soft mode ☐ Standard ☐ Sport ☐

- SIMP-001 ☐ **£1,995.00** Simplex smart attendant control unit [Simplex link](#)
- MED1200 ☐ **£895.00** Upcharge for dual attendant control with bracket
- MED7660 ☐ **£160.00** Upcharge for parallel swivel arm bracket
- MED1356 ☐ **£110.00** Upcharge for attendant joystick holder

Joystick mounting selection - Choose joystick mounting MED No. from this link:

Attendant control push handle diameter 19-22mm ☐ 25-26mm ☐ 28-32mm ☐

- MED1347 ☐ Joystick mounted to frame
- MED-JOY-A ☐ Joystick mounted to armrest (check on adapter website first for correct armrest type MED No.)

[Adapter link](#)

Tyre: Please choose tyre style

- MS2 ☐ **N.C.O** Flex Insert solid tyre
- PNE-16 ☐ **N.C.O** Pneumatic tyre

Anti-tipper: (Tippskydd med stolsresning) - Choose from below and complete TB No.

- RL-KIPPS ☐ **N.C.O** Wheelchair has its own anti-tippers left and right supplied by wheelchair manufacturer
- XOF ☐ **£120.00** Anti-tip bracket on the wheel
- TBB2000 ☐ **£352.00** Anti-tippers with jack up function mounted to frame (Choose TB No.)

T.B No. [Adapter link](#)

Battery, charger and bag options

Battery placement: ☐ Backrest mounted ☐ Underseat (not li-ion) or some models

- MED2001 ☐ **N.C.O** NiMH Battery
- MED-L ☐ **£995.00** Upcharge to lithium-ion battery with charger
- MED12V ☐ **£230.00** Upcharge -Converter for charging in boat, caravan etc - output 230v ac 50hz
- MCG2018 ☒ **Included** 150w UK-Plug

Additional options and accessories

- MIMH-S ☐ **£1,095.00** Additional spare NiMH battery
- LITH-S ☐ **£1,495.00** Additional spare Li-ion battery
- MEM7611 ☐ **£250.00** Programming cable and bluetooth module (Dealers only)

Notes: