

E-Drive Blackline, Electric drive unit, specifications:

E-Drive, onboard programmable or programmable with Bluetooth Windows laptop/tablet:

- Speed 0-6km/h, braking distance, turning speed, acceleration, engine power
- 3 programs pre-programmed, Soft, Standard, Sport
- Program and joystick individually programmable
- Motors maximum 2 x 120W
- Wheel self-shutdown time 10, 60 minutes or not shutdown
- Weight per wheel 24" only 6.4 kg
- Max. user weight 125 kg, optional 150 kg "Heavy Duty"
- Normal charge about 2-3 hours, Refresh charge 7-8 hours
- Maximum camber 3 degrees
- 3-in-1: drive yourself with pushrims, drive yourself with joystick or drive with attendant control
- Standard USB-Anschluss im Joystick zum Laden von z.B. Handy, 5V, 500 mA maximal



Institute: TÜV Süd, Hanover
Certification: EN12184:2014

Crash test according to ISO 7176-19



CE Conforms to the EU directive for:
MDR 2017/745 class 1



[Reset form](#)

[Save form](#)

Company name:

P.O No.

Date:

E-Drive Plus - Choose model below

- MED ☒ **£5,495.00** E-Drive Blackline wheels with QR axle, Aluminium push rims, nickel-hydride battery, charger, battery bag
- HD ☐ **£400.00** Heavy duty wheel upcharge, up to 150kg user weight

E-Drive mounting information - complete wheelchair info and fill in MED number

Wheelchair manufacturer

Wheelchair model

MED3

Wheel mounting MED No.

Date of manufacture

[Adapters link](#)

Wheel selection - Choose wheel size below

- | | | | | | |
|-----|-----------------------|---|-------|-----------------------|--|
| 316 | ☐ | £200 16" Flex insert, no push rim, (No HD version) | 324-M | ☐ | £325 24" With off-road bike tyres |
| 320 | ☐ | £200 20" | 325 | ☐ | £200 25" |
| 322 | ☐ | £200 22" | 326 | ☐ | £200 26" |
| 324 | ☐ | N.C.O 24" | | | |

Joystick control selection -

- JOY-0 ☐ Without standard joystick
- JOY-R ☐ Mount joystick on right
- JOY-L ☐ Mount joystick on left

Joystick programme

Soft mode ☐ Standard ☐ Sport ☐

- SIM001 ☐ **£1,995.00** Upcharge for Simplex smart attendant control [Simplex link](#)
- MED1200 ☐ **£895.00** Upcharge for dual attendant control with bracket
- MED7660 ☐ **£160.00** Upcharge for parallel swivel arm bracket
- MED1356 ☐ **£110.00** Upcharge for attendant joystick holder

Joystick mounting selection - Choose joystick mounting

- Attendant control push handle diameter 19-22mm ☐ 25-26mm ☐ 28-32mm ☐
- MED1347 ☐ MED1347 Joystick mounted to frame
- MED**** ☐ Joystick mounted to armrest (check on adapter website first for correct armrest type MED No.)
- MED No.

[Adapters link](#)

Anti-tipper: (Tippskydd med stolsresning) - Choose from below and complete TB No. if applicable from this link:

RL-KIPPS	☐	N.C.O	Wheelchair has its own anti-tippers left and right supplied by wheelchair manufacturer
XOF	☐	£120.00	Anti-tip bracket on wheel (16", 20", 22", 24")
TBB2000	☐	£375.00	Anti-tippers with jack up function mounted to frame (Choose TB No.)

[Adapters link](#)

Choose Push rim option

PR-A	☐	N.C.O	Aluminium push rim
PR-SS	☐	£64.00	Upcharge Stainless Steel push rim
PR-PC	☐	£165.00	Upcharge Poly push rim
PR-FC	☐	£95.00	Upcharge foam compact push rim
PR-MG	☐	£415.00	Upcharge Max Grip push rim
A-COVER	☐	£85.00	Upcharge push rim removable cover for aluminium or stainless
SG	☐	£125.00	Upcharge for transparent spoke protectors

Choose tyre style

MED-PNE	☐	N.C.O	1" standard tyre
MED7692	☐	£95.00	Upcharge Schwalbe Marathon Plus tyres
DM2412	☐	£95.00	Upcharge Coarse tread 24 x1 3/8" tyre
MS2	☐	£95.00	Flex Insert (20" x 1" - 26" x 1") - standard on 16"
MEM7691	☐	£95.00	PU Tyre (20", 22", 24", 26")

Battery, charger and bag options

Battery placement:



Backrest mounted



Underseat (not with li-ion) or some models

MED2001	☐	N.C.O	Standard NiMH battery
MEDL	☐	£995.00	Upcharge to lithium-ion battery with charger
MED12V	☐	£230.00	Upcharge -Coverter for charging in boat, caravan or camper - output 230v ac 50hz 150w
MCG2018	☒	Included	UK-Plug

Additional options and accessories

NIMHS	☐	£1,095.00	Additional spare NiMH battery
LITHS	☐	£1,495.00	Additional spare Li-ion battery
MEM7611	☐	£250.00	Programming cable and bluetooth module (Dealers only)

Notes: