



eFOLDi[®]

eFOLDi MK 1.5

eFOLDi MK 1.5 is the lightest, most compact, folding mobility scooter in its class, weighing only 17Kg and no larger than a small suitcase when folded. The eFOLDi MK1.5 is fully UK road legal with a top speed of 8mph and a range of 14 miles.

eFOLDi fits easily into the smallest car and is fully airline compliant - the perfect travel companion.

Enjoy freedom, friendship and fun!



eFOLDi®

Features

Speed	4mph (6kmh) / 8mph (12kmh)
Range	14 miles / 22 Km
Carry capacity	120Kg / 265 lbs / 19 Stone
Climbing ability	1:6 Gradient (9.5 degrees)
Folded dimensions	66 x 56 x 40 cm / 26 x 22 x 15.6 inches
Tyre type	Pneumatic
Battery standard	12AH Li-ion
Motor	240W 24V Hub Motor
Ignition	Key operated

Version: MK1.5 05052021

www.efoldi.com