

The

# Denwick

Tilt-in-Space Single-motor  
Tailored Lateral

Made  
in the UK



# Denwick

## Tailored Lateral Back

Choose our furniture with confidence -  
All riser recliner chairs are tested by FIRA -  
The leading authority in furniture safety by  
design.



Riser Recliners available from STOCK in the following fabrics



VN  
Dove Grey



VN  
Slate Grey



VN  
Cocoa



VN  
Oatmeal



Colours may vary slightly from images shown which are for illustrative purposes only.

### Main Features:

- 3-button power on handset
- Tailored Lateral or Waterfall Back design
- Pocket sprung seat base
- Scroll arms with wings
- Twin Pockets
- Arm covers and Antimacassar included
- 5 Year Guarantee on frames, mechanisms & all electrics
- Quality Craftsmanship, **BRITISH MADE**

The Denwick Cosi Chair is a Tilt-in-Space Single-motor Riser Recliner inviting you to relax with its exceptional level of comfort & support. Maintaining a consistent angle between the seat & backrest for stability & to reduce shear/friction, it also offers the perfect neutral sitting position & level of back support for those needing high leg elevation. Featuring a **NEW** Tailored Lateral back design, it boasts style & structure coupled with excellent lumbar back support. Available in four new contemporary colours.

### Denwick Size Guide:

	One size
Seat height	48.3cm / 19"
Seat width	50.8cm / 20"
Seat depth	50.8cm / 20"
Arm height	20.3cm / 8"
Back height	71.1cm / 28"
Overall width	82.5cm / 32½"
Overall height	111.8cm / 44"
Weight limit	20 stone / 127kg

Local Cosi Chair Retailer: