The

# Denwick

Tilt-in-Space Single-motor Waterfall





## **Denwick** Waterfall Back

Choose our furniture with confidence -All riser recliner chairs are tested by FIRA -The leading authority in furniture safety by desian.




### Riser Recliners available from STOCK in the following fabrics













Colours may vary slightly from images shown which are for illustrative purposes only.

The Denwick Cosi Chair is a Tilt-in-Space Single-motor Riser Recliner inviting you to relax with its exceptional level of comfort & support. Maintaining a consistent angle between the seat & backrest for stability & to reduce shear/friction, it also offers the perfect neutral sitting position & level of back support & is ideal for when legs need elevating above the hips. Featuring a deep fibre filled Waterfall backrest, available in four new stylish & contemporary colours.

#### Main Features:

- 3-button power on handset
- Waterfall or Tailored Lateral Back design
- Pocket sprung seat base
- Scroll arms with wings
- Twin Pockets
- Arm covers and Antimacassar included
- 5 Year Guarantee on frames, mechanisms & all electrics
- Quality Craftsmanship, BRITISH MADE



#### **Denwick Size Guide:**

	One size
Seat height	48.3cm / 19"
Seat width	50.8cm / 20"
Seat depth	50.8cm / 20"
Arm height	20.3cm / 8"
Back height	71.1cm / 28"
Overall width	82.5cm / 32½"
Overall height	111.8cm / 44"
Weight limit	20 stone / 127kg

Local Cosi Chair Retailer:







