

388 S Specification:

Type Class	B/2
Overall length	125cm / 49"
Overall width	54cm / 21"
Floor to seat height (cm)	56-64
Seat width	45.7cm / 18"
Seat depth	40.6cm/16"
Maximum carrying capacity	136kgs/21 stone
Battery voltage and capacity	12V/33Ahx2
Maximum speed	6.4km/h / 4mph
Ground clearance	10cm / 4"
Maximum kerb height	5cm / 2"
Range* - standard battery pack	up to 32km / 20 miles
The maximum safe slope	8°
Turn-around width	1.7m / 5'6"
Colour - 388 S	Red/Silver
+D 6 11 1 10 :	

^{*}Range on full charge and flat ground

Main Features:

- Outstanding reliability
- Flip-up width adjustable armrests
- · Comfortable high back, height adjustable, sliding swivel seat
- Narrow compact chassis

UPGRADE

Upgrade:

• The Rascal 388 XL, offers all the great features of the 388 Deluxe, plus a maximum speed of 6 mph and a ground clearance of 13cm

Rascal 388 S Profile:

The Rascal 388 S is a well established, great value for money mid-size scooter to suit your everyday mobility needs. At 21" wide, it is the narrowest width of any scooter in its class & is equally at home inside a busy supermarket as it is on the pavement. With a comfortable high back seat with headrest, smooth controls and adjustable tiller, this high performance, reliable scooter is your ideal mobility partner.

Local Rascal Mobility Retailer:









