

ELEVATION



The ***Elevation by PDG*** allows its user to embrace mobility related activities of daily living by quickly and easily adapting their seating position to their mobility needs through ***on-the-fly seat height and back angle adjustment.***



Hand
Propulsion



Foot
Propulsion



Self-Tilt



Desk/Table
Accessible



Transit Tie
Downs

Elevation by PDG Award-winning ultralight chair

visit **elevationbypdg.com**
for more information



DYNAMIC TILT
35° Range



WEIGHT CAPACITY
250 lbs



SEAT WIDTH
12" to 19"



SEAT DEPTH
14" to 20"



SEAT HEIGHT
18" to 21"



BACK SEAT ANGLE
84° to 120°



OVERALL WIDTH
+9.0"



CHAIR WEIGHT
22 lbs



ELEVATION



Exclusive features

- ✓ Dynamic seat height adjustment (anterior/posterior tilt)
- ✓ Gas strut suspension
- ✓ Fully adjustable back angle
- ✓ Patented stability system

Available options *See order form*

- ✓ Anti-shear side guard system
- ✓ Fold down back
- ✓ Compatible with SmartDrive and Alber e-motion
- ✓ Numerous wheel and caster options




User-controlled dynamic seat height adjustment
from 15° of seat dump to 20° of seat lift



On-the-Fly back angle
adjustment



35° of lift/tilt range



Optional SmartDrive mount

