

E | F55

Retail Order Form



ORDER FORM

Billing Address Or Customer-ID			Shipping Address	Identical To Billing Address
C.-ID:			C.-ID:	
Name:			Name:	
Street:			Street:	
Town/PC:			Town/PC:	
E-Mail:			E-Mail:	
Phone:			Phone:	
Order Information			User Weight	
Order	Calculation	Built-4-Me®	Max. user weight: 110 kg	
Purchase Order No.:			User Weight:	kg
Art. No.: 47600000				
Order Date:	Day	Month	Year	

Technical Data Overview

Max. user weight	110 kg	Max. speed forward	6 kph; optional 10 kph or 15 kph
Front frame width	260 – 480 mm ⁽¹⁾	Wheel size	8,5" air-chambered solid-tyre / 14" pneumatic tyre
Battery weight	2 kg	Battery type	Li-Ion, 36V, 8.1 Ah
Total weight w/o battery	11 kg	Brake type	Twin-brake system: mechanical disc-brake and separate electric brake
Total weight with battery	13 kg	Operating Voltage	36 Volt
Colour	Black	Range of action	25 km max.
Colour accents (included in delivery)	Black, red, blue, orange ⁽²⁾	Climbing power	up to 11% (6°)

Base		£
	Pull-device F55 400 W (max.) motor, 8.1 Ah Li-ion battery, forward- and reverse gear, cruise-control, dual-brake-system (disc-brake, e-brake), mudguard, front-light and colour-kit included.	0
Drive wheels		
EMF070085	Drive wheel, 8,5", solid tyre	0
EMF070086	Drive wheel, 14", pneumatic tyre	0
Ordering/Bundling		£
EMF090213	Without wheelchair order	2,605
EMF090217	With wheelchair order, not mounted to the wheelchair (1)	2,605
Docking		£
EMF010027	Docking for fixed-front wheelchairs	0
EMF010028	Docking for swing-away wheelchairs	0
Docking-Clamps		
EMF010042	Docking-clamps for Ki Mobility Rogue (3) (4)	0
EMF010041	Docking-clamps for oval tubes (Progeo Joker / Ego / Tekna) (3) (4)	0
EMF010029	Docking-clamps for round tube frames (adaptable to 19.5 mm, 23 mm, 25 mm, 28.6 mm & 30 mm tubes) (3) (4)	0
EMF010030	Docking-clamps for oval tubes (Quickie Xenon ² FF, Helium) (3)	0
EMF010031	Docking-clamps for Quickie Xenon ² SA (3)	0
Adjustability		£
	Complete easy-adapt kit included to fit F55 to front frame widths from 260 – 480 mm (5). The F55 is fully adjustable to all seat heights and frame angles.	0
Speed Setting		£
EMF000025	6 kph	0
EMF000027	10 kph (6)	0
EMF000028	15 kph (6)	0
Display		£
	Display: Battery status, Odometer, trip counter, avg. speed, real time power consumption	0
Accessories		£
	LED front light	0
	Mudguard for front wheel	0
	Colour-kit included: mudguard and sidebars in 4 colours	0
EMF010032	Kickstand: mounts to the sub-frame	158
EMF010033	Padding for sub-frame (enclosed)	0
EMF160010	Additional Li-Ion battery, Standard: 36V, 8.1 Ah, 291 Wh (approx. 25 km range)	424
EMF110155	Additional charger	112

(1) Wheelchair order on separate order form. Separate shipment of wheelchair and F55.

(2) Wheelchair-order on separate order form. Shipment of readily mounted wheelchair/F55-unit

(3) Always check installation space for the docking clamps before ordering. Depending on configuration and option, a particular configuration of a chair might not offer enough space.

(4) Always check combination agreements for non-Sunrise wheelchairs.

(5) Frame width measured from center frame tube to center frame tube at the mounting position.

Use only the measured dimension and not the seat width to check compatibility.

(6) Check your local laws and regulations for max. permitted speed in your country. Speed-settings > 6 km/h can only be ordered via Fax with a signed waiver (next page). Max. speed-settings can only be altered by Sunrise Medical.

Accessories



LED-light, Fender (accent colour red) and disc-brake (pictured on EMF070085 Drive wheel 8,5", solid tyre)



// EMF010032 Kickstand: mounted to the sub-frame



Colour-kit (mudguard) included in delivery



Colour-kit (sidebars) included in delivery (accent colour blue)

Letter of Acknowledgement

We herewith advise, that using the Empulse F55 with max. speed-settings higher than 6 kph is only permitted on public roads in certain countries and otherwise may only be used on private property.

In the absence of approval, in accordance with the road traffic laws, the Empulse F55 with speed-settings > 6 kph may not be driven on public roads, bicycle-lanes and pavements. You need to inform yourself about the country specific legislation regarding usage and your legal obligation to insure the Empulse F55 speed settings higher than 6 kph (10/15 kph). When you're travelling with your F55, always check the specific legislation in the respective country.

When placing an order for a max. speed-setting > 6 kph, please confirm with your signature below, that you have understood the above information and that you only intend the usage of your F55 on private property.

I herewith confirm, that I am aware of the information above, that I will follow strictly the country specific legislation.

Date

Customer Signature

Day Month Year

Comments

--

F55

£

Base Price

00,000.00

Total Of Additional Equipment

Packaging And Transport

Exclusive VAT

VAT

Total Price

Date

Signature

Day Month Year