

Wheel Add-On



03/2025

ORDER FORM

Errors and omissions excepted.

Billing Address Or Customer-ID		Shipping Address		☐ Identical To Billing Address
C.-ID:		C.-ID:		
Name:		Name:		
Street:		Street:		
Town/PC:		Town/PC:		
E-Mail:		E-Mail:		
Phone:		Phone:		
Order Information		User Weight		
☐ Order ☐ Calculation ☐ Built-4-Me®		Max. user weight: 100 kg		
Purchase Order No.:		User Weight: kg		
Art. No.: 99200001		Reimbursement Code: xxx		
Order Date: Day Month Year		Base Price £4,995		

Technical Data Overview

Max. user weight	100 kg (220 lbs)	Motor	150W
Rear Wheel Sizes	24" (540mm)	Operating Voltage	24V
Total Device Weight	12 kg (26.4 lbs)	Maximum Motor Noise	60db
Per Battery Weight (4.0Ah)	0.6 kg (1.3 lbs)	Battery Type	Lithium Ion (21.6V DC and 86.4Wh each battery)
Max. Range Per Charge	15 km (9.3 mi)	Device Charger Type	28V DC 5.7A 160W
Maximum Speed	6 kph (3.7 mph)	Docking Charger	28V DC 1.5A IP44
Maximum Slope	6°	Control Box	3.7V 1Ah 3.7Wh
Bluetooth Standard	BLE (125Kbps Coded PHY)		

Base			£
☒	M90 unit Incl. Control box, Joystick +Two-Way-Toggle, Batteries, long and short mounting hardware, Quick Release Axles, Cable Charger)		4995

Sunrise Medical Wheelchair Models - Specific mounting hardware will be provided for the chair selected

Specific wheelchair model the M90 Wheels will be used on					
☐	EMM000035	Parix ²	☐	EMM000184	Neon ²
☐	EMM000036	Basix ²	☐	EMM000052	Ibis
☐	EMM000037	Rubix ²	☐	EMM000075	SR45
☐	EMM000038	Helix ²	☐	EMM000065	Xenon ²
☐	EMM000055	Easy Life Eco	☐	EMM000076	IRIS
☐	EMM000056	Easy Life i	☐	EMM000103	Argon ²
☐	EMM000102	Easy Life T	☐	EMM000174	QS5 X
☐	EMM000111	Easy Life RT	☐	EMM000101	Helium or Nitrum
☐	EMM000059	Style X / Style X Ultra	☐	EMM000182	Zippie IRIS
☐	EMM000108	Easy Life	☐	EMM000183	SIMBA
☐	EMM000109	Easy Life R	☐	EMM000149	Youngster 3
☐	EMM000129	Easy 200/300			
RGK - Empulse M90 is compatible with RGK models. Contact your local RGK representative for more information.					

Handrims			£
☒	EMM070096	Handrims, aluminium, anodised, black	0
☐	EMM070207	Handrims, Supergrip R	250
☐	EMM070212	Handrims, Ellipse, 3R, black	125

Tyres			£
☒	EMM070104	Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution	0
☐	EMM070102	Solid tyre, puncture proof, tread	0

Options and Accessories - Batteries are interchangeable and can be used with either wheel			£
☐	EMM090227	Leg Mount Controller Includes thigh strap and hardware	50
☐	EMM160042	Additional Li Ion Battery Set 2 additional batteries (21.6V DC and 86.4Wh each battery)	550
☐	EMM260243	Docking charger, 3 A - 36 V For charging battery external	130

Comments

Comments

M90

£

Base Price	00,000.00
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	

Date

Signature

Day

Month

Year

Send E-Mail

**Sunrise Medical GmbH**

Kahlbachring 2-4

D-69254 Malsch/HD

Phone: +49 (0) 72 53 / 9 80-0

Fax: +49 (0) 72 53 / 9 80-222

E-Mail: kundenservice@sunrisemedical.de

www.SunriseMedical.eu