

**AVAILABLE FOR TRIAL
AT OUR PARTNERS**



"I truly realized the relief that Skeed provides for the hands, wrists, arms... You just activate the lever with your fingers and it brakes! It's an everyday advantage, even indoors and at low speed, it's incomparable!"

Clément
Skeed user



CONTACT US

✉ hello@eppur.eu ☎ +33. 6. 59. 95.02

WWW.EPPUR.EU



**BOOK YOUR
FREE TRIAL !**



Do not dispose of in public areas.

This medical device is a regulated health product which, under this regulation, bears the CE marking.

EPPUR
-SKEED

GROUND BRAKING TECHNOLOGY



MADE IN FRANCE

EPPUR, a simplified joint-stock company with a share capital of €28,901, registered with the Registry of Commerce and Companies of Lille Métropole under number B 881 789 408.

STEER AND BRAKE WITH YOUR FINGERTIPS

Skeed is a pair of wheels that integrates a unique braking system at its center. A fixed activation lever directly connected to the hub allows you to slow down and immobilize your wheelchair, even without grip!



HANDLEBAR
TETRA OPTION



HANDLEBAR
PARA OPTION



PROGRESSIVE BRAKE

The more you pull, the more you brake! Activate the brake lever backwards to slow down and steer your wheelchair with unprecedented precision.



PARKING BRAKE

Push the lever forward to immobilize your wheelchair. A powerful brake ensures total stability, regardless of your tire pressure.

Voir details et options sur la fiche de mesure.

LESS EFFORT, + MORE SAFETY SKEED SAVES YOUR HANDS AND ARMS

5X LESS EFFORT AND ZERO FRICTION

With Skeed, no more burns and risk of injury ! The lever eliminates all friction on your hands and divides the braking effort by 5. Your arms and hands are thus preserved !

MAXIMUM LIGHTNESS

A light wheel, to facilitate transfers and reduce the overall weight of your wheelchair.

Naked rim, hub & lever : 1.5 Kg
Standard configuration : 2.15 Kg



100% TWO-WHEEL MANEUVERABILITY PRESERVED

Thanks to its fixed handrim, Skeed guarantees you perfect two-wheel control. Overcome obstacles and maneuver with agility.

COMPATIBLE WITH MOTORIZATION

Skeed wheels are particularly useful for mastering 5th rear propulsion wheels safely, without straining your shoulders on every turn or slope.