


Date	UK (English)	Folding frame wheelchair
User	Company stamp	
End user weight (kg)		£ 2290,00
Language user manual	Responsible	Retail price - VAT not included

Wheelchair measures	
Seat width (LS)	☐ 33 ☐ 36 ☐ 39 ☐ 42 ☐ 45 ☐ 48
Seat depth (PS)	☐ 35 ☐ 37,5 ☐ 40 ☐ 42,5 ☐ 45 ☐ 47,5
Footplate distance (DP)	 (cm) ± 0.5 cm - adjustment steps 1.5cm (Real distance between seat and footplate) - IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm.)
Setting	☒ Standard ☐ Prudential ☐ Active ☐ Extreme
Backrest Height (TS)	Std.Value 36(min 30-max 54) adjustment steps 1.5cm
Crossbrace (black)	☒ Single ☐ Double (Not seat width 33-36)
Front height	(STD 49) - Rear height (STD 44)

Rear frame	
Aluminium frame	☒ Aluminium standard frame - Same colour of front frame  The real seat depth, from the front part of the seat cloth to the backrest tube will be 1 cm more than the declared one. camber ☒ 0° ☐ 1,5° ☐ 3°
	☐ Aluminium frame DYNAMIC 2.0 - (COLOUR: BLACK) 09955 315,00  The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one. With rake adjustable backrest ($\pm 9^\circ$) Exclude the installation of : Brake handles for companion (Drum Brake) camber ☒ 0° ☐ 1,5° ☐ 2° ☐ 4°
Titanium frame TMAX	Front frame + Titanium frame T MAX 06858 687,00  The real seat depth, from the front part of the seat cloth to the backrest tube will be 1 cm more than the declared one. <u>Brushed Titanium frame - not varnished</u> <u>Weight reduction approx. 450gr</u> camber ☒ 0° ☐ 1,5° ☐ 3°

Front frame	
☐ Short frame - Ideal with Seat depth (PS) 35 - 37,5 ☐ Medium frame - Ideal with Seat depth (PS) 40 - 42,5 ☐ Long frame - Ideal with Seat depth (PS) 45 - 47,5	Frame Angle ☒ 110° ☐ 100°  








Frame colour (aluminium)	
☒ Opaque black ☐ Bright silver ☐ Bright fuchsia	☐ Bright anthracite ☐ Bright orange ☐ Bright blue ☐ Bright lime ☐ Bright blue china ☐ Pearl white ☐ Bright red

Backrest	
☒ Adjustable velcro straps and padded transpiring upholstery	

Handle tube	
☒ Tubes with standard handle ☐ Folding handles	 02958 252,00

Colour set: Rear wheels and forks - To be matched only with: Rear wheels Ultra-light Progeo - Standard!			
☒ Polished black	22" 24" 25" 26"		
☐ Bright silver	22" 24"	00833	126,00
☐ Bright anthracite	24"	00833	126,00
☐ Bright orange	22" 24"	00833	126,00
☐ Bright fuchsia	22" 24"	00833	126,00
☐ Bright blue	22" 24"	00833	126,00
☐ Bright red	22" 24"	00833	126,00
☐ Bright lime	22" 24"	00833	126,00
☐ Pearl white	22" 24"	00833	126,00

Date	Order	User	Address
Casters			
☆ 5" Soft			
☐ 4" Soft			
☐ 3" Sport			
☐ 6" Soft (Available only with 110° angle)			
☐ 6" Pneumatic (Available only with 110° angle)			
☐ 7" Soft (Available only with 110° angle)			
☐ 3" Active Soft (wide)	(Aluminium alloy rim)	06824	126,00
☐ 4" Active Sport (narrow)	(Aluminium alloy rim)	01398	126,00
☐ 4" Active Soft (wide)	(Aluminium alloy rim)	03943	126,00
☐ 5" Active Soft (wide)	(Aluminium alloy rim)	07518	126,00
☐ 5" Pneumatic		10844	63,00
Rear wheels			
☆ 24" Ultra-light Progeo - Standard!			
☐ 22" Ultra-light Progeo - Standard!			
☐ 25" Ultra-light Progeo - Standard! (colour: black)			
☐ 26" Progeo - Standard! (colour: black)			
☐ 24" ULTRA-LIGHT PROGEO KEY SPORT	Weight reduction approx. 450gr.	04383	876,00
☐ 24" Ultra-light "Spinergy LX"		01685	876,00
☐ 24" Ultra-light "Spinergy LXL"	☐ Black hub and white spokes ☐ Hub and spokes black	07846	876,00
☐ 24" Ultra-light "Spinergy Flex rim"	Integrated aluminium/rubber handrim	02270	1.374,00
☐ 25" Ultra-light "Spinergy LX"		02101	876,00
☐ 25" Ultra-light "Spinergy LXL"	☐ Black hub and white spokes	07847	876,00
☐ 26" Ultra-light "Spinergy LX"		02102	876,00
☐ 24" Off-road (aluminium handrim black)		01392	224,00
☐ 22" with brake handles for companion (drum brake) (aluminium handrim)		03385	360,00
Available only with camber 0°			
Not available with Frame Dynamic 2.0°			
☐ 24" with brake handles for companion (drum brake) (aluminium handrim)		03384	360,00
Available only with camber 0°			
Not available with Frame Dynamic 2.0°			
☐ One arm drive - (24" rear wheels Aluminium handrim High pressure Tyre) - camber 0°			
☐ Right ☐ Left		02103	401,00
Rear wheel tyre			
☆ High pressure tyre (22"-24"-25"-26")			
☐ High pressure profile tyre (24")			
☐ Marathon Plus	(24"-25") with anti-puncture protection	01669	63,00
☐ Anti-puncture - Colour Black	(24")	11083	126,00
☐ Anti-puncture - Colour Grey	(22")	04181	98,00
Handrim			
Wheel distance ☆ 3 cm ☐ 4 cm			
☆ Aluminium (22"-24"-25"-26")			
☐ Aluminium black SOFT-TOUCH	(24"-25") (24" handrim, on 25" wheel is installed with spacers)	09170	126,00
☐ Chromed Steel	(24"-25") (24" handrim, on 25" wheel is installed with spacers)	01684	148,00
☐ Titanium	(24"-25") (24" handrim, on 25" wheel is installed with spacers)	00016	275,00
☐ Ultra Grepp	(24"-25") (24" handrim, on 25" wheel is installed with spacers)	00017	315,00
☐ Max Grepp	(24"-25") (24" handrim, on 25" wheel is installed with spacers)	00834	349,00
☐ Ergo Para	(24")	00632	349,00
☐ TOP FRICTION	(24"-25") (24" handrim, on 25" wheel is installed with spacers)	11023	399,00
Oval countour with a rubber insert			
☐ NOVA H	(24"-25") (24" handrim, on 25" wheel is installed with spacers)	11025	263,00
Ergonomic black profile			
☐ NOVA D	(24"-25") (24" handrim, on 25" wheel is installed with spacers)	11026	286,00
Ergonomic profile with rubber coating (ideal for quadriplegics)			
☐ SURGE LT	(25") Oval countour with a rubber insert	06427	360,00
☐ Quad	(24")	00693	200,00
☐ Silicon cover for standard aluminium handrim	Available for measures (22"- 24" - 26")	04050	160,00

Date	Order	User	Address
Wheel locks			
☒ Bent lever push to lock ☐ Straight lever push to lock ☐ Straight lever pull to lock ☐ Sport brake ☐ Brake with extension lever push to lock ☐ Brake with extension lever pull to lock		     	 06428 126,00 01949 63,00 03372 63,00
Footplate - See reference tables N°3			
☒ Plastic footrest, rake adjustable (Minimum footplate distance 33cm.) ☒ 2/3 Internal ☐ 2/3 External			11525 no list
☐ Automatic Closure footplate, rake and depth adjustable (Minimum footplate distance 30cm.) ☐ Internal ☒ 2/3 Internal ☐ External ☐ Minimum footplate distance 25cm. Available only with Multi-function and rotative system - 03945 189,00			11526 no list
☐ One piece footplate - Flip-up! rake and depth adjustable (Minimum footplate distance 30cm.) ☐ Internal ☒ 2/3 Internal ☐ External ☐ Minimum footplate distance 25cm. Available only with Multi-function and rotative system - 03945 189,00			11527 no list
☐ One piece aluminium footplate - Flip-up! rake and depth adjustable (Minimum footplate distance 30cm.) ☐ Internal ☒ 2/3 Internal ☐ 2/3 External ☐ External			11528 no list
☐ One piece Carbon fibre footplate rake adjustable, swing away, flip-up. (Minimum footplate distance 30cm.) ☒ Internal ☐ External			03291 275,00
☐ Footplate COMPACT, automatic Closure, rake and depth adjustable Minimum footplate distance 35cm. ☒ Internal ☐ External			09957 252,00
Mudguard			
Mudguard - FIXED			
☒ ABS mudguard ☐ ABS Straight sideguard ☐ Carbon fibre mudguard ☐ Carbon fibre Straight sideguard			03946 315,00 08918 315,00
Mudguard - REMOVABLE			
☐ Removable ABS mudguard ☐ Removable ABS fibre straight sideguard ☐ Removable carbon fibre mudguard ☐ Removable Carbon fibre straight sideguard			08921 126,00 08921 126,00 08919 435,00 08920 435,00
Armrests			
☐ Tubular armrests (Swing away, removable and height adjustable armrest) ☐ Standard Support (Not with frame Dynamic) ☐ Steel support with safety locking system			04828 177,00 00694 252,00
Removable armrest with ABS mudguard			
☐ Desk armrests PROGEO (Removable and height adjustable) Not available with Frame Dynamic 2.0 + camber 4°			09956 224,00
☐ Removable armrest with ABS straight sideguard ☐ Removable armrest with carbon-fiber straight sideguard			08881 224,00 08882 453,00

Date	Order	User	Address
Accessories			
☐	Pad for front frame protectors		07058 52,00
☐	Backrest-seat connection cloth		10241 17,00
☐	Removable storage bag (with quick coupling)		09811 126,00
Optional 1			
☐	Transit wheels		00836 148,00
☐	Anti-tipper with rapid hookup	(inside rotation to reduce the encumbrance) ☐ Right ☐ Left	00023 126,00 00024 126,00
☐	Tipping aid	☐ Right ☐ Left	00025 74,00 01046 74,00
☐	Crutches support	☐ Right ☐ Left	00838 74,00 00837 74,00
Optional 2			
☐	Tetraclip	(Quick release axle for quadriplegic)	03947 177,00
☐	Rear wheels plate extension	(Not with frame Dynamic 2.0)	01950 177,00
☐	Handybike adaptor	(Not with frame Dynamic 2.0)	00506 74,00
☐	Push handles for companion	(Removable and height adjustable)	01048 189,00
☐	Harness		04103 126,00
☐	Safety pelvic belt		03972 63,00
☐	Tray table (removable)		00508 275,00
☐	Lumbar support	(to insert in the backrest upholstery)	00013 40,00
☐	Super Head rest	(Folding with the wheelchair)	04049 453,00
☐	Swing away lateral support	☐ Right ☐ Left	02097 252,00 02098 252,00
	Not available with Posture backrest PHYSIO		
☐	Spoke guards Black	(22"-24")	01043 148,00
☐	Spoke guards Transparent	(22"-24")	00835 148,00
☐	Spoke guards Progeo design	(22"-24")	00019 177,00
☐	Spoke guards Dragon	(22"-24")	02106 177,00
Additional Kit			
☐	24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)	03378 498,00
☐	24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)	09357 876,00

Date Order User Address

PHYSIO - Posture backrest - Std. Backrest upholstery not provided

Removable with quick release
Suitable for tubes with diameter 18÷26mm
Depth adjustment (4cm)
Angle adjustment (±18°)

☐ **AIR** Carbon shell

not available for seat width 48

BACKREST HEIGHT		28	34	40			
Shell's shape		Lateral retaining wing A=high - R=regular					
7.5 cm		A	R	A	R	A	R
							
		☐	☐	☐	☐	☐	☐

1.047,00

☐ **POSTURE** Aluminum shell

BACKREST HEIGHT	34	40	52	
Shell's shape	Lateral retaining wing B=regular			
	R	R	R	
7,5 cm				
	☐	☐	☐	

1.013,00

Configuration DEEP with side support PAD (height adjustable)

☐ Aluminum lateral with padded, breathable upholstery.



Small (h 13cm)



☐ Right
10272

29,00



☐ Left
10272

29,00

Large (h 20cm)



☐ Right
10273

40,00




☐ Left
10273

40,00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Aluminum lateral with padded, breathable upholstery.



☐ Large (18x12)

05647 183,00

☐ Small (14x9)

05643 183,00

☐ Lateral support extension

Quantity (number pcs).....

10271 69,00

☐ Lateral support extension

Quantity (number pcs).....

☐ Small (14x9)

05642 183,00

☐ Large (18x12)

05646 183,00

☐ Headrest adjustable - Double ball mounting system
Available only with backrest Physio
Not available for backrest Physio AIR carbon shell

☐ Large (34x14)
☐ Small (28x14)



05648 252,00

Permobil LimitedUnit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts, AL7 1ANUnited Kingdom - Phone: +44 1484 722888

customer_support_uk@permobil.com - https://www.permobil.com/en-gb


Date	Order	User	Address
------------	-------------	------------	---------------

Note

Wheelchair : Progeo EXELLE - Rev. 2023.1	Prices - VAT excluded to be added to the price list	Send this page only if any optional have been chosen	
		☐ standard	(supplied if no requests are indicated)
		☐ on request	(without extra charge)
		☐ optional	(cost to be added at the basic price)

REFERENCE TABLES

1	Front height	
With closed frame (95°) max 5" caster		
Casters	Front height	
	min.	max
3"	41	52
4"	42	53
5"	44	54
6"	48	55
7"	49	56
The size of the fork will be defined only by Progeo in accordance with wheelchair final configuration.		

2	Rear height				
Rear frame		Rear wheels			
		22"	24"	25"	26"
STANDARD	 Straight plate	min. 37	min. 39	min. 40	min 41
		max 45	max 47	max 48	max 49
	 Shaped plate	min. 36	min. 38	min. 39	min. 40
		max 49	max 51	max 52	max 53
DYNAMIC 2.0		min. 36	min. 38	min. 39	min. 40
		max 49	max 51	max 52	max 53

3	Footplate option	Footplate distance	
		Frame Angle 110°	Frame Angle 100°
	Plastic footrest	min. 33	min. 31
	Automatic Closure footplate	min. 30	min. 28
	Automatic Closure footplate , with multifunction rotative system	min. 25	min. 23
	One piece footplate - Flip-up!	min. 30	min. 28
	One piece footplate - Flip-up! , with multifunction rotative system	min. 25	min. 23
	One piece aluminium footplate - Flip-up!	min. 30	min. 28
	One piece Carbon fibre footplate	min. 30	min. 28
Calculated with the best configuration and with the footplates positioned on inside.			
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm.)			

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all coulors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)