




GEPARD

HD

CARBON FIBRE ROLLATOR

GEPARD HD ROLLATOR



The Mobilex heavy duty Gepard rollator has a frame made of carbon fibre and aluminium, making it one of the lightest rollators on the market today starting from just 5kgs*. Although very light, it offers a high degree of stability coupled with extreme ease of movement.

The push handles are height-adjustable with an easy to read numerical marking and a simple thumb turn lock. It features a textile seat and padded back strap, and a locking system that stops it reopening once folded.

The large 8" smooth-running wheels are equipped with 2 brakes for optimum speed control that are also lockable for parking. The detachable bag has a load capacity of 5kgs and can also be used as a convenient shopping bag.

This heavy duty rollator is tested and approved for a user weight of up to 200kgs, a 50kg increase from the standard version.

Specifications	MG62HD
Overall Width	68cm
Overall Length	66cm
Product Weight	6.1kgs
Width Folded	22.5cm
Seat Height	62cm
Seat Width	56cm
Height of Push Handles	74cm – 102cm
Max. User Weight	200kgs
Wheels	200mm x 35mm
Max. Bag Load	5kgs

*Weight of Gepard MG55

OPTIONAL ACCESSORIES



WATERPROOF BAG



PADDED SEAT



SOFT WHEELS



ONE HANDED BRAKE



SPEED REDUCTION WHEELS

Sold and maintained by: