




GEPARD

CARBON FIBRE ROLLATOR

GEPARD ROLLATOR



The Mobilex Gepard has a frame made of carbon fibre, making it one of the lightest rollators on the market at just 5kgs*. Although very light, it offers a high degree of stability coupled with ease of movement.

The push handles are height-adjustable with easy to read numerical markings and a simple thumb turn lock. It features a textile seat, padded back strap, walking stick holder and a unique locking system that stops it reopening once folded.

The large 200mm (8") smooth-running wheels have soft tyres for enhanced comfort, and the rear wheels are equipped with brakes for optimum speed and movement control. These are also lockable for parking. The mesh bag has a load capacity of 5kgs and can be removed and used as a convenient shopping bag.

Specifications	MG55	MG62
Overall Width	60cm	60cm
Overall Length	66cm	66cm
Product Weight	5.5kgs	5.8kgs
Width Folded	22.5cm	22.5cm
Seat Height	55cm	62cm
Seat Width	47cm	47cm
Height of Push Handles	66cm – 86cm	74cm – 102cm
Max. User Weight	150kgs	150kgs
Wheels	200mm x 35mm	200mm x 35mm
Max. Bag Load	5kgs	5kgs

*Weight of Gepard MG55

OPTIONAL ACCESSORIES



WATERPROOF BAG



LARGE BACKSTRAP



PADDED SEAT



SOFT WHEELS



ONE HANDED BRAKE



NON-SLIP TRAY

Sold and maintained by: