




JAGUAR

INDOOR ROLLATOR

JAGUAR ROLLATOR



The Mobilex Jaguar indoor rollator has a light frame made of aluminium with the ability to be neatly folded for easy storage.

Overall the structure is very light, but offers high stability with high rolling ease. The large running wheels are equipped with a brake on each of the rear wheels which can be locked using the wide single pull bracket running the width of the frame.

The rollator is only approved for indoor use and has been tested and approved for a user weight of up to 100kg. The Jaguar comes complete with bag and detachable tray.

Specifications	MJAG
Overall Width	54.5cm
Overall Length (Depth) Opened	66cm
Length (Depth) Folded	27cm
Height of the Push Handle (Adjustable)	83.5cm – 94.5cm
Size of the Rollator Folded (HxWxD)	83.5 x 54.5 x 27cm
Weight	6.3kgs
Max. User Weight	100kgs
Brakes	2 with parking function
Front and Rear Wheel Size	175mm x 27mm
Max. Bag Load	5kgs

Sold and maintained by: