

# FUZE JR



The *Fuze JR* is a 50° tilt chair for mature pediatrics or petite adults. It features a unique center of gravity pivot and provides *effortless and maintenance-free performance*. The JR is a durable chair and comes with a growing frame.



Hand  
Propulsion



Foot  
Propulsion



Self-Tilt



Desk/Table  
Accessible



Transit Tie  
Downs

## *Fuze JR* A durable pediatric tilt chair

### *Available options* See order form

- ✓ Transit tie downs WC-19 & ISO-7176-19
- ✓ 30° manual dynamic recline
- ✓ Adjustable width frame

visit [pdgmobility.com/jr](https://pdgmobility.com/jr)  
for more information



DYNAMIC TILT  
Up to 50°



WEIGHT CAPACITY  
250 lbs



SEAT WIDTH  
12" to 18"



SEAT DEPTH  
13" to 18"



SEAT HEIGHT  
13" to 20"



BACK SEAT ANGLE  
84° to 120°



OVERALL WIDTH  
+7.5" to 8.5"



CHAIR WEIGHT  
42 lbs