

# EASY-OUT



*Easy-Out* by PDG is a new generation of plug-in front rigging. Its design features a low profile and small on-chair footprint. Removing and refitting is simple and quick, plus once the Easy-Outs come off, nothing remains on the chair.

## Exclusive features

- ✓ Effortless insertion and removal with a single lever
- ✓ When the hangers are removed, no hardware or brackets remain on the chair
- ✓ Tool-free leg length adjustment with a single lever
- ✓ The Angle Adjustable option offers angles from 0° (horizontal) to 110° (contracture)

## Easy-Out Fixed 70

- ✓ 70° fixed angle

## Easy-Out Universal Angle Adjustable

- ✓ Angles from 0° (horizontal) to 110° (contracture)
- ✓ Angle can be changed with a simple allen key



Side-view of chair with and without the Easy-Outs fitted



110° range of Angle Adjustable model



Tool-free leg length adjustment

Available on all  
PDG chairs  
(except Stellar HD and Elevation)

visit [pdgmobility.com/easy-out](https://pdgmobility.com/easy-out) for more information