

Sleek Style Luxurious Features



GO GO
ELITE TRAVELLER LX



Every day is a new adventure with the Go-Go Elite Traveller LX®, an exceptionally stable and manoeuvrable travel scooter. With easy disassembly and a convenient storage basket, the Go-Go Elite Traveller LX keeps up with your active life!

PRIDE[®]
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST*

www.pride-mobility.co.uk

Specifications


Weight Capacity	19 st 9 lbs (125 kg)
Maximum Speed ¹	Up to 4 mph
Ground Clearance ³	38 mm (1.5")
Turning Radius ³	1121 mm (44.25")
Overall Length ³	1003 mm (39.5")
Overall Width ³	495 mm (19.5")
Max. Obstacle Climb	50 mm (2")
Drive Tyres	Front: 178 mm (7") - solid Rear: 203 mm (8") - solid
Suspension	Limited
Total Weight Without Batteries ^{3, 4}	35.75 kg (79 lbs)
Heaviest Piece When Disassembled ³	15.25 kg (33.5 lbs) - front section
Drivetrain	24 volt DC motor
Dual Breaking System	Regenerative & electromechanical
Battery Requirements ³	2 x 12 Ah deep cycle (standard) 2 x 18 Ah deep cycle (optional)
Battery Weight ⁴	12 Ah - 8.75 kg (19.5 lbs) 18 Ah - 13.5 kg (29.5 lbs)
Battery Charger	2 amp off-board
Colour	Red & blue

Features

- Feather-touch disassembly with one hand
- Front and rear LED lighting
- Wrap around delta tiller
- Seat post suspension
- Includes 2 sets of interchangeable shrouds
- Exclusive & stylish low profile tyres
- Convenient front basket

Also available in:

- Red
- Blue



¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

² Includes standard seat weight.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.