



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

JANUARY 2022
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 3NG Comfort

Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 19
- Self Propulsion or Attendant Propulsion
- Choice of 4 seat widths (16" - 20") 405mm to 505mm
- Choice of 2 seat depths (16" or 18") 400mm or 450mm
- Choice of 3 seat height (16" - 18") 410mm to 460mm
- Black nylon seat upholstery
- Reclining (Ratchet) backrest 0° - 30° with Laguna shell 20° (510mm)
- Comfort Black multi stretch PU coated fabric Passad 2 backrest
- Elevating legrest hangers 0° to 75° including calf pad
- Flip up, angle and depth adjustable footplates
- Removable, height and depth adjustable comfort armrests
- Black multi stretch PU coated fabric trunk supports and headrest
- Choice of 8 colours (with 7 special colour options)
- Push or Pull to lock brakes (configuration dependant)
- Comfort transportation kit and buckle belt (adjustable both sides)
- Adjustable heel straps
- User tool kit



Standard Base price: £1,044

Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
CEG0005	ACTION 3NG COMFORT	☐
CEGSPEC	SPECIAL	POA
PROPULSION*		
CEG0100	SELF PROPULSION (NC 1062, 1074)	£1,044
CEG0110	TRANSIT PROPULSION (NC 0290, 0300, 1068, 1150, 1180, 1260, 1275, 1390)	£1,044

SEAT WIDTH*				
CEG0190		16" / 405 MM SEAT WIDTH		○
CEG0200		17" / 430 MM SEAT WIDTH		◆
CEG0210		18" / 455 MM SEAT WIDTH		○
CEG0230		20" / 505 MM SEAT WIDTH		○
SEAT DEPTH*				
CEG0240		16" / 400 MM SEAT DEPTH		○
CEG0250		18" / 450 MM SEAT DEPTH		◆
SEAT HEIGHT*				
CEG0280		18" / 460 MM SEAT HEIGHT		◆
CEG0290		17" / 435 MM SEAT HEIGHT	(NC 1010, 1062)	○
CEG0300		16" / 410 MM SEAT HEIGHT	(NC 1010, 1062, 1370, 1440, 1441)	○
SEAT UPHOLSTERY*				
CEG0340		NYLON WITH HOOK AND LOOP STRIPS		◆
CEG0345		TENSION ADJUSTABLE SEAT UPHOLSTERY & DEPTH ADJUSTABLE -50MM		○
CEG0372		CUSHION COMFORT MISTRAL 2 LAYER (INCLUDING SEAT UPHOLSTERY 3 VELCROS)		○
BACKREST*				
CEG0441		RECLINING (RATCHET) 0/30° LAGUNA SHELL 20"		◆
BACKREST UPHOLSTERY*				
CEG0490		COMFORT BACKREST LAGUNA BLACK MULTI STRETCH PU COATED 20" / 510 MM HEIGHT		○
CEG0491		COMFORT BACKREST PASSAD 2 BLACK MULTI STRETCH PU COATED 20" / 510 MM HEIGHT		◆
LEGREST HANGER RH*				
CEG0520		ELEVATING LEGREST RH SWING IN/SWING OUT	(NC 66666)	◆
CEG0545		AMPUTEE LEGREST RH	(NC 0610, 1111)	£124
LEGREST HANGER LH*				
CEG0570		ELEVATING LEGREST LH SWING IN/SWING OUT	(NC 55555)	◆
CEG0595		AMPUTEE LEGREST LH	(NC 0630, 2222)	£124
FOOTREST RH*				
CEG0610		FLIP UP DEPTH AND ANGLE ADJ RH		◆
SWO1111		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)	(NC 0545)	◆
CEE66666		NO FOOTREST		○
FOOTREST LH*				
CEG0630		FLIP UP DEPTH AND ANGLE ADJ LH		◆
SWO2222		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)	(NC 0595)	◆
CEE55555		NO FOOTREST		○
ARMREST RH*				
CEG0691		REMOVABLE RH (N°4) COMFORT ARM		◆
CEG0732		REMOVABLE RH (N°4) HEMI ARM	(NC 1481)	£55
ARMREST LH*				
CEG0791		REMOVABLE LH (N°4) COMFORT ARM		◆
CEG0832		REMOVABLE LH N°4 HEMI ARM	(NC 1481)	£55

OTHER SEATING OPTIONS			
CEG0881		TRUNK SUPPORT COMFORT BLACK MULTI STRETCH PU COATED FABRIC (PAIR) ONLY FOR LAGUNA SHELL	◆
CEG0906		HEADREST BLACK MULTI STRETCH PU COATED FABRIC	◆
CEG0937		BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED	£109
CEG0938		BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED	£109
COLOUR*			
CEG0961		SLATE GREY (RAL7024 FINE TEXTURE)	○
CEG0964		PEARL WHITE	○
CEG0965		ELECTRIC GREEN	○
CEG0966		MATT BLACK	○
CEG0967		BLUE POP	○
CEG0968		HYPNOTIC PURPLE	○
CEG0969		RED MATT	○
CEG0972		WATER BLUE	○
CEG0980		ANTHRACITE BLACK	○
CEG0991		SPECIAL COLOUR ON DEMAND ONLY	£161
SWO3333		DARK RED - SPECIAL	£161
SWO4444		PEARL GREY - SPECIAL	£161
SWO5555		INTENSE YELLOW	£161
SWO6666		VIBRANT ORANGE	£161
SWO7777		POLAR BLUE	£161
SWO8888		CANDY PINK	£161
CASTOR HOUSING*			
CEG0995		STANDARD	◆
CEG0996		LIGHT WEIGHT	£28
FRONT CASTORS*			
CEG1010		8" X 1.1/4" BLACK SOLID TYRE	◆
CEG1050		6" X 1.1/4" BLACK SOLID TYRE	○
REAR WHEEL SIZE*			
CEG1062		12" (NC 1071, 1072, 1084, 1088, 1140, 1150, 1180, 1260, 1390)	○
CEG1068		24"	◆
REAR WHEELS*			
CEG1071		STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES) (NC 1280)	◆
CEG1072		SELF-PROPEL REINFORCED REAR WHEEL (NC 1280)	£22
CEG1073		DRUM BRAKE SELF-PROPEL 22 AND 24" or TRANSIT 12" (NC 1201, 1270, 1275)	○
CEG1074		TRANSIT BLACK (NC 1088, 1260, 1280)	○
TYRES*			
CEG1082		SOLID BLACK	◆
CEG1084		PNEUMATIC PROFILE BLACK	£20
CEG1088		PNEUMATIC LOW PROFILE BLACK	○
HANDRIM			
CEG1140		ALUMINIUM PRE-ANODIZED HANDRIM	◆
CEG1150		HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC	£31

REAR WHEEL AXLE*			
CEG1170		FIXED AXLE	○
CEG1180		Q/R AXLE	◆
REAR WHEEL POSITION*			
CEG1200		PASSIVE REAR WHEELS POSITION	◆
CEG1201		DUAL AMPUTEE REAR WHEEL POSITION (NC 1275, 1280)	£69
PUMP			
CEG1230		PUMP AND SUPPORT	£6
CHASSIS OPTION			
CEG1241		FRONT WHEELS EXTENSION	£21
CEG1260		SPOKE GUARDS PAIR 24" BLK	£73
BRAKES*			
CEG1270		PUSH TO LOCK	◆
CEG1275		PULL TO LOCK BRAKES	○
CEG1280		ATTENDANT HUB BRAKES	£220
COMFORT TRANSPORTATION KIT*			
CEG1346		COMFORT TRANSPORTATION KIT	◆
SAFETY OPTIONS			
CEG1361		BODYPOINT 2-PT POSTURE BELT M (NC 1362)	£76
CEG1362		BODYPOINT 2-PT POSTURE BELT L	£76
CEG1370		REAR ANTI-TIPPERS PAIR (NC 1375, 1440, 1441)	○
CEG1375		REAR ANTI-TIPPERS PAIR WITH TIPPING AID (NC 1440, 1441)	£67
CEG1390		BRAKE EXTENSION PAIR	£14
OTHER OPTIONS			
CEG1430		CANE HOLDER	£22
CEG1440		TIPPING AID RH	£14
CEG1441		TIPPING AID LH	£14
CEG1481		COMFORT ACRYLIC TRAY	£97
SWO10001		ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618	£102

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:	
CEG0290	17" / 435 MM SEAT HEIGHT
CEG0300	16" / 410 MM SEAT HEIGHT
CEG0545	AMPUTEE LEGREST RH
CEG0595	AMPUTEE LEGREST LH
CEG0610	FLIP UP DEPTH AND ANGLE ADJ RH
SWO1111	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEE66666	NO FOOTREST
CEG0630	FLIP UP DEPTH AND ANGLE ADJ LH
SWO2222	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEE55555	NO FOOTREST
CEG1010	8" X 1.1/4" BLACK SOLID TYRE
CEG1062	12"
CEG1068	24"
CEG1071	STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES)
CEG1072	SELF-PROPEL REINFORCED REAR WHEEL
CEG1074	TRANSIT BLACK
CEG1084	PNEUMATIC PROFILE BLACK
CEG1088	PNEUMATIC LOW PROFILE BLACK
CEG1140	ALUMINIUM PRE-ANODIZED HANDRIM
CEG1150	HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC
CEG1180	Q/R AXLE
CEG1201	DUAL AMPUTEE REAR WHEEL POSITION
CEG1260	SPOKE GUARDS PAIR 24" BLK
CEG1270	PUSH TO LOCK
CEG1275	PULL TO LOCK BRAKES
CEG1280	ATTENDANT HUB BRAKES
CEG1362	BODYPOINT 2-PT POSTURE BELT L
CEG1370	REAR ANTI-TIPPERS PAIR
CEG1375	REAR ANTI-TIPPERS PAIR WITH TIPPING AID
CEG1390	BRAKE EXTENSION PAIR
CEG1440	TIPPING AID RH
CEG1441	TIPPING AID LH
CEG1481	COMFORT ACRYLIC TRAY



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Fax:

Description of special request (please include relevant option code or part number):