



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Power Retail Prescription Form

05 JANUARY 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Esprit Action Junior

Standard Features

- Max User Weight 80kg (12.6 Stone)
- Crash tested ISO7176 - 19*
- Single crossbar (Width Adjustable)
- Seat widths 305mm to 380mm 12" - 15"
- Seat depth 325mm (13") to 375mm (15") adjustable (set at 350mm)
- Seat height 460mm (18")
- Nylon seat upholstery
- FIXED 7° BACKREST 430mm backrest
- Standard nylon sling type backrest upholstery
- 80° swing in / swing out legrest hangers
- Flip up angle and depth adjustable footplates
- Height adjustable, removable arm armrests
- Choice of 3 standard colours (with 10 Special colour options)
- *
- * Battery MUST be removed when used as a seat in a vehicle

Standard Base price: £3,706



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
SPD0050	ESPRIT ACTION JUNIOR (SINGLE CROSSBAR) (NC 2300)	£3,706
CHAIR TYPE*		
SPDSPEC	SPECIAL	POA
SPD0100	BASE MODEL ESPRIT ACTION	◆
CHASSIS*		
SPD0200	SINGLE CROSSBAR (SW 305 MM TO 380 MM / SD 325 MM TO 375 MM) USER MAX WEIGHT = 80Kg	◆

BRACKET TYPE*			
SPD0300		STANDARD SWING AWAY BRACKET	○
SPD0310		ALBER E35 SWIVEL ARM	£360
USER REMOTE*			
SPD0400		USER REMOTE CONTROL RH	○
SPD0410		USER REMOTE CONTROL LH	○
DUAL REMOTE			
SPD0440		INTUITIVE DUAL REMOTE CONTROL RH	£853
SPD0450		INTUITIVE DUAL REMOTE CONTROL LH	£853
SPD0460		INTUITIVE SINGLE DUAL REMOTE CONTROL (ONE HANDLE) (ADJUSTABLE RH TO LH SIDE)	£608
SPEED*			
SPD0500		6 KM/H	◆
BATTERY*			
SPD0600		BATTERY 7,65AH/36V (NON SPILLABLE LEAD ACID BATTERY)	◆
POWER OPTIONS			
SPD0610		ADDITIONAL BATTERY 7,65AH/36V	£271
SPD0630		REMOTE PROTECTOR	£62
SEAT WIDTH*			
SPD1000		305 MM SEAT WIDTH	○
SPD1010		330 MM SEAT WIDTH	○
SPD1020		355 MM SEAT WIDTH	○
SPD1030		380 MM SEAT WIDTH	○
SEAT DEPTH*			
SPD1200		325 TO 375 MM ADJUSTABLE SEAT DEPTH SET AT 350 MM	○
SEAT HEIGHT*			
SPD1300		460 MM SEAT HEIGHT	○
SPD1310		485 MM SEAT HEIGHT	○
SEAT UPHOLSTERY TYPE*			
SPD1400		NYLON WITH LOOP BANDS STRIPS SEAT UPHOLSTERY, WIDTH & DEPTH ADJUSTABLE (SINGLE CROSSBAR ONLY)	◆
BACKREST TYPE*			
SPD2000		FIXED 7° BACKREST 2 HEIGHTS 400/430 MM (DELIVERED 430 MM)	◆
SPD2020		FOLDING BACKREST 430 MM HEIGHT	£32
SPD2040		ANGLE ADJUSTABLE AND FOLD DOWN (-15° TO + 15°) FIXED 7° HEIGHT 425-450 MM / 17"-18" SET AT 425 MM / 17" (2 POSITIONS; 25 MM / 1" STEPS THROUGH UPHOLSTERY)	£97
SPD2050		RECLINING/FOLDING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£64
SPD2060		RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£124
SPD2080		RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£51
BACKREST UPHOLSTERY TYPE*			
SPD2300		LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER	◆
SPD2320		CONTOUR ADJUSTABLE (FULL STRAPS) BI MATERIAL COVER	○

HANGERS RH*			
SPD3000		80° SHORT, UL FROM 200 TO 330MM	◆
HANGERS LH*			
SPD3100		80° SHORT, UL FROM 200 TO 330MM	◆
FOOTPLATES RH*			
SPD3200		FLIP UP RH FOOTPLATE INCLUDING HEEL STRAP	○
SPD3210		FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP	£42
SPD3220		FLIP UP ONE PIECE FOOTPLATE, DEPTH & ANGLE ADJUSTABLE (NC 3310)	£76
FOOTPLATES LH*			
SPD3300		FLIP UP LH FOOTPLATE, INCLUDING HEEL STRAP	○
SPD3310		FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP	£42
SPD9925		No LH footplate required due to FLIP UP ONE PIECE BEING SELECTED (NC 3210)	○
ARMRESTS RH*			
SPD4000		REMOVABLE RH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	○
SPD4010		REMOVABLE RH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (+35 MM BACKWARD) & WIDTH (+25MM OUTWARDS)	◆
ARMRESTS LH*			
SPD4100		REMOVABLE LH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	○
SPD4110		REMOVABLE LH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (35 MM BACKWARD) & WIDTH (+25MM OUTWARDS)	◆
SEAT CUSHION*			
SPD5000		ESSENTIAL CUSHION 375 MM DEPTH	○
SEATING OPTIONS			
SPD5110		HEADREST BLACK PU TR26	£153
SPD5120		HEADREST WITH CHEEK SUPPORT BLACK PU TR26 (INCLUDING BACKBRACE EXCEPT FOR LAGUNA)	£177
SPD5140		LATERAL POSITIONING PADS (PAIR)	£84
SPD5150		ADJUSTABLE CALF STRAP	○
BACK BRACE			
SPD5160		BACK BRACE WIDTH ADJUSTABLE	○

COLOURS*			
SPD6000	☐	ANTHRACITE BLACK	○
SPD6015	☐	DARK PURPLE	○
SPD6020	☐	WATER BLUE	○
SWO6529	☐	POP BLUE - 3 week lead time	£164
SWO6259	☐	SLATE GREY - 3 week lead time	£164
SWO0806	☐	Electric Green - 3 week lead time	£164
SWO6033	☐	Red Matt - 3 week lead time	£164
SWO0607	☐	Dark Red - 3 week lead time	£164
SWO0515	☐	Pearl Grey - 3 week lead time	£164
SWO0647	☐	Intense Yellow - 3 week lead time	£164
SWO0648	☐	Vibrant Orange - 3 week lead time	£164
SWO0651	☐	Polar Blue - 3 week lead time	£164
SWO0652	☐	Candy Pink - 3 week lead time	£164
SPD6050	☐	SPECIAL COLOUR (ON DEMAND ONLY)	POA
ACCENT COLOUR*			
SPD6100	☐	BLACK	○
SPD6120	☐	LIME	○
FRONT WHEELS*			
SPD7000	☐	200 X 50 MM (8" X 2") BLACK SOLID TYRE	£47
SPD7010	☐	200 X 50 MM (8" X 2") BLACK PNEUMATIC TYRE	◆
REAR WHEELS*			
SPD7100	☐	355 X 70 MM (14" X 27") BLACK SOLID WITH PUNCTURE PROOF TYRE	£82
SPD7110	☐	355 X 70 MM (14" X 27") BLACK PNEUMATIC TYRE	◆
ANTI TIPPER*			
SPD8000	☐	REAR ANTI-TIPPERS, FIXED, PAIR	◆
SPD8010	☐	REAR ANTI-TIPPERS REMOVABLE WITH CRUTCHES (E25), PAIR	£276
LAP BELT*			
SPD8020	☐	POSTURE BELT	◆
TRANSPORTATION BRACKET*			
SPD8030	☐	TRANSPORTATION BRACKET	◆
OTHER OPTIONS			
SPD8100	☐	CANE HOLDER	£20
SPD8200	☐	ACRYLIC TRAY FLIP UP	£308
TOTAL PRICE			

Non Conforming Options:

SPD2300	LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER
---------	--

SPD3210	FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP
---------	---

SPD3310	FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP
---------	---



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):