



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

JANUARY 2022
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Azalea Minor

Standard Features

- Max User Weight 75KG/12 Stone
- Crash Tested ISO 7176-Part 19
- Transit Propulsion with Attendant Hub Brakes and Black Solid Tyres
- Seat Width 340mm (+50mm or +100mm Width Adjustment)
- Choice of 2 Seat Heights 400mm or 450mm
- Depth Adjustable Seat 380mm-450mm
- Backrest Plate with Braced Push Handles
- Assistant Operated Gas Piston Seat Tilt +1° - +25° and Backrest Recline 30° with Tilt Lock
- Angle Adjustable Legrests with Calf Support
- Height Adjustable Wide Pad Armrests
- Pearl Silver Frame Colour
- User Manoeuvred Tyre Brakes
- Pelvic Belt
- Anti Tippers

Standard Base price: £2,162



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
BNC0014	Azalea Minor	£2,162
BNC8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	POA
PROPULSION*		
BNC0020	Transit 16" wheels with hub brakes (NC 1055, 1062, 1074, 1100, 1120)	<input checked="" type="radio"/>
BNC0025	Self Propel 24"/22" wheels with user brakes (NC 0075, 0605, 0635, 1047, 1076, 1305)	<input type="radio"/>
BNC0030	Self Propel 24"/22" wheels with assistant hub brakes (NC 0075, 0605, 0635, 1047, 1074, 1305)	£133

SEAT DIMENSIONS*				
BNC0037		Minor (- 50 mm In seat depth and - 30 mm In chassis length)		◆
SEAT WIDTH*				
BNC0038		340 mm/13" (Minor)		◆
SEAT WIDTH EXTENSION*				
BNC0065		+ 50 mm seat width pre-adjustment on armrests and legrests		○
BNC0075		+ 100 mm seat width pre-adjustment on armrests and legrests	(NC 0110)	○
BNC0085		Without seat width extension, set at 0mm adjustment but can be adjusted wider if required		○
SEAT DEPTH*				
BNC0095		380-450 mm / 15" - 18" (Azalea Minor)		◆
SEAT HEIGHT*				
BNC0110		400 mm / 16"	(NC 1062)	○
BNC0120		450 mm / 18"	(NC 1055)	○
SEAT ANGLE*				
BNC0000		Standard Seat Angle +1° to +25°		◆
BNC0127		Negative Seat Angle -1 To +23°		○
SEAT TYPE*				
BNC0130		Depth Adjustable		◆
SEAT CUSHION*				
BNC0165		Cushion Flosshape seat upholstery	(NC 1111)	£54
BNC0180		Without Seat Cushion, Seat Plate only	(NC 0201, 0202)	◆
SEAT UPHOLSTERY*				
BNC0201		Black Multi-stretch Black PU coated fabric		◆
BNC0202		Black Spacer TR36		○
BNC1111		No upholstery, select with options 0170 Multifunctional and 0180 Without seat cushion		○
BACKREST*				
BNC0217		Flex 2 Backrest Tension adjustable (Azalea Minor)	(NC 0925)	£131
BNC0218		Flex 2 Backrest Fixed wooden plate (Azalea Minor)	(NC 0240, 0905)	◆
BACKREST CUSHION*				
BNC0240		Cover for tension adjustable upholstery		○
BNC0250		Cushion, Laguna (Lateral Support)	(NC 0905)	○
BNC0260		Cushion, Mistral 2 (Waist Support)	(NC 0322, 0905)	£137
BACKREST UPHOLSTERY*				
BNC0311		Black Plush TR35		○
BNC0321		Black Multi-stretch Black PU coated fabric		◆
BNC0322		Black Spacer TR36		○
PUSH HANDLES*				
BNC0340		Push Handles braced (Azalea Minor)		◆
BNC0345		Push Bar, angle adjustable		£97
SEAT TILT AND BACKREST REGULATION*				
BNC0350		Assistant operated gas piston seat tilt and backrest recline		◆
TILT LOCK*				
BNC1263		Tilt Lock (locking of pre-set angle)		◆

LEGREST HANGER ATTACHMENT*			
BNC0375		Standard position for legrest mounting	◆
SWO0189		Standard position with 0-30° Legrest Abduction function	£95
BNC0378		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest only (NC 0390, 0410, 0440, 0460)	○
LEGREST HANGER LH*			
BNC0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH (NC 0475, 0480, 0490, 0942, 0945, 5555)	◆
BNC0390		Legrest Hanger fixed 80 deg, swing in/swing out LH (NC 0928, 0929, 0930)	○
BNC0410		Legrest Amputee LH (NC 0475, 0476, 0477, 0480, 0490, 0510, 0520, 0928, 0929, 0930, 0941, 0942, 0945)	£130
LEGREST HANGER RH*			
BNC0430		Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH (NC 0475, 0540, 0550, 0942, 0945, 5555)	◆
BNC0440		Legrest Hanger fixed 80 deg, swing in/swing out RH (NC 0931, 0932, 0940)	○
BNC0460		Legrest Amputee RH (NC 0475, 0476, 0477, 0540, 0550, 0570, 0580, 0931, 0932, 0940, 0941, 0942, 0945)	£130
FOOTREST LH*			
BNC0475		Footrest one piece, depth and angle adjustable short tube (NC 0477, 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940)	£112
BNC0476		Footrest one piece, depth and angle adjustable long tube (NC 0477, 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940)	£112
BNC0477		Footbox, angle adjustable, padded with short tube (NC 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940, 0941, 0942, 0945)	£516
BNC0480		Footrest lower part, fixed with short tube and heel strap, LH (NC 0930, 0941)	○
BNC0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.	○
BNC0510		Footrest lower part, fixed with long tube and heel strap, LH. (NC 0930, 0941)	○
BNC0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNC3333		Without footrest, only select with option 0410 Amputee legrest	○
FOOTREST RH*			
BNC0540		Footrest lower part, fixed with short tube and heel strap, RH (NC 0940, 0941)	○
BNC0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.	○
BNC0570		Footrest lower part, fixed with long tube and heel strap, RH (NC 0940, 0941)	○
BNC0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNC4444		Without footrest, only select with options 0460 Amputee legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○

ARMREST LH*				
BNC0600		Armrest Ducktail Pad LH, height adjustable with autolock	(NC 1232)	○
BNC0605		Armrest Wide Pad LH, height adjustable with autolock. Compatible with 16" wheels only	(NC 1232)	◆
BNC0607		Armrests Reinforced Wide Pad, LH, height adjustable with autolock	(NC 1230, 1231, 1232, 1260)	£85
BNC0610		Armrest Padded LH, height adjustable with autolock	(NC 1230, 1231, 1232, 1260)	£61
BNC0615		Armrest Hemiplegic LH, height adjustable with autolock	(NC 0695, 1230, 1231, 1232, 1260)	£170
BNC0620		Without Armrest LH Select when choosing BNC1232 Half Tray RH, Narrow armrest is included	(NC 0695, 1230, 1231, 1260)	○
ARMREST RH*				
BNC0630		Armrest Ducktail Pad RH, height adjustable with autolock	(NC 1233)	○
BNC0635		Armrest Wide Pad RH, height adjustable with autolock. Compatible with 16" wheels only	(NC 1233)	◆
BNC0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock	(NC 1230, 1231, 1233, 1260)	£85
BNC0640		Armrest Padded RH, height adjustable with autolock.	(NC 1230, 1231, 1233, 1260)	£61
BNC0645		Armrest Hemiplegic, RH, height adjustable with autolock	(NC 0695, 1230, 1231, 1233, 1240, 1260)	£170
BNC0650		Without Armrest RH Select when choosing BNC1233 Half Tray LH, Narrow armrest is included	(NC 0695, 1230, 1231, 1240, 1260)	○
ARMREST OPTIONS. Side support cover 0680, must be selected if 0690 Foam for side support is required				
BNC0680		Side Support Cover, removable with tools (wipeable cover for the sideguard)	(NC 0695)	£48
BNC0690		Foam for side support cover (Internal cushioned pad, which is only possible in combination with 0680)	(NC 0695)	○
BNC0695		Pad for sideguard, removable with hook and loop fastening		£48
SEAT OPTIONS				
BNC0710		Incontinence Cover		£38
BNC0720		Pommel Abduction Cushion		£112
BNC0725		Tilt Scale for seat tilt		£42
BACKREST OPTIONS				
BNC0730		Headrest, Sideway Adjustable	(NC 0750, 0770)	◆
BNC0750		Headrest with cheek support, sideway adjustable	(NC 0740, 0770, 0780)	£89
SWO0190		Junior Headrest, sideway adjustable. Pad size 230mm x 100mm	(NC 0730, 0740, 0750, 0770, 0780)	£89
BNC0740		Cover for headrest		£29
BNC0770		Neckrest		○
BNC0780		Cover For Neckrest		£21
BNC0790		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)		£46
BNC0900		Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included		○
BNC0905		Lateral Positioning Pads. Attaches with hook and loop between the strap and the cover of the tension adjustable upholstery		£58
BNC0910		Lumbar Support Cushion		£29
BNC0915		Tilt Scale for backrest recline		£42
BNC0920		Wedges (One Pair)		£10
BNC0925		Privacy Attachment (Material connecting backrest and seat)		£32

TRUNK SUPPORT LH				
BNC0820		Trunk Support Multi Adjustable with extended height arm LH	(NC 0900)	○
BNC0880		Trunk Support Swingaway LH	(NC 0900)	£100
BNC0895		Trunk Support Multi Adjustable Swingaway LH	(NC 0900)	£107
TRUNK SUPPORT RH				
BNC0830		Trunk Support Multi Adjustable with extended height arm RH	(NC 0900)	○
BNC0890		Trunk Support Swingaway RH	(NC 0900)	£100
BNC0896		Trunk Support Multi Adjustable Swingaway RH	(NC 0900)	£107
ELEVATING LEGREST OPTIONS*				
BNC0928		Soft foam padded PU fabric cover, for LH elevating legrest calf pad		◆
BNC0929		Large PU fabric coated knee pad for LH elevating legrest		◆
BNC0931		Soft foam padded PU fabric cover, for RH elevating legrest calf pad		◆
BNC0932		Large PU fabric coated knee pad for RH elevating legrest		◆
BNC5555		Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)		○
LEGREST OPTIONS				
BNC0930		Footplate Extender LH	(NC 0941)	£32
BNC0940		Footplate Extender RH	(NC 0941)	£32
BNC0941		Cover Footplate, Black Plush Fabric		£51
BNC0942		Calf Strap Padded		£69
BNC0945		Calf Strap (for One Piece Footrest)		○
CHASSIS COLOUR*				
BNC0944		Silver Pearl		◆
BNC0946		Azurite (Dark Blue)		○
BNC0947		Sand (Beige)		○
BNC0948		Electric Green		○
BNC0949		Happy Red		○
CASTOR ATTACHMENT*				
BNC0960		Fixed castor attachment		◆
FRONT CASTORS*				
BNC0982		140x45 mm (5 1/2" x 1 3/4") Medium Hard Black tyre		£19
BNC0992		150x30 mm (6" x 1 1/4") Medium Hard Black tyre		◆
REAR WHEEL SIZE*				
BNC1047		16" / 405 mm Transit rear wheel with fixed attachment	(NC 1074, 1082, 1084, 1100, 1120, 1160)	◆
BNC1055		22" / 560 mm Self Propel rear wheel with QR attachment		○
BNC1062		24" / 600 mm Self Propel rear wheel with QR attachment		○
REAR WHEEL TYPE*				
BNC1074		Standard self propel with colourless hub, clear anodized rim and Stainless steel spokes		○
BNC1076		Hub Brake 16" Transit or 22" or 24" Self Propel Transit with Black Plastic Hub and Rim Self Propel with colourless hub, clear anodized rim and Stainless steel spokes		○

REAR WHEEL TYRE*			
BNC1082	<input type="checkbox"/>	Pneumatic Profile Black 1 3/8" tyres	○
BNC1084	<input type="checkbox"/>	Pneumatic Puncture proof / flatfree 1 3/8" Black tyres	£39
BNC1092	<input type="checkbox"/>	Solid Black tyres	◆
HANDRIM			
BNC1100	<input type="checkbox"/>	Aluminium handrim	◆
BNC1120	<input type="checkbox"/>	Plastic coated handrim	£46
SPOKEGUARD			
BNC1160	<input type="checkbox"/>	Frost spokeguard	£158
BRAKE*			
BNC1170	<input type="checkbox"/>	User manoeuvred tyre brakes	◆
PELVIC BELT*			
BNC1200	<input type="checkbox"/>	Pelvic belt	◆
SAFETY OPTIONS			
BNC1220	<input type="checkbox"/>	Reflector Kit (wheels and chassis)	£21
TRANSPORTATION BRACKETS*			
BNC1225	<input type="checkbox"/>	Transportation Brackets, mounted to the front of the wheelchair	◆
OTHER OPTIONS			
BNC1230	<input type="checkbox"/>	Fixed Table Tray without elbow pads includes mounting kit (NC 1237, 1260)	£131
BNC1231	<input type="checkbox"/>	Fixed Table Tray with elbow Pads includes mounting kit (NC 1237, 1260)	£131
BNC1232	<input type="checkbox"/>	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1235, 1260)	£98
BNC1233	<input type="checkbox"/>	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1235, 1260)	£98
BNC1240	<input type="checkbox"/>	Swing away table without elbow pads (NC 1237, 1260)	£144
BNC1235	<input type="checkbox"/>	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)	£48
BNC1237	<input type="checkbox"/>	Half table tray padded cushion	£48
BNC1260	<input type="checkbox"/>	Mounting kit for table (if ordered separately)	£8
BNC1270	<input type="checkbox"/>	Pump	£21
BNC1280	<input type="checkbox"/>	Cane Holder	£47
BNC1290	<input type="checkbox"/>	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	£52
BNC1300	<input type="checkbox"/>	Harness attachment , must be selected in combination with Head/Neckrest	£81
BNC1305	<input type="checkbox"/>	Front castor locks, includes front tiedown brackets. Allows castors to remain parallel for easy access into confined spaces such as lifts or vehicles	£174
BNC1310	<input type="checkbox"/>	Attachment for positioning belt	£30

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:	
BNC0075	+ 100 mm seat width pre-adjustment on armrests and legrests
BNC0110	400 mm / 16"
BNC0201	Black Multi-stretch Black PU coated fabric
BNC0202	Black Spacer TR36
BNC1111	No upholstery, select with options 0170 Multifunctional and 0180 Without seat cushion
BNC0240	Cover for tension adjustable upholstery
BNC0322	Black Spacer TR36
BNC0390	Legrest Hanger fixed 80 deg, swing in/swing out LH
BNC0410	Legrest Amputee LH
BNC0440	Legrest Hanger fixed 80 deg, swing in/swing out RH
BNC0460	Legrest Amputee RH
BNC0475	Footrest one piece, depth and angle adjustable short tube
BNC0476	Footrest one piece, depth and angle adjustable long tube
BNC0477	Footbox, angle adjustable, padded with short tube
BNC0480	Footrest lower part, fixed with short tube and heel strap, LH
BNC0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.
BNC0510	Footrest lower part, fixed with long tube and heel strap, LH.
BNC0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNC0540	Footrest lower part, fixed with short tube and heel strap, RH
BNC0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.
BNC0570	Footrest lower part, fixed with long tube and heel strap, RH
BNC0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNC0605	Armrest Wide Pad LH, height adjustable with autolock. Compatible with 16" wheels only
BNC0635	Armrest Wide Pad RH, height adjustable with autolock. Compatible with 16" wheels only
BNC0695	Pad for sideguard, removable with hook and loop fastening
BNC0730	Headrest, Sideway Adjustable
BNC0750	Headrest with cheek support, sideway adjustable
BNC0740	Cover for headrest
BNC0770	Neckrest
BNC0780	Cover For Neckrest
BNC0900	Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included
BNC0905	Lateral Positioning Pads. Attaches with hook and loop between the strap and the cover of the tension adjustable upholstery
BNC0925	Privacy Attachment (Material connecting backrest and seat)
BNC0928	Soft foam padded PU fabric cover, for LH elevating legrest calf pad
BNC0929	Large PU fabric coated knee pad for LH elevating legrest
BNC0931	Soft foam padded PU fabric cover, for RH elevating legrest calf pad
BNC0932	Large PU fabric coated knee pad for RH elevating legrest
BNC5555	Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)
BNC0930	Footplate Extender LH
BNC0940	Footplate Extender RH
BNC0941	Cover Footplate, Black Plush Fabric
BNC0942	Calf Strap Padded

BNC0945	Calf Strap (for One Piece Footrest)
BNC1047	16" / 405 mm Transit rear wheel with fixed attachment
BNC1055	22" / 560 mm Self Propel rear wheel with QR attachment
BNC1062	24" / 600 mm Self Propel rear wheel with QR attachment
BNC1074	Standard self propel with colourless hub, clear anodized rim and Stainless steel spokes
BNC1076	Hub Brake 16" Transit or 22" or 24" Self Propel Transit with Black Plastic Hub and Rim Self Propel with colourless hub, clear anodized rim and Stainless steel spokes
BNC1082	Pneumatic Profile Black 1 3/8" tyres
BNC1084	Pneumatic Puncture proof / flatfree 1 3/8" Black tyres
BNC1100	Aluminium handrim
BNC1120	Plastic coated handrim
BNC1160	Frost spokeguard
BNC1230	Fixed Table Tray without elbow pads includes mounting kit
BNC1231	Fixed Table Tray with elbow Pads includes mounting kit
BNC1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BNC1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BNC1240	Swing away table without elbow pads
BNC1235	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)
BNC1237	Half table tray padded cushion
BNC1260	Mounting kit for table (if ordered separately)
BNC1305	Front castor locks, includes front tiedown brackets. Allows castors to remain parallel for easy access into confined spaces such as lifts or vehicles



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Fax:

Description of special request (please include relevant option code or part number):