



TDX2 Low Rider

The new **TDX SP2 Low Rider** features a drop down seat unit to boast unrivalled low seat to ground heights of either 403mm or 420mm. This allows practical driving positions combined with the ability to easily access WAVs, fit under tables and align with a variety of transfer surfaces.

TDX SP2 is the power chair for individuals who require full functionality, performance and style combined! The most stable and high performance midwheel drive base on the market also has a compact turning radius of 560mm - making the **TDX SP2 Low Rider** equally suitable to outdoor and indoor environments.



Super Low Seat to Ground Height

The **TDX SP2 Low Rider** is available with either a 403mm or 420mm seat to ground height, utilising a drop down seat unit and 30° centre of gravity shift seat tilt within the base specification. The modern aesthetics of the **TDX2** can be complimented with integrated LED lights and a choice of seven exciting shroud colours.

Technical data

| | | | | | | |
|--------------------------------------|---------------------------|---|--------------------------------|----------------------------------|--|---------------------------|
| | | | | | | |
| Invacare TDX SP2 low rider | 330 - 580 mm ¹ | 410 - 510 mm 510 - 560 mm ² | 403 / 420mm ^{3,10} | 480, 540, 560 mm ⁶ | 275 - 340 mm | 290 - 460 mm ⁷ |
| | | | | | | |
| Invacare TDX SP2 low rider | 90° - 120° | 1010 - 1170 mm ¹⁰ | 1160 mm ⁸ | 980 mm | 650 mm ⁹ | 160 kg ⁴ |
| | | | | | | |
| Invacare TDX SP2 low rider | 1120 mm | 75 mm | 9° / 16 % | 2 x 340 W | 32 km (6 km/h) ¹² 39 km (8 km/h) ¹² | 6 and 8 km/h |
| | | | | | | |
| Invacare TDX SP2 low rider | | | | | | 154 kg ¹⁰ |

- 1/ 50 mm by armrest adjustment
- 2/ Extended profile rail package for larger seat depth
- 3/ ground to seat plate without cushion
- 4/ with seat lifter or powered tilt
- 5/ with lifter
- 6/ measured from seat plate
- 7/ measured from footplate to seat plate
- 8/ depending on seat position
- 9/ width of wheelbase
- 10/ depending on configuration
- 11/ Depending on seat to ground height
- 12/ according to ISO 7176-4:2008

- 13/ Invacare **Modulite™** seat components are cross compatible across the Invacare power wheelchair range and partly with the Invacare Rea™ family.
- 14/ Compliant with the directive on the Restriction of the use of certain Hazardous Substances in electrical and electronic equipment - RoHS (2011/65/EC).
- 15/ Compliant with the regulation concerning the Registration, Evaluation, Authorization and Restriction of Chemicals.

Invacare Limited
Pencoed Technology Park -
Pencoed - Bridgend CF35 5AQ -
United Kingdom
Tel: +44 1656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: ireland@invacare.com
www.invacare.ie

© 2016 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

TDX SP2 LR - UK - 07/2016 - 1626896

Invacare® **Modulite™**



Bariatric use



G-Trac™ option available



Shroud colours



Wet black

Electric red

Jade green



Electric blue

Electric silver

Arctic white

Transparent purple



Follow us on Twitter:
[@Think_Mobility](https://twitter.com/Think_Mobility)

Follow us on YouTube:
Invacare Ltd - Think Mobility

More detailed information, technical data, references for transport or description of certain options can be found in the wheelchair's user manual.

