

# **RETAIL PRESCRIPTION FORM**

# **JAY UNION RP**

(reduced profile - cushion height 9cm)

Invoice to:		Deliver to:		
Name:		Name:		
Street:		Street:		
Town:	Postal code:	Town:	Postal code:	
Tel:	Fax:	Tel:	Fax:	
		<del>.</del>		_
Order date:		Marked for:		
Purchase order:				
		Sunrise order no:		

### Jay® UNION (Reduced Profile) Features

Medium to High skin protection cushion

Deep pre-contoured foam cushion with viscoelastic top layer

Integral JAY Flow Fluid pad in seat well Water Resistant (Dartex) inner cover

Microclimatic outer cover as standard

Reduced profile (height) is recommended for medium to active manual chairs to enable a lower seat to ground height

2 Year guarantee

Price £250

E-mail: orders@sunmed.co.uk

Fax: 0845 605 66 89

Maximum user weight 136kg (SW 34cm-50cm)

227kg (SW 52-60cm)

## Application;

Clients with moderate positioning needs and at moderate to high risk of skin breakdown (medium to high risk skin protection).

Lateral pelvic wedges provides lateral stability and alignment of the pelvis. Contoured base encourages lateral and medial support for the thighs.

Reduced profile option to lower the overall height by 2,5cm (1") for a lower seat-to-floor height.

					C	LIEN	T M	EAS	URE	ME	NT					
				width C	JU <sup>36</sup>	N38 ,	P&B	6 July 68	2 JU	PR JUN	ko jut	D JUST	PR	pR July	\$ 100 \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \	
			C	cm	36	38	40	42	44	48	50	52	54	56	60	
		JU 36RP	3	36	$\Diamond$	$\Diamond$	$\Diamond$	$\Diamond$	$\Diamond$							Basic Price
		JU 38RP	3	38	$\Diamond$	0	<b>\( \)</b>	$\Diamond$	$\Diamond$	$\Diamond$	$\Diamond$					£250
		JU 40RP	4	40	0	0	0	0	0	$\Diamond$	$\Diamond$					
	_	JU 42RP	4	42	0	0	0	0	0	0	$\Diamond$					Product Code
	Ĕ	JU 44RP	4	44	$\Diamond$	$\Diamond$	0	0	0	0	$\Diamond$					JU RP
	Ä	JU 48RP	4	48		$\Diamond$	0	0	0	$\Diamond$	$\Diamond$					
	Ę	JU 50RP		50		$\Diamond$	<b>\( \)</b>	$\Diamond$	$\Diamond$	$\Diamond$	$\Diamond$					Lead Times
	SEAT DEPTH	JU 52RP		52												O 5 day lead time
	•	JU54RP		54												
		JU 56RP		56												
		JU 60RP	6	60												Not available
				•			•								••••••	
Example		JU	Width 44	4cm [	Depth	44cm		JU	1444	RP =	Con	plete	Cus	hion	(Reduced	Profile height)

SPARE CUSHION COVER			Quantity	Unit Price	it Price Item Code	
CIJU RPEU		Spare Incontinence Resistant Inner Cover		£51	CIJU	RPEU
CJURPMEU		Spare Microclimatic Outer Cover		£65	CJU	RPMEU
JAY UNION RP					£250	
Option Price					£	
VAT					£	
Total Price					£	



# **JAY UNION RP CUSHION**

Stability and positioning with skin integrity protection

### **Product Features:**

- · Medium to High Skin protection risk
- · Contoured foam base with anatomical seat well and femoral loading
- MemGel® soft memory foam top layer for even pressure distribution, heat dissipation and comfort
- · JAY Fluid Tripad for ischial immersion and envelopment
- · Incontinence resistant inner cover, Microclimatic outer cover
- Reduced height profile (1" lower height) is ideal for medium to active manual chair users allowing reduced seat to floor heights for easier transfers and to aid optimal propulsion
- Maximum User Weight 136kg upto SW 36-50cm, 227kg SW 52-60cm



<b>JAY UNION CUS</b>	HION	
JU3636RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 36x36	£250
JU3638RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 36x38	£250
JU3640RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 36x40	£250
JU3642RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 36x42	£250
JU3644RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 36x44	£250
JU3836RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 38x36	£250
JU3838RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 38x38	£250
JU3840RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 38x40	£250
JU3842RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 38x42	£250
JU3844RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 38x44	£250
JU3848RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 38x48	£250
JU3850RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 38x50	£250
JU4036RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 40x36	£250
JU4038RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 40x38	£250
JU4040RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 40x40	£250
JU4042RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 40x42	£250
JU4044RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 40x44	£250
JU4048RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 40x48	£250
JU4050RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 40x50	£250
JU4236RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 42x36	£250
JU4238RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 42x38	£250
JU4240RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 42x40	£250
JU4242RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 42x42	£250
JU4244RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 42x44	£250
JU4248RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 42x48	£250
JU4250RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 42x50	£250
JU4436RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 44x36	£250
JU4438RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 44x38	£250
JU4440RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 44x40	£250
JU4442RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 44x42	£250
JU4444RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 44x44	£250

JU4448RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 44x48	£250
JU4450RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 44x50	£250
JU4838RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 48x38	£250
JU4840RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 48x40	£250
JU4842RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 48x42	£250
JU4844RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 48x44	£250
JU4848RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 48x46	£250
JU4850RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 48x50	£250
JU5036RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 50x38	£250
JU5040RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 50x40	£250
JU5042RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 50x42	£250
JU5044RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 50x44	£250
JU5048RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 50x48	£250
JU5050RP	JAY UNION RP Cushion (Reduced Profile Height - 9cm) - 50x50	£250