

Elevate Your Life To The Next Level!



Jazzy
air₂



At the touch of a switch, the Jazzy Air® 2 elevates 12" in just 11 seconds. With a sleek, more modern look, this power chair safely drives at speeds up to 4 mph while elevated, allowing users to enjoy face-to-face social engagement.



Specifications

Weight Capacity	21 st 6 lbs (136 kg)
Maximum Speed ¹	Up to 4 mph
Ground Clearance ³	25 mm (1")
Turning Radius ³	528 mm (20.75")
Overall Length ³	1067 mm (41.93")
Overall Width ³	590 mm (23.19")
Max. Obstacle Climb	55 mm (2.2")
Drive Tyres	Drive: 254 mm (10") Flat Free Castors: 152 mm (6") Flat Free
Suspension	Active-Trac Suspension with Mid-Wheel 6 Technology
Total Weight Without Batteries ^{3, 4}	81.53 kg (180 lbs)
Heaviest Piece When Disassembled ³	81.53 kg (180 lbs)
Drivetrain	Two motor, mid-wheel-6 drive
Dual Breaking System	"Intelligent Braking" electronic regenerative, disc park brake
Battery Requirements ³	2 x 12v, deep cycle U-1 batteries
Battery Weight ⁴	11.5 kg (25.4 lbs) each
Battery Charger	5 Amp off-board
Colour	Citrine yellow, pearly white, ruby red, sapphire blue, black pearl, garnet red, black onyx, pink topaz, silver, tanzanite

Features

- Fast Rise: 12" of power adjustable seat height in 11 seconds
- Active-Trac Suspension with Mid-Wheel 6 Technology
- Bright LED lighting package
- Storage Compartment
- Angle, height and depth-adjustable armrest
- High-back or solid seat pan, reclining sport seat
- Integrated foot platform

Also available in:

- Citrine Yellow
- Pearly White
- Ruby Red
- Sapphire Blue
- Black Pearl (matte)
- Garnet Red (matte)
- Black Onyx
- Pink Topaz
- Silver
- Tanzanite (matte)



¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

² Includes standard seat weight.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.