

JOKER JUNIOR 2.0



Seat width	24 - 27 - 30 - 33 - 36 cm
Seat depth	27 - 30 - 32,5 - 35 - 37,5 cm
Seat height	Front from 36 to 53 cm (adjustable) Rear from 36 to 47 cm (adjustable)
Footplate distance	From 8 to 38 cm (adjustable)
Backrest height	From 29 to 51 cm (adjustable)
Front frame angle	95°
Total length	75 cm approx. (std. conf.)
Total width	Seat depth + 22,5 cm (with 0° camber)
Backrest angle	From 74° to 94° to the seat (adjustable)
Camber	0° - 2° - 4°
Weight without rear wheels	Starting from 7,3 kg approx.
Wheelchair weight	Starting from 9,9 kg approx.
Maximum load	75 kg