

JUNIOR LIGHT



Seat width	27 - 30 - 33 - 36 cm
Seat depth	30 - 32,5 - 35 - 37,5 cm
Seat height	Front from 41 a 50 cm (adjustable) Rear from 38 a 47 cm (adjustable)
Footplate distance	From 18 to 40 cm (adjustable)
Backrest height	From 27 to 47 cm (adjustable)
Front frame angle	100°
Total length	80 cm approx. (conf. std)
Total width	Seat width + 18 cm (with 0° camber)
Encumbrance	28 cm approx. (with 0° camber)
Backrest angle	90° to the seat
Weight without rear wheels	Starting from 9,7 kg approx.
Wheelchair weight	Starting from 13,0 kg approx.
Maximum load	80 kg