# Base Price - £3795.00 Valid from 01.04.22 to 31.03.23 Please ensure you are using the latest order form ☐ Provide estimate only □ Order ☐ Attachment for custom fabrication enclosed Billing address/customer number **Delivery address** Sales person Name Date Signature

By placing the order, you agree to comply with the instructions for use/service ma $nuals \ for \ accessories \ and \ products \ of \ other \ manufacturers \ when \ combining \ them$ with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products, and you confirm to have read, understood and accepted our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

#### Standard equipment

### Juvo B4

490E75=1\_AS01\_C

Power wheelchair [product name]

Max. load: 160 kg (standard equipment)







Additional options installed on the power wheelchair can reduce the maximum load. Please observe the maximum gross weight on the nameplate.

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders. Please use the spare parts catalogue for spare parts orders.

Basic data (volu	ntary	information)				Mandatory informatio	n*
Gender		Inc	dication			User height in cm	
Age		Pa	ying party _			User weight in kg	
	<u> </u>					1 Seat depth	
						2 Trunk depth	
	1		Swing			3 Arm support height	
		1 / / )	8	, \		4 Armpit height	
		7	(4)	<b>)</b>		5 Shoulder height	
6	5	4 3				6 Head height	
,	$\downarrow$					7 Forearm length to the m	niddle of the palm
		1	10			8 Trunk width	
		9				9 Lower leg length	
						10 Seat width	
		1	3				
* This information	will h	elp to support user-oriented	l settings. Otto	obock guarar	ntees the confi	idential and purely internal use of these	data.
	catio	n (mandatoy information)					
EA001		Category A (indoor)					£3795.00
EA002	Ш	Category B (in/outdoor)	(electric light	ting mandato	ory in Europe)		£3795.00
Max. user weig	ht (ma	andatoy information; user +	personal lugg	age, informa	tion required	for chassis suspension)	
		75 kg					
		100 kg (for category A up	to 12.5"/8")				
		113 kg (not possible in cor	nbination wit	h category A	(indoor) exclu	isions see 2.1)	
		125 kg (not possible in co	mbination wit	h category A	(indoor) exclu	usions see 2.1)	
		140 kg (not possible in co	mbination wit	h category A	(indoor) exclu	usions see 2.1)	
		150 kg (not possible in co	mbination wit	h category A	(indoor) exclu	usions see 2.1)	
		160 kg (not possible in co	mbination wit	h category A	(indoor) exclu	usions see 2.1)	
Seating system	ıs 💽	(mandatoy information)					
EC001		Standard seat Junior (ma	ax. load: 75 kg	)			£128.00
EC002		Standard seat Small					Inclusive
EC003		Standard seat, Large					Inclusive
EC004		Standard seat, XL					£165.00
EC005		VAS (Vario Adjust Seat)	(possible only	in combinati	ion with max.	load <= 125 kg) exclusions see 3.3)	£380.00
Settings for Sta	ındarı	d seat Junior exclusions see	e 3.2) 4.0) 5.0)				
Seat depth SD		□34	☐36	38	☐ 40		
Seat width SW	in cm	□ 34	□36		40		
		BH 47 cm (adjustable to 42					
Settings for Sta	ındarı	d seat Small exclusions see	3.2) 4.0) 4.1) 4	.4) 5.0)			
Seat depth SD		□38	<b>4</b> 0	42	44	46	
Seat width SW		□ 38*	40	42			
			* Possible o		ination with El	L022	
Back support h	eight	BH 55 cm (adjustable to 50	cm or 45 cm)				

Settings for Star	Settings for Standard seat Large exclusions see 4.0) 5.0)												
Seat depth SD is	n cm	42   44   46   48   50											
Seat width SW i	n cm	□44 □46 □48											
Back support he	eight	BH 55 cm (adjustable to 50 cm or 45 cm)											
Settings for Star	ndar	seat XL exclusions see 4.0) 5.0)											
Seat depth SD i	n cm	42   44   46   48   50											
Seat width SW i	n cm	□ 50 □ 52 □ 54 □ 56											
Back support he	eight	BH 55 cm (adjustable to 50 cm or 45 cm)											
Settings for VAS	5 excl	usions see 4.0) 5.0)											
Seat depth SD is	n cm	□ 38     □ 40     □ 42     □ 44     □ 46     □ 48       □ 50     □ 52     □ 54     □ 56     □ 58											
Seat width SW i	n cm	□ 38     □ 40     □ 42     □ 44     □ 46     □ 48       □ 50     □ 52     □ 54											
Back support he	eight	BH 55 cm (adjustable to 53 cm, 51 cm or 45 cm)											
Lower leg lengt	th LL												
Seat cushions (	<b>①</b> (n	nandatory information)											
EC011		Without cushion	Inclusive										
EC012		Seat cushion, black cover (seat cushion thickness 5 cm)	Inclusive										
EC013		Seat cushion with incontinence cover (seat cushion thickness 5 cm)	£86.00										
EC034		Seat cushion Livo (max. user weight up to 125 kg; seat cushion thickness 5 cm) exclusions see 3.3)	£71.00										
EC035		Seat cushion Evoflair (max. user weightup to 200 kg; seat cushion thickness 8 cm)	£256.00										
EC014		Terra seat cushion (max. user weight up to 150 kg; seat cushion thickness 6 cm) exclusions see 7.1) 7.2)	£174.00										
EC015		Terra Aquos seat cushion (max. user weight up to 150 kg; seat cushion thickness 6 cm) exclusions see 7.1) 7.2)	£197.00										
EC016		Terra Flair seat cushion (max. user weight up to 150 kg; seat cushion thickness 6 cm) exclusions see 7.1) 7.3)	£395.00										
EC018		Cloud seat cushion (max. user weight up to 250 kg; seat cushion thickness 11 cm) exclusions see 7.5)	£483.00										
EC020		Flat contoured (seat cushion thickness 7 cm) exclusions see 7.8) to 7.11)	£213.00										
EC021		Deep contoured (seat cushion thickness 9 cm) exclusions see 7.8) to 7.11)	£213.00										
Please refer to ou	r seat	cushion brochure for information on additional seat cushions.											
Cover for seat c	ushic	n, contoured (possible only in combination with EC020 or EC021) exclusions see 8.0) 8.1)											
EC022		Fabric cover	Inclusive										
EC023		Synthetic leather cover	£260.00										
EC024		Fabric cover with incontinence barrier	£283.00										
Back support ve	ersior	n/with contour 📵 (mandatory information)											
ED032		Without back support upholstery	Inclusive										
ED033		Back support upholstery, adaptable	Inclusive										
ED034		Back support pad flat contoured (SD reduction by approx. 2 cm) exclusions see 7.8) 7.9)	£304.00										
ED035		Back support pad deep contoured (SD reduction by approx. 2 cm) exclusions see 7.8) 7.9)	£304.00										
ED050		Back support pad flat contoured AD (adaptable, SD reduction by approx. 2 cm) exclusions see 9.1)	£355.00										
ED051		Back support pad deep contoured AD (adaptable, SD reduction by approx. 2 cm) exclusions see 9.1)	£355.00										
ED036		Baxx aluminium flat top back support (max. user weight 113 kg, installation accessories 4-point-PRO,	£588.00										
		SD reduction by approx. 2 cm, back support width = chosen seat width) exclusions see 9.1) to 9.3)											

Cover for back	supp	ort pad, contoured (possible only in combin. with ED034, ED035, ED050 or ED051) exclusions see 10.0) 10.1)	
ED037		Fabric cover	Inclusive
ED038		Synthetic leather cover	£42.00
Accessories for	back	support pad, contoured 📵	
ED039		Side panels (only for back support pad deep contoured) exclusions see 10.0)	£80.00
ED040		Lumbar support exclusions see 11.1)	£60.00
Accessories for	r arm s	support exclusions see 12)	
EL021		Padded covers for arm supports exclusions see 12.1)	£79.00
EL022		Padded covers for clothing guard	£60.00
Back support h	neight	for Baxx exclusions see 14.0) 14.1)	
ED041		High (18"/46 cm)	Inclusive
ED042		Extra high (20"/51 cm)	Inclusive
Back support o	onto	ur for Baxx exclusions see 14.0) 14.1)	
ED043		Active (2.5"/6 cm)	Inclusive
ED044		Deep (4"/10 cm)	£24.00
ED045		Extra deep (6"/15 cm)	£47.00
Colour for Bax	x exclı	usions see 14.0) 14.1)	
ED005		Sparkle light orange	Inclusive
ED009		Vinho sparkle	Inclusive
ED010		Marine blue	Inclusive
ED013		Jet black RAL 9005 matte	Inclusive
ED014		Candy red	Inclusive
ED021		Sparkle granny-smith	Inclusive
ED023		White RAL 9016	Inclusive
Seat inclination	າ (man	idatoy information)	
		-3° (not in combination with EC005) exclusions see 17.1)	
		0°	
		3° (not in combination with EC005) exclusions see 17.1)	
		6°	
		9°	
Power seat adj	ustme	ent (seat inclinition 0°)	
EC025		Power seat tilt 45° with centre of gravity shifting (max. load: VAS 100 kg, standard/contour seat 140 kg) exclusions see 3.1) 18.1) 24.9)	£729.00
EC026		Power seat height adjustment 35 cm (max. load: VAS 100 kg, standard/contour seat 140 kg; not possible in combination with EU060 or EU051) exclusions see 3.1) 18.1) 24.9)	£2012.00
EC027		Power seat height adjustment 35 cm with 45° seat tilt, with centre of gravity shifting (max. load: VAS 100 kg, standard/contour seat 140 kg; not possible in combination with EU060 or EU051) exclusions see 3.1) 18.1) 24.9)	£2351.00
EC028		Power seat tilt 20° (max. load: standard/contour seat 140 kg, not in combination with EC005) exclusions see 3.1) 18.1) 24.9)	£508.00
Caster fork (ma	ndato	y information)	
EF001		Suspension (not possible in combination with EF007 and EF008) exclusions see 2.1)	£217.00
EF002		No suspension	Inclusive

Castar																						
Proper	Caster wh	neels (ı	manda	toy info	rmatio	n)																
	EF005			9" pneumatic, puncture-resistant exclusions see 2.1) Inclusive																		
	EF006			9" pun	·																	
Part	EF007			10" pn	euma	tic, pu	ncture	-proo	f (not p	ossible	in con	nbinati	on with	EF001	) exclus	ions se	e 2.1) 2	20.1)			Ir	nclusive
Privacy   Priv	EF008			10" pu	ncture	e-proc	f (not p	oossible	e in cor	nbinati	on witl	n EF001	) exclu	sions s	ee 2.1) 2	20.1)					lr	nclusive
	EF004			8" pne	umati	c, pun	cture-	proof	(possib	le only	in com	binatio	n with	EG003	) exclu	sions se	ee 20.2)	)				£62.00
	Drive whe	eels (m	nandat	oy infor	mation	)																
Type	EG001			14" pn	pneumatic, puncture-resistant exclusions see 2.1) Inclusive																	
Type color   Typ	EG002			14" pu	ncture	e-proc	f exclu	sions s	ee 2.1)												Ir	nclusive
Tright   Sign	EG003			12" pu	puncture-proof (possible only in combination with EF004) exclusions see 20.2) £139.00																	
Tread	Tyre colou	ur (mai	ndatoy	informa	mation)																	
Tread   Cols	EG011			Grey																		
Figure   Second   Content   Conten	EG012			Black (																		
Figure   Second   Content   Conten	Tread 💽	l (man	datoy	informa	tion)																	
Casta	EG013			Ribs (o	nly wit	h tyre	colour	grey) ex	clusio	ns see 2	23.1)										Ir	nclusive
Second   S	EG014			Lugs (r	not pos	sible w	/ith 12"	and 8"	wheels	s) exclu	sions s	ee 22.1)									Ir	nclusive
Seat height   Seath reliable   Seath r	EG015			Caster	whee	ls with	ribs, o	drive v	vheels	with I	ugs (n	ot poss	ible wit	th EG00	)3 and	EF004)	exclusi	ons			lı	nclusive
Notices and set the light standard set the l				see 22.	1) 23.1)																	
Seath height hei	Seat heigl	hts (at	0° sea	t inclina	tion, w	ithout	seat cu	shion, s	suspen	sion fu	lly exte	nded/v	vithout	t user; r	mandat	oy info	rmatio	n)				
Regint   R	Wheels			1.1)	1.1)	1.2)	1.3)	1.1)	1.2)	1.3)	1.1)	1.1)	1.2)	1.3)	1.2)	1.1)	1.2)	1.3)	1.3)	1.1)	1.2)	1.3)
Notices and the contribution of the contributi	Seat	cm		38		40			42,5		44		45		46		47,5		48,5		50	
Notices and exclusions   24.37	height			11.5		45.6					47.0		47.6		40		40.6		40		40.5	
Standard seat   Standard sea	Notices and e		ns	+	24.37)		6) 24.0)	24.		4.29)		24.3		1.27)		24.3		1.24)		24.3		1.23)
Standard seat   Standard sea			-												24.0)							
Standard seat				24.4)		24.14) 24.3	)	24.11	) 24.9) 24.2)	24.3)	24.4)				24.7)	24	.18) 24.0) 24	1.2)	24.8)		24.3) 24.2)	
Standard seat   Standard sea															24.6)				24.6)			
Standard seat																						
Standard seat         24.21)         24.27																						
Standard seat   2,2   2   2   3   3   4   5   5   5   5   5   5   5   5   5				24.2)																		
VAS         Image: Control of the				24.21)					I					I								
VAS  VAS  VAS  VAS  VAS  VAS  Standard seat 20° seat tilt  VAS  VAS  45° seat tilt  VAS, Lift seat 35 cm + 45° seat tilt  VAS, Lift seat 35 cm +  VAS  VAS  VAS  VAS  VAS  VAS  VAS  VA	Standard seat	t		_																		
Standard seat 20° seat tilt         Company of the company of th	VAS																					
20° seat tilt  Standard seat 45° seat tilt  VAS, Lift seat 35 cm +  Output  Ou					2.2)			2.2)		2.2)												
Standard seat 45° seat tilt         Image: Control of the contro		t																				
45° seat tilt       2.2)       1		+																				
VAS 45° seat tilt         45° seat tilt         1		-																				
Standard seat, Lift seat 35 cm + 45° seat tilt  VAS, Lift seat 35 cm + 45° seat tilt											_											
35 cm + 45° seat tilt  VAS, Lift seat 35 cm +	45° seat tilt										2.1)				Ц				Ц			
VAS, Lift seat 35 cm +			at																			
											2.1)								_			

Wheels			Weight		
1.1)	8"	/12.5"	2.1)	Max. 75 kg	
1.2)	9"	caster wheels	2.2)	Max. 100 kg	
1.3)	9"	caster wheels with suspension or 10" wheels			
	Co	ombination possible			
	Co	ombination not possible			
Back support a	ngle a	adjustment (mandatory information)			
ED046		Mechanical in increments (Standard seat 0/10/20/30°, -9/1/1	1/21°; VAS -5/0/1	0/20/30°)	Inclusive
ED047		Power (30°)			£621.00
Mounting kit fo	or hea	nd/neck supports			
EL001		Gooseneck			£114.00
EL002		Multi-axis, offset			£106.00
EL003		Multi-axis, straight			£102.00
EL031		UniLink holder for head support, multi-axis (suitable only possible in combination with EL006, EL007, EL013 to EL016)	for Advanced lir	ne head supports; not	£212.00
EL039		UniLink extension element, 38 mm (may be used to compere)	nsate the distan	ce at kyphotic head postu-	£73.00
EL038		UniLink extension element, 50 mm (may be used to compere)	nsate the distan	ce at kyphotic head postu-	£54.00
EL040		UniLink extension element, 76 mm (may be used to compere)	nsate the distan	ce at kyphotic head postu-	£73.00
Head/neck sup	ports	exclusions see 27.0) 27.1)			
Care line					
EL006		Head support, small			£80.00
EL007		Head support, large			£80.00
EL013		Head/neck support, small			£76.00
EL014		Head/neck support, large			£76.00
EL015		Head/neck support, X-shape, small			£91.00
EL016		Head/neck support, X-shape, large			£91.00
Advanced line					
EL033		Comfort Plus, size S			£141.00
EL034		Comfort Plus, size M			£164.00
EL035		Comfort Plus, size L			£164.00
EL036		Combo (3-pad), size S			£282.00
EL037		Combo (3-pad), size L			£293.00
Belts • exclu	sions s	see 28.0)			
EP001		Standard lap belt			Inclusive
EP002		2-point lap belt, padded, metal buckle			£134.00
EP003		2-point lap belt, padded, plastic buckle			£121.00
EP020		2-point lep belt with metal buckle and rewinder (buckle of	on side of joystic	k) exclusions see 7.8)	£206.00

EPO05	Chest/shoulde	r harn	esses 📵	
EP007	EP005		Dynamic harness, size S	£174.00
EP008   Sternum harness, size S	EP006		Dynamic harness, size M	£174.00
EP009   Sternum harness, size M	EP007		Dynamic harness, size L	£174.00
EP010	EP008		Sternum harness, size S	£167.00
EP011	EP009		Sternum harness, size M	£167.00
EP012	EP010		Sternum harness, size L	£167.00
EP013	EP011		Slimline sternum harness, size S	£162.00
EP014   Shoulder harness, size S   £180.00   EP015   Shoulder harness, size M   £180.00   EP016   Shoulder harness, size M   £180.00   EP016   Shoulder harness, size L   £180.00   VAS accessories exclusions see 9.1) ED048   Push handles   £60.00   Arm supports	EP012		Slimline sternum harness, size M	£162.00
EPO15   Shoulder harness, size M	EP013		Slimline sternum harness, size L	£162.00
EP016	EP014		Shoulder harness, size S	£180.00
WS accessories exclusions see 9.1)  ED048	EP015		Shoulder harness, size M	£180.00
ED048	EP016		Shoulder harness, size L	£180.00
Arm supports (mandatory information)  EE023	VAS accessorie	s exclu	usions see 9.1)	
EE023	ED048		Push handles	£60.00
Flip-up (effective seat width +/- 1 cm per side; on Standard seat Small, SW 42 cm; Standard seat Large, SW 48 cm; Standard seat XL, SW 54 cm; VAS plus 2 cm; max. user weight: 125 kg; not in combination with EC001) exclusions see 3.3) 7.8) 67.1) to 67.4     Leg support version	Arm supports	(r	nandatoy information)	
SW 48 cm; Standard seat XL, SW 54 cm; VAS plus 2 cm; max. user weight: 125 kg; not in combination with EC001) exclusions see 3.3) 7.8) 67.1) to 67.4)    Leg support version	EE023		Standard	Inclusive
EB011	EE024		SW 48 cm; Standard seat XL, SW 54 cm; VAS plus 2 cm; max. user weight: 125 kg; not in combination with	£233.00
EB001	Leg support ve	ersion	(mandatoy information)	
EB002   Mechanically elevating with length adjustment (max. user weight 140 kg; not available for LL < 21   £540.00   cm; not in combination with EB012 / EB009) exclusions see 3.1) 31.1) to 31.4)   Power elevating with length adjustment and separate control (max. load 140 kg; not available   £1101.00   for LL < 21 cm; not in combination with EB012 / EB009) exclusions see 3.1) 31.1) to 31.4)   Power elevating with length adjustment and separate control (max. load 140 kg; not available   £1101.00   for LL < 21 cm; not in combination with EB012 / EB009) exclusions see 3.1) 31.1) to 31.4)   EB005   Without leg supports exclusions see 31.5)   Inclusive   Foot support version   exclusions see 32.0) 32.1)   EB012   Individual, plastic foot plates (max. load 140 kg; not available on elevating leg supports; not in combination with EB002 / EB003 / EB004) exclusions see 3.1) 31.6) 32.2) to 32.4)   Individual aluminium foot plates, new design (max. user weight 160 kg; not in combination with EB002 / EB001 / EB001) exclusions see 32.6) to 32.8) 32.10)   EB009   Single-panel, aluminium foot plate (max. user weight 160 kg; not in combination with EB002 / £136.00   EB003) exclusions see 32.6) to 32.9)   Leg support accessories   EB010   Lateral leg pad (for mechanical and power elevating leg support adjustment; pair) exclusions see 33.1)   £198.00   Batteries (mandatoy information)   EU036   AGM 56 Ah (C5) / 60,4 Ah (C20)   Inclusive   EU032   Gel 63 Ah (C5) / 74 Ah (C20)   £551.00	EB011		Standard (basic design) exclusions see 7.8) 31.1) 31.6)	Inclusive
EB003	EB001		Standard (new design) exclusions see 31.1)	£194.00
Foot support version    Continuous points   Foot support version   F	EB002			£540.00
EB005	EB003			£1101.00
Foot support version  exclusions see 32.0) 32.1)  EB012	EB004			£1101.00
Individual, plastic foot plates (max. load 140 kg; not available on elevating leg supports; not in combination with EB002 / EB003 / EB004) exclusions see 3.1) 31.6) 32.2) to 32.4)   EB007	EB005		Without leg supports exclusions see 31.5)	Inclusive
nation with EB002 / EB003 / EB003 / EB004) exclusions see 3.1) 31.6) 32.2) to 32.4)   EB007	Foot support v	ersior	exclusions see 32.0) 32.1)	
EB007	EB012			Inclusive
Leg support accessories  EB010	EB007			£233.00
EB010	EB009			£136.00
Batteries (mandatoy information)         EU036       ☐ AGM 56 Ah (C5) / 60,4 Ah (C20)       Inclusive         EU032       ☐ Gel 63 Ah (C5) / 74 Ah (C20)       £551.00	Leg support ac	cesso	ries	
EU036       AGM 56 Ah (C5) / 60,4 Ah (C20)       Inclusive         EU032       Gel 63 Ah (C5) / 74 Ah (C20)       £551.00				£198.00
EU032 Gel 63 Ah (C5) / 74 Ah (C20) £551.00	Batteries (mand	latoy ii	nformation)	
	EU036		AGM 56 Ah (C5) / 60,4 Ah (C20)	Inclusive
EU032 Gel 63 Ah (C5) / 74 Ah (C20) Basic £551.00	EU032		Gel 63 Ah (C5) / 74 Ah (C20)	£551.00
	EU032		Gel 63 Ah (C5) / 74 Ah (C20) Basic	£551.00

Battery charge	r 🗿	(mandatoy information)	
EU040		Without battery charger	Inclusive
EU041		Battery charger 8 A, fanless (suitable only for batteries with 56 Ah (C5)) exclusions see 35.3) 35.4)	Inclusive
EU044		Battery charger 10 A, fanless (IP54) (suitable for batteries 56 Ah (C5) and 63 Ah (C5)) exclusions see 35.2)	£242.00
EU046		Battery charger 10 A, with fan (suitable for batteries 56 Ah (C5) and 63 Ah (C5))	£252.00
EU042		Battery charger 12 A, incl. external charging receptacle (magnetic, installation side right) exclusions see 35.4)	£428.00
Additional cha	rging	receptacle 📵	
EU043		External magnetic charging receptacle (installation side right; possible only with EU042) exclusions see 36.1)	£173.00
Wheelchair col	our (c	olour for caster fork and metal strip; frame colour: KTL black; mandatoy information)	
005		Sparkle light orange	Inclusive
009		Vinho sparkle	Inclusive
010		Marine blue	Inclusive
013		Jet black RAL 9005 matte	Inclusive
014		Candy red	Inclusive
021		Sparkle granny-smith	Inclusive
023		White RAL 9016	Inclusive
Chassis accesso	ories	<b>⊙</b>	
EA010		Splash guard for drive wheels exclusions see 2.1) 39.5) 39.6)	£100.00
EA011		Splash guard for caster wheels exclusions see 2.1)	£58.00
Caster wheel sv	wivel	lock	
EF020		Mechanical	£399.00
Chassis safety a	access	sories 💽	
EF021		Curb climbing assist exclusions see 2.1)	£322.00
			2322:00
ES002		Kit for vehicle transport (not in combination with ED036; max. load 140 kg) exclusions see 3.1) 41.1) 41.3)	Inclusive
ES002 ES013		· · · · · · · · · · · · · · · · · · ·	
		41.3) Ottobock installation kit for Dahl docking system (not possible in combination with EA001 and EF021; max. load 140 kg) exclusions see 2.1) 3.1) 24.5) 41.1) 41.2)	Inclusive
ES013		41.3) Ottobock installation kit for Dahl docking system (not possible in combination with EA001 and EF021; max. load 140 kg) exclusions see 2.1) 3.1) 24.5) 41.1) 41.2)	Inclusive
ES013 General access	ories	41.3) Ottobock installation kit for Dahl docking system (not possible in combination with EA001 and EF021; max. load 140 kg) exclusions see 2.1) 3.1) 24.5) 41.1) 41.2)	Inclusive £100.00
ES013  General access ES003	ories	Ottobock installation kit for Dahl docking system (not possible in combination with EA001 and EF021; max. load 140 kg) exclusions see 2.1) 3.1) 24.5) 41.1) 41.2)  Pocket for mobile phone	£100.00 £49.00
ES013  General access ES003 ES004	ories	41.3) Ottobock installation kit for Dahl docking system (not possible in combination with EA001 and EF021; max. load 140 kg) exclusions see 2.1) 3.1) 24.5) 41.1) 41.2)  Pocket for mobile phone  Luggage carrier exclusions see 2.1) 3.1)	£100.00 £49.00 £223.00
ES013  General access ES003 ES004 ES005	ories	Ottobock installation kit for Dahl docking system (not possible in combination with EA001 and EF021; max. load 140 kg) exclusions see 2.1) 3.1) 24.5) 41.1) 41.2)  Pocket for mobile phone  Luggage carrier exclusions see 2.1) 3.1)  Crutch holder  Rearview mirror (installation side across from control panel side; not possible in combination with	£100.00 £49.00 £223.00 £114.00
ES013  General access ES003 ES004 ES005 ES006	ories	41.3) Ottobock installation kit for Dahl docking system (not possible in combination with EA001 and EF021; max. load 140 kg) exclusions see 2.1) 3.1) 24.5) 41.1) 41.2)  Pocket for mobile phone  Luggage carrier exclusions see 2.1) 3.1) Crutch holder  Rearview mirror (installation side across from control panel side; not possible in combination with ES012) exclusions see 42.1)	£100.00 £49.00 £223.00 £114.00
ES013  General access ES003 ES004 ES005 ES006 ES007	ories	Ottobock installation kit for Dahl docking system (not possible in combination with EA001 and EF021; max. load 140 kg) exclusions see 2.1) 3.1) 24.5) 41.1) 41.2)  Pocket for mobile phone Luggage carrier exclusions see 2.1) 3.1) Crutch holder  Rearview mirror (installation side across from control panel side; not possible in combination with ES012) exclusions see 42.1) Tool kit	£100.00 £49.00 £223.00 £114.00 £58.00
ES013  General access ES003 ES004 ES005 ES006 ES007 ES010	ories	41.3) Ottobock installation kit for Dahl docking system (not possible in combination with EA001 and EF021; max. load 140 kg) exclusions see 2.1) 3.1) 24.5) 41.1) 41.2)  Pocket for mobile phone  Luggage carrier exclusions see 2.1) 3.1) Crutch holder  Rearview mirror (installation side across from control panel side; not possible in combination with E5012) exclusions see 42.1)  Tool kit  Beverage holder	£100.00 £49.00 £223.00 £114.00 £170.00 £58.00
ES013  General access ES003 ES004 ES005 ES006 ES007 ES010 ES011	ories	41.3) Ottobock installation kit for Dahl docking system (not possible in combination with EA001 and EF021; max. load 140 kg) exclusions see 2.1) 3.1) 24.5) 41.1) 41.2)  Pocket for mobile phone Luggage carrier exclusions see 2.1) 3.1) Crutch holder Rearview mirror (installation side across from control panel side; not possible in combination with ES012) exclusions see 42.1) Tool kit Beverage holder Airman pump	£100.00 £49.00 £223.00 £114.00 £58.00 £25.00 £241.00
ES013  General access ES003 ES004 ES005 ES006 ES007 ES010 ES011 ET001	ories	Ottobock installation kit for Dahl docking system (not possible in combination with EA001 and EF021; max. load 140 kg) exclusions see 2.1) 3.1) 24.5) 41.1) 41.2)  Pocket for mobile phone  Luggage carrier exclusions see 2.1) 3.1)  Crutch holder  Rearview mirror (installation side across from control panel side; not possible in combination with ES012) exclusions see 42.1)  Tool kit  Beverage holder  Airman pump  External horn  Tray (swivel mechanism across from control panel side; not possible in comb. with ES006) exclusions see	£100.00 £49.00 £223.00 £114.00 £170.00 £58.00 £25.00 £241.00

Receptacle for	exter	nal power supply, right 📵	
ET002		12 V exclusions see 43.1) 43.4)	£300.00
ET003		24 V exclusions see 43.2)	£120.00
ET004		USB exclusions see 43.3)	£135.00
Lighting  e	kclusio	ins see 44.0)	
ES021		Electric lighting	£459.00
Control device	•	(mandatoy information)	
EU059		VR2 (70 A), without power functions and lighting (recommended for load up to 140 kg) exclusions see 46.1) 46.2) 46.4) 46.7)	Inclusive
EU060		VR2 (70 A), max. 2 power functions and lighting possible (recommended for load up to 140 kg; not possible in combination with EC026 or EC027) exclusions see 46.2) 46.3) 46.5)	Inclusive
EU050		VR2 (90 A), without power functions and lighting (recommended for load up to 160 kg) exclusions see 2.1) 46.1) 46.4) 46.5) 46.7)	Inclusive
EU051		VR2 (90 A), max. 2 power functions and lighting possible (recommended for load up to 160 kg; not possible in combination with EC026 or EC027) exclusions see 2.1) 46.3) 46.7)	£170.00
EU052		TEN° (90 A), 1 power function, choose option EU056 for more power functions, lighting possible (recommended for load up to 160 kg)	£793.00
EU058		R-net LED-L (90 A), 1 power function, choose option EU056 for more power functions, lighting possible (recommended for load up to 160 kg)	£549.00
EU064		R-net Colour (90 A), 1 power function, choose option EU056 for more power functions, lighting not possible (recommended for load up to 160 kg) exclusions see 46.4)	£823.00
TEN° control p	anel (	possible only in combination with EU052) exclusions see 47.0) 47.1)	
EU065		Standard colour shimano matte	Inclusive
EU066		Colour selectable	Inclusive
Colour for TEN	° cont	trol panel (possible only in combination with EU066) exclusions see 48.0) 48.1)	
EU005		Sparkle light orange	Inclusive
EU009		Vinho sparkle	Inclusive
EU010		Marine blue	Inclusive
EU013		Jet black RAL 9005 matte	Inclusive
EU014		Candy red	Inclusive
EU021		Sparkle granny-smith	Inclusive
EU023		White RAL 9016	Inclusive
Control panel	butto	n (possible only in combination with EU052) exclusions see 47.0)	
EU067		Piko button, diameter 50 mm (red/black)	£356.00
EU068		Piko button, diameter 30 mm (red/black)	£356.00
Control panel	holde	r (mandatoy information)	
EU070		Control panel holder, standard (not possible in combination with TEN° (EU052)) exclusions see 50.1)	Inclusive
EU071		Control panel holder, swing-away in parallel, 2 cm height-adjustable	£113.00
EU072		Control panel holder, swing-away in parallel, 9.5 cm height-adjustable	£136.00
Control panel	holde	r mounting side (mandatoy information)	
		Right	
		Left	,

Joystick top	)		
EU080		Fork for tetraplegics, horizontal (adjustable for hand width 70 - 90 mm)	£113.00
EU081		Stick S80 (special joystick handle 80 mm)	£100.00
EU082		Soft ball, black cover	£170.00
EU083		Ball top, small (diameter 43 mm)	£114.00
EU084		Ball top, large (diameter 48 mm)	£105.00
EU085		Soft ball, coloured (diameter 30 mm, 4 pieces)	£206.00
Modules for TE	N°/R-	net LED-L/R-net colour control device 📵 exclusions see 53.0)	
		TEN° LCD module (when using 1-switch scan drive, programming by service required; please order but-	
EU054		ton separately; adapter switch box not required)	£1622.00
Holder for TEN	° LCD	module exclusions see 54.0) 54.1)	
EU073		Standard (installation across from control panel side; not possible in combination with ES012) exclusions see 42.1)	£83.00
EU074		With gooseneck (installation across from control panel side; not possible in combination with ES012) exclusions see 42.1)	£256.00
EU075		With gooseneck for tray (included)	£192.00
Actuator modu	ıle for	TEN°/R-net LED-L/R-net colour control device exclusions see 53.0) 55.0) 55.1)	
EU056		Actuator module (required for power seat functions on R-net LED-L, R-net colour and TEN°, when more	£420.00
L0030	Ш	than one power seat function is needed)	1420.00
Control device	acces	ssories 📵	
EU057		Attendant control	£706.00
Attendant con	trol in	stallation side exclusions see 57.0) 57.1)	
		Right exclusions see 57.2)	
		Left	
Track stabiliser	exclu	sions see 53.0) 59.0)	
EF023		Power (gyro-module) (not in combination with EU059, EU060, EU050 or EU051; depending on terrain, the specified driving distance range can be reduced when using the gyro)	£432.00
Drive unit (man	datoy	information)	
ET016		Performance (usable up to a gross weight of 320 kg)	Inclusive
ET011		High performance (only in combination with EU052, EU058 or EU064) exclusions see 2.1) 53.0)	£559.00
Speed (mandate	ov info	prmation)	
ET012		6 km/h	Inclusive
ET013		7.2 km/h exclusions see 2.1) 61.1)	Inclusive
ET014		10 km/h exclusions see 2.1) 61.1) 61.2)	Inclusive
Drive unit onti	on ev	:lusions see 62.0) 62.1)	
ET008		Individual type approval for speed > 6 km/h	£505.00
Drive unit war	antv		
EW002		From 24 to 48 months	£80.00
			253.00
		ons for use (additional IFU for non-national languages)	
ES028		Qualified personnel	Inclusive
ES029	Ш	User	Inclusive

#### Seating systems 🝮



EC001 - EC004 Standard seat



EC005 VAS

#### Seat cushion 5



EC034 Livo



EC035 Evoflair

#### Seat cushions/back support version with contour 5



EC020/ED034 Flat contoured



EC021/ED035 Deep contoured

#### Cover for seat cushion/back support pad contoured 5



EC022 Fabric cover



EC023 Synthetic leather cover

#### Back support version 5



ED050/ED051 Back support pad flat/deep contoured



ED036 Baxx aluminium flat top back support

#### Accessories for back support pad, contoured 5



ED039 Side panels



ED040 Lumbar support

#### Power seat adjustment 🝮



EC025 45° seat tilt



EC026 Seat height adjustment 35 cm



ED047 Power back support angle adjustment

#### Tread 5



EG013 Ribs



EG014 Lugs



EG015 Caster wheel with ribs

#### Mounting Kit for head/neck supports 5



EL001 Gooseneck



EL002 Multi-axis, offset



EL003 Multi-axis, straight



EL031 UniLink multi axis, including seat back mount

#### Head/neck supports Care line 5



EL006/EL007 Head support small/large



EL013/EL014 Combination head/neck support small/large



EL015/EL016 Head/neck support, X-shape, small/ large

### Head/neck supports Advanced line ◆



EL033/EL034/EL035 Comfort Plus



EL036/EL037 Combo

#### Arm supports 👈



EE024 Flip-up arm supports

#### Belts ◆



EP001 Standard lap belt



EP002 2-point lap belt, padded, metal buckle



EP003 Lap belt, padded, plastic buckle



EP020 Lap belt with rewinder

#### Chest/shoulder harnesses 5



EP005 - EP007 Dynamic harness S/M/L



EP008 - EP010 Sternum harness S/M/L



EP011 - EP013 Slimline sternum harness S/M/L



EP014 - EP016 Shoulder harness S/M/L

#### Leg support version 5



EB011 Standard (basic design)



EB001 Standard (new design)



EB002 Mechanically elevating with length adjustment



EB003 Power elevating, separate control

#### Foot plate version 🝮



EB004 Power elevating, simultaneous control



EB012 Individual, plastic foot plates



EB007 Individual aluminium foot plates

#### Add.charging receptacle 5



EU043 External magnetic charging receptacle

#### Battery charger 👈



EU041 Battery charger 8 A, fanless



EU044 Battery charger 10 A, fanless



EU046 Battery charger 10 A, with fan



EU042 Battery charger 12 A, incl. external charging receptacle

#### Chassis accessories 5



EA010 Splash guard for drive wheels

#### Chassis safety accessories 👈



EF021 Curb climbing assist



ES013 Dahl mounting kit

#### General accessories 5



ES003 Pocket for mobile phone



ES004 Luggage carrier



ES005 Crutch holder



ES006 Rearview mirror

ES010 Beverage holder



ES023 Backpack hooks

#### Receptacle f. ext. power 🔿



ET002 Receptacle f. 12 V ET003 Receptacle f. 24 V ET004 USB

#### Lighting 5



ES021 Electric lighting

#### Control device 👈



EU050/EU059 VR2



EU051/EU060 VR2



EU052 TEN°



EU058 R-net LED-L

### Joystick top 👈



EU064 R-net Colour



EU080 Fork for tetraplegics, horizontal



EU081 Stick S80



EU082 Soft ball, black cover

EU083 Ball top, small



EU084 Ball top, large



EU085 Soft ball, coloured

### Modul for TEN° ◆



EU054 TEN° LCD module

### Control device accessories 5



EU057 Attendant control



EU095 Hand warmer (single hood for handmodule)

#### Tables for notes and exclusions

#### Seat depth

Cook						Seat de	pth SD						
Seat	34 cm	36 cm	38 cm	40 cm	42 cm	44 cm	46 cm	48 cm	50 cm	52 cm	54 cm	56 cm	58 cm
Standard seat Junior	х	Х	Х	Х									
Standard seat, small			Х	Х	Х	Х	Х						
Standard seat, large					Х	Х	Х	Х	Х				
Standard seat, XL					Х	Х	Х	Х	Х				
VAS			Х	Х	Х	Х	Х	Х	Х	Х	Х	Х	Х

X = combination possible

### Seat width

6 .						Seat wid	dth (SW)					
Seat	34 cm	36 cm	38 cm	40 cm	42 cm	44 cm	46 cm	48 cm	50 cm	52 cm	54 cm	56 cm
Standard seat Junior	Х	Х	Х	Х								
Standard seat, small			Х	Х	Х							
Standard seat, large						Х	Х	Х				
Standard seat, XL									Х	Х	Х	Х
VAS			Х	Х	Х	Х	Х	Х	Х	Х	Х	

X = combination possible

					Liv	o seat cush	nion					
SD						Seat wi	dth [cm]					
[cm]	34	36	38	40	42	44	46	48	50	52	54	56
34												
36												
38				Х								
40			Х		Х							
42				Х	Х	Х	Х	Х		Х		
44					Х	Х	Х	Х	Х	Х		
46					Х	Х	Х	Х	Х	Х		
48					Х	Х	Х				•	
50						Х	Х					
52					Х	Х	Х			Х		
54												
56												
58												



EC034 Livo

	Evoflair seat cushion											
SD	Seat width [cm]											
[cm]	34	36	38	40	42	44	46	48	50	52	54	56
34												
36												
38												
40			Х	Х	Х							
42				х	Х		Х	Х	Х			
44					Х		Х	Х		Х		
46					Х		Х	Х		Х		
48									,			
50												
52					Х		Х			Х		
54												
56												
58												



EC035 Evoflair

	Terra, Terra Aquos seat cushion											
SD		Seat width [cm]										
[cm]	34	36	38	40	42	44	46	48	50	52	54	56
34		X										
36		X									_	
38		Х	Х	X	Х	Х	Χ	Х	Χ	Х		
40				Χ	Х	Х	Χ					
42				Х	X	Х	Χ					
44				Х	Х	Х	Χ	Х	Χ	Х		
46				Χ	Х	Х	Χ	Х	Χ	Х		
48				Х	Х	Х	Χ	Х	Х	Х		
50												
52												
54												
56												
58												

Terra Flair seat cushion												
SD	Seat width [cm]											
[cm]	34	36	38	40	42	44	46	48	50	52	54	56
34												
36												
38			Х	Х	Х	Х	Χ	Х	Х	Х		
40				Х	Х	Х	Х				-	
42				Х	Х	Х	X				_	
44				Х	X	Х	Χ	Х	X	Х		
46				X	Х	Х	Х	Х	Х	Х		
48				X	Х	Х	Χ	Х	X	Х		
50												
52												
54												
56												
58												







EC014 Terra EC015 Terra Aquos

EC016 Terra Flair

	Cloud seat cushion												
SD	Seat width [cm]												
[cm]	34	36	38	40	42	44	46	48	50	52	54	56	
34	Χ	Х											
36	Х	Х	Х										
38	Χ	Х	Х	Х	Х		Х	Х	Х	Х	Х	Х	
40	Х	Х	Х	Х	Х	Х	Х	Х	Х	Х	Х	Х	
42			Х	Х	Х	Х	Х	Х	Х	Х	Х	Х	
44				Х	Х	Х	Х	Х	Х	Х	Х	Х	
46				Х	Х		Х	Х	Х	Х	Х	Х	
48				Х	Х		Х	Х	Х	Х	Х	Х	
50						,	Х			Х	Х	Х	
52							Х			Х	Х	Х	
54								1					
56													
58													



EC018 Cloud

	Seat cushion flat contoured, deep contoured								
CD.[]	SW [cm]								
SD [cm]	38	40	42	44	46	48	50	52	54
38	Х	Х					•		
40	Х	Х							
42	Х	Х	Х	Х	Х	Х	Х	Х	Х
44	Х	Х	Х	Х	Х	Х	Х	Х	Х
46	Х	Х	Х	Х	Х	Х	Х	Х	Х
48			Х	Х	Х	Х	Х	Х	Х
50			Х	Х	Х	Х	Х	Х	Х
52			Х	Х	Х	Х	Х	Х	Х
54			Х	Х	Х	Х	Х	Х	Х
56			Х	Х	Х	Х	Х	Х	Х
58				'	1		'	1	1

 $\ensuremath{\mathsf{X}}\xspace$  possible dimension in combination with VAS

#### Notes and exclusions

- 2 1) not possible in conjunction with product classification category A (indoor)
- 3.1) possible only in combination with maximum user weight <= 140 kg
- 3.2) possible only in combination with maximum user weight <= 75 kg
- 3.3) possible only in combination with maximum user weight <= 125 kg
- see table "seat depth" 4.0)
- 4.1) not possible in conjunction with the combination of maximum user weight > 75 kg and Standard seat Junior or VAS
- 4.3)  $not possible in conjunction with the combination of maximum user weight > 100 \, kg and Standard seat Small or Standard seat Large or Standard$ seat XL or VAS
- 4.4) possible only in combination with maximum user weight <= 113 kg möglich
- 5.0) see table "seat width"
- possible only in combination with maximum load capacity <= 150 kg 7.1)
- see table "seat cushion Terra, Terra Aquos" 7.2)
- see table "seat cushion Terra-Flair" 7.3)
- 7.5) see table "seat cushion Cloud"
- possible only in combination with maximum user weight <= 125 kg 7.6)
- 7.8) cannot be combined with Standard seat, Junior
- Not possible in conjunction with the combination of VAS, seat width 38 cm or 40 cm and seat depth > 46 cm, see table "seat cushion with flat con-7.9) tours, with deep contours" (490E75=1-konturiert)
- Not possible in conjunction with the combination of VAS, seat width > 40 cm and seat depth 38 cm or 40 cm, see table "seat cushion with flat con-7.10) tours, with deep contours" (490E75=1-konturiert)
- Not possible in conjunction with the combination of VAS and seat depth 58 cm, see table "seat cushion with flat contours, with deep contours" 7.11) (490E75=1-konturiert)
- 7.12) cannot be combined with seat depth 58 cm
- 8.0) possible only in combination with seat cushion with flat or deep contours
- 8.1) mandatory information in combination with seat cushion with flat or deep contours
- possible only in combination with VAS 9.1)
- possible only in combination with maximum user weight <= 113 kg 9.2)
- 9.3) cannot be combined with seath width 54 cm
- 10.0) possible only in combination with back cushion with flat or deep contours or back cushion with flat or deep contours, for adaptable back upholstery
- 10.1) mandatory information in combination with back cushion with flat or deep contours or back cushion with flat or deep contours, for adaptable back upholstery
- possible only in combination with back cushion with deep contours 11.1)
- pandatory information in combination with Standard seat Small and SW 38 cm 12)
- 12.1) . not possible in conjunction with Standard seat Small and SW 38 cm and unselected upholstery covers for arm supports
- 14.0) possible only in combination with Baxx aluminium flat top back support
- 14.1) mandatory information in combination with Baxx aluminium flat top back support
- cannot be combined with VAS 17.1)
- 18.1) possible only with seat inclination 0°
- 18.2) cannot be combined with seat inclination -3° / 6° / 9°
- 20.1) not possible in conjunction with caster wheel swing arm spring mounted
- 20.2) not possible in conjunction with product classification category B (indoor/outdoor)
- 22.1) not possible in conjunction with caster wheels 8" or drive wheels 12'
- possible only in conjunction with tyre colour grey 23.1)
- seat height >= lower leg length +3 cm 24.0)
- 24.1) possible only in combination with product classification category A (indoor)
- 24.2) cannot be combined with seat tilt 45° or seat height adjustment 35cm or seat height adjustment 35 cm with seat tilt 45°
- 24.3) not possible in conjunction with VAS and caster wheel swing arm unsprung and caster wheels 9'
- 24.4) not possible in combination with VAS and maximum user weight > 75 kg
- 24.5) cannot be combined with seat tilt 20°
- possible only in combination with seat tilt 45° or seat height adjustment 35 cm or seat height adjustment 35 cm with seat tilt 45° 24.6)
- 24 7) possible only in combination with caster wheel swing arm unsprung and caster wheels 9'
- possible only in combination with caster wheel swing arm spring mounted or caster wheels 10' 24.8)
- not possible in conjunction with the combination of maximum user weight > 100 kg and VAS 24.9)
- 24.10) not possible in conjunction with the combination of seat height adjustment up to 35 cm with seat tilt 45° and VAS
- not possible in conjunction with the combination of VAS and SW > 42 cm and seat inclination 6°, 9° and caster wheel swing arm spring mounted or 24.11)
- 24.12) not possible in conjunction with the combination of VAS and caster wheel swing arm spring mounted or caster wheels 10"
- 24.13) cannot be combined with caster wheels 10" or caster wheel swing arm spring mounted
- 24.14) cannot be combined seat inclination -3°
- 24.15) not possible in conjunction with the combination of seat inclination -3° and caster wheels 10" or caster wheel swing arm spring mounted
- 24.16) not possible in conjunction with Standard seat Large/ XL and seat tilt 20° and caster wheels 8"
- 24.17)  $not possible in conjunction with Standard seat Large/\ XL\ and\ seat\ tilt\ 20^{\circ}\ and\ caster\ wheel\ swing\ arm\ unsprung\ and\ caster\ wheels\ 9^{\circ}\ and\ caster\ wheel\ swing\ arm\ unsprung\ and\ caster\ wheels\ 9^{\circ}\ and\ caster\ wheel\ swing\ arm\ unsprung\ and\ caster\ wheels\ 9^{\circ}\ and\ caster\ wheel\ swing\ arm\ unsprung\ and\ caster\ wheels\ 9^{\circ}\ and\ caster\ wheel\ swing\ arm\ unsprung\ arm\ unsprun$
- not possible in conjunction with the combination of VAS and seat inclination 9° and caster wheel swing arm unsprung and caster wheels 9" 24.18) not possible in conjunction with the combination of VAS and seat inclination 9° and caster wheel swing arm spring mounted or caster wheels 10" 24.19)
- not possible in conjunction with VAS and seat inclination 9° and caster wheel swing arm unsprung and caster wheels 10" or caster wheel swing 24.20) arm sprung and caster wheels 9"
- not possible in conjunction with VAS and caster wheels 8" 24.21)
- not possible in conjunction with Standard seat Junior/Small/Large/XL and seat inclination 6°, 9° and caster wheel swing arm unsprung and caster 24.22) wheels 10" or caster wheel swing arm sprung and caster wheels 9"
- 24.23) not possible in conjunction with Standard seat Junior/Small/Large/XL and seat inclination 3°, 6°, 9° and caster wheel swing arm unsprung and caster wheels 9"
- 24.24) not possible in conjunction with Standard seat Junior/Small/Large/XL and seat inclination 9° and caster wheel swing arm unsprung and caster wheels 9"
- not possible in conjunction with Standard seat Junior/Small/Large/XL and seat inclination 9° and caster wheels 8" 24.26)
- not possible in conjunction with Standard seat Large/XL and seat tilt 20° and caster wheel swing arm unsprung and caster wheels 10" or caster 24.27) wheel swing arm sprung and caster wheels 9"
- 24.28) not possible in conjunction with caster wheel swing arm spring mounted or caster wheels 10"
- 24.29) not possible in conjunction with the combination of VAS and SW > 42 cm and caster wheels 8"
- 24.31) not possible in conjunction with the combination of VAS and seat inclination 9° and caster wheels 8"
- not possible in conjunction with Standard seat Large/XL 24.32)

- Notes and exclusions 24 33) not possible in conjunction with Standard seat Large/XL and caster wheel swing arm unsprung and caster wheels 9" 24.34) not possible in combination with seat height adjustment up to 35 cm with seat tilt 45° and maximum user weight > 75 kg 24.35) not possible in combination with seat tilt 20° and seat inclination 3° and caster wheel swing arm unsprung and caster wheels 9" not possible in combination with seat tilt 20° and seat inclination 3° and caster wheel swing arm unsprung and caster wheels 8" 24.36) 24.37) see table "seat hight" 27.0) possible only in combination with head/neck support mounting kit 27.1) mandatory information in combination with head/neck support mounting kit 28.0) mandatory information in combination with category B not possible in conjunction with the combination of standard seat, small, seat depth < 42 cm and seat width > 38 cm 28.1) possible only in combination with selected lower lea length 31.1) not possible in conjunction with the combination of Standard seat or VAS and lower leg length < 21 cm 31.2) not possible in conjunction with the combination of maximum user weight > 75 kg and Standard seat and lower leg length < 28 cm 31.3) not possible in conjunction with the combination of maximum user weight > 75 kg and VAS and lower leg length < 27 cm 31.4) 31.5) not possible in conjunction with selected lower leg length not possible in conjunction with the combination of seat height SH < 45 cm and lower leg length LLL <= 30 cm 31.6) 32.0) cannot be combined with leg support adjustment, without 32.1) mandatory information in combination with leg support adjustment, standard, mechanically elevating or electrically elevating 32.2) not possible in conjunction with the combination of Standard seat and lower leg length  $\,<\!25\,cm$ not possible in conjunction with the combination of VAS and lower leg length < 24 cm 32.3) cannot be combined with leg support adjustment, mechanically elevating or electrically elevating 32 4) cannot be combined with leg support adjustment standard (basic design) 32.6) not possible in conjunction with the combination of Standard seat or VAS and lower leg length < 17 cm and leg support adjustment, standard 32.7) not possible in conjunction with the combination of maximum user weight > 75 kg and Standard seat or VAS and lower leg length < 24 cm and leg 32.8) support adjustment, standard 32.9) cannot be combined with leg support adjustment, mechanically elevating or electrically elevating, separate controllable 32.10) not possible in conjunction with the combination of maximum user weight > 140 kg and leg support adjustment, electrically elevating, simultaneously controllable 331) possible only in combination with leg support adjustment, mechanically or electrically elevating 35.2) not possible in conjunction with the combination of Argentina, Australia, Brazil, Hong Kong, Japan, Canada, Columbia, Mexico, Thailand, United Kingdom or United States of America as the country of destination cannot be combined with AGM batteries 63 Ah (C5) or Gel batteries 63 Ah (C5) 35.3) not possible in conjunction with the combination of Argentina, Brazil, Canada, Columbia, Mexico or USA as the country of destination 35.4)possible only in conjunction with battery charger, without 36.1) 39.5) not possible in conjunction with the combination of SH 40 cm and caster wheel swing arm unsprung and caster wheeles 9" and VAS or Standard 39.6) not possible in conjunction with the combination of SH 42,5 cm and caster wheel swing arm spring mounted or caster wheels 10" and VAS or eat systems Standard seat Large/XL 41.1) cannot be combined with Baxx aluminium flat top back support 41.2) cannot be combined with chassis safety accessories curb climbing assist 41.3) not possible in conjunction with the combination of Electic seat adjustment Seat tilt 20° Electrical, Belts 2-point Lap belt with rewinder, metal buckle and SD <44 cm cannot be combined with therapy tray 42.1) cannot be combined with rearview mirror 42.2) 42.3) cannot be combined with Germany as the country of destination 43.1) cannot be combined with socket for external power supply, right 24 V and socket for external power supply, right USB 43.2) cannot be combined with socket for external power supply, right 12 V and socket for external power supply, right USB 43.3) cannot be combined with socket for external power supply, right 12 V and socket for external power supply, right 24 V 43.4) cannot be combined with batteries/battery capacity AGM 63 Ah (C5) / 75 Ah (C20) mandatory information in conjunction with the combination of product classification category B and all countries of destination except Other coun-44.0) tries, Argentina, Australia, Brazil, France, Hong Kong, Japan, Canada, Columbia, Korea, Mexico, Russia, Switzerland, Thailand, United States of America or United Kingdom Not possible in conjunction with the combination of Germany as the country of destination and product classification category B 441) cannot be combined with electric functions 46.1) 46.2) possible only in combination with total weight <= 260 kg 46.3) possible only in combination with max. 3 electric functions cannot be combined with lighting acc. to German StVZO 46.4) not possible in conjunction with Germany as the country of destination and more than 2 electric functions 46.5) 46.7) cannot be combined with seat height adjustment 35 cm or seat height adjustment 35 cm with seat tilt 45° 47.0) possible only in combination with control unit, TEN° 47.1) mandatory information in combination with control unit, TEN° possible only in combination with control panel, TEN°, colour selectable 48.0) 48 1) mandatory information in combination with control panel, TEN°, colour selectable 50.1) cannot be combined with control unit TEN° possible only in combination with control unit, TEN° or R-net LED-L (90 A) or R-net Colour (90 A) 53.0) possible only in combination with TEN° LCD module 54.0) 54.1) mandatory information in combination with TEN° LCD module  $mandatory information in conjunction with the combination of control unit, TEN^{\circ} (90 \, A) or \, R-net \, LED-L (90 \, A) or \, R-net \, Colour (90 \, A) and > 1 \, electronic conjunction (90 \, A) and (90$ 55.0) tric function 55.1)  $mandatory information in conjunction with the combination of lighting acc. to German StVZO and control unit, TEN^{\circ} or R-net LED-L (90 A)$ 57.0) possible only in combination with attendant control 57.1) mandatory information in combination with attendant control 57.2) not possible in conjunction with the combination of Standard seat Junior and actuator module for control unit, TEN°
- possible only in combination with control unit, VR2 58.1)
- not possible in conjunction with Australia. Canada or the United States of America as the country of destination 59.0) 59.1) not possible in conjunction with the combination of batteries Gel 63 Ah (C5) / 74 Ah (C20) and general accessories external horn
- 59.2) not possible in conjunction with the combination of batteries AGM 63 Ah (C5) / 75 Ah (C20) and general accessories external horn
- 61.1)
- for Germany as the country of destination, this is possible only in conjunction with the combination of rear warning sign and lighting acc. to Ger-
- 61.2) not possible in conjunction with drive unit high performance
- possible only in conjunction with the combination of Germany as the country of destination and speed > 6 km/h 62.0)

#### Notes and exclusions

62.1)	mandatory information in conjunction with the combination	n of Germany as the country of destination and speed > 6 km/h

- 63.0) not possible in conjunction with Australia as the country of destination
- 64.0) possible only in combination with additional instructions for use, qualified personnel
- 64.1) mandatory information in combination with additional instructions for use, qualified personnel
- 65.0) possible only in combination with additional instructions for use, user
- mandatory information in combination with additional instructions for use, user 65.1)
- 67.1) not possible in conjunction with the combination of Standard seat Small and seat width 38 cm or 40 cm  $\,$
- 67.2) 67.3) not possible in conjunction with the combination of Standard seat Large and seat width 44 cm or 46 cm not possible in conjunction with the combination of Standard seat XL and seat width 50 cm or 52 cm not possible in conjunction with 2-point Lap belt with rewinder, metal buckle
- 67.4)