Base Price -EA004 - FWD - £5462.00 EA005 - RWD - £4984.00 EA006 - MWD - £7215.00 Valid from 01.01.24 to 31.12.24 Please ensure you are using the latest order form ☐ Provide estimate only □ Order ☐ Attachment for custom fabrication enclosed Billing address/customer number **Delivery address** Sales person Name Date Signature By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of

the respective Ottobock products, and you confirm to have read, understood and accepted our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

Standard equipment

Juvo B5 Power Wheelchair

EA004 - 490E75=0_AS01_C

EA005 - 490E75=0_AS01_C

EA006 - 490E75=0_AS05_C

Max. load: 160 kg (standard equipment)







Additional options installed on the power wheelchair can reduce the maximum load. Please observe the maximum gross weight on the nameplate.

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders. Please use the spare parts catalogue for spare parts orders.

Basic data (volur	ntary	nformation)	Mandatory information*
Gender		Indication	User height in cm
Age		Paying party	User weight in kg
			1 Seat depth
			2 Trunk depth
	1		3 Arm support height
			4 Armpit height
		7	5 Shoulder height
6	5	3 2	6 Head height
			7 Forearm length to the middle of the palm
		1 10	8 Trunk width
		9	9 Lower leg length
		9	10 Seat width
* This information	will h	elp to support user-oriented settings. Ottobock guarantees the confidential a	nd purely internal use of these data.
Product classifi	catio	1	
EA002		Category B (in/outdoor) (electric lighting mandatory in Europe)	Inclusive
Drive type			
EA004		Front-wheel drive (FWD)	£5462.00
EA005		Rear-wheel drive (RWD)	£4984.00
EA006		Mid-wheel drive (MWD) exclusions see 2.1)	£7215.00
Max. user weig	ht (us	er + personal luggage, information required to adapt chassis suspension)	
		70 kg	Inclusive
		75 - 100 kg	Inclusive
		100 - 140 kg	Inclusive
		140 - 160 kg	Inclusive
Chassis size			
EA007		Size 1 (overall width approx. 60 cm; max. user weight 140 kg) exclusions see	4.1) Inclusive
Seating system	s 💽		
EC001		Standard seat Junior (max. user weight 75 kg)	£147.00
LCOOT	Ш	(SD 3440 cm, SW 3440 cm, BH 47 cm to 37 cm, adjustable in 5 cm increme	ents) exclusions see 5.2)
EC002		Standard seat Small (SD 3846 cm, SW 3842 cm, BH 55 cm to 45 cm, adjuments)	stable in 5 cm incre- Inclusive
LCOOZ	Ш	exclusions see 4.1)	metasive
	_	Standard seat Large (SD 4250 cm, SW 4448 cm, BH 55 cm to 45 cm, adju	stable in 5 cm incre-
EC003		ments) exclusions see 5.1)	Inclusive
	_	Standard seat XL (SD 4250 cm, SW 5056 cm, BH 55 cm to 45 cm, adjustal	ole in 5 cm increments)
EC004		exclusions see 5.1)	£189.00
EC005		VAS (Vario Adjust Seat) (SD 3858 cm, SW 3854 cm, BH 55 cm, 53 cm, 5	(cm, 45 cm) £435.00
EC006		Recaro n-Joy* (manually adjustable back support, cushion thickness appro	x. 16.5 cm) exclusions see 5.1) £1759.00
EC007		Recaro n-Joy-Plus* (adaptable on site, manually adjustable back support, scushion thickness approx. 16.5 cm) exclusions see 5.1)	5D 38 56 cm, £4186.00

Back support for i	mod	lular Recaro	seating s	system (exclus	ions see 6	.0) and 6.1)							
ED030	ן כ	LX (manually	adjustab	le)										£1489.0	00
ED031		LT (adaptable lateral supports, manually adjustable) £1547.00													
Seat bottom for n	nodi	ular Recaro	back sup	port syst	em 📵	exclusions	s see 6.0) a	nd 6.1)							
EC008	J	F (flat padde	d, seat bot	tom sewi	n to latera	l guides, c	ushion thi	ckness ap	prox. 14 cı	m)				£958.0)0
EC009	ם [X (contoured	lateral gu	ides with	seat exte	nsion, cusl	hion thick	ness appr	ox. 16.5 cn	n)				£1126.0)0
EC010	⊒ [W (flat padde	ed, with se	at extens	ion, cushi	on thickne	ess approx	. 14.5 cm)						£965.0)0
* All Recaro seats are equipped with a head support.															
Settings for Stanc	dard	seat Junior	exclusion	s see 8.0) 1	to 8.12) an	d 9.0) to 9.	.2)								
Seat depth SD in o	cm		□34		36	□38	□40)							
Seat width SW in	cm		□34		36	□38	□40)							
Back support heig	ght E	3H 47 cm (ad	ljustable t	o 42 cm o	r 37 cm)										
Settings for Stanc	dard	seat Small e	xclusions	see 8.0) to	o 8.12) and	l 9.0) to 9.2	2)								
Seat depth SD in o	cm		□38		40	42	□ 44]46						
Seat width SW in	cm		□38		40	42									
Back support heig	ght E	3H 55 cm (ac	ljustable t	o 50 cm o	r 45 cm)										
Settings for Stanc	dard	seat Large	exclusions	see 8.0) t	o 8.16) and	d 9.0) to 9.2	2)								
Seat depth SD in	cm		42		44	☐ 46	□ 48	3 []50						
Seat width SW in	cm		44		46	□ 48									
Back support heig	ght E	3H 55 cm (ac	ljustable t	o 50 cm o	r 45 cm)										
Settings for Stanc	dard	seat XL excl	usions see	8.0) to 8.	16) and 9.0)) to 9.2)									
Seat depth SD in	cm		42		44	☐46	□ 48]50						
Seat width SW in	cm		□ 50		52	☐ 54	☐ 56	i							
Back support heig	ght E	3H 55 cm (ac	ljustable t	o 50 cm o	r 45 cm)										
Settings for VAS e	exclu	sions see 8.0)	to 8.2) an	d 9.0) to 9	.2)										
Seat depth SD in o	cm	□38	□ 40 □ 50		42	□ 44 □ 54	□ 46 □ 56]48]58						
Seat width SW in	cm.	□38			52 42	□ 34 □ 44] 48						_
Scat Width SW III	CIII	30	☐ 50			☐ 54		, <u> </u>	140						
Back support heig	ght E	3H 55 cm (ac	ljustable t	o 53 cm, 5	1 cm or 4	5 cm)									
Overview to seating	na si	vstems, seat	denths a	and seat	widths										
	95.	, , , , , , , , , , , , , , , , , , , ,	. асрано (arra seat											
Seat depth															П
Seat	Seat depth SD [cm]														
Seat		34	36	38	40	42	44	46	48	50	52	54	56	58	
Standard seat Jun	nior														
Standard seat Sm	all														
Standard seat Lar	ge														
Standard seat XL															
VAS															

Recaro n-Joy Plus

Seat width

Seat			Seat width SW [cm]											
Seat			34	36	38	40	42	44	46	48	50	52	54	56
Standard seat .	Junior	r												
Standard seat S	Small													
Standard seat I	Large													
Standard seat 2	ΧL													
VAS														
Lower leg leng	th LL						nout seat co	ushion; > 50 9.3)	0 cm not p	ossible in c	combinatio	n with Rec	aro)	
Seat cushion (p	ossibl	le only i	in combir	nation with	Standard	seat or VAS	S, mandato	ry informat	tion) exclu	sions see 9	.0) 9.1)			
EC011		т -	out cush				-	•						Inclusive
EC012		Seat	cushion	, black co	ver (seat o	cushion th	nickness 5	cm)						Inclusive
EC013		Seat	cushion	with inco	ntinence	cover (sea	at cushion	thickness 5	cm)					£99.00
EC034		Seat 11.6)	cushion	Livo (max	. user weig	ht up to 12	25 kg; seat	cushion thi	ckness 5 cr	n) exclusio	ns see 11.4	.)		£82.00
EC035		Seat	cushion	Evoflair (ı	nax. user v	veightup to	o 200 kg; se	eat cushion	thickness	8 cm) exclu	usions see	11.7)		£293.00
EC014		Terra	seat cus	shion (up	to 150 kg; s	eat cushio	n thicknes	s 6 cm) excl	usions see	11.1) 11.2)				£200.00
EC015		Terra	Aquos s	seat cushi	on (up to	150 kg; sea	t cushion t	hickness 6	cm) exclusi	ons see 11.	.1) 11.2)		-	£226.00
EC016		Terra	a Flair sea	at cushior	1 (up to 150) kg; seat c	ushion thic	kness 6 cm) exclusion	ns see 11.1)	11.3)			£452.00
EC018		Clou	d seat cu	d seat cushion (seat cushion thickness 11 cm) exclusions see 11.5) £553.00										
EC020		Cont	our, flat	our, flat contoured (seat cushion thickness 7 cm) exclusions see 11.8) 11.10) to 11.12) £244.00										
EC021		Cont	our, dee	p contou	red (seat c	ushion thic	ckness 9 cn	n) exclusior	ns see 11.8)	11.10) to 1	1.12)			£244.00
Please refer to ou	ur seat	t cushio	n brochu	re for info	mation on	additiona	l seat cushi	ons.						
Cover for seat of exclusions see 12			ntoured	(poss	ible only in	combinat	ion with fla	it or deep c	ontoured :	seat cushio	on, mandat	ory inform	ation)	
EC022		Fabri	ic cover											£290.00
EC023		Syntl	hetic lea	ther cove	r									£298.00
EC024		Fabri	ic cover	with inco	ntinence	cover								£324.00
Back support v		n/with	contoui	r 📵 (pos	sible only i	in combina	ation with S	Standard se	at or VAS,	mandatory	information	on)		
ED032		With	out back	k support	upholste	ry exclusio	ons see 9.0)							Inclusive
ED033		Back	support	tupholste	ry, adapt	able exclu	sions see 9	0.0)						Inclusive
ED034			our back sions see		pad, flat o	contoure	d (seat dep	th reductio	on by appro	ox. 2 cm)				£348.00
ED035			tour back sions see		pad, dee	p contour	red (seat d	epth reduc	tion by app	orox. 2 cm)				£348.00
ED050			tour back sions see		pad, flat	contoure	d AD (adar	otable, seat	depth red	uction by a	approx. 2 c	m)		£407.00
ED051			tour back sions see		pad, dee	p contoui	red AD (ad	aptable, se	at depth re	eduction b	y approx. 2	2 cm)		£407.00
ED036		sories	s 2-point-		seat depth			kg, withou a. 2 cm, not			stallation a	cces-		£673.00

Cover for back	SUDDA	ort pad, contoured (possible only in combination with flat or deep contoured back support pad,	
		n) exclusions see 14.0) 14.1)	
ED037		Fabric cover	Inclusive
ED038		Synthetic leather cover	£49.00
Accessories for	back	support pad, contoured ((possible only in combination with flat or deep contoured back support pad,	
		n) exclusions see 14.0) 14.1)	
ED039		Side panels (possible only in combination with deep contoured back support pad) exclusions see 15.1)	£92.00
ED040		Lumbar support	£69.00
Back support h	eight	for Baxx (possible only in combination with Baxx aluminium flat top back support; mandatory information) exclusions	
ED041		High (18"/46 cm)	Inclusive
ED042		Extra high (20"/51 cm)	Inclusive
Back support c	onto	ur for Baxx (possible only in combination with Baxx aluminium flat top back support; mandatory information) exclusions	
ED043		Active (2.5"/6 cm)	Inclusive
ED044		Deep (4"/10 cm)	£28.00
ED045		Extra deep (6"/15 cm)	£54.00
Colour for Baxx	x (poss	sible only in combination with Baxx aluminium flat top back support; mandatory information) exclusions see 16.0) 16.1)	
ED002		Anthrazit metallic	Inclusive
ED003		Jet black RAL 9005	Inclusive
ED004		Ice blue	Inclusive
ED005		Sparkle light orange	Inclusive
ED009		Vinho sparkle	Inclusive
ED012		Sky blue	Inclusive
ED013		Jet black RAL 9005 matte	Inclusive
ED014		Candy red	Inclusive
ED016		Silver metallic	Inclusive
ED017		Yellow RAL 1023	Inclusive
ED018		Signal red RAL 3001	Inclusive
ED021		Sparkle granny-smith	Inclusive
ED022		Creame RAL 9001	Inclusive
ED023		White RAL 9016	Inclusive
ED027		Shimano matte	Inclusive
ED029		Pink RAL 4010	Inclusive
Seat inclination	າ (pre-	tilt)	
		-3° (not possible in combination with Recaro n-Joy, n-Joy-Plus, modular seating system; front-wheel drive; mid-wheel dri combination of rear-wheel drive and VAS) exclusions see 10.1) 19.1) 19.2) 19.3)	ve;
		0°	
		3° (not possible in combination with Recaro n-Joy, n-Joy-Plus or modular seating system) exclusions see 10.1)	
		6° (not possible in combination with Recaro n-Joy, n-Joy-Plus or modular seating system) exclusions see 5.1) 10.1)	
		9° (not possible in combination with Recaro n-Joy, n-Joy-Plus or modular seating system) exclusions see 5.1) 10.1)	

		at tilt, without seat c additional approx. 1							+ 3 cm	I	I	
FWD, RWD, MV	VD		cm	41	43	45	47	49	51	53	55	57
Sizes 1 & 2			Inch	16.5	17.5	18	18.5	19.5	20.5	21	22	23
Notices and excl	usions			21.0) bis 21.6) 21.9) 21.10)	21.0) 21.5) bis 21.7) 21.9) 21.10) 21.14) bis 21.17)	21.0) 21.5) 21.6) 21.15) bis 21.17)	21.0) 21.15) 21.17)	8.12) 21.0) 21.6) 21.15)	8.12) 10.2) 21.0) 21.6)	5.1) 19.1) 21.0) 21.6) 21.12)	5.1) 19.1) 21.0) 21.6) 21.11) bis 21.13)	10.1) 19.1) 19.3) 21.0) 21.6) 21.8) 21.12) 21.13)
Standard seat												
Vario Adjust Se	eat (VA	NS)										
		Plus and seat modu seat cushion thickn										
Standard seat;	tilt 45	0										
Vario Adjust se	at (VA	S); tilt 45°										
Recaro (all); tilt	:45°											
Standard seat; seat height adj												
Vario Adjust se seat height adj	•	**										
Recaro (all); tilt seat height adj												
		without load from floo er load the specified s					, at 0° pre-t	tilt, without	seat height	adjustmer	nt and wit-	
Power seat adj	ustme	ent 📵 (possible onl	y in comb	ination wit	h 0° seat inc	clination)						
EC025		Power seat tilt 45° 20.29) 20.30)	with cer	ntre of grav	vity shiftin	g exclusior	ns see 20.0)	20.11) to 20	.13) 20.24)			£834.00
EC026		Power seat height 42.3)	adjustm	ent 35 cm	exclusions	see 20.0) 20	0.10) bis 20.	.14) 20.25) 2	0.29) 20.30)		£2302.00
EC027		Power seat height exclusions see 20.0)	-				th centre o	of gravity s	shifting			£2690.00
EC033		Without exclusions	see 20.15) to 20.23) 2	0.26) bis 20).28) 20.31) 2	20.32)					Inclusive
Back support a	ngle a	adjustment (mandate	ory inforn	nation in co	mbination	with Standa	ard seat or \	VAS) exclus	ions see 9.1)		
ED046		Mechanical in incr						S: -5/0/10/2	0/30°9; pos	sible		Inclusive
ED047		Power (Standard sea	at 30°; VA	S -5° to 30°)	exclusions	see 22.1)						£711.00
Mounting kit fo	or hea	d/neck supports (n	ot possib	le in combi	nation with	Recaro sea	t and seati	ng system)	exclusions s	see		
EL001		Gooseneck										£131.00
EL002		Multi-axis, offset										£122.00
EL003		Multi-axis, straigh	t									£179.00
EL031		UniLink holder for sible in combination		• •		,	for Advanc	ed line head	d supports;	not pos-		£243.00

UniLink access	ories	(may be used to compensate the distance at kyphotic head posture) exclusions see 83.0)	
EL039		UniLink extension, 38 mm	£84.00
EL038		UniLink extension, 50 mm	£62.00
EL040		UniLink extension, 75 mm	£84.00
Head/neck sup	ports	(possible only in comb. with mounting kit for head/neck supports, mandatory information) exclusions see 24.0) 24.1)	
Care line			
EL006		Head support, small exclusions see 24.2)	£92.00
EL007		Head support, large exclusions see 24.2)	£92.00
EL013		Head/neck support, small exclusions see 24.2)	£87.00
EL014		Head/neck support, large exclusions see 24.2)	£87.00
EL015		Head/neck support, X-shape, small exclusions see 24.2)	£105.00
EL016		Head/neck support, X-shape, large exclusions see 24.2)	£105.00
Advanced line			
EL033		Comfort Plus, size S	£162.00
EL034		Comfort Plus, size M	£188.00
EL035		Comfort Plus, size L	£188.00
EL036		Combo (3-pad), size S	£323.00
EL037		Combo (3-pad), size L	£336.00
Belts (man	datory	r information in combination with category B) exclusions see 25.0)	
EP001		Standard lap belt	Inclusive
EP002		2-point lap belt, padded, metal buckle (not possible in combination with Recaro seat and seating system) exclusions see 10.1)	£156.00
EP003		2-point lap belt, padded, plastic buckle (not possible in combination with Recaro seat and seating system) exclusions see 10.1)	£141.00
EP020		2-point lap belt with metal buckle and rewinder (buckle on side of joystick) exclusion see 10.1) 11.8) 25.2)	£236.00
Chest/shoulde	r harn	esses (not possible in combination with Recaro) exclusions see 10.1)	
EP005		Dynamic harness, size S	£200.00
EP006		Dynamic harness, size M	£200.00
EP007		Dynamic harness, size L	£200.00
EP008		Sternum harness, size S	£192.00
EP009		Sternum harness, size M	£192.00
EP010		Sternum harness, size L	£192.00
EP011		Slimline sternum harness, size S	£186.00
EP012		Slimline sternum harness, size M	£186.00
EP013		Slimline sternum harness, size L	£186.00
EP014		Shoulder harness, size S	£206.00
EP015		Shoulder harness, size M	£206.00
EP016		Shoulder harness, size L	£206.00
VAS accessorie	s (pos	sible only in combination with VAS)	
ED048		Push handles	£69.00

Arm supports	•		
EE023		Standard	Inclusive
EE024		Flip-up (effective seat width (+/- 1 cm per side) on Standard seat Small, SW 42 cm; Standard seat Large, SW 48 cm; Standard seat XL, SW 56 cm; VAS plus 2 cm) (max. user weight 125 kg; not in combination with Standard seat Junior and Recaro) exclusions see 10.1) 11.6) 11.8) 79.1) to 79.3)	£262.00
Padding for arm	suppo	ort and clothing guard can be ordered from the spare parts catalogue.	
Leg support ve	ersion	⊙	
EB001		Standard (LL is required, exclusions see 20.29) 28.1)	Inclusive
EB002		Mechanically elevating with length adjustment (max. load 140 kg; with length adjustment, not available for LL < 21 cm; LL is required) exclusions see 4.1) 28.1) to 28.5) 28.45)	£618.00
EB003		Power elevating with length adjustment and separate control (max. load 140 kg; not available for LL < 21 cm; LL is required) exclusions see 4.1) 28.1) to 28.7)	£1260.00
EB004		Power elevating with length adjustment and simultaneous control (not available for LL < 21 cm; LL is required) exclusions see 28.1) to 28.7)	£1260.00
EB013		Centrally mounted leg support (LL is required) exclusions see tables from page 22 (3); 28.1) 28.9) to 28.44) 28.46) to 28.48)	£895.00
EB005		Without leg support (not possible in combination with selected lower leg length) exclusions see 28.8)	Inclusive
Foot plate vers	ion [exclusions see 29.0) 29.1)	
EB006		Individual, plastic foot plate	Inclusive
EB007		Individual, aluminium foot plate exclusions see 5.1) 29.11)	£292.00
EB009		Single-panel, aluminium foot plate exclusions see 29.10)	£195.00
Leg support ac	cesso	ries	
EB010		Lateral leg pad for power and mechanical elevating leg supports (pair) (possible only in combination with mechanically or power elevating leg support adjustment) exclusions see 30.1)	£227.00
Batteries			
EU035		Gel 39 Ah (C5) / 50 Ah (C20) exclusions see 31.2) 31.3)	Inclusive
EU036		AGM 53 Ah (C5) / 60,4 Ah (C20) exclusions see 5.1) 31.3)	Inclusive
Battery charge	r 💽		
EU041		Battery charger 8 A, without fan (only for batteries with 39 Ah (C5)) exclusions see 32.1) 32.3)	Inclusive
EU042		Battery charger 12 A, incl. external charging receptacle (magnetic, installation side right) exclusions see 32.3)	Inclusive
Additional cha	rging	receptacle 📵	
EU043		External charging receptacle (magnetic, installation side right; only for battery charger EU040; recommended with existing 12A battery charger) exclusions see 75.1)	£198.00
Caster wheel s	wing	arm 📵	
EF001		Suspension	£156.00
EF002		No suspension exclusions see 5.1) 19.3)	Inclusive
Caster wheels			
EF003		6" PU Komfort (only for MWD) exclusions see 19.1) 64.6)	Inclusive
EF009		6" PU Performance (only for MWD, harder rubber mixture)	Inclusive
EF005		9" pneumatic, puncture-resistant exclusions see 19.3)	Inclusive
EF006		9" puncture-proof exclusions see 19.3)	Inclusive
EF007		10" pneumatic, puncture-resistant exclusions see 19.3)	Inclusive
EF008		10" puncture-proof exclusions see 19.3)	Inclusive

Drive wheels			
EG001		14" pneumatic, puncture-resistant	Inclusive
EG002		14" puncture-proof	Inclusive
Tyre colour			
EG011		Grey	Inclusive
EG012		Black	Inclusive
Tread 📵			
EG013		Ribs (only with tyre colour grey) exclusions see 37.1)	Inclusive
EG014		Lugs exclusions see 19.3)	Inclusive
EG015		Caster wheels with ribs, drive wheels with lugs exclusions see 37.2)	Inclusive
Wheelchair col	lour	·	
002		Anthrazit metallic	Inclusive
003		Jet black RAL 9005	Inclusive
004		Ice blue	Inclusive
005		Sparkle light orange	Inclusive
009		Vinho sparkle	Inclusive
012		Sky blue	Inclusive
013		Jet black RAL 9005 matte	Inclusive
014		Candy red	Inclusive
016		Silver metallic	Inclusive
017		Yellow RAL 1023	Inclusive
018		Signal red RAL 3001	Inclusive
021		Sparkle granny-smith	Inclusive
022		Cream RAL 9001	Inclusive
023		White RAL 9016	Inclusive
027		Shimano matte	Inclusive
037		Pink RAL 4010	Inclusive
Chassis accesso	ories		
EA010		Splash guard for drive wheels (included on FWD; not possible in conjunction with the combination of	£98.00
EA011		RWD, VAS and SH 43 cm) exclusions see 39.1) to 39.3) Splash guard for caster wheels	£38.00
EAUTI	Ш	Spiasii guaru for caster wheels	£38.00
Caster wheel s	wivel		
EF020		Mechanical (possible only in conjunction with the combination of splash guard for drive wheels and splash guard for caster wheels) exclusions see 19.3) 40.1)	£457.00
Chassis safety	acces	sories 💽	
EF021		Curb climbing assist (for RWD only) exclusions see 19.1) 19.3) 41.2)	£369.00
ES002		Kit for vehicle transport (not in combination with Baxx ED036; max. load 140 kg) exclusions see 4.1) 41.1)	Inclusive
ES013		Dahl docking system installation kit exclusions see 4.1) 41.1) 41.3)	£121.00
General access	ories	<u> </u>	
ES003		Pocket for mobile phone	£57.00
ES004		Luggage carrier exclusions see 5.1)	£259.00
	1		

ES005		Crutch holders	£131.00
ES006		Rearview mirror (installation side across from control panel side or swivel mechanism for mid-tray control or TEN° tray module) exclusions see 42.1)	£195.00
ES007		Tool kit	£67.00
ES010		Beverage holder	£29.00
ES011		Airman pump	£276.00
ET001		External horn	£536.00
ES012		Tray (swivel mechanism across from control panel side; not possible in combination with rearview mirror) exclusions see 42.2) 42.3)	£457.00
ES023		Backpack hook exclusion see 10.1)	£84.00
Socket for exte	rnal p	ower supply, right (front right next to safety switch, mandatory information) exclusions see 43.0) 43.1)	
ET002		12 V	£344.00
ET003		24 V	£138.00
ET004		USB	£155.00
Socket for exte	rnal r	power supply, left (front left next to safety switch, mandatory information) exclusions see 44.0) 44.1)	
ET005		12 V (not possible in combination with socket for external power supply, right, 12 V or 24 V) exclusions see 44.2)	£344.00
ET006		24 V (not possible in combination with socket for external power supply, right, 12 V or 24 V) exclusions see 44.2)	£138.00
ET007		USB	£155.00
Lighting 💽 e	xclusio	ons see 63.0)	
ES021		Electric lighting according to German Motor Vehicle Safety Standards (StVZO)	£735.00
Control device	(•)		
FLIOFO		VR2 (90 A), without power functions and lighting	la alcator
EU050		(not recommended for >120 kg user weight) exclusions see 45.1) 45.4) to 45.6)	Inclusive
EU051		VR2 (90 A), max. 2 power functions and lighting possible (not recommended for >120 kg user weight) exclusions see 45.3) to 45.5) 45.7)	Inclusive
EU061		VR2 (90 A) power module with tray module, manual flip over (max. 2 power funktions and lighting are possible) exclusions see 42.1) 45.3) to 45.5) 45.7)	£1519.00
EU053		R-net/TEN° (120 A) (mandatory for 14 km/h, recommended for >150 kg user weight)	£1126.00
Input for R-net	/TEN°	Power module exclusions see 46.0)	
EU065		TEN° hand module, standard colour Shimano matte	Inclusive
EU066		TEN° hand module, colour selectable	£1107.00
EU063		TEN° (120 A), tray module, electrical flip over (incl. 1 Piko-Button for rotation) exclusions see 42.1) 42.3) Note: The supplied horn can be perceived as too quiet; Ottobock recommends the use of an external horn, if no TEN° Display module is used.	£2499.00
493U75=ST188		Jog Wheel with Toggle	£41.80
Colour for TEN	° hand	d module (possible only in combination with control panel TEN°, mandatory information) exclusions see 47.0) 47.1)	
EU002		Anthrazit metallic	Inclusive
EU003		Jet black RAL 9005	Inclusive
EU004		Ice blue	Inclusive
EU005		Sparkle light orange	Inclusive
EU009		Vinho sparkle	Inclusive
EU012		Sky blue	Inclusive
EU013		Jet black RAL 9005 matte	Inclusive
EU014		Candy red	Inclusive

EU017	EU016		Silver metallic	Inclusive
EU021 Sparkle granny-smith Inclusive EU022 Cream RAL 9001 Inclusive EU023 White RAL 9016 Inclusive EU029 Pink RAL 4010 Inclusive EU029 Pink RAL 4010 Inclusive EU029 Pink RAL 4010 Inclusive EU116 Required force 250 g Inclusive EU117 Required force 120 g Inclusive EU117 Required force 120 g Inclusive EXternal buttons	EU017		Yellow RAL 1023	Inclusive
EU022	EU018		Signal red RAL 3001	Inclusive
EU023 White RAL 9016 Inclusive EU029 Pink RAL 4010 Inclusive Joystick for TEN* tray module exclusions see 80.0 80.1) EU116 Required force 250 g Inclusive EU117 Required force 120 g Inclusive EU117 Required force 120 g Inclusive EU117 Required force 120 g Inclusive EXERT ALL 4010 PROBLEM Set 10 Inclusive EXERT A	EU021		Sparkle granny-smith	Inclusive
EU029 Pink RAL 4010 Inclusive Diopytick for TEN** Temodule exclusions see 80.0) 80.1) EU116	EU022		Cream RAL 9001	Inclusive
Button module (for power seat functions) incl. 1 Piko button diameter 50 mm (for installation on TEN* control panel, Piko button red for on/off) Button module (for power seat functions) incl. 1 Piko button diameter 50 mm (for installation on TEN* control panel, Piko button red for on/off) Control panel holder, swing-away in parallel, height-adjustable Control panel holder, swing-away in parallel, neight-adjustable Control panel holder, swing-away in parallel, neight-adjustable Control panel holder, swing-away in parallel, neight-adjustable Control panel by the dealer) Control panel for TEN* tray module (for individual placement by the dealer) Control panel	EU023		White RAL 9016	Inclusive
EU116 Carrier Required force 250 g Inclusive Inclusive Required force 120 g Inclusive Required for 120 g Inclusive Required force 120 g	EU029		Pink RAL 4010	Inclusive
Etternal buttors Connection on TEN* control panel or tray module; possible only in combination with TEN* control device, electrical flip-over) exclusions see 48.0; EU067	Joystick for TEN	√l° tray	y module exclusions see 80.0) 80.1)	
External buttons (connection on TEN° control panel or tray module; possible only in combination with TEN° control device, electrical filip-over) exclusions see 48.0 EU067	EU116		Required force 250 g	Inclusive
EU067	EU117		Required force 120 g	Inclusive
EU068 2 Piko buttons, diameter 30 mm (red/black) £408.00 EU097 Button module (for power seat functions, installation opposite control panel side) £1451.00 EU098 Button module (for power seat functions) incl. 1 Piko button diameter 30 mm (for installation on TEN* control panel, Piko button red for on/off) EU099 Button module (for power seat functions) incl. 1 Piko button diameter 50 mm £1656.00 (for installation on TEN* control panel, Piko button red for on/off) Holder for hand the exclusions see 49.0) 49.2) EU070 Control panel holder, standard (possible only in combination with VR2) exclusions see 49.3) Inclusive EU071 Control panel holder, swing-away in parallel, height-adjustable £130.00 EU072 Control panel holder, swing-away in parallel, 9.5 cm height-adjustable £156.00 Mounting material tray for TEN* tray module exclusions see 80.0) 80.1) E5014 Tray with attachment, TEN* tray module centrally mounted £1016.00 Tray (without cut-out) with attachment, incl. cutting template for TEN* tray module £722.00 (for individual placement by the dealer) Attachment only, incl. cutting template for TEN* tray module (for individual placement by the dealer) Mounting side for tolder for holder for holder for holder for tray module (for individual placement by the dealer) £545.00 Mounting side for tolder for holder for tray module (for individual placement by the dealer) £722.00 (for individual placement by the dealer) Attachment only, incl. cutting template for TEN* tray module (for individual placement by the dealer) £99.00 Mounting side for tolder for half module or tray (for VR2, for mid-tray-told (MTC) or TEN* tray module; mandatory information) exclusions see 50.0) 50.1) Right		ns 졑	(connection on TEN° control panel or tray module; possible only in combination with TEN° control device, electrical	flip-over) exclu-
EU097 Button module (for power seat functions, installation opposite control panel side) £1451.00 EU098 Button module (for power seat functions) incl. 1 Piko button diameter 30 mm £1656.00 (for installation on TEN° control panel, Piko button red for on/off) EU099 Button module (for power seat functions) incl. 1 Piko button diameter 50 mm £1656.00 (for installation on TEN° control panel, Piko button red for on/off) EU097 Control panel holder, standard (possible only in combination with VR2) exclusions see 49.3) Inclusive EU071 Control panel holder, swing-away in parallel, height-adjustable £130.00 EU072 Control panel holder, swing-away in parallel, 9.5 cm height-adjustable £156.00 Mounting material tray for TEN° tray module exclusions see 80.0) 80.1) ES014 Tray with attachment, TEN° tray module centrally mounted £1016.00 ES015 Tray (without cut-out) with attachment, incl. cutting template for TEN° tray module £722.00 (for individual placement by the dealer) ES016 Attachment only, incl. cutting template for TEN° tray module (for individual placement by the dealer) MOUNTING SIDE for holder for hand module or tray (for VR2, for mid-tray control (MTC) or TEN° tray module; mandatory information) exclusions see 50.0) 50.1) Right Left	EU067		2 Piko buttons, diameter 50 mm (red/black)	£408.00
Button module (for power seat functions) incl. 1 Piko button diameter 30 mm (for installation on TEN° control panel, Piko button red for on/off) Button module (for power seat functions) incl. 1 Piko button diameter 50 mm (for installation on TEN° control panel, Piko button red for on/off) Holder for hand module exclusions see 49.0) 49.2) EU070	EU068		2 Piko buttons, diameter 30 mm (red/black)	£408.00
International Control panel, Piko button red for on/off)	EU097		Button module (for power seat functions, installation opposite control panel side)	£1451.00
Holder for hand motule exclusions see 49.0) 49.2) EU070	EU098			£1656.00
EU070	EU099			£1656.00
EU071	Holder for hand	d mod	dule exclusions see 49.0) 49.2)	
EU072	EU070		Control panel holder, standard (possible only in combination with VR2) exclusions see 49.3)	Inclusive
Mounting material tray for TEN° tray module exclusions see 80.0) 80.1) ES014	EU071		Control panel holder, swing-away in parallel, height-adjustable	£130.00
ES014	EU072		Control panel holder, swing-away in parallel, 9.5 cm height-adjustable	£156.00
ES015 Tray (without cut-out) with attachment, incl. cutting template for TEN° tray module (for individual placement by the dealer) Attachment only, incl. cutting template for TEN° tray module (for individual placement by the dealer) ES017 Cutting template for TEN° tray module (for individual placement by the dealer) Mounting side for holder for hand module or tray (for VR2, for mid-tray control (MTC) or TEN° tray module; mandatory information) exclusions see 50.0) 50.1) Right Left Left	Mounting mate	erial t	ray for TEN° tray module exclusions see 80.0) 80.1)	
(for individual placement by the dealer) Attachment only, incl. cutting template for TEN° tray module (for individual placement by the dealer) ES017 Cutting template for TEN° tray module (for individual placement by the dealer) Mounting side for holder for hand module or tray (for VR2, for mid-tray control (MTC) or TEN° tray module; mandatory information) exclusions see 50.0) 50.1) Right Left Left	ES014		Tray with attachment, TEN° tray module centrally mounted	£1016.00
ES016 ES017 Cutting template for TEN° tray module (for individual placement by the dealer) Mounting side for holder for hand module or tray (for VR2, for mid-tray control (MTC) or TEN° tray module; mandatory information) exclusions see 50.0) 50.1) Right Left Left	ES015			£722.00
Mounting side for holder for hand module or tray (for VR2, for mid-tray control (MTC) or TEN° tray module; mandatory information) exclusions see 50.0) 50.1) Right Left Logstick top	ES016			£545.00
(for VR2, for mid-tray control (MTC) or TEN° tray module; mandatory information) exclusions see 50.0) 50.1) Right Left Joystick top	ES017		Cutting template for TEN° tray module (for individual placement by the dealer)	£99.00
Right Left Joystick top				
Joystick top 💽	, , , , , , ,			
			Left	
	Joystick top	ל		
			Fork for tetraplegics, horizontal (adjustable for hand width 70 - 90 mm) exclusions see 51.1) 51.2)	£130.00
EU081 Stick S80 (special joystick handle 80 mm) exclusions see 51.1) 51.2) £115.00	EU081		Stick S80 (special joystick handle 80 mm) exclusions see 51.1) 51.2)	£115.00
EU082	EU082		Soft ball, black cover exclusions see 51.1) 51.2)	£195.00
EU083 Ball top, small (diameter 43 mm) exclusions see 51.1) 51.2) £131.00	EU083		Ball top, small (diameter 43 mm) exclusions see 51.1) 51.2)	£131.00
EU084	EU084		Ball top, large (diameter 48 mm) exclusions see 51.1) 51.2)	£121.00
EU085	EU085		Soft ball, coloured (diameter 30 mm, 4 pieces) exclusions see 51.1)	£236.00

Module for cor	ntrol c	levice R-net/TEN° (possible only in combination with TEN° control device) exclusions see 46.0) 52.0)	
		TEN° LCD module (when using 1-switch scan drive, programming by service required; please order but-	
EU054		ton separately; adapter switch box not required) exclusions see 52.0) 52.1)	£1856.00
EU055		Input/output module (required e.g. for Easywave radio control) exclusions see 52.2) 52.3)	£1193.00
			21175.00
Holder for TEN	° LCD	module (possible only in combination with TEN° LCD module, mandatory information) exclusions see 53.0) 53.1)	
EU073		Standard (installation across from control panel side; not possible in combination with tray) exclusions see 42.1) 51.2)	£95.00
EU074		With gooseneck (installation across from control panel side; not possible in combination with tray) exclusions see 42.1) 51.2)	£293.00
EU075		With gooseneck for tray (included)	£220.00
Special control	ls for ⁻	TEN° LCD module (possible only in combination with TEN° LCD module) exclusions see 52.0) 53.0)	
EU100		mo-Vis micro-joystick (actuating force 10 g)	£4319.00
EU101		mo-Vis multi-joystick (actuating force 50 g)	£3945.00
EU102		mo-Vis all-round joystick light (actuating force 120 g)	£3269.00
EU128		mo-Vis All-round Joystick (actuating force 250 g)	£2081.00
EU129		mo-Vis All-round Heavy Duty Joystick (actuating force 650 g)	£4014.00
EU103		Sip and puff control (possible only in combination with Recaro seat and seating system) exclusions see 25.1)	£1298.00
EU104		Button control (3 or 4 button switch control application, on-/off and mode included (at TEN° LCD module))	£765.00
Number of but	tons	for button control (possible only in comb. with button control, mandatory information) exclusions see 55.0) 55.1)	
		3 buttons	
		4 buttons	
Putton 1 for hi	ıtton	. Control (mandatory information with number of buttons for button control = 3, 4; only in combination with number of but	•
		l = 3, 4) exclusions see 72.0) 72.1)	
EU086		Piko button, diameter 50 mm, red	£195.00
EU087		Piko button, diameter 50 mm, yellow	£195.00
EU088		Piko button, diameter 50 mm, black	£195.00
EU089		Piko button, diameter 50 mm, silver	£195.00
EU090		Piko button, diameter 30 mm, red	£195.00
EU091		Piko button, diameter 30 mm, yellow	£195.00
EU092		Piko button, diameter 30 mm, black	£195.00
EU093		Piko button, diameter 30 mm, silver	£195.00
Dutton 2 for by			
		control (mandatory information with number of buttons for button control = 3, 4; only in combination with number of but = 3, 4) exclusions see 72.0) 72.1)	[-
EU186			£195.00
		Piko button, diameter 50 mm, red	2175.00
EU187		Piko button, diameter 50 mm, red Piko button, diameter 50 mm, yellow	£195.00
EU187 EU188	1		
		Piko button, diameter 50 mm, yellow	£195.00
EU188		Piko button, diameter 50 mm, yellow Piko button, diameter 50 mm, black	£195.00 £195.00
EU188 EU189		Piko button, diameter 50 mm, yellow Piko button, diameter 50 mm, black Piko button, diameter 50 mm, silver	£195.00 £195.00 £195.00
EU188 EU189 EU190		Piko button, diameter 50 mm, yellow Piko button, diameter 50 mm, black Piko button, diameter 50 mm, silver Piko button, diameter 30 mm, red	£195.00 £195.00 £195.00 £195.00

		control (mandatory information with number of buttons for button control = 3, 4; only in combination with number of = 3, 4) exclusions see 72.0) 72.1)	of but-
EU286		Piko button, diameter 50 mm, red	£195.00
EU287		Piko button, diameter 50 mm, yellow	£195.00
EU288		Piko button, diameter 50 mm, black	£195.00
EU289		Piko button, diameter 50 mm, silver	£195.00
EU290		Piko button, diameter 30 mm, red	£195.00
EU291		Piko button, diameter 30 mm, yellow	£195.00
EU292		Piko button, diameter 30 mm, black	£195.00
EU293		Piko button, diameter 30 mm, silver	£195.00
		control (mandatory information with number of buttons for button control = 4; only in combination with number of exclusions see 74.0) 74.1)	outtons
EU386		Piko button, diameter 50 mm, red	£195.00
EU387		Piko button, diameter 50 mm, yellow	£195.00
EU388		Piko button, diameter 50 mm, black	£195.00
EU389		Piko button, diameter 50 mm, silver	£195.00
EU390		Piko button, diameter 30 mm, red	£195.00
EU391		Piko button, diameter 30 mm, yellow	£195.00
EU392		Piko button, diameter 30 mm, black	£195.00
EU393		Piko button, diameter 30 mm, silver	£195.00
		nout TEN° LCD module (direct connection to bus system; possible only in combination with TEN° control devices (SN° LCD module) exclusions see 46.0) 57.0)	e; not possible
EU105		mo-Vis micro-joystick (actuating force 10 g)	£4479.00
EU106		mo-Vis multi-joystick (actuating force 50 g)	£4218.00
EU107		mo-Vis all-round joystick light (actuating force 120 g)	£3532.00
EU108		mo-Vis all-round joystick (actuating force 250 g)	£2465.00
EU109		mo-Vis all-round heavy duty joystick (actuating force 650 g)	£4408.00
Holder for mo-	Vis jo	ystick	
EU076		Power swivel arm, installation side left on the back support (to hold joysticks, sip & puff and button controls) exclusions see 42.3) 58.2) 58.4) 58.7)	£2867.00
EU077		Power swivel arm, installation side right on the back support (to hold joysticks, sip & puff and button controls) exclusions see 42.3) 58.2) 58.4) 58.7)	£2867.00
EU078		Holder for tray (must be installed by the customer) exclusions see 42.3)	£633.00
EU110		Holder for arm support incl. 1 finger button (for on/off, installation on control panel side) exclusions see 58.6)	£320.00
EU111		Holder for arm support incl. 2 finger buttons (for on/off and mode, installation on control panel side) exclusions see 58.5)	£499.00
EU112		Holder for arm support incl. 1 piko button (diameter 30 mm red for on/off, installation on control panel side) exclusions see 58.6)	£347.00
EU113		Holder for arm support incl. 2 piko buttons (diameter 30 mm red/black for on/off and mode, installation on control panel side) exclusions see 58.5)	£553.00

Actuator modu	ıle for	R-net/TEN° control device exclusions see 46.0) 59.0) 59.1)	
EU056		Actuator module (required for lighting and/or power seat functions on TEN° (120 A))	Inclusive
Control device	acces	ssories 📵	
EU057		Attendant control exclusions see 60.2)	£808.00
EU094		Easywave radio module (possible only in combination with input/output module) exclusions see 60.1)	£368.00
		stallation side (possible only in combination with attendant control; mandatory information)	
exclusions see 61	1.0) 61.		
		Right exclusions see 61.2) Left	
		datory information in conjunction with the combination of front-wheel drive and VR2 control device or front-wheel drive e) exclusions see 64.0)	
EF022		Mechanical (not possible in comb. with drive type MWD; possible only in comb. with chassis acc. splash guard for caster wheels; not possible in comb. of drive type FWD and speed > 7,2 km/h) exclusions see 19.3) 40.2) 64.3)	Inclusive
EF023		Power (gyro-module, included for FWD 10 km/h; not in combination with FWD and VR2; depending on terrain, the specified driving distance range can be reduced when using a gyro) exclusions see 2.1) 46.0)	£638.00
EF024		Mechanical incl. ASM (stability control; only in combination with TEN°, only for Australia, Canada or USA, not for RWD) exclusions see 19.3) 46.0) 64.5)	£455.00
Memory modu		power seat functions (possible only in combination with TEN° (120A) control device or TEN° tray module, electrical e 46.0)	
EU122		Back support angle adjustment, pre-programmed, 20° exclusions see 77.1) 77.2)	£550.00
EU123		Seat tilt, pre-programmed, 20° exclusions see 77.2) 77.3)	£564.00
EU124		Back support angle adjustment/seat tilt, pre-programmed 35°/15° exclusions see 77.4) to 77.6)	£1175.00
Option for elec	tric se	eat adjustment	
EU127		Speed limiter (sum angle >35°) exclusions see 77.3) 82.0) 82.1)	£252.00
Drive unit			
ET010		Standard (max. user weight 140 kg; possible only in comb. with overall weight <= 260 kg; not possible in comb. with TEN° control device (120A) or TEN° tray module, electrical flip-over)) exclusions see 4.1) 45.2) 65.1)	Inclusive
ET016		Performance (recommended for user weight > 100 kg; widening size 1 to 620 mm)	Inclusive
ET011		High performance (recommended for user weight > 100 kg; mandatory 14 km/h; (widening size 1 to 620 mm)	£640.00
Speed			
ET012		6 km/h	Inclusive
ET013		7.2 km/h exclusions see 66.1)	Inclusive
ET014		10 km/h (not in combination with FWD and VR2) exclusions see 5.1) 45.4) 64.2) 66.1) 66.5) 66.6)	Inclusive
ET015		14 km/h (only in comb. with RWD and TEN° (120A)) exclusions see 19.1) 19.3) 45.2) 46.0) 66.1) 66.3) 66.4)	Inclusive
Drive unit warr	anty	extension (only for drive motors, not for actuators)	
EW002		From 24 months to 48 months (the Mobility Solutions GmbH terms and exclusions apply) exclusions see 5.1)	£92.00
Additional inst	ructio	ons for use (additional IFU for non-national languages) exclusions see 68.0)	
ES028		Qualified personnel	£23.00
ES029		User	£23.00


Seating system 5



EC001 - EC004 Standard seat



EC005 VAS



EC006 Recaro n-Joy

Back support version 5

ED050/ED051 Contour back support pad AD

Seat cushionbottom for modular Recaro back support system 🤝



EC008 F



FC009 X



EC010 W

Back support version 5



ED036 Baxx aluminium flat top back support

Back support for modular Recaro seating system 5



ED031 Recaro LT



ED030 Recaro LX

Seat cushion/back support version contoured 5



EC020/ED034 Flat contoured



EC021/ED035 Deep contoured

Cover f. seat cushion/back support version contoured 5



EC022/ED037 Fabric cover



EC023/ED038 Synthetic leather cover

Accessories for back support pad, contoured 5



ED039 Siede panels



ED040 Lumbar support

Power seat adjustment 5



EC025 45° seat tilt



EC026 Seat height adjustment 35 cm

Belts 👈



EP002 Lap belt, metal buckle



EP003 Lap belt, plastic buckle

Belts ጛ



EP020 Belt with rewinder

Chest/shoulder harnesses 5



EP005 - EP007 Dynamic harness



EP008 - EP010 Sternum harness



EP011 - EP013 Slimline sternum harness

Chest/shoulder harn. \circlearrowleft



EP014 - EP017 Shoulder harness

Mounting kits for head/neck supports 5



EL001 Gooseneck



EL002 Mehrachsig, gekröpft



EL003 Multi-axis, straight

Mounting kits 5



EL031 UniLink holder for head support, multi-axis

Head/neck supports Advanced line 5



EL033/EL034/EL035 Comfort Plus



EL036/EL037 Combo

Head/neck supports Care line 5



EL006/EL007 Head support small/large



EL013/EL014 Combination head/neck support small/large



EL015/EL016 Head/neck support, X-shape, small/ large



Arm supports 5

EE024 Flip-up arm support

Leg support version 🝮



EB003 Power elevating, separate control



EB013 Centrally mounted leg support

Foot plate version 5



EB007 Individual aluminium foot plate

Add. charg. receptacle 5



EU043 External magnetic charging receptacle

Battery charger 👈



EU041 Battery charger 8 A, fanless



EU044 Battery charger 10 A, fanless (IP54)



EU046 Battery charger 10 A, with fan



EU042 Battery charger 12 A, incl. external charging receptacle

Tread 5



EG013 Ribs



EG014 Lugs



EG015 Caster wheel with ribs

Caster wheel swivel lock 5



EF020 Caster wheel swivel lock

Chassis safety accessories 5



EF021 Curb climbing assist



ES002 Kit for vehicle transport



ES013 Dahl docking system installation kit

General accessories 5



ES003 Pocket for mobile phone

General accessories 5



ES004 Luggage carrier



ES005 Crutch holder



ES006 Rearview mirror



ES010 Beverage holder

General accessories 5



ES023 Backpack holder

Socket for external power supply 5



ET002/005 12 V



ET004/007 USB

Lighting 5



ES020/021 Electric lighting for footpaths

Control device 5



EU050/EU051 VR2



EU053 TEN°



EU063 TEN° tray module, electrical flip over

External buttons 👈



EU098 Button module

Module for TEN° ★



EU054 TEN° LCD module

Control device accessorie 🝮



EU057 Attendant control

Caster wheel swing arm 5



EF001 Caster wheel swing arm, spring-mounted

Joystick top 👈



EU080 Fork for tetraplegics, horizontal



EU081 Stick S80



EU082 Soft ball, black cover



EU083 Ball top, small

Joystick top 👈



EU084 Ball top, large



EU085 Soft ball, coloured

Holder for mo-Vis joystick 👈



EU110-113 Holder for arm support



EU111 mo-Vis Joystick incl. 2 finger buttons

Special controls with TEN° LCD module 5/without TEN° LCD module 5



EU100/105 mo-Vis microjoystick



EU101/106 mo-Vis multijoystick




EU102/107 mo-Vis allround joystick light



EU103 Sip and puff control

 $Overview\ possible\ combinations\ centrally\ mounted\ leg\ support\ -\ Standard\ seat\ junior\ and\ Standard\ seat\ Small\ \ \ref{Standard}$

		SH [cm]	45		47			49			51		5	3	55		57
SD [cm]	Chassis size	Equipment	MW	FW	MW	RW	FW	MW	RW	FW	MW	RW	MW	RW	MW	RW	RW
		without tilt															
	1	with seat lift and tilt					$\overline{}$										
24		with tilt															
34		without tilt															
	2	with seat lift and tilt															
		with tilt															
		without tilt															
	36	with seat lift and tilt															
26		with tilt															
36		without tilt															
	2	with seat lift and tilt					$\overline{}$										
		with tilt															
		without tilt															
	1	with seat lift and tilt															
		with tilt															
38		without tilt															
	2	with seat lift and tilt															
		with tilt															
		without tilt															
	1	with seat lift and tilt															
		with tilt										Ĭ		ĺ			
40	2	without tilt															
		with seat lift and tilt															
		with tilt												ĺ			
		without tilt															
	1	with seat lift and tilt															
		with tilt															
42		without tilt															
	2	with seat lift and tilt															
		with tilt															
		without tilt															
	1	with seat lift and tilt					$\overline{/}$										
4.6		with tilt															
44		without tilt															
	2	with seat lift and tilt															
		with tilt															
		without tilt															
	1	with seat lift and tilt															
, ,		with tilt															
46		without tilt															
	2	with seat lift and tilt															
		with tilt															



Overview possible combinations centrally mounted leg support - Standard seat Large and Standard seat XL 🕏

		SH [cm]	4	43 45					7 SS larg	ie		47 SS XI		4	9 SS larg	ne .		49 SS XI	
SD [cm]	Chassis size	Equipment	FW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW
	1	Α																	
	1	В																	
		С																	
42		Α																	
	2	В																	
		С																	
		А																	
	1	В																	
	·	С																	
44		А																	
	2	В																	
		С																	
		А																	
	1	В																	
16		С																	
46	2	А																	
		В																	
		С																	
		Α																	
	1	В																	
48		С																	
40		А																	
	2	В																	
		С																	
		А																	
	1	В			/														
50		С	/_																
		А	/_													/_			
	2	В	/_		/_														
		С																	

Combination possible



- A without tilt
- B with seat lift and tilt
- C with tilt

Continuation overview possible combinations centrally mounted leg support - Standard seat Large and Standard seat XL 🕏

		SH [cm]	5	1 SS lar	ge		51 SS XI	_	5	53	5	5	57
SD [cm]	Chassis size	Equipment	FW	MW	RW	FW	MW	RW	MW	RW	MW	RW	RW
		without tilt											
	1	with seat lift and tilt											
42		with tilt											
42		without tilt											
	2	with seat lift and tilt											
		with tilt											
		without tilt											
	1	with seat lift and tilt											
		with tilt											
44		without tilt											
	2	with seat lift and tilt											
		with tilt			ĺ			ĺ					
		without tilt											
	1	with seat lift and tilt											
		with tilt											
46		without tilt											
	2	with seat lift and tilt											
		with tilt											
		without tilt											
	1	with seat lift and tilt											
		with tilt						ĺ					
48		without tilt											
	2	with seat lift and tilt											
		with tilt			ĺ								
		without tilt											
	1	with seat lift and tilt		/	/								
		with tilt											
50		without tilt											
	2	with seat lift and tilt		/									
		with tilt											

Combination possible

Combination not possible

Overview possible combinations centrally mounted leg support - VAS \circlearrowleft

		SH [cm]	4	13		45			47			49			51		5	3	55		5	7
SD [cm]	Chassis size	Equip- ment	FW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW	MW	RW	MW	RW	MW	RW
		А		/																		
	1	В															/					
		С																				
38		А																				
	2	В																				
		С																				
		А																				
	1	В																				
40		С																				
40		А																				
	2	В																				
		С																				
		Α																				
	1	В																				
42		С																				
72	2	А												\angle								
		В			\angle																\angle	
		С																				
	1	А			\angle						\angle			\angle								
		В			\angle																\angle	
44		С																				
		Α	\angle		\angle						<u>/</u>			<u>/</u>								
	2	В		\angle	\angle	\angle	/_														\angle	
		С																				
		А			\angle						/			/_								
	1	В			\angle																//	
46		С																				
		A			\angle									/								
	2	В			\angle																\angle	
		C																				
		A	/		/																	
	1	В			\angle																\angle	
48		C			/																	
		A			/						/											
	2	В	/																		/	
		С																				

Combination possible

Combination not possible

- A without tilt
- B with seat lift and tilt
- C with tilt

Continuation overview possible combinations centrally mounted leg support - VAS \circlearrowleft

		SH [cm]	4	3		45			47			49			51		53		5	5	5	7
SD [cm]	Chassis size	Equip- ment	FW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW	MW	RW	MW	RW	MW	RW
		А																				
	1	В																				
50		С																				
50		А																				
	2	В																				
		С																				
		А																				
	1	В																				
52		С																				
32		А																				
	2	В																				
		С																				
		А																				
	1 2	В	\angle		\angle	\angle															\angle	
54		С	<u>/</u>		\angle																	
	2	А	<u>/</u>		\angle																	
		В	<u>/</u>	/_	\angle	/_															<u>/</u> ,	
	ļ	С																				
		А																				
	1	В	4	\angle	4	/_															/_	
56		С																				
		А																				
	2	В	/_			/															/_	
		С																				
		A																				
	1	В	/_			/															/_	
58		С																				
		A																				
	2	В	/	/	/	/															/	
		С																				

Combination possible

Combination not possible

- A without tilt
- B with seat lift and tilt
- C with tilt