

2020/21

ottobock.

v2

Base Price -

EA004 - FWD - £4514.00

EA005 - RWD - £4119.00

EA006 - MWD - £5963.00

Valid 01.04.20 to 31.03.21

Please ensure you are using the latest order form

☐ Provide estimate only

☐ Order

☐ Attachment for custom fabrication enclosed

Billing address/customer number

Delivery address

Sales person

Name

.....
Date Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products, and you confirm to have read, understood and accepted our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

Standard equipment

-
-
-

Juvo B5 Power Wheelchair

☐ EA004 - 490E75=0_AS01_C

☐ EA005 - 490E75=0_AS01_C

☐ EA006 - 490E75=0_AS05_C

Max. load: 140 kg

(standard equipment)



Additional options installed on the power wheelchair can reduce the maximum load.

Please observe the maximum gross weight on the nameplate.

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders.

Please use the spare parts catalogue for spare parts orders.


Basic data (voluntary information)

Gender _____

Indication _____

Age _____

Paying party _____



Mandatory information*

User height in cm _____

User weight in kg _____

1 Seat depth _____

2 Trunk depth _____

3 Arm support height _____

4 Armpit height _____

5 Shoulder height _____

6 Head height _____

7 Forearm length to the middle of the palm _____

8 Trunk width _____

9 Lower leg length _____

10 Seat width _____

* This information will help to support user-oriented settings. Ottobock guarantees the confidential and purely internal use of these data.

Product classification			
EA001	<input type="checkbox"/>	Category A (indoor)	
EA002	<input type="checkbox"/>	Category B (in/outdoor)	
Drive type			
EA004	<input type="checkbox"/>	Front-wheel drive (FWD)	£4514.00
EA005	<input type="checkbox"/>	Rear-wheel drive (RWD)	£4119.00
EA006	<input type="checkbox"/>	Mid-wheel drive (MWD) exclusions see 2.1)	£5963.00
Max. user weight (user + personal luggage, information required to adapt chassis suspension)			
	<input type="checkbox"/>	70 kg	
	<input type="checkbox"/>	75 - 100 kg	
	<input type="checkbox"/>	100 - 140 kg	
Chassis size			
EA007	<input type="checkbox"/>	Size 1 (overall width approx. 60 cm; max. user weight 140 kg) exclusions see 4.1)	Inclusive
Seating systems			
EC001	<input type="checkbox"/>	Standard seat Junior (max. user weight 75 kg) (SD 34...40 cm, SW 34...40 cm, BH 47 cm to 37 cm, adjustable in 5 cm increments) exclusions see 5.2)	£117.00
EC002	<input type="checkbox"/>	Standard seat Small (SD 38...46 cm, SW 38...42 cm, BH 55 cm to 45 cm, adjustable in 5 cm increments) exclusions see 4.1)	Inclusive
EC003	<input type="checkbox"/>	Standard seat Large (SD 42...50 cm, SW 44...48 cm, BH 55 cm to 45 cm, adjustable in 5 cm increments) exclusions see 5.1)	Inclusive
EC004	<input type="checkbox"/>	Standard seat XL (SD 42...50 cm, SW 50...56 cm, BH 55 cm to 45 cm, adjustable in 5 cm increments) exclusions see 5.1)	£117.00
EC005	<input type="checkbox"/>	VAS (Vario Adjust Seat) (SD 38...58 cm, SW 38...54 cm, BH 55 cm, 53 cm, 51 cm, 45 cm)	£348.00
EC006	<input type="checkbox"/>	Recaro n-Joy* (manually adjustable back support, cushion thickness approx. 16.5 cm) exclusions see 5.1)	£1410.00
EC007	<input type="checkbox"/>	Recaro n-Joy-Plus* (adaptable on site, manually adjustable back support, SD 38 ... 56 cm, cushion thickness approx. 16.5 cm) exclusions see 5.1)	£2296.00

Back support for modular Recaro seating system exclusions see 6.0) and 6.1)			
ED030	<input type="checkbox"/>	LX (manually adjustable)	£538.00
ED031	<input type="checkbox"/>	LT (adaptable lateral supports, manually adjustable)	£601.00

Seat bottom for modular Recaro back support system exclusions see 6.0) and 6.1)			
EC008	<input type="checkbox"/>	F (flat padded, seat bottom sewn to lateral guides, cushion thickness approx. 14 cm)	£767.00
EC009	<input type="checkbox"/>	X (contoured lateral guides with seat extension, cushion thickness approx. 16.5 cm)	£902.00
EC010	<input type="checkbox"/>	W (flat padded, with seat extension, cushion thickness approx. 14.5 cm)	£773.00
* All Recaro seats are equipped with a head support.			

Settings for Standard seat Junior exclusions see 8.0) to 8.12) and 9.0) to 9.2)				
Seat depth SD in cm	<input type="checkbox"/> 34	<input type="checkbox"/> 36	<input type="checkbox"/> 38	<input type="checkbox"/> 40
Seat width SW in cm	<input type="checkbox"/> 34	<input type="checkbox"/> 36	<input type="checkbox"/> 38	<input type="checkbox"/> 40
Back support height BH 47 cm (adjustable to 42 cm or 37 cm)				

Settings for Standard seat Small exclusions see 8.0) to 8.12) and 9.0) to 9.2)					
Seat depth SD in cm	<input type="checkbox"/> 38	<input type="checkbox"/> 40	<input type="checkbox"/> 42	<input type="checkbox"/> 44	<input type="checkbox"/> 46
Seat width SW in cm	<input type="checkbox"/> 38	<input type="checkbox"/> 40	<input type="checkbox"/> 42		
Back support height BH 55 cm (adjustable to 50 cm or 45 cm)					

Settings for Standard seat Large exclusions see 8.0) to 8.16) and 9.0) to 9.2)					
Seat depth SD in cm	<input type="checkbox"/> 42	<input type="checkbox"/> 44	<input type="checkbox"/> 46	<input type="checkbox"/> 48	<input type="checkbox"/> 50
Seat width SW in cm	<input type="checkbox"/> 44	<input type="checkbox"/> 46	<input type="checkbox"/> 48		
Back support height BH 55 cm (adjustable to 50 cm or 45 cm)					

Settings for Standard seat XL exclusions see 8.0) to 8.16) and 9.0) to 9.2)					
Seat depth SD in cm	<input type="checkbox"/> 42	<input type="checkbox"/> 44	<input type="checkbox"/> 46	<input type="checkbox"/> 48	<input type="checkbox"/> 50
Seat width SW in cm	<input type="checkbox"/> 50	<input type="checkbox"/> 52	<input type="checkbox"/> 54	<input type="checkbox"/> 56	
Back support height BH 55 cm (adjustable to 50 cm or 45 cm)					

Settings for VAS exclusions see 8.0) to 8.2) and 9.0) to 9.2)								
Seat depth SD in cm	<input type="checkbox"/> 38	<input type="checkbox"/> 40	<input type="checkbox"/> 42	<input type="checkbox"/> 44	<input type="checkbox"/> 48	<input type="checkbox"/> 50	<input type="checkbox"/> 52	<input type="checkbox"/> 54
		<input type="checkbox"/> 56	<input type="checkbox"/> 58					
Seat width SW in cm	<input type="checkbox"/> 38	<input type="checkbox"/> 40	<input type="checkbox"/> 42	<input type="checkbox"/> 44	<input type="checkbox"/> 48	<input type="checkbox"/> 50	<input type="checkbox"/> 52	<input type="checkbox"/> 54
Back support height BH 55 cm (adjustable to 53 cm, 51 cm or 45 cm)								

Overview to seating systems, seat depths and seat widths

Seat depth

Seat	Seat depth SD [cm]												
	34	36	38	40	42	44	46	48	50	52	54	56	58
Standard seat Junior	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>									
Standard seat Small			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>						
Standard seat Large					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
Standard seat XL					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
VAS			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Recaro n-Joy Plus			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Seat width

Seat	Seat width SW [cm]											
	34	36	38	40	42	44	46	48	50	52	54	56
Standard seat Junior	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>								
Standard seat Small			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>							
Standard seat Large						<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
Standard seat XL									<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
VAS			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Lower leg length LL cm (15 ... 54 cm) (without seat cushion; > 50 cm not possible in combination with Recaro)
exclusions see 5.1) 8.12) 10.1) 10.2) 19.1) 19.3)

Seat cushion (possible only in combination with Standard seat or VAS, mandatory information) exclusions see 9.0) 9.1)

EC011	<input type="checkbox"/>	Without cushion	Inclusive
EC012	<input type="checkbox"/>	Seat cushion, black cover (seat cushion thickness 5 cm)	Inclusive
EC013	<input type="checkbox"/>	Seat cushion with incontinence cover (seat cushion thickness 5 cm)	£78.00
EC014	<input type="checkbox"/>	Terra seat cushion (up to 150 kg; seat cushion thickness 6 cm) exclusions see 11.1) 11.2)	£159.00
EC015	<input type="checkbox"/>	Terra Aquos seat cushion (up to 150 kg; seat cushion thickness 6 cm) exclusions see 11.1) 11.2)	£180.00
EC016	<input type="checkbox"/>	Terra Flair seat cushion (up to 150 kg; seat cushion thickness 6 cm) exclusions see 11.1) 11.3)	£330.00
EC018	<input type="checkbox"/>	Cloud seat cushion (seat cushion thickness 11 cm) exclusions see 11.5	£442.00
EC020	<input type="checkbox"/>	Contour, flat contoured (seat cushion thickness 7 cm) exclusions see 11.8) to 11.12)	£195.00
EC021	<input type="checkbox"/>	Contour, deep contoured (seat cushion thickness 9 cm) exclusions see 11.8) to 11.12)	£195.00
Please refer to our seat cushion brochure for information on additional seat cushions.			

Cover for seat cushion, contoured (possible only in combination with flat or deep contoured seat cushion, mandatory information) exclusions see 12.0) 12.1)

EC022	<input type="checkbox"/>	Fabric cover	£232.00
EC023	<input type="checkbox"/>	Synthetic leather cover	£37.00
EC024	<input type="checkbox"/>	Fabric cover with incontinence cover	£24.00

Back support version/with contour (possible only in combination with Standard seat or VAS, mandatory information) exclusions see 9.1)

ED032	<input type="checkbox"/>	Without back support upholstery exclusions see 9.0)	Inclusive
ED033	<input type="checkbox"/>	Back support upholstery, adaptable exclusions see 9.0)	Inclusive
ED034	<input type="checkbox"/>	Contour back support pad, flat contoured (seat depth reduction by approx. 2 cm) exclusions see 9.0) 11.8) 11.9)	£278.00
ED035	<input type="checkbox"/>	Contour back support pad, deep contoured (seat depth reduction by approx. 2 cm) exclusions see 9.0) 11.8) 11.9)	£278.00
ED050	<input type="checkbox"/>	Contour back support pad, flat contoured AD (adaptable, seat depth reduction by approx. 2 cm) exclusions see 11.9) 13.2)	£325.00
ED051	<input type="checkbox"/>	Contour back support pad, deep contoured AD (adaptable, seat depth reduction by approx. 2 cm) exclusions see 11.9) 13.2)	£325.00
ED036	<input type="checkbox"/>	Baxx aluminium flat top back support (max. load 113 kg, without shoulder recess, installation accessories 2-point-PRO, long, seat depth reduction by approx. 2 cm, not for Standard seat) exclusions see 13.1) 13.2) 13.3)	£539.00

Cover for back support pad, contoured (possible only in combination with flat or deep contoured back support pad, mandatory information) exclusions see 14.0) 14.1)			
ED037	<input type="checkbox"/>	Fabric cover	Inclusive
ED038	<input type="checkbox"/>	Synthetic leather cover	£37.00
Accessories for back support pad, contoured (possible only in combination with flat or deep contoured back support pad, mandatory information) exclusions see 14.0) 14.1)			
ED039	<input type="checkbox"/>	Side panels (possible only in combination with deep contoured back support pad) exclusions see 15.1)	£72.00
ED040	<input type="checkbox"/>	Lumbar support	£54.00
Back support height for Baxx (possible only in combination with Baxx aluminium flat top back support; mandatory information) exclusions see 16.0) 16.1)			
ED041	<input type="checkbox"/>	High (18"/46 cm)	Inclusive
ED042	<input type="checkbox"/>	Extra high (20"/51 cm)	Inclusive
Back support contour for Baxx (possible only in combination with Baxx aluminium flat top back support; mandatory information) exclusions see 16.0) 16.1)			
ED043	<input type="checkbox"/>	Active (2.5"/6 cm)	Inclusive
ED044	<input type="checkbox"/>	Deep (4"/10 cm)	£21.00
ED045	<input type="checkbox"/>	Extra deep (6"/15 cm)	£42.00
Colour for Baxx (possible only in combination with Baxx aluminium flat top back support; mandatory information) exclusions see 16.0) 16.1)			
ED002	<input type="checkbox"/>	Anthracite metallic	Inclusive
ED003	<input type="checkbox"/>	Jet black RAL 9005	Inclusive
ED004	<input type="checkbox"/>	Ice blue	Inclusive
ED005	<input type="checkbox"/>	Sparkle light orange	Inclusive
ED009	<input type="checkbox"/>	Vinho sparkle	Inclusive
ED012	<input type="checkbox"/>	Sky blue	Inclusive
ED013	<input type="checkbox"/>	Jet black RAL 9005 matte	Inclusive
ED014	<input type="checkbox"/>	Candy red	Inclusive
ED016	<input type="checkbox"/>	Silver metallic	Inclusive
ED017	<input type="checkbox"/>	Yellow RAL 1023	Inclusive
ED018	<input type="checkbox"/>	Signal red RAL 3001	Inclusive
ED021	<input type="checkbox"/>	Sparkle granny-smith	Inclusive
ED022	<input type="checkbox"/>	Creame RAL 9001	Inclusive
ED023	<input type="checkbox"/>	White RAL 9016	Inclusive
ED027	<input type="checkbox"/>	Shimano matte	Inclusive
ED029	<input type="checkbox"/>	Pink RAL 4010	Inclusive
Seat inclination (pre-tilt)			
	<input type="checkbox"/>	-3° (not possible in combination with Recaro n-Joy, n-Joy-Plus, modular seating system; front-wheel drive; mid-wheel drive; combination of rear-wheel drive and VAS) exclusions see 10.1) 19.1) 19.2) 19.3)	
	<input type="checkbox"/>	0°	
	<input type="checkbox"/>	3° (not possible in combination with Recaro n-Joy, n-Joy-Plus or modular seating system) exclusions see 10.1)	
	<input type="checkbox"/>	6° (not possible in combination with Recaro n-Joy, n-Joy-Plus or modular seating system) exclusions see 5.1) 10.1)	
	<input type="checkbox"/>	9° (not possible in combination with Recaro n-Joy, n-Joy-Plus or modular seating system) exclusions see 5.1) 10.1)	

Seat height (at 0° seat tilt, without seat cushion, suspension fully extended/without user) SH ≥ LL + 3 cm
 Recaro seat/system: additional approx. 13 cm for padding, recommended max. SH 49 cm

FWD, RWD, MWD Sizes 1 & 2	cm	41	43	45	47	49	51	53	55	57
	Inch	16.5	17.5	18	18.5	19.5	20.5	21	22	23
Notices and exclusions		21.0) bis 21.6) 21.9) 21.10)	21.0) 21.5) bis 21.7) 21.9) 21.10) 21.14) bis 21.17)	21.0) 21.5) 21.6) 21.15) bis 21.17)	21.0) 21.15) 21.17)	8.12) 21.0) 21.6) 21.15)	8.12) 10.2) 21.0) 21.6)	5.1) 19.1) 21.0) 21.6) 21.12)	5.1) 19.1) 21.0) 21.6) 21.11) bis 21.13)	10.1) 19.1) 19.3) 21.0) 21.6) 21.8) 21.12) 21.13)
Standard seat	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Vario Adjust Seat (VAS)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Recaro (n-Joy, n-Joy Plus and seat modules) + approx. 13 - 16 cm seat cushion thickness	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Standard seat; tilt 45°	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Vario Adjust seat (VAS); tilt 45°	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Recaro (all); tilt 45°	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Standard seat; tilt 45° & seat height adjustment	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Vario Adjust seat (VAS); tilt 45° & seat height adjustment	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Recaro (all); tilt 45° & seat height adjustment	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Seat height measured without load from floor to top of seat plate/seat frame (for Recaro), at 0° pre-tilt, without seat height adjustment and without seat cushion. Under load the specified seat height will be min. 1 cm lower.										

Power seat adjustment (possible only in combination with 0° seat inclination)

EC025	<input type="checkbox"/>	Power seat tilt 45° with centre of gravity shifting exclusions see 20.0) 20.10) to 20.13) 20.24) 20.29) 20.30)	£668.00
EC026	<input type="checkbox"/>	Power seat height adjustment 35 cm exclusions see 20.0) 20.10) bis 20.14) 20.25) 20.29) 20.30)	£1846.00
EC027	<input type="checkbox"/>	Power seat height adjustment 35 cm with 45° seat tilt with centre of gravity shifting exclusions see 20.0) 20.10) bis 20.14) 20.25) 20.29) 20.30)	£2158.00
EC033	<input type="checkbox"/>	Without exclusions see 20.15) to 20.23) 20.26) bis 20.28) 20.31) 20.32)	Inclusive

Back support angle adjustment (mandatory information in combination with Standard seat or VAS) exclusions see 9.1)

ED046	<input type="checkbox"/>	Mechanical in increments (Standard seat: 0/10/20/30° or -9/1/11/21°; VAS: -5/0/10/20/30°9; possible only in combination with Standard seat or VAS) exclusions see 9.0)	Inclusive
ED047	<input type="checkbox"/>	Power (Standard seat 30°; VAS -5° to 30°) exclusions see 22.1)	£569.00

Mounting kit for head/neck supports (not possible in combination with Recaro seat and seating system) exclusions see 10.1)

EL001	<input type="checkbox"/>	Gooseneck	£103.00
EL002	<input type="checkbox"/>	Multi-axis, offset	£97.00
EL003	<input type="checkbox"/>	Multi-axis, straight	£93.00
EL005	<input type="checkbox"/>	Rock-n-Lock, flip down	£296.00

Head/neck supports (possible only in comb. with mounting kit for head/neck supports, mandatory information) exclusions see 24.0) 24.1)

EL006	<input type="checkbox"/>	Head support, small	£72.00
EL007	<input type="checkbox"/>	Head support, large	£72.00
EL009	<input type="checkbox"/>	Head support, concave, size M	£141.00
EL011	<input type="checkbox"/>	Head support with ABS exterior shell, small	£114.00
EL012	<input type="checkbox"/>	Head support with ABS exterior shell, large	£144.00
EL013	<input type="checkbox"/>	Head/neck support, small	£68.00
EL014	<input type="checkbox"/>	Head/neck support, large	£68.00

EL015	<input type="checkbox"/>	Head/neck support, X-shape, small	£83.00
EL016	<input type="checkbox"/>	Head/neck support, X-shape, large	£83.00
EL017	<input type="checkbox"/>	Head support, concave, adjustable, one size only	£244.00
EL020	<input type="checkbox"/>	Head support, 3-piece, angle-adjustable, size L	£266.00

Belts (mandatory information in combination with category B) exclusions see 25.0)			
EP001	<input type="checkbox"/>	Standard lap belt	Inclusive
EP002	<input type="checkbox"/>	2-point lap belt, padded, metal buckle (not possible in combination with Recaro seat and seating system) exclusions see 10.1)	£124.00
EP003	<input type="checkbox"/>	2-point lap belt, padded, plastic buckle (not possible in combination with Recaro seat and seating system) exclusions see 10.1)	£112.00

Chest/shoulder harnesses (not possible in combination with Recaro) exclusions see 10.1)			
EP005	<input type="checkbox"/>	Dynamic harness, size S	£159.00
EP006	<input type="checkbox"/>	Dynamic harness, size M	£159.00
EP007	<input type="checkbox"/>	Dynamic harness, size L	£159.00
EP008	<input type="checkbox"/>	Sternum harness, size S	£152.00
EP009	<input type="checkbox"/>	Sternum harness, size M	£152.00
EP010	<input type="checkbox"/>	Sternum harness, size L	£152.00
EP011	<input type="checkbox"/>	Slimline sternum harness, size S	£148.00
EP012	<input type="checkbox"/>	Slimline sternum harness, size M	£148.00
EP013	<input type="checkbox"/>	Slimline sternum harness, size L	£148.00
EP014	<input type="checkbox"/>	Shoulder harness, size S	£165.00
EP015	<input type="checkbox"/>	Shoulder harness, size M	£165.00
EP016	<input type="checkbox"/>	Shoulder harness, size L	£165.00

VAS accessories (possible only in combination with VAS)			
ED048	<input type="checkbox"/>	Push handles	£54.00

Arm supports			
EE023	<input type="checkbox"/>	Standard	Inclusive
EE024	<input type="checkbox"/>	Flip-up (effective seat width (+/- 1 cm per side) on Standard seat Small, SW 42 cm; Standard seat Large, SW 48 cm; Standard seat XL, SW 56 cm; VAS plus 2 cm) (max. user weight 125 kg; not in combination with Standard seat Junior and Recaro) exclusions see 10.1) 11.6) 11.8) 79.1) to 79.3)	£82.00

Padding for arm support and clothing guard can be ordered from the spare parts catalogue.

Leg support version			
EB001	<input type="checkbox"/>	Standard (LL is required, exclusions see 28.1)	Inclusive
EB002	<input type="checkbox"/>	Mechanically elevating with length adjustment (max. load 140 kg; with length adjustment, not available for LL < 21 cm; LL is required) exclusions see 4.1) 28.1) to 28.5)	£495.00
EB003	<input type="checkbox"/>	Power elevating with length adjustment and separate control (max. load 140 kg; not available for LL < 21 cm; LL is required) exclusions see 4.1) 28.1) to 28.7)	£1010.00
EB004	<input type="checkbox"/>	Power elevating with length adjustment and simultaneous control (not available for LL < 21 cm; LL is required) exclusions see 28.1) to 28.7)	£1010.00
EB013	<input type="checkbox"/>	Centrally mounted leg support (LL is required) exclusions see tables from page 22 28.1) 28.9) to 28.44) for MWD	£446.00 Inclusive
EB005	<input type="checkbox"/>	Without leg support (not possible in combination with selected lower leg length) exclusions see 28.8)	Inclusive

Foot plate version exclusions see 29.0) 29.1)		
EB012	<input type="checkbox"/>	Individual, plastic foot plate (max. user weight 140 kg; not available on elevating leg supports) exclusions see 4.1) 29.2) to 29.4) 29.6) Inclusive
EB007	<input type="checkbox"/>	Individual, aluminium foot plate exclusions see 5.1) 29.7) to 29.9) 29.11) £112.00
EB009	<input type="checkbox"/>	Single-panel, aluminium foot plate exclusions see 29.7) to 29.10) £124.00
Leg support accessories		
EB010	<input type="checkbox"/>	Lateral leg pad for power and mechanical elevating leg supports (pair) (possible only in combination with mechanically or power elevating leg support adjustment) exclusions see 30.1) £114.00
Batteries		
EU035	<input type="checkbox"/>	Gel 39 Ah (C5) / 50 Ah (C20) exclusions see 31.2) 31.3) Inclusive
EU036	<input type="checkbox"/>	AGM 56 Ah (C5) / 60,4 Ah (C20) (exclusions see 5.1) 31.3) Inclusive
Battery charger		
EU041	<input type="checkbox"/>	Battery charger 8 A, without fan (only for batteries with 39 Ah or 53 Ah (C5)) exclusions see 32.1) Inclusive
EU042	<input type="checkbox"/>	Battery charger 12 A incl. external charging receptacle (magnetic, installation side right) Inclusive
Additional charging receptacle		
EU043	<input type="checkbox"/>	External magnetic charging receptacle (installation side right; only for battery charger EU040; recommended with existing 12A battery charger) exclusions see 75.1)0 £158.00
Caster wheel swing arm		
EF001	<input type="checkbox"/>	Suspension £124.00
EF002	<input type="checkbox"/>	No suspension exclusions see 5.1) 19.3) Inclusive
Caster wheels		
EF003	<input type="checkbox"/>	6" PU (only for MWD) exclusions see 19.1) 64.6) Inclusive
EF009	<input type="checkbox"/>	6" PU Performance (only for MWD, firmer rubber mixture) Inclusive
EF005	<input type="checkbox"/>	9" pneumatic, puncture-resistant exclusions see 19.3) Inclusive
EF006	<input type="checkbox"/>	9" puncture-proof exclusions see 19.3) Inclusive
EF007	<input type="checkbox"/>	10" pneumatic, puncture-resistant exclusions see 19.3) Inclusive
EF008	<input type="checkbox"/>	10" puncture-proof exclusions see 19.3) Inclusive
Drive wheels		
EG001	<input type="checkbox"/>	14" pneumatic, puncture-resistant Inclusive
EG002	<input type="checkbox"/>	14" puncture-proof Inclusive
Tyre colour		
EG011	<input type="checkbox"/>	Grey Inclusive
EG012	<input type="checkbox"/>	Black Inclusive
Tread		
EG013	<input type="checkbox"/>	Ribs (only with tyre colour grey) exclusions see 37.1) Inclusive
EG014	<input type="checkbox"/>	Lugs exclusions see 19.3) Inclusive
EG015	<input type="checkbox"/>	Caster wheels with ribs, drive wheels with lugs exclusions see 37.2) Inclusive

Wheelchair colour			
002	<input type="checkbox"/>	Anthracite metallic	Inclusive
003	<input type="checkbox"/>	Jet black RAL 9005	Inclusive
004	<input type="checkbox"/>	Ice blue	Inclusive
005	<input type="checkbox"/>	Sparkle light orange	Inclusive
009	<input type="checkbox"/>	Vinho sparkle	Inclusive
012	<input type="checkbox"/>	Sky blue	Inclusive
013	<input type="checkbox"/>	Jet black RAL 9005 matte	Inclusive
014	<input type="checkbox"/>	Candy red	Inclusive
016	<input type="checkbox"/>	Silver metallic	Inclusive
017	<input type="checkbox"/>	Yellow RAL 1023	Inclusive
018	<input type="checkbox"/>	Signal red RAL 3001	Inclusive
021	<input type="checkbox"/>	Sparkle granny-smith	Inclusive
022	<input type="checkbox"/>	Cream RAL 9001	Inclusive
023	<input type="checkbox"/>	White RAL 9016	Inclusive
027	<input type="checkbox"/>	Shimano matte	Inclusive
037	<input type="checkbox"/>	Pink RAL 4010	Inclusive
Chassis accessories			
EA010	<input type="checkbox"/>	Splash guard for drive wheels (included on FWD; not possible in conjunction with the combination of RWD, VAS and SH 43 cm) exclusions see 39.1) to 39.3) for FWD and MWD	£56.00 Inclusive
EA011	<input type="checkbox"/>	Splash guard for caster wheels for FWD and MWD	£30.00 Inclusive
ES001	<input type="checkbox"/>	Rear marker plate in accordance with German Motor Vehicle Safety Standards (StVZO)	£94.00
Caster wheel swivel lock			
EF020	<input type="checkbox"/>	Mechanical (possible only in conjunction with the combination of splash guard for drive wheels and splash guard for caster wheels) exclusions see 19.3) 40.1)	£366.00
Chassis safety accessories			
EF021	<input type="checkbox"/>	Curb climbing assist (for RWD only) exclusions see 19.1) 19.3) 41.2)	£295.00
ES002	<input type="checkbox"/>	Kit for vehicle transport (not in combination with Baxx ED036; max. load 140 kg) exclusions see 4.1) 41.1)	Inclusive
ES013	<input type="checkbox"/>	Dahl docking system installation kit exclusions see 4.1) 41.1) 41.3)	£96.00
General accessories			
ES003	<input type="checkbox"/>	Pocket for mobile phone	£22.00
ES004	<input type="checkbox"/>	Luggage carrier exclusions see 5.1)	£206.00
ES005	<input type="checkbox"/>	Crutch holders	£66.00
ES006	<input type="checkbox"/>	Rearview mirror (installation side across from control panel side or swivel mechanism for mid-tray control or TEN° tray module) exclusions see 42.1)	£100.00
ES007	<input type="checkbox"/>	Tool kit	£32.00
ES010	<input type="checkbox"/>	Beverage holder	£22.00
ES011	<input type="checkbox"/>	Airman pump	£124.00
ET001	<input type="checkbox"/>	External horn	£177.00
ES012	<input type="checkbox"/>	Tray (swivel mechanism across from control panel side; not possible in combination with rearview mirror) exclusions see 42.2) 42.3)	£366.00

Socket for external power supply, right (front right next to safety switch, mandatory information) exclusions see 43.0) 43.1)			
ET002	<input type="checkbox"/>	12 V	£274.00
ET003	<input type="checkbox"/>	24 V	£109.00
ET004	<input type="checkbox"/>	USB	£123.00

Socket for external power supply, left (front left next to safety switch, mandatory information) exclusions see 44.0) 44.1)			
ET005	<input type="checkbox"/>	12 V (not possible in combination with socket for external power supply, right, 12 V or 24 V) exclusions see 44.2)	£274.00
ET006	<input type="checkbox"/>	24 V (not possible in combination with socket for external power supply, right, 12 V or 24 V) exclusions see 44.2)	£109.00
ET007	<input type="checkbox"/>	USB	£123.00

Lighting exclusions see 63.0)			
ES021	<input type="checkbox"/>	Electric lighting according to German Motor Vehicle Safety Standards (StVZO)	£421.00

Control device			
EU050	<input type="checkbox"/>	VR2 (90 A), without power functions and lighting (not recommended for >120 kg user weight) exclusions see 45.1) 45.4) to 45.6)	Inclusive
EU051	<input type="checkbox"/>	VR2 (90 A), max. 2 power functions and lighting possible (not recommended for >120 kg user weight) exclusions see 45.3) to 45.5) 45.7)	Inclusive
EU061	<input type="checkbox"/>	VR2 (90 A) mid-tray control, manually reversible (max. 2 power funktions and lighting are possible) exclusions see 42.1) 45.3) to 45.5) 45.7)	£1217.00
EU053	<input type="checkbox"/>	TEN° (120 A) (mandatory for 14 km/h, recommended for >150 kg user weight)	£902.00
EU063	<input type="checkbox"/>	TEN° (120A) tray module, electrical flip over (incl. 1 Piko button for rotation) exclusions see 42.1) 42.3) Note: The supplied horn can be perceived as too quiet; Ottobock recommends the use of an external horn, if no TEN° Display module is used.	£2057.00

TEN° control panel (possible only in combination with control device TEN°, mandatory information) exclusions see 46.0) 46.1)			
EU065	<input type="checkbox"/>	Standard colour shimano matte	Inclusive
EU066	<input type="checkbox"/>	Colour selectable	£24.00
493U75=ST188	<input type="checkbox"/>	Jog wheel with Toggle	£35.40

Colour for TEN° control panel (possible only in combination with control panel TEN°, mandatory information) exclusions see 47.0) 47.1)			
EU002	<input type="checkbox"/>	Anthracyte metallic	Inclusive
EU003	<input type="checkbox"/>	Jet black RAL 9005	Inclusive
EU004	<input type="checkbox"/>	Ice blue	Inclusive
EU005	<input type="checkbox"/>	Sparkle light orange	Inclusive
EU009	<input type="checkbox"/>	Vinho sparkle	Inclusive
EU012	<input type="checkbox"/>	Sky blue	Inclusive
EU013	<input type="checkbox"/>	Jet black RAL 9005 matte	Inclusive
EU014	<input type="checkbox"/>	Candy red	Inclusive
EU016	<input type="checkbox"/>	Silver metallic	Inclusive
EU017	<input type="checkbox"/>	Yellow RAL 1023	Inclusive
EU018	<input type="checkbox"/>	Signal red RAL 3001	Inclusive
EU021	<input type="checkbox"/>	Sparkle granny-smith	Inclusive
EU022	<input type="checkbox"/>	Cream RAL 9001	Inclusive
EU023	<input type="checkbox"/>	White RAL 9016	Inclusive
EU029	<input type="checkbox"/>	Pink RAL 4010	Inclusive

Joystick for TEN° tray module exclusions see 80.0) 80.1)			
EU116	<input type="checkbox"/>	Required force 250 g	Inclusive
EU117	<input type="checkbox"/>	Required force 120 g	Inclusive
External buttons (connection on TEN° control panel or tray module; possible only in combination with TEN° control device, electrical flip-over) exclusions see 46.0)			
EU067	<input type="checkbox"/>	2 Piko buttons, diameter 50 mm (red/black)	£145.00
EU068	<input type="checkbox"/>	2 Piko buttons, diameter 30 mm (red/black)	£145.00
EU097	<input type="checkbox"/>	Button module (for power seat functions, installation opposite control panel side)	£1164.00
EU098	<input type="checkbox"/>	Button module (for power seat functions) incl. 1 Piko button diameter 30 mm (for installation on TEN° control panel, Piko button red for on/off)	£1328.00
EU099	<input type="checkbox"/>	Button module (for power seat functions) incl. 1 Piko button diameter 50 mm (for installation on TEN° control panel, Piko button red for on/off)	£1328.00
Control panel holder exclusions see 49.0) 49.2)			
EU070	<input type="checkbox"/>	Control panel holder, standard (possible only in combination with VR2) exclusions see 65.1)	Inclusive
EU071	<input type="checkbox"/>	Control panel holder, swing-away in parallel, height-adjustable	£102.00
EU072	<input type="checkbox"/>	Control panel holder, swing-away in parallel, 9.5 cm height-adjustable	£124.00
Mounting material tray for TEN° tray module exclusions see 80.0) 80.1)			
ES014	<input type="checkbox"/>	Tray with attachment, TEN° tray module centrally mounted	£814.00
ES015	<input type="checkbox"/>	Tray (without cut-out) with attachment, incl. cutting template for TEN° tray module (for individual placement by the dealer)	£578.00
ES016	<input type="checkbox"/>	Attachment only, incl. cutting template for TEN° tray module (for individual placement by the dealer)	£436.00
ES017	<input type="checkbox"/>	Cutting template for TEN° tray module (for individual placement by the dealer)	£78.00
Mounting side, control panel holder or tray (for VR2, for mid-tray control (MTC) or TEN° tray module; mandatory information) exclusions see 50.0) 50.1)			
	<input type="checkbox"/>	Right	
	<input type="checkbox"/>	Left	
Joystick top			
EU080	<input type="checkbox"/>	Fork for tetraplegics, horizontal (adjustable for hand width 70 - 90 mm) exclusions see 51.1) 51.2)	£102.00
EU081	<input type="checkbox"/>	Stick S80 (special joystick handle 80 mm) exclusions see 51.1) 51.2)	£46.00
EU082	<input type="checkbox"/>	Soft ball, black cover exclusions see 51.1) 51.2)	£83.00
EU083	<input type="checkbox"/>	Ball top, small (diameter 43 mm) exclusions see 51.1) 51.2)	£78.00
EU084	<input type="checkbox"/>	Ball top, large (diameter 48 mm) exclusions see 51.1) 51.2)	£96.00
EU085	<input type="checkbox"/>	Soft ball, coloured (diameter 30 mm, 4 pieces) exclusions see 51.1)	£188.00
Module for TEN° control device (possible only in combination with TEN° control device) exclusions see 48.0)			
EU054	<input type="checkbox"/>	TEN° LCD module (when using 1-switch scan drive, programming by service required; please order button separately; adapter switch box not required) exclusions see 52.1)	£1488.00
EU055	<input type="checkbox"/>	Input/output module (required e.g. for Easywave radio control) exclusions see 52.2)	£423.00
Holder for TEN° LCD module (possible only in combination with TEN° LCD module, mandatory information) exclusions see 53.0) 53.1)			
EU073	<input type="checkbox"/>	Standard (installation across from control panel side; not possible in combination with tray) exclusions see 42.1) 51.2)	£75.00
EU074	<input type="checkbox"/>	With gooseneck (installation across from control panel side; not possible in combination with tray) exclusions see 42.1) 51.2)	£234.00
EU075	<input type="checkbox"/>	With gooseneck for tray (included)	£175.00

Special controls with TEN° LCD module (possible only in combination with TEN° LCD module) exclusions see 53.0)			
EU100	<input type="checkbox"/>	mo-Vis micro-joystick (actuating force 10 g)	£2094.00
EU101	<input type="checkbox"/>	mo-Vis multi-joystick (actuating force 50 g)	£1912.00
EU102	<input type="checkbox"/>	mo-Vis all-round joystick light (actuating force 120 g)	£1624.00
EU103	<input type="checkbox"/>	Sip and puff control (possible only in combination with Recaro seat and seating system) exclusions see 25.1)	£1040.00
EU104	<input type="checkbox"/>	Button control (3 or 4 button switch control application, on-/off and mode included (at TEN° LCD module))	£440.00

Number of buttons for button control (possible only in comb. with button control, mandatory information) exclusions see 55.0) 55.1)			
	<input type="checkbox"/>	3 buttons	
	<input type="checkbox"/>	4 buttons	

Button 1 for button control (mandatory information with number of buttons for button control = 3, 4; only in combination with number of buttons for button control = 3, 4) exclusions see 72.0) 72.1)			
EU086	<input type="checkbox"/>	Piko button, diameter 50 mm, red	£94.00
EU087	<input type="checkbox"/>	Piko button, diameter 50 mm, yellow	£94.00
EU088	<input type="checkbox"/>	Piko button, diameter 50 mm, black	£94.00
EU089	<input type="checkbox"/>	Piko button, diameter 50 mm, silver	£94.00
EU090	<input type="checkbox"/>	Piko button, diameter 30 mm, red	£94.00
EU091	<input type="checkbox"/>	Piko button, diameter 30 mm, yellow	£94.00
EU092	<input type="checkbox"/>	Piko button, diameter 30 mm, black	£94.00
EU093	<input type="checkbox"/>	Piko button, diameter 30 mm, silver	£94.00

Button 2 for button control (mandatory information with number of buttons for button control = 3, 4; only in combination with number of buttons for button control = 3, 4) exclusions see 72.0) 72.1)			
EU186	<input type="checkbox"/>	Piko button, diameter 50 mm, red	£94.00
EU187	<input type="checkbox"/>	Piko button, diameter 50 mm, yellow	£94.00
EU188	<input type="checkbox"/>	Piko button, diameter 50 mm, black	£94.00
EU189	<input type="checkbox"/>	Piko button, diameter 50 mm, silver	£94.00
EU190	<input type="checkbox"/>	Piko button, diameter 30 mm, red	£94.00
EU191	<input type="checkbox"/>	Piko button, diameter 30 mm, yellow	£94.00
EU192	<input type="checkbox"/>	Piko button, diameter 30 mm, black	£94.00
EU193	<input type="checkbox"/>	Piko button, diameter 30 mm, silver	£94.00

Button 3 for button control (mandatory information with number of buttons for button control = 3, 4; only in combination with number of buttons for button control = 3, 4) exclusions see 72.0) 72.1)			
EU286	<input type="checkbox"/>	Piko button, diameter 50 mm, red	£94.00
EU287	<input type="checkbox"/>	Piko button, diameter 50 mm, yellow	£94.00
EU288	<input type="checkbox"/>	Piko button, diameter 50 mm, black	£94.00
EU289	<input type="checkbox"/>	Piko button, diameter 50 mm, silver	£94.00
EU290	<input type="checkbox"/>	Piko button, diameter 30 mm, red	£94.00
EU291	<input type="checkbox"/>	Piko button, diameter 30 mm, yellow	£94.00
EU292	<input type="checkbox"/>	Piko button, diameter 30 mm, black	£94.00
EU293	<input type="checkbox"/>	Piko button, diameter 30 mm, silver	£94.00

Button 4 for button control (mandatory information with number of buttons for button control = 4; only in combination with number of buttons for button control = 4) exclusions see 74.0) 74.1)			
EU386	<input type="checkbox"/>	Piko button, diameter 50 mm, red	£94.00
EU387	<input type="checkbox"/>	Piko button, diameter 50 mm, yellow	£94.00
EU388	<input type="checkbox"/>	Piko button, diameter 50 mm, black	£94.00
EU389	<input type="checkbox"/>	Piko button, diameter 50 mm, silver	£94.00
EU390	<input type="checkbox"/>	Piko button, diameter 30 mm, red	£94.00
EU391	<input type="checkbox"/>	Piko button, diameter 30 mm, yellow	£94.00
EU392	<input type="checkbox"/>	Piko button, diameter 30 mm, black	£94.00
EU393	<input type="checkbox"/>	Piko button, diameter 30 mm, silver	£94.00
Special controls without TEN° LCD module (direct connection to bus system; possible only in combination with TEN° control device; not possible in combination with TEN° LCD module) exclusions see 48.0) 57.0)			
EU105	<input type="checkbox"/>	mo-Vis micro-joystick (actuating force 10 g)	£2200.00
EU106	<input type="checkbox"/>	mo-Vis multi-joystick (actuating force 50 g)	£2008.00
EU107	<input type="checkbox"/>	mo-Vis all-round joystick light (actuating force 120 g)	£1769.00
EU108	<input type="checkbox"/>	mo-Vis all-round joystick (actuating force 250 g)	£1231.00
EU109	<input type="checkbox"/>	mo-Vis all-round heavy duty joystick (actuating force 650 g)	£2196.00
Holder for mo-Vis joystick exclusions see 58.0) 58.1)			
EU076	<input type="checkbox"/>	Power swivel arm, installation side left on the back support (to hold joysticks, sip & puff and button controls) exclusions see 42.3) 58.2) to 58.4)	£1424.00
EU077	<input type="checkbox"/>	Power swivel arm, installation side right on the back support (to hold joysticks, sip & puff and button controls) exclusions see 42.3) 58.2) to 58.4)	£1424.00
EU078	<input type="checkbox"/>	Holder for tray (must be installed by the customer) exclusions see 42.3)	£200.00
EU110	<input type="checkbox"/>	Holder for arm support incl. 1 finger button (for on/off, installation on control panel side) exclusions see 58.5)	£255.00
EU111	<input type="checkbox"/>	Holder for arm support incl. 2 finger buttons (for on/off and mode, installation on control panel side) exclusions see 58.5)	£400.00
EU112	<input type="checkbox"/>	Holder for arm support incl. 1 piko button (diameter 30 mm red for on/off, installation on control panel side) exclusions see 58.5)	£277.00
EU113	<input type="checkbox"/>	Holder for arm support incl. 2 piko buttons (diameter 30 mm red/black for on/off and mode, installation on control panel side) exclusions see 58.5)	£442.00
Actuator module for TEN° control device exclusions see 48.0) 59.0) 59.1)			
EU056	<input type="checkbox"/>	Actuator module (required for lighting and/or power seat functions on TEN° (120 A))	Inclusive
Control device accessories			
EU057	<input type="checkbox"/>	Attendant control	£647.00
EU094	<input type="checkbox"/>	Easywave radio module (possible only in combination with input/output module) exclusions see 60.1)	£294.00
Attendant control installation side (possible only in combination with attendant control; mandatory information) exclusions see 61.0) 61.1)			
	<input type="checkbox"/>	Right exclusions see 61.2)	
	<input type="checkbox"/>	Left	

Track stabiliser (mandatory information in conjunction with the combination of front-wheel drive and VR2 control device or front-wheel drive and TEN° control device) exclusions see 64.0)		
EF022	<input type="checkbox"/>	Mechanical (not possible in combination with drive type MWD; possible only in combination with chassis accessories splash guard for caster wheels; not possible in combination of drive type FWD and speed > 7,2 km/h) exclusions see 19.3) 40.2) 64.3) Inclusive
EF023	<input type="checkbox"/>	Power (gyro-module, included for FWD 10 km/h; not in combination with FWD and VR2; depending on terrain, the specified driving distance range can be reduced when using a gyro) exclusions see 48.0) 64.4) £510.00
Memory module for power seat functions (possible only in combination with TEN° (120A) control device or TEN° tray module, electrical flip-over) exclusions see 48.0)		
EU122	<input type="checkbox"/>	Back support angle adjustment, pre-programmed, 20° exclusions see 77.1) 77.2) £439.00
EU123	<input type="checkbox"/>	Seat tilt, pre-programmed, 20° exclusions see 77.2) 77.3) £452.00
EU124	<input type="checkbox"/>	Back support angle adjustment/seat tilt, pre-programmed 35°/15° exclusions see 77.4) to 77.6) £942.00
Option for electric seat adjustment		
EU127	<input type="checkbox"/>	Speed limiter (sum angle >35°) exclusions see 77.3) 82.0) 82.1) £155.00
Drive unit		
ET010	<input type="checkbox"/>	Standard (max. user weight 140 kg; possible only in combination with overall weight ≤ 260 kg; not possible in combination with TEN° control device (120A) or TEN° tray module, electrical flip-over)) exclusions see 4.1) 45.2) 65.1) Inclusive
ET016	<input type="checkbox"/>	Performance (recommended for user weight > 100 kg; widening size 1 to 620 mm) Inclusive
ET011	<input type="checkbox"/>	High performance (recommended for user weight > 100 kg; mandatory 14 km/h; (widening size 1 to 620 mm) £512.00
Speed		
ET012	<input type="checkbox"/>	6 km/h Inclusive
ET013	<input type="checkbox"/>	7.2 km/h exclusions see 66.1) Inclusive
ET014	<input type="checkbox"/>	10 km/h (not in combination with FWD and VR2) exclusions see 5.1) 45.4) 64.2) 66.1) 66.5) 66.6) (must order ES021 lighting with this option) Inclusive
ET015	<input type="checkbox"/>	14 km/h (only in combination with RWD and TEN° (120A)) exclusions see 19.1) 19.3) 45.2) 48.0) 66.1) 66.3) 66.4) (must order ES021 lighting with this option) Inclusive
Drive unit warranty extension (only for drive motors, not for actuators)		
EW002	<input type="checkbox"/>	From 24 months to 48 months (the Mobility Solutions GmbH terms and exclusions apply) exclusions see 5.1) £72.00
Additional instructions for use (additional IFU for non-national languages) exclusions see 68.0)		
ES028	<input type="checkbox"/>	Qualified personnel £17.00
ES029	<input type="checkbox"/>	User £17.00

Seating system



EC001 - EC004 Standard seat



EC005 VAS



EC006 Recaro n-Joy

Back support for modular Recaro seating system



ED031 Recaro LT



ED030 Recaro LX

Seat bottom for modular Recaro back support system



EC008 F



EC009 X



EC010 W

Seat cushion/back support version contoured



EC020/ED034 Flat contoured



EC021/ED035 Deep contoured

Cover f. seat cushion/back support version contoured



EC022/ED037 Fabric cover



EC023/ED038 Synthetic leather cover

Back support version



ED050/ED051 Contour back support pad AD



ED036 Baxx aluminium flat top back support

Accessories for back support pad, contoured



ED039 Siede panels




ED040 Lumbar support

Power seat adjustment



EC025 45° seat tilt



EC026 Seat height adjustment 35 cm

Mounting kit for head/neck supports



EL005 Rock-n-Lock flip down

Head/neck supports



EL006 Head support



EL011/EL012 Head support



EL013/EL014 Head support with ABS exterior shell



EL015/EL016 X-shaped head/neck support

Belts



EL017 Head support, concave



EL020 Head support, angle-adjustable



EP002 Lap belt, metal buckle




EP003 Lap belt, plastic buckle


Chest/shoulder harnesses




EP005 - EP007 Dynamic harness



EP008 - EP010 Sternum harness



EP011 - EP013 Slimline sternum harness



EP014 - EP017 Shoulder harness

Arm supports



EE024 Flip-up arm support

Leg support version



EB003 Power elevating, separate control



EB013 Centrally mounted leg support

Foot plate version



EB007 Individual aluminium foot plate

Battery charger



EU041 Battery Charger



EU044 Battery Charger



EU042 Battery Charger

Add. charg. receptacle



EU043 External magnetic charging receptacle

Caster wheel swing arm



EF001 Caster wheel swing arm, spring-mounted

Tread



EG013 Ribs



EG014 Lugs



EG015 Caster wheel with ribs

Caster wheel swivel lock



EF020 Caster wheel swivel lock

Chassis safety accessories



EF021 Curb climbing assist



ES002 Kit for vehicle transport



ES013 Dahl docking system installation kit

General accessories




ES003 Pocket for mobile phone



ES004 Luggage carrier



ES005 Crutch holder



ES006 Rearview mirror

Socket for external power supply



ES010 Beverage holder



ET002/005 12 V



ET004/007 USB

Lighting



ES020/021 Electric lighting for footpaths


Control device



EU050/EU051 VR2



EU053 TEN°



EU063 TEN° tray module, electrical flip over



External buttons



EU098 Button module

Joystick top



EU080 Fork for tetraplegics, horizontal



EU081 Stick S80




EU082 Soft ball, black cover

Module for TEN° ↻



EU083 Ball top, small



EU084 Ball top, large



EU085 Soft ball, coloured



EU054 TEN° LCD module

Special controls with TEN° LCD module / without TEN° LCD module



EU100/105 mo-Vis micro-joystick



EU101/106 mo-Vis multi-joystick



EU102/107 mo-Vis all-round joystick light



EU103 Sip and puff control

Holder for mo-Vis joystick



EU110-113 Holder for arm support



EU111 mo-Vis Joystick incl. 2 finger buttons

Control device accessories



EU057 Attendant control




EU095 Hand warmer (single hood for hand module)




EU096 Hand warmer (double hood for mid-tray control)

Overview possible combinations centrally mounted leg support - Standard seat junior and Standard seat Small



SD [cm]	Chassis size	Equipment	SH [cm]	45				47				49				51				53				55				57			
				MW	FW	MW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW	MW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW	FW
34	1	without tilt																													
		with seat lift and tilt																													
		with tilt																													
	2	without tilt																													
		with seat lift and tilt																													
		with tilt																													
36	1	without tilt																													
		with seat lift and tilt																													
		with tilt																													
	2	without tilt																													
		with seat lift and tilt																													
		with tilt																													
38	1	without tilt																													
		with seat lift and tilt																													
		with tilt																													
	2	without tilt																													
		with seat lift and tilt																													
		with tilt																													
40	1	without tilt																													
		with seat lift and tilt																													
		with tilt																													
	2	without tilt																													
		with seat lift and tilt																													
		with tilt																													
42	1	without tilt																													
		with seat lift and tilt																													
		with tilt																													
	2	without tilt																													
		with seat lift and tilt																													
		with tilt																													
44	1	without tilt																													
		with seat lift and tilt																													
		with tilt																													
	2	without tilt																													
		with seat lift and tilt																													
		with tilt																													
46	1	without tilt																													
		with seat lift and tilt																													
		with tilt																													
	2	without tilt																													
		with seat lift and tilt																													
		with tilt																													

 Combination possible

 Combination not possible

Overview possible combinations centrally mounted leg support - Standard seat Large and Standard seat XL


SD [cm]	Chassis size	SH [cm]		43		45			47 SS large			47 SS XL			49 SS large			49 SS XL		
		Equipment	FW	RW	FW	MW	RW		FW	MW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW
42	1	A																		
		B																		
		C																		
	2	A																		
		B																		
		C																		
44	1	A																		
		B																		
		C																		
	2	A																		
		B																		
		C																		
46	1	A																		
		B																		
		C																		
	2	A																		
		B																		
		C																		
48	1	A																		
		B																		
		C																		
	2	A																		
		B																		
		C																		
50	1	A																		
		B																		
		C																		
	2	A																		
		B																		
		C																		


-  Combination possible
 Combination not possible

- A without tilt
 B with seat lift and tilt
 C with tilt

Continuation overview possible combinations centrally mounted leg support - Standard seat Large and Standard seat XL


SD [cm]	Chassis size	SH [cm]	Equipment	51 SS large			51 SS XL			53		55		57
				FW	MW	RW	FW	MW	RW	MW	RW	MW	RW	RW
42	1		without tilt											
			with seat lift and tilt											
			with tilt											
	2		without tilt											
			with seat lift and tilt											
			with tilt											
44	1		without tilt											
			with seat lift and tilt											
			with tilt											
	2		without tilt											
			with seat lift and tilt											
			with tilt											
46	1		without tilt											
			with seat lift and tilt											
			with tilt											
	2		without tilt											
			with seat lift and tilt											
			with tilt											
48	1		without tilt											
			with seat lift and tilt											
			with tilt											
	2		without tilt											
			with seat lift and tilt											
			with tilt											
50	1		without tilt											
			with seat lift and tilt											
			with tilt											
	2		without tilt											
			with seat lift and tilt											
			with tilt											

 Combination possible

 Combination not possible

Overview possible combinations centrally mounted leg support - VAS



SD [cm]	SH [cm]		43		45			47			49			51			53		55		57	
	Chassis size	Equip- ment	FW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW	MW	RW	MW	RW	MW	RW
38	1	A																				
		B																				
		C																				
	2	A																				
		B																				
		C																				
40	1	A																				
		B																				
		C																				
	2	A																				
		B																				
		C																				
42	1	A																				
		B																				
		C																				
	2	A																				
		B																				
		C																				
44	1	A																				
		B																				
		C																				
	2	A																				
		B																				
		C																				
46	1	A																				
		B																				
		C																				
	2	A																				
		B																				
		C																				
48	1	A																				
		B																				
		C																				
	2	A																				
		B																				
		C																				



- A without tilt
- B with seat lift and tilt
- C with tilt

Continuation overview possible combinations centrally mounted leg support - VAS


SD [cm]	Chassis size	SH [cm]	43		45			47			49			51			53		55		57	
			FW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW	FW	MW	RW	MW	RW	MW	RW	MW	RW
50	1	A																				
		B																				
		C																				
	2	A																				
		B																				
		C																				
52	1	A																				
		B																				
		C																				
	2	A																				
		B																				
		C																				
54	1	A																				
		B																				
		C																				
	2	A																				
		B																				
		C																				
56	1	A																				
		B																				
		C																				
	2	A																				
		B																				
		C																				
58	1	A																				
		B																				
		C																				
	2	A																				
		B																				
		C																				

 Combination possible
 Combination not possible

- A without tilt
- B with seat lift and tilt
- C with tilt

Terra, Terra Aquos seat cushion												
SD [cm]	Seat width [cm]											
	34	36	38	40	42	44	46	48	50	52	54	56
34		X										
36		X										
38		X	X	X	X	X	X	X	X	X		
40				X	X	X	X					
42				X	X	X	X					
44				X	X	X	X	X	X	X		
46				X	X	X	X	X	X	X		
48				X	X	X	X	X	X	X		
50												
52												
54												
56												
58												

Terra Flair seat cushion																							
SD [cm]	Seat width [cm]																						
	34	36	38	40	42	44	46	48	50	52	54	56											
34																							
36																							
38													X	X	X	X	X	X	X	X	X		
40														X	X	X	X						
42														X	X	X	X						
44														X	X	X	X	X	X	X			
46														X	X	X	X	X	X	X			
48														X	X	X	X	X	X	X			
50																							
52																							
54																							
56																							
58																							



EC014 Terra



EC015 Terra Aquos



EC016 Terra Flair

Cloud seat cushion												
SD [cm]	Seat width [cm]											
	34	36	38	40	42	44	46	48	50	52	54	56
34	X	X										
36	X	X	X									
38	X	X	X	X	X		X	X	X	X	X	X
40	X	X	X	X	X	X	X	X	X	X	X	X
42			X	X	X	X	X	X	X	X	X	X
44				X	X	X	X	X	X	X	X	X
46				X	X		X	X	X	X	X	X
48				X	X		X	X	X	X	X	X
50							X			X	X	X
52							X			X	X	X
54												
56												
58												



EC018 Cloud

Notes and exclusions

- 2.1) Not possible in combination with Australia, Canada or the USA as the country of destination
- 4.1) Possible only in combination with max. user weight <= 140 kg
- 5.1) Possible only in combination with maximum user weight <= 160 kg
- 5.2) Possible only in combination with maximum user weight <= 100 kg
- 6.0) Possible only in combination with modular Recaro seating system
- 6.1) Mandatory information only in combination with modular Recaro seating system
- 8.0) Possible only in combination with Standard seat, VAS or Recaro n-Joy Plus
- 8.1) Mandatory information in combination with Standard seat, VAS or Recaro n-Joy-Plus
- 8.2) See table "Seat depth"
- 8.3) Not possible in conjunction with the combination of max. user weight > 75 kg and mid-wheel drive
- 8.4) Not possible in conjunction with the combination of max. user weight > 100 kg, chassis size 1 and Standard seat Small or VAS
- 8.5) Not possible in conjunction with the combination of max. user weight > 100 kg, rear-wheel drive and VAS
- 8.6) Not possible in conjunction with the combination of max. user weight > 100 kg, rear-wheel drive, chassis size 1 and VAS
- 8.7) Not possible in conjunction with the combination of max. user weight > 100 kg, rear-wheel drive, chassis size 2 and VAS
- 8.8) Not possible in conjunction with the combination of max. user weight > 100 kg, mid-wheel drive, chassis size 2 and Standard seat Small or VAS
- 8.9) Not possible in conjunction with the combination of max. user weight > 140 kg, mid-wheel drive and VAS
- 8.10) Not possible in conjunction with the combination of max. user weight > 140 kg, rear-wheel drive and VAS
- 8.11) Not possible in conjunction with the combination of maximum user weight > 140 kg, drive type front-wheel drive (FWD) and VAS
- 8.12) Not possible in conjunction with the combination of maximum user weight > 160 kg and drive type front-wheel drive (FWD)
- 8.13) Not possible in conjunction with the combination of maximum user weight > 160 kg and drive type mid-wheel drive (MWD)
- 8.14) Not possible in conjunction with the combination of maximum user weight > 180 kg and drive type mid-wheel drive (MWD)
- 8.15) Not possible in conjunction with the combination of maximum user weight > 180 kg and drive type front-wheel drive (FWD)
- 8.16) Not possible in conjunction with the combination of maximum user weight > 140 kg, drive type rear-wheel drive (RWD) or mid-wheel drive (MWD) and standard seat, large or standard seat, XL
- 9.0) Possible only in combination with Standard seat or VAS
- 9.1) Mandatory information in combination with Standard seat or VAS
- 9.2) See table "Seat width"
- 10.1) Not possible in combination with Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system
- 10.2) Not possible in conjunction with the combination of front-wheel drive and Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system
- 11.1) Possible only in combination with max. user weight <= 150 kg
- 11.2) See table "Terra, Terra Aquos seat cushion"
- 11.3) See table "Terra-Flair seat cushion"
- 11.4) See table "Z-Flo seat cushion"
- 11.5) See table "Cloud seat cushion"
- 11.6) Possible only in combination with max. user weight <= 125 kg
- 11.8) Not possible in combination with Standard seat Junior
- 11.9) Not possible in conjunction with the combination of VAS and SW 38 cm, 42 cm, 46 cm, 48 cm or 52 cm
- 11.10) Not possible in conjunction with the combination of VAS, SW 40 cm and SD 48 cm, 50 cm, 52 cm, 54 cm or 56 cm
- 11.11) Not possible in conjunction with the combination of VAS, SW 44 cm or 50 cm or 54 cm and SD 38 cm or 40 cm
- 11.12) Not possible in combination with SD 58 cm
- 12.0) Possible only in combination with seat cushion with flat or deep contours
- 12.1) Mandatory information in combination with seat cushion with flat or deep contours
- 13.1) Possible only in combination with max. user weight <= 113 kg
- 13.2) Possible only in combination with VAS
- 13.3) Not possible in combination with SW 54 cm
- 14.0) Possible only in combination with back version back cushion with flat or deep contours
- 14.1) Mandatory information in combination with back version back cushion with flat or deep contours
- 15.1) Possible only in combination with back version back cushion with deep contours
- 16.0) Possible only in combination with back version ADI aluminium flat top back
- 16.1) Mandatory information in combination with back version ADI aluminium flat top back
- 19.1) Not possible in combination with front-wheel drive
- 19.2) Not possible in conjunction with the combination of rear-wheel drive and VAS
- 19.3) Not possible in combination with mid-wheel drive
- 20.0) Possible only in combination with seat inclination 0°
- 20.10) Not possible in conjunction with the combination of max. user weight > 140 kg, front-wheel drive, VAS and SD 38 cm
- 20.11) Not possible in conjunction with the combination of max. user weight > 100 kg, mid-wheel drive, chassis size 1, Standard seat Small and SD 40 cm
- 20.12) Not possible in conjunction with the combination of max. user weight > 100 kg, rear-wheel drive, chassis size 1, Standard seat Large or Standard seat XL and SD 42 cm
- 20.13) Not possible in conjunction with the combination of max. user weight > 100 kg, mid-wheel drive, chassis size 1, Standard seat XL and SD 42 cm
- 20.14) Not possible in combination with LL >= 45 cm
- 20.15) Not possible in conjunction with the combination of max. user weight > 75 kg, front-wheel drive and SD 34 cm
- 20.16) Not possible in conjunction with the combination of max. user weight > 100 kg, front-wheel drive, chassis size 2, Standard seat Small or VAS and SD 38 cm
- 20.17) Not possible in conjunction with the combination of max. user weight > 100 kg, front-wheel drive, Standard seat Small or VAS and SD 40 cm
- 20.18) Not possible in conjunction with the combination of max. user weight > 100 kg, mid-wheel drive, chassis size 1, VAS and SD 40 cm
- 20.19) Not possible in conjunction with the combination of max. user weight > 100 kg, front-wheel drive, chassis size 1, Standard seat Large or Standard seat XL and SD 42 cm
- 20.20) Not possible in conjunction with the combination of max. user weight > 140 kg, front-wheel drive, VAS and SD 42 cm
- 20.21) Not possible in conjunction with the combination of max. user weight > 100 kg, rear-wheel drive, chassis size 1, VAS and SD 44 cm
- 20.22) Not possible in conjunction with the combination of max. user weight > 75 kg, mid-wheel drive, chassis size 2, Standard seat Small and SD 46 cm
- 20.23) Not possible in conjunction with the combination of max. user weight <= 100 kg, mid-wheel drive, chassis size 2, Standard seat XL and SD 50 cm
- 20.24) Not possible in conjunction with the combination of France as the country of destination, drive type Rear-wheel Drive (RWD), maximum user weight >= 150 kg and lower leg length >= 49 cm
- 20.25) Not possible in conjunction with the combination of maximum user weight > 100 kg and chassis, size 1
- 20.26) Not possible in conjunction with the combination of maximum user weight > 160 kg, drive type front-wheel drive (FWD) and seat depth 50 cm
- 20.27) Not possible in conjunction with the combination of maximum user weight > 160 kg, drive type mid-wheel drive (MWD) and seat depth 44 cm
- 20.28) Not possible in conjunction with the combination of maximum user weight > 180 kg, drive type mid-wheel drive (MWD) and seat depth 46 cm or 48 cm

- 20.29) Not possible in conjunction with the combination of maximum user weight > 70 kg, VAS and lower leg length 23 cm
- 20.30) Not possible in conjunction with the combination of maximum user weight > 140 kg, drive type rear-wheel drive (RWD), standard seat, large or standard seat, XL and seat depth 44 cm
- 20.31) Not possible in conjunction with the combination of maximum user weight > 140 kg, drive type front-wheel drive (FWD), standard seat, large or standard seat, XL and seat depth 42 cm
- 30.32) Not possible in conjunction with the combination of maximum user weight > 140 kg, drive type mid-wheel drive (MWD), standard seat, large and seat depth 44 cm
- 21.0) SH >= LL +3 cm
- 21.1) Not possible in combination with VAS, Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system
- 21.2) Not possible in conjunction with the combination of front-wheel drive, chassis size 1 and Standard seat Large
- 21.3) Not possible in conjunction with the combination of front-wheel drive and Standard seat XL
- 21.4) Not possible in conjunction with the combination of rear-wheel drive and Standard seat
- 21.5) Not possible in combination with 45° seat tilt
- 21.6) Not possible in combination with 35 cm seat height adjustment or 35 cm seat height adjustment combined with 45° seat tilt
- 21.7) Not possible in conjunction with the combination of rear-wheel drive and Standard seat Large or Standard seat XL and -3° seat inclination
- 21.8) Not possible in combination with Standard seat when 45° seat tilt is not selected
- 21.9) Not possible in conjunction with the combination of mid-wheel drive and Standard seat Large or Standard seat XL
- 21.10) Not possible in conjunction with the combination of mid-wheel drive, Standard seat Small and SW 42 cm
- 21.11) Not possible in conjunction with the combination of mid-wheel drive and Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system
- 21.12) Not possible in conjunction with the combination of France as the country of destination, drive type Rear-wheel Drive (RWD), maximum user weight >= 150 kg and electric seat adjustment seat tilt 45°
- 21.13) Cannot be combined with seat inclination -3°
- 22.1) Not possible in combination with Recaro n-Joy or Recaro n-Joy Plus
- 24.0) Possible only in combination with mounting kit for head/neck supports
- 24.1) Mandatory information in combination with the mounting kit for head/neck supports
- 25.0) Mandatory information in combination with category B
- 25.1) Possible only in combination with Recaro n-Joy or Recaro n-Joy-Plus or modular Recaro seating system
- 28.1) Possible only in combination with selected LL
- 28.2) Not possible in conjunction with the combination of Standard seat or VAS and LL < 21 cm
- 28.3) Not possible in conjunction with the combination of max. user weight > 70 kg and Standard seat and LL < 28 cm
- 28.4) Not possible in conjunction with the combination of max. user weight > 70 kg and VAS and LL < 27 cm
- 28.5) Not possible in conjunction with the combination of max. user weight > 70 kg and Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system and LL < 21 cm
- 28.6) Not possible in conjunction with the combination of rear-wheel drive, chassis size 2, Standard seat Small and SH 43 cm
- 28.7) Not possible in conjunction with the combination of rear-wheel drive, chassis size 2, Standard seat Small, -3° seat inclination and SH 43 cm
- 28.8) Not possible in combination with selected lower leg length
- 28.9) Not possible in conjunction with the combination of drive type Front-wheel Drive (FWD), Standard seat Junior or Standard seat Small and electric seat adjustment, without
- 28.10) Not possible in conjunction with the combination of drive type Front-wheel Drive (FWD), Standard seat Large or Standard seat XL, seat depth <= 46 cm and electric seat adjustment, without
- 28.11) Not possible in conjunction with the combination of drive type Front-wheel Drive (FWD), chassis, size 2, Standard seat Large or Standard seat XL, seat depth >= 48 cm, seat height 43 cm and electric seat adjustment, without
- 28.12) Not possible in conjunction with the combination of drive type Front-wheel Drive (FWD), chassis, size 2, Standard seat Large, seat depth 42 cm and electric seat adjustment seat height adjustment up to 35 cm or seat height adjustment up to 35 cm with seat tilt 45°
- 28.13) Not possible in conjunction with the combination of drive type Front-wheel Drive (FWD), chassis, size 2, Standard seat Large, seat depth 42 cm, seat height 47 cm or 49 cm and electric seat adjustment seat tilt 45°
- 28.14) Not possible in conjunction with the combination of drive type Front-wheel Drive (FWD), Standard seat and lower leg length < 28 cm
- 28.15) Not possible in conjunction with the combination of drive type Front-wheel Drive (FWD), VAS (Vario Adjust Seat), seat depth <= 46 cm and electric seat adjustment, without
- 28.16) Not possible in conjunction with the combination of drive type Front-wheel Drive (FWD), chassis, size 1, VAS (Vario Adjust Seat), seat depth 48 cm, seat height <= 47 cm and electric seat adjustment, without
- 28.17) Not possible in conjunction with the combination of drive type Front-wheel Drive (FWD), chassis, size 1, VAS (Vario Adjust Seat), seat depth 50 cm, seat height 43 cm and electric seat adjustment, without
- 28.18) Not possible in conjunction with the combination of drive type Front-wheel Drive (FWD), chassis, size 2, VAS (Vario Adjust Seat), seat depth 48 cm, seat height <= 49 cm and electric seat adjustment, without
- 28.19) Not possible in conjunction with the combination of drive type Front-wheel Drive (FWD), chassis, size 2, VAS (Vario Adjust Seat), seat depth 50 cm or 52 cm or 54 cm, seat height <= 45 cm and electric seat adjustment, without
- 28.20) Not possible in conjunction with the combination of drive type Front-wheel Drive (FWD), VAS (Vario Adjust Seat), Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system, lower leg length < 28 cm or > 43 cm and electric seat adjustment, without
- 28.21) Not possible in conjunction with the combination of drive type Front-wheel Drive (FWD), VAS (Vario Adjust Seat), Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system, lower leg length < 28 cm and electric seat adjustment seat tilt 45°, seat height adjustment up to 35 cm or seat height adjustment up to 35 cm with seat tilt 45°
- 28.22) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), chassis, size 1, Standard seat Junior or Standard seat Small, seat height <= 51 cm and electric seat adjustment, without
- 28.23) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), chassis, size 1, Standard seat Junior or Standard seat Small, seat depth <= 40 cm, seat height 53 cm or 55 cm and electric seat adjustment, without
- 28.24) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), chassis, size 1, Standard seat Junior, seat depth 34 cm or 36 cm, seat height 49 cm or 51 cm or 53 cm and electric seat adjustment seat tilt 45°
- 28.25) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), chassis, size 2, Standard seat Junior or Standard seat Small or Standard seat XL and electric seat adjustment, without
- 28.26) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), chassis, size 2, Standard seat Junior and electric seat adjustment seat height adjustment up to 35 cm or seat height adjustment up to 35 cm with seat tilt 45°
- 28.27) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), chassis, size 2, Standard seat Small, seat depth <= 42 cm and electric seat adjustment seat height adjustment up to 35 cm or seat height adjustment up to 35 cm with seat tilt 45°
- 28.28) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), chassis, size 2, Standard seat Junior, seat height 47 cm and electric seat adjustment seat tilt 45°
- 28.29) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), chassis, size 2, Standard seat Small, seat depth <= 42 cm, seat height 47 cm and electric seat adjustment seat tilt 45°
- 28.30) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), Standard seat, lower leg length < 28 cm or > 48 cm and electric seat adjustment, without
- 28.31) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), Standard seat, lower leg length < 28 cm and electric seat adjustment seat height adjustment up to 35 cm or seat height adjustment up to 35 cm with seat tilt 45°

- 28.32) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), Standard seat, lower leg length < 28 cm or > 49 cm and electric seat adjustment seat tilt 45°
- 28.33) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), chassis, size 1, VAS (Vario Adjust Seat), seat depth <= 42 cm, seat height 43 cm and electric seat adjustment, without
- 28.34) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), chassis, size 2, VAS (Vario Adjust Seat), seat depth <= 40 cm, seat height 43 cm and electric seat adjustment, without
- 28.35) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), VAS (Vario Adjust Seat), Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system, lower leg length < 28 cm and electric seat adjustment, without
- 28.36) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), VAS (Vario Adjust Seat), Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system, lower leg length < 28 cm and electric seat adjustment seat height adjustment up to 35 cm or seat height adjustment up to 35 cm with seat tilt 45°
- 28.37) Not possible in conjunction with the combination of drive type Rear-wheel Drive (RWD), VAS (Vario Adjust Seat), Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system, lower leg length < 28 cm or > 48 cm and electric seat adjustment seat tilt 45°
- 28.38) Not possible in conjunction with the combination of drive type Mid-wheel Drive (MWD), Standard seat Junior or Standard seat Small, seat height 41 cm or 43 cm and electric seat adjustment, without
- 28.39) Not possible in conjunction with the combination of drive type Mid-wheel Drive (MWD), Standard seat, lower leg length < 28 cm or > 48 cm and electric seat adjustment, without
- 28.40) Not possible in conjunction with the combination of drive type Mid-wheel Drive (MWD), Standard seat, lower leg length < 27 cm or > 40 cm and electric seat adjustment seat height adjustment up to 35 cm or seat height adjustment up to 35 cm with seat tilt 45°
- 28.41) Not possible in conjunction with the combination of drive type Mid-wheel Drive (MWD), Standard seat, lower leg length < 27 cm or > 48 cm and electric seat adjustment seat tilt 45°
- 28.42) Not possible in conjunction with the combination of drive type Mid-wheel Drive (MWD), VAS (Vario Adjust Seat), Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system, seat height 43 cm and electric seat adjustment, without
- 28.43) Not possible in conjunction with the combination of drive type Mid-wheel Drive (MWD), VAS (Vario Adjust Seat), Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system, lower leg length < 28 cm or > 49 cm and electric seat adjustment set tilt 45° or without
- 28.44) Not possible in conjunction with the combination of drive type Mid-wheel Drive (MWD), VAS (Vario Adjust Seat), Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system, lower leg length < 28 cm or > 41 cm and electric seat adjustment seat height adjustment up to 35 cm or seat height adjustment up to 35 cm with seat tilt 45°
- 29.0) Not possible in combination with leg support adjustment, none or centrally mounted
- 29.1) Mandatory information in combination with standard, mechanically elevating or power elevating leg support adjustment
- 29.2) Not possible in conjunction with the combination of Standard seat and LL < 25 cm
- 29.3) Not possible in conjunction with the combination of VAS and LL < 24 cm
- 29.4) Not possible in conjunction with the combination of Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system and LL < 18 cm
- 29.6) Not possible in combination with mechanically elevating or power elevating leg support adjustment
- 29.7) Not possible in conjunction with the combination of Standard seat or VAS and LL < 17 cm and standard leg support adjustment
- 29.8) Not possible in conjunction with the combination of max. user weight > 70 kg and Standard seat or VAS and LL < 24 cm and standard leg support adjustment
- 29.9) Not possible in conjunction with the combination of max. user weight > 70 kg and Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system and LL < 17 cm and standard leg support adjustment
- 29.10) Not possible in combination with mechanically elevating or power elevating leg support adjustment, separate control
- 29.11) Not possible in conjunction with the combination of max. user weight > 140 kg and power elevating leg support adjustment, simultaneous control
- 30.1) Possible only in combination with mechanically or power elevating leg support adjustment
- 31.1) Possible only in combination with chassis size 2
- 31.2) Possible only in combination with chassis, size 1
- 31.3) Not possible in conjunction with the combination of maximum user weight > 100 kg and electric seat adjustment seat height adjustment up to 35 cm or seat height adjustment up to 35 cm with seat tilt 45°
- 32.1) Possible only in combination with AGM 53 Ah (C5) or Gel 39 Ah (C5) batteries
- 32.2) Not possible in combination with Argentina, Australia, Brazil, Hong Kong, Japan, Canada, Colombia, Mexico, Thailand, the UK or USA as the country of destination
- 37.1) Possible only in combination with grey tyres
- 37.2) Not possible in conjunction with the combination of front-wheel drive or rear-wheel drive and black tyres
- 39.1) Mandatory information in combination with front-wheel drive or mid-wheel drive (derivation)
- 39.2) Not possible in conjunction with the combination of rear-wheel drive, VAS and SH 43 cm
- 39.3) Not possible in conjunction with the combination of rear-wheel drive, Recaro n-Joy, Recaro n-Joy-Plus or modular Recaro seating system and SH 43 cm
- 40.1) Possible only in conjunction with the combination of splash guard for drive wheels and splash guard for caster wheels
- 40.2) Possible only in combination with splash guard for caster wheels
- 41.1) Not possible in combination with ADI aluminium flat top back support
- 41.2) Not possible in combination with leg support adjustment, with central mounting
- 41.3) Not possible in combination with curb climbing assist
- 42.1) Not possible in combination with general accessories therapy tray
- 42.2) Not possible in combination with rearview mirror
- 42.3) Possible only in combination with maximum user weight <= 180 kg
- 43.0) Possible only in combination with external power supply, installation side right
- 43.1) Mandatory information in combination with external power supply, installation side right
- 44.0) Possible only in combination with external power supply, installation side left
- 44.1) Mandatory information in combination with external power supply, installation side left
- 44.2) Not possible in combination with socket for external power supply, right, 12 V or 24 V
- 45.1) Not possible in combination with power functions
- 45.2) Possible only in combination with overall weight <= 260 kg
- 45.3) Possible only in combination with max. 2 power functions
- 45.4) Possible only in combination with overall weight <= 320 kg
- 45.5) Possible only in combination with rear-wheel drive, mid-wheel drive or the combination of front-wheel drive and splash guard for caster wheels
- 45.6) Not possible in combination with lighting in accordance with Germany Motor Vehicle Safety Standards (StVZO)
- 45.7) Not possible in combination with modular Recaro seating system and power back support angle adjustment
- 46.0) Possible only in combination with TEN° control device
- 46.1) Mandatory information in combination with TEN° control device
- 47.0) Possible only in combination with TEN° control panel, colour selectable
- 47.1) Mandatory information in combination with TEN° control panel, colour selectable
- 48.0) Possible only in combination with control unit, TEN° (120A) or TEN° tray module, electrical flip over
- 49.0) Possible only in combination with VR2 control device, without power functions and lighting, VR2, max. 2 power functions and lighting, or TEN° (120 A)

- 49.1) Possible only in combination with VR2 control device
- 49.2) Mandatory information in combination with VR2 control device, without power functions and lighting, VR2, max. 2 power functions and lighting or TEN° (120 A)
- 50.0) Possible only in combination with control unit, VR2 or TEN° (120A) or in conjunction with the combination of control unit, TEN° tray module, electrical flip over and mounting material tray for tray module - tray with attachment, tray module centrally mounted or tray with attachment, incl. template or only attachment, incl. template
- 50.1) Mandatory information in combination with control unit, VR2 or TEN° (120A) or in conjunction with the combination of control unit TEN° tray module, electrical flip over and mounting material tray for tray module - tray with attachment, tray module centrally mounted or tray with attachment, incl. template or only attachment, incl. template
- 51.1) Not possible in combination with VR2, mid-tray control, manually reversible
- 51.2) Not possible in combination with control unit, TEN° tray module, electrical flip over
- 52.1) Not possible in conjunction with the combination of Standard seat and modules for control unit TEN° Input/Output module
- 52.2) Not possible in conjunction with the combination of Standard seat and modules for control unit TEN° - TEN° LCD module
- 53.0) Possible only in combination with TEN° LCD module
- 53.1) Mandatory information in combination with TEN° LCD module
- 55.0) Possible only in combination with modules for control unit TEN° - TEN° LCD-Modul
- 55.1) Mandatory information in combination with special controls with TEN° LCD modul button control
- 57.0) Not possible in combination with modules for control unit, TEN° - TEN° LCD module
- 58.0) Possible only in combination with special controls with TEN° LCD module – micro-joystick, multi-joystick, all-round joystick light or special controls – micro-joystick, multi-joystick, all-round joystick light, all-round joystick, all-round joystick heavy duty
- 58.1) Mandatory information in combination with special controls with TEN° LCD module – micro-joystick, multi-joystick, all-round joystick light or special controls – micro-joystick, multi-joystick, all-round joystick light, all-round joystick, all-round joystick heavy duty
- 58.2) Not possible in combination with special controls – all-round joystick or all-round joystick heavy duty
- 58.3) Not possible in conjunction with the combination of Standard seat and installation kit for head/neck supports Rock-n-Lock or Rock-n-Lock, flip down
- 58.4) Not possible in combination with external buttons – 2 piko buttons diameter 50 mm or 2 piko buttons diameter 30 mm
- 58.5) Not possible in combination with special controls with TEN° LCD module – all-round joystick light or special controls – all-round joystick light, all-round joystick or all-round joystick heavy duty
- 59.0) Mandatory information in combination with >= 1 power function
- 59.1) Mandatory information in conjunction with the combination of lighting according to German Motor Vehicle Safety Standards (StVZO) and TEN° control device or TEN° tray module, electrical flip over
- 60.1) Possible only in combination with input/output module
- 61.0) Possible only in combination with attendant control
- 61.1) Mandatory information in combination with attendant control
- 61.2) Not possible in conjunction with the combination of Standard seat Junior and actuator module for TEN° control device
- 63.0) Mandatory information **in conjunction with the combination of** all countries of destination except other countries, Argentina, Australia, Brazil, France, Hong Kong, Japan, Canada, Colombia, Korea, Mexico, Russian Federation, Switzerland, Thailand, the UK or USA **and product classification category B (in-/outdoor)**
- 63.1) Not possible in conjunction with the combination of Germany as the country of destination and category B
- 64.0) Mandatory information in conjunction with the combination of front-wheel drive and VR2 control device or front-wheel drive and TEN° control device
- 64.2) Not possible in conjunction with the combination of front-wheel drive and VR2 control device
- 64.3) Not possible in conjunction with the combination of front-wheel drive and speed > 7.2 km/h
- 64.4) Not possible in combination with Australia, Canada or the USA as the country of destination
- 64.5) Only possible in combination with Australia, Canada or the USA as the country of destination
- 64.6) Not possible in combination with rear-wheel drive
- 65.1) Not possible in combination with TEN° control device (120 A) or TEN° Tray module, electrical flip over
- 66.1) For Germany as the country of destination, possible only in conjunction with the combination of rear marker plate and lighting in accordance with Germany Motor Vehicle Safety Standards (StVZO)
- 66.3) Possible only in combination with drive unit high performance
- 66.4) Cannot be combined with batteries, 39Ah Gel (C5) or AGM 53Ah (C5))
- 66.5) Not possible in conjunction with the combination of maximum user weight > 100 kg and batteries, Gel 39Ah (C5) or AGM 53Ah (C5)
- 66.6) Cannot be combined with drive unit high performance
- 67.0) Possible only in conjunction with the combination of Germany as the country of destination and speed > 6 km/h
- 67.1) Mandatory information in conjunction with the combination of Germany as the country of destination and speed > 6 km/h
- 68.0) Not possible in conjunction with Australia as the country of destination
- 69.0) Possible only in combination with instructions for use (qualified personnel)
- 69.1) Mandatory information in combination with instructions for use (qualified personnel)
- 70.0) Possible only in combination with instructions for use (user)
- 70.1) Mandatory information in combination with instructions for use (user)
- 71.1) Theoretical overall weight (without special requests) = sum of the chosen characteristic values
- 71.2) Cannot be combined with theoretical total weight > 380 kg (the evaluated options exceed the permissible total weight of 380 kg)
- 72.0) Mandatory information with number of buttons for button control = 3, 4
- 72.1) Possible only in combination with number of buttons for button control = 3, 4
- 74.0) Mandatory information with number of buttons for button control = 4
- 74.1) Possible only in combination with number of buttons for button control = 4
- 75.1) Possible only in combination with battery charger
- 77.1) Possible only in combination with power back support angle adjustment
- 77.2) Not possible in combination with > 5 power functions
- 77.3) Possible only in combination with 45° seat tilt or 35 cm seat height adjustment combined with 45° seat tilt
- 77.4) Possible only in conjunction with the combination of power back support angle adjustment and 45° seat tilt or power back support angle adjustment and 35 cm seat height adjustment combined with 45° seat tilt
- 77.5) Not possible in combination with > 4 power functions
- 77.6) Not possible in combination with mechanical track stabiliser incl. ASM
- 79.1) Not possible in conjunction with the combination of Standard seat Small and SW 38 cm or SW 40 cm
- 79.2) Not possible in conjunction with the combination of Standard seat Large and SW 44 cm or SW 46 cm
- 79.3) Not possible in conjunction with the combination of Standard seat XL and SW 50 cm or SW 52 cm
- 80.0) Possible only in combination with control unit, TEN° tray module, electrical flip over
- 80.1) Mandatory information in combination with control unit, TEN° tray module, electrical flip over
- 82.0) Not possible in conjunction with France as the country of destination
- 82.1) Mandatory information in conjunction with the combination of maximum user weight > 160 kg and electric seat adjustment seat tilt 45° or seat height adjustment up to 35 cm with seat tilt 45° except with France as the country of destination