



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Company: _____

Account Number: _____ Date: _____

☐ Quote ☐ PO Number: _____

RTS: _____

Marked For: _____

Select if chair quote has: ☐ Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** _____

☐ **Little Wave Arc XP - Rigid Base** **£2995**

165 lb weight capacity. Price includes one time free growth kit.
Tilt range: 0° - 45°. Maximum Seat Width of 18".

☐ **Little Wave Arc XP - Folding Base** **£2995**

165 lb weight capacity. Price includes one time free growth kit.
Tilt range: 0° - 45°. Maximum Seat Width of 18".

Effective 18/07/2022

Frame

Frame Width

Select box corresponding to desired Frame Width

Frame Width Range	Frame Width		
☐ Small 10" - 12"	☐ 10	☐ 11	☐ 12
☐ Medium 12" - 14"	☐ 12	☐ 13	☐ 14
☐ Large 14" - 16"	☐ 14	☐ 15	☐ 16
☐ X-Large 16" - 18"	☐ 16	☐ 17	☐ 18

Frame Depth

Select Frame Depth within desired Frame Depth Range

Frame Depth Range	Frame Depth			
☐ X-Short 12" - 15"	☐ 12 N/C	☐ 13 N/C	☐ 14 N/C	☐ 15 N/C
☐ Short 14" - 17"	☐ 14 N/C	☐ 15 N/C	☐ 16 N/C	☐ 17 N/C
☐ Medium 16" - 19"	☐ 16 N/C	☐ 17 N/C	☐ 18 N/C	☐ 19 N/C
☐ Long 18" - 20"	☐ 18 N/C	☐ 19 N/C	☐ 20 N/C	

Transit Option

- ☐ Transit System _____ £247
☐ Omit Transit _____ N/C

Offset Hanger Option

- ☐ Offset Hanger _____ £32
Not available with X-Short 15" Frame Depth

Tilt Mechanism

- ☐ Hand Tilt - Right _____ N/C
☐ Hand Tilt - Left _____ N/C
☐ Hand Tilt - Dual _____ N/C
☐ Foot Tilt - Right _____ N/C
☐ Foot Tilt - Left _____ N/C

Seat Pan

- ☐ Growing Solid Seat Pan _____ N/C
☐ Omit Seat Pan _____ N/C

Vent and Battery Tray

Requires Rigid Base and Anti-Tips. Not available on frame widths less than 16" or Attendant Foot Lock or Foot Tilt.

- ☐ Vent Tray _____ £263
☐ Battery Box Holder _____ £263
Not available with X-Short Frame Depth Range
☐ Vent Tray and Battery Box Holder _____ £525
Not available with X-Short Frame Depth Range or Foot Tilt or Drum Brake or Hemi Wheel Lock. Must remove Ventilator and Battery for transit. Includes extended frame.

Backrest

Back Post / Height Range / Back Height

Stroller Handle Back Post (Folding)

- ☐ Short 20" £95
☐ Tall 24" £95

Fixed Height with Folding Adjustable Handle (Folding)

Not available with Seating Dynamics Rocker Back

- ☐ Short 22" N/C
☐ Tall 26" N/C

Reclining Back with Folding Adjustable Handle (Folding)

Not available with Seating Dynamics Rocker Back or O2 Holder. Includes extended frame. Requires Anti-tips.

- ☐ Short 22" £761
☐ Tall 26" £761

Stroller Handle Back Post Option

Available only with Stroller Handle Back Post

- ☐ Removable Stroller Handle Extension £110

Adjustable Height with Folding Adjustable Handle (Folding)

- ☐ Short 18"-25" N/C

Not available with Seating Dynamics Rocker Back

- ☐ 18" ☐ 19" ☐ 20" ☐ 21"
☐ 22" ☐ 23" ☐ 24" ☐ 25"

- ☐ Tall 20"-27" N/C

20"-22" Heights not available with Seating Dynamics Rocker Back

- ☐ 20" ☐ 21" ☐ 22" ☐ 23"
☐ 24" ☐ 25" ☐ 26" ☐ 27"

Other Backrest Options

- ☐ Non Folding Back N/C
Standard with Seating Dynamics Rocker Back

- ☐ Seating Dynamics Rocker Back £1680
Not available with Flip Up Armrests. Available only with Non Folding Back.

- ☐ Adjustable Height Rigidizer Bar £63

- ☐ Adjustable Height Folding Rigidizer Bar £231
Not available with Reclining Back

Footrest

Hanger / Front Frame Bend

Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount	☐	☐	☐	☐		Std
Front Mount		☐	☐	☐		N/C
Pediatric Pro Elevating Leg Rests (Peds PRO ELR)					☐	£247
Pro Elevating Leg Rests (PRO ELR)					☐	£247
Contracture Footrest ^{3, 4, 5} Bilateral Mount					☐	£336
Contracture Footrest ^{3, 5} Center Mount					☐	£336
Omit Hangers					☐	-£32

¹ Not available with 7" or 8" casters with Offset Hangers or 8" casters with Aluminum Angle Adjustable Footplates

² Not available with Frame Widths less than 14" without Offset Hangers or Frame Widths less than 12" with Offset Hangers

³ Aluminum Angle Adjustable footplates included. Requires Rigid Base. Not available with Transit. See image to right.

⁴ Not available with X-Short 12", Short 14", Medium 16", or Long 18" Frame Depths

⁵ Not available with X-Short 15" Frame Depth

PRO ELR Cover

Available only on PRO ELR


- ☐ Gel Pad Covers (Pair) £264

3" Extension


Available only with Front Mount Hangers. Not available with 90° hangers in combination with Offset Hangers and 8" casters.

- ☐ 3" Extension £42

- ☐ 3" Extension (Provided as Parts) £42



Bilateral Mount



Center Mount

Footrest continued

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length									
Extension Tube Hangers (60°, 70°, 80°)	☐ Short	☐ 10	☐ 10.5	☐ 11	☐ 11.5	☐ 12	☐ 12.5	☐ 13	☐ 13.5		
	☐ Medium	☐ 11.5	☐ 12	☐ 12.5	☐ 13	☐ 13.5	☐ 14	☐ 14.5	☐ 15	☐ 15.5	
	☐ Long	☐ 14	☐ 14.5	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	
Front Mount Hangers ¹		Standard Selections									Requires 3" Extension for £40
	☐ Short (70°, 80°, 90°)	☐ 4 ⁵	☐ 4.5 ⁵	☐ 5	☐ 5.5	☐ 6					☐ 8.5 ☐ 9 ☐ 9.5
		☐ 6.5	☐ 7	☐ 7.5	☐ 8						☐ 10 ☐ 10.5 ☐ 11
	☐ Medium (70°, 80°, 90°)	☐ 4 ⁵	☐ 4.5 ⁵	☐ 5	☐ 5.5	☐ 6	☐ 6.5	☐ 7			☐ 10.5 ☐ 11 ☐ 11.5
		☐ 7.5	☐ 8	☐ 8.5	☐ 9	☐ 9.5	☐ 10				☐ 12 ☐ 12.5 ☐ 13
	☐ Long (70°, 80°)	☐ 5	☐ 5.5	☐ 6	☐ 6.5	☐ 7	☐ 7.5	☐ 8	☐ 8.5		
		☐ 9	☐ 9.5	☐ 10	☐ 10.5	☐ 11	☐ 11.5	☐ 12			☐ 12.5 ☐ 13 ☐ 13.5
Contracture Footrest	☐ Short	☐ 10.5	☐ 11	☐ 11.5	☐ 12	☐ 12.5	☐ 13	☐ 13.5	☐ 14	☐ 14.5	
	☐ Medium	☐ 13.5	☐ 14	☐ 14.5	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	
	☐ Long	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5	
Peds ² PRO ELR	☐ Low Pivot	☐ 6 ³	☐ 6.5 ³	☐ 7 ⁴	☐ 7.5 ⁴	☐ 8 ⁴	☐ 8.5	☐ 9	☐ 9.5	☐ 10	☐ 10.5 ☐ 11 ☐ 11.5
	☐ High Pivot	☐ 6 ³	☐ 6.5 ³	☐ 7 ⁴	☐ 7.5 ⁴	☐ 8 ⁴	☐ 8.5	☐ 9	☐ 9.5	☐ 10	
PRO ELR Low Pivot Setting	☐ Short	10.5 - 15.5									
	☐ Long	13 - 19									
PRO ELR High Pivot Setting	☐ Short	8.5 - 17									
	☐ Long	11.5 - 17.5									

¹ See Minimum Front Seat Height chart to the right

² Comes standard with 3" Extension

³ Not available on chairs 13" wide or more

⁴ Not available on chairs 15" wide or more

⁵ Available only with 90°

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	14.5"	14.5"	14.5"	14.5"	14.5"	15.5"
Minimum Front Seat Height with 3" Extension	14.5"	14.5"	15.5"	16"	17"	18"

Footrest continued

Footplate


No selection needed with Contracture Footrest

- ☐ Composite _____ Std
Not available with Front Mount Hanger, Frame Widths less than 14" without Offset Hangers, Frame Widths less than 12" with Offset Hanger
- ☐ Composite Angle Adjustable _____ £158
Not available with Front Mount Hanger or 10"-12" Frame Widths. 13" and 14" Frame Widths require Offset Hanger.
- ☐ Aluminum Angle Adjustable _____ N/C
- ☐ Aluminum Locking Multi-Angle Adjustable _____ £263
Available only with Front Mount Hanger
- ☐ One Piece Flip Up Angle Adjustable _____ £152
Not available with Front Mount Hanger, PRO ELR, Frame Widths less than 14" without Offset Hanger or Frame Widths less than 12" with Offset Hanger

Hanger Release

No selection needed with Contracture Footrest. If no selection is made 4-Way - In line will be used.

- ☐ Classic _____ Std
- ☐ 4-Way - In line _____ N/C
See reference A below. Allows chair to be folded without removing footrest hanger. Not available with PRO ELR or Peds PRO ELR.
- ☐ 4-Way - Out _____ N/C
See reference B below. Not available with PRO ELR or Peds PRO ELR.



Heel Loop

Not available with One Piece Flip Up Angle Adjustable Footplate

- ☐ Heel Loops _____ £37
Not available with Frame Widths less than 14" without Offset Hangers or Frame Widths less than 12" with Offset Hangers
- ☐ Adjustable Heel Loops _____ £42
- ☐ Omit Heel Loops _____ N/C

Calf Strap

- ☐ Velcro Adjustable _____ £42
- ☐ Padded Velcro Adjustable _____ £58
- ☐ Bodypoint Padded AeroMesh _____ £89

Residual Limb Support

Residual Limb Support

- ☐ Right and Left _____ £315
- ☐ Left _____ £158 ☐ Right _____ £158
- ☐ Omit Left Hanger _ -£16 ☐ Omit Right Hanger _ -£16

Cover

- ☐ Smooth ☐ Air Mesh ☐ Smooth ☐ Air Mesh

Pad Size

- | | | | |
|---------------------------------|----------------------------------|---------------------------------|----------------------------------|
| ☐ 5 x 8 | ☐ 9 x 14 | ☐ 5 x 8 | ☐ 9 x 14 |
| ☐ 7 x 10 | ☐ 11 x 10 | ☐ 7 x 10 | ☐ 11 x 10 |
| ☐ 9 x 10 | ☐ 11 x 14 | ☐ 9 x 10 | ☐ 11 x 14 |

Seat Height, Casters and Rear Wheels

Seat Height / Caster Size / Rear Wheel Size

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size			
	5"	6"	7"	8"
Seat Height				
12" with Pneumatic or Pneumatic with Airless Insert Tire	☐ 14.5			
	☐ 15			
	☐ 15.5	☐ 15.5		
	☐ 16	☐ 16		
	☐ 16.5	☐ 16.5		
12" with Full Poly Tire	☐ 14.5			
	☐ 15			
	☐ 15.5	☐ 15.5		
12" Drum Brake	☐ 16	☐ 16		
	☐ 14.5			
	☐ 15			
	☐ 15.5	☐ 15.5		
16" with Pneumatic, Pneumatic with Airless Insert or Full Poly Tire	☐ 16	☐ 16		
	☐ 16.5	☐ 16.5		
	☐ 17	☐ 17	☐ 17	
	☐ 17.5	☐ 17.5	☐ 17.5	
	☐ 18	☐ 18	☐ 18	☐ 18
	☐ 18.5	☐ 18.5	☐ 18.5	☐ 18.5
	☐ 14.5			
	☐ 15			
	☐ 15.5	☐ 15.5		

Rear Wheel Size	Caster Size			
	5"	6"	7"	8"
Seat Height				
16" with Low Poly Tire	☐ 14.5			
	☐ 15			
	☐ 15.5	☐ 15.5		
	☐ 16	☐ 16		
	☐ 16.5	☐ 16.5		
	☐ 17	☐ 17	☐ 17	
16" Drum Brake	☐ 17.5	☐ 17.5	☐ 17.5	
	☐ 18	☐ 18	☐ 18	☐ 18
	☐ 15			
	☐ 15.5	☐ 15.5		
	☐ 16	☐ 16		
20"	☐ 16.5	☐ 16.5		
	☐ 17	☐ 17	☐ 17	
	☐ 17.5	☐ 17.5	☐ 17.5	
	☐ 18	☐ 18	☐ 18	☐ 18
	☐ 18.5	☐ 18.5	☐ 18.5	☐ 18.5
	☐ 19	☐ 19	☐ 19	☐ 19
	☐ 14.5			
	☐ 15			
	☐ 15.5	☐ 15.5		

Seat Height, Casters and Rear Wheels continued

Caster					
Caster Type	Caster Size				
	Price	5"	6"	7"	8"
	.75" Lighted Rollerblade	£74	☐		
	1" Poly	Std	☐	☐	☐
	1" Poly Aluminum	£53	☐		
	1" Pneumatic	N/C		☐	☐
	1.5" Poly	£53	☐	☐	
	1.5" Soft Roll Aluminum	£53	☐	☐	
	2" Poly	£65		☐	

Axle

No selection needed with Drum Brake

- ☐ Fixed _____ N/C
☐ Quick Release _____ N/C

Wheel Locks

Wheel Lock

- ☐ Push to Lock _____ Std
☐ Pull to Lock _____ N/C
Not available with 12" or 16" Full Poly tires.
☐ Attendant Foot Lock _____ £105
Not available with 20" wheels. Mounted to the right side unless Foot Tilt - Right is selected.
☐ Drum Brake _____ £462
Available only with 12" Mag, 16" Mag, or MAXX Performance Spoke wheels. Not available with 16" Low Poly tire or Recline Back or Vent and Battery Tray or Hand Tilt. Axle is Quick Release.
☐ Omit Wheel Locks _____ N/C

Armrests

Armrest

Height Adjustable T-Arm Tall (11.5"- 15") **New!**

Not available with 12" Frame Depth

- ☐ Desk _____ £236
☐ Full _____ £236

Height Adjustable Pediatric

Best fit with Frame Depths less than 14"

- ☐ Desk (6.5"- 9") _____ N/C
☐ Desk Tall (8"-10.5") _____ N/C

Height Adjustable T-Arm (8.5"- 13.5")

Not available with 12" Frame Depth

- ☐ Desk _____ £236
☐ Full _____ £236

Height Adjustable T-Arm Low (6.5"- 10.5")

Not available with 12" Frame Depth. If requested configuration is unavailable, Standard Height T-Arm will be used.

- ☐ Desk _____ £236
☐ Full _____ £236

Rear Wheel

Rear Wheel Type	Rear Wheel Size			
	Price	12"	16"	20" ¹
	Mag (Black)	Std	☐	☐
	Colored Mag	£158		☐
	MAXX Performance Spoke	£105		☐

¹ Not available with Attendant Foot Lock

Colored Mag

Mag Color: ☐ Florescent Red ☐ White ☐ Blue ☐ Pink

Rear Tire Type

- ☐ Full Poly _____ Std
☐ Low Poly _____ N/C
Available only with 16" wheel size
☐ Pneumatic _____ N/C
☐ Pneumatic with Airless Insert _____ £105

Handrim

Handrims not available on 12" or 16" wheel sizes

Standard Handrims

- ☐ Aluminum Anodized _____ Std
☐ Plastic Coated _____ £84
☐ Omit Handrims _____ N/C

Dual - Drum Brake _____ £567
Available only with 12" Mag, 16" Mag, or MAXX Performance Spoke wheels. Not available with 16" Low Poly tire or Recline Back or Vent and Battery Tray or Hand Tilt. Axle is Quick Release.

- ☐ Push to Lock ☐ Pull to Lock

Hemi Wheel Lock _____ £263
Not available with 20" wheels. Extension Handle included.
☐ Left ☐ Right

Wheel Lock Options

- ☐ Extension Handles _____ £47
Available only with Push to Lock or Pull to Lock Wheel Locks.

Angle Adjustable Locking Extendable Flip Up (9.5"-15")

Not available with Seating Dynamics Rocker Back. Available only with 12" and 16" Rear Wheels.

- ☐ Desk _____ £236
☐ Full _____ £236

- ☐ Tubular Flip Up _____ Std
Not available with Seating Dynamics Rocker Back, Reclining Back or Arm Pads
☐ Omit Armrests _____ N/C

Arm Pad

- ☐ Standard Arm Pad _____ N/C
☐ Waterfall Arm Pad _____ £26
Not available with Angle Adjustable Locking Extendable Flip Up Armrest
☐ Foam Grip _____ N/C
Available only with Angle Adjustable Locking Extendable Flip Up Armrest

Accessories

Removable Side Guard

Not available with Height Adjustable Armrests

- ☐ Adult Composite Side Guards (14" - 16") £147
☐ Pediatric Composite Side Guards (10" - 14") £147

Anti-Tip

- ☐ Rear Anti-Tips £184
☐ Omit Rear Anti-Tips N/C

Pelvic Positioning Belt

- ☐ 1" Auto Buckle £45
☐ 1" Padded Auto Buckle £131
☐ 1 1/2" Auto Buckle £47
☐ 1 1/2" Padded Auto Buckle £131
☐ 2" Airline Buckle £90

Bodypoint Evoflex - Pelvic Stabilizer £163
Size: ☐ Small (16-24") ☐ Medium (19-29") ☐ Large (25-37")

Belt Mounting Kit

- ☐ Bodypoint 1" Band Clamp Belt Mounting Kit £89
☐ Seat Belt Mounting Hardware £14

Spoke Guards

- ☐ Clear Spoke Guards £131
Available only with MAXX Performance Spoke


Designer Spoke Guards £158
Available only with MAXX Performance Spoke

Design:

- ☐ Blue Monster ☐ Fairy Garden ☐ Havin' a Hoot
☐ Robo Fun ☐ Rockin' Reptiles ☐ Super Puppies
☐ Tutu Cute ☐ Tuxedo Time

Other Accessories

- ☐ Tool Kit £26
☐ Removable Underseat Pouch £53
☐ Pediatric Ki Mobility Backpack N/C
☐ Ki Mobility Backpack £74
☐ Black Neoprene Impact Guards £42
☐ Tilt Stop Posterior (Pair) £26
☐ Tilt Stop Anterior (Pair) £26
☐ IV Pole £278
Must remove for transit
☐ O2 Holder £263
Mounted to the left side. Must remove for transit.
☐ Canopy with Hardware £236
See image below



Color

Frame

Select desired base frame and seat frame / hanger colors. If no Seat Frame color is selected, Base Frame color will be used.

Gloss Finish		Base Frame	Seat Frame	Matte Finish		Base Frame	Seat Frame
Black	N/C	☐	☐	Charcoal Gray	N/C	☐	☐
Blue Skies	N/C	☐	☐	Emerald Green	N/C	☐	☐
Bubblegum Pink	N/C	☐	☐	Iridescent Burgundy	N/C	☐	☐
Candy Blue	N/C	☐	☐	Iridescent Navy	N/C	☐	☐
Candy Purple	N/C	☐	☐	Military Green	N/C	☐	☐
Candy Red	N/C	☐	☐	Retro Red	N/C	☐	☐
Copper	N/C	☐	☐	Royal Blue	N/C	☐	☐
Granny Smith	N/C	☐	☐	Shadow Black	N/C	☐	☐
Metallic Orange	N/C	☐	☐	Smooth White	N/C	☐	☐
Silver	N/C	☐	☐				
Teal	N/C	☐	☐				
Yellow	N/C	☐	☐				

Axiom G, S and P Cushions

Select size in chart at right. Weight capacity: 350 lbs.

Axiom G, S and P Cushions

- ☐ Axiom G £79
☐ Axiom S £110
☐ Axiom P £155

Axiom G, S or P Cushion Quantity: ☐ 1 ☐ 2

If no selection is made, quantity 1 will be selected

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Depth	14"	☐								
	15"	☐	☐	☐						
	16"	☐	☐	☐	☐	☐	☐	☐		
	17"	☐	☐	☐	☐	☐	☐	☐		
	18"	☐	☐	☐	☐	☐	☐	☐		☐
	19"		☐	☐	☐	☐	☐	☐		☐
	20"		☐	☐	☐	☐	☐	☐		☐
	21"			☐		☐	☐	☐		☐
	22"			☐		☐	☐	☐		☐

Axiom SP Cushions

Select size in chart at right. Weight capacity: 350 lbs.

Axiom SP Cushions

- ☐ Axiom SP - Visco £340
☐ Axiom SP - Fluid £315

Axiom SP Cushion Quantity: ☐ 1 ☐ 2

If no selection is made, quantity 1 will be selected

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Depth	14"	☐	☐	☐						
	15"	☐	☐	☐	☐	☐				
	16"	☐	☐	☐	☐	☐	☐	☐		☐
	17"	☐	☐	☐	☐	☐	☐	☐		☐
	18"	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19"	☐	☐	☐	☐	☐	☐	☐	☐	☐
	20"	☐	☐	☐	☐	☐	☐	☐	☐	☐
	21"		☐	☐	☐	☐	☐	☐	☐	☐
	22"		☐	☐	☐	☐	☐	☐	☐	☐

Axiom Cushion Accessories

Outer Covers

- ☐ Axiom G Extra Outer Cover £59
☐ Axiom S Extra Outer Cover £67
☐ Axiom P Extra Outer Cover £76
☐ Axiom SP Extra Visco Outer Cover £84
☐ Axiom SP Extra Fluid Outer Cover £84

Outer Cover Quantity: ☐ 1 ☐ 2

If no selection is made, quantity 1 will be selected

Inner Covers

- ☐ Axiom SP Extra Visco Inner Cover £84
☐ Axiom SP Extra Fluid Inner Cover £84

Inner Cover Quantity: ☐ 1 ☐ 2

If no selection is made, quantity 1 will be selected

Solid Seat Insert

Select size in chart below. Weight capacity: 350 lbs.

Solid Seat Insert £46

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Depth	14"	☐	☐	☐						
	15"	☐	☐	☐	☐					
	16"	☐	☐	☐	☐	☐	☐	☐		☐
	17"	☐	☐	☐	☐	☐	☐	☐		☐
	18"	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19"	☐	☐	☐	☐	☐	☐	☐	☐	☐
	20"	☐	☐	☐	☐	☐	☐	☐	☐	☐
	21"		☐	☐	☐	☐	☐	☐	☐	☐
	22"		☐	☐	☐	☐	☐	☐	☐	☐

Solid Seat Insert Quantity: ☐ 1 ☐ 2


If no selection is made, quantity 1 will be selected

Axiom Back


Select back that is the same or up to 2" narrower than the width of the frame. Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports. Weight capacity: 300 lbs.

Axiom Posterior Back (PB)

		Width					
		12"	14"	16"	18"	20"	22"
		Retail £499	£499	£499	£499	£499	£499
Height	Low	8"	☐ AXPB1208	☐ AXPB1408	☐ AXPB1608	☐ AXPB1808	☐ AXPB2008
		10"		☐ AXPB1410	☐ AXPB1610	☐ AXPB1810	☐ AXPB2010
		12"	☐ AXPB1212	☐ AXPB1412	☐ AXPB1612	☐ AXPB1812	☐ AXPB2012
	Mid	14"		☐ AXPB1414	☐ AXPB1614	☐ AXPB1814	☐ AXPB2014
		16"	☐ AXPB1216	☐ AXPB1416	☐ AXPB1616	☐ AXPB1816	☐ AXPB2016
		18"		☐ AXPB1418	☐ AXPB1618	☐ AXPB1818	☐ AXPB2018



Axiom Posterior Back - Low




Axiom Posterior Back - Mid

Axiom Posterior Back Measurement Reference				
		Back Height	A	B
Height	Low	8"	10"	6.5"
		10"	12"	8.5"
		12"	14"	10.5"
	Mid	14"	16"	8.5"
		16"	18"	9.5"
		18"	20"	10.5"


2.25"

Axiom Posterior Lateral Back (PL)


		Width					
		12"	14"	16"	18"	20"	22"
		Retail £530	£530	£530	£530	£530	£530
Height	Low	8"	☐ AXPL1208	☐ AXPL1408	☐ AXPL1608	☐ AXPL1808	☐ AXPL2008
		10"		☐ AXPL1410	☐ AXPL1610	☐ AXPL1810	☐ AXPL2010
		12"	☐ AXPL1212	☐ AXPL1412	☐ AXPL1612	☐ AXPL1812	☐ AXPL2012
	Mid	14"		☐ AXPL1414	☐ AXPL1614	☐ AXPL1814	☐ AXPL2014
		16"	☐ AXPL1216	☐ AXPL1416	☐ AXPL1616	☐ AXPL1816	☐ AXPL2016
		18"		☐ AXPL1418	☐ AXPL1618	☐ AXPL1818	☐ AXPL2018
	Tall	20"	☐ AXPL1220	☐ AXPL1420	☐ AXPL1620	☐ AXPL1820	☐ AXPL2020
		22"		☐ AXPL1422	☐ AXPL1622	☐ AXPL1822	☐ AXPL2022



Axiom Posterior Lateral Back - Low



Axiom Posterior Lateral Back - Mid



Axiom Posterior Lateral Back - Tall

Axiom Posterior Lateral Back Measurement Reference				
		Back Height	A	B
Height	Low	8"	10"	6.5"
		10"	12"	8.5"
		12"	14"	10.5"
	Mid	14"	16"	8.5"
		16"	18"	9.5"
		18"	20"	10.5"
	Tall	20"	22"	11.5"
		22"	24"	12.5"


3.65"

Axiom Back continued


Select back that is the same or up to 2" narrower than the width of the frame. Hardware is included with back. Weight capacity: 300 lbs.

Axiom Deep Lateral Back (DL)

			Width					
			12"	14"	16"	18"	20"	22"
			Retail	£610	£610	£610	£610	£610
Height	Mid	14"	☐ AXDL1214	☐ AXDL1414	☐ AXDL1614	☐ AXDL1814	☐ AXDL2014	☐ AXDL2214
		16"		☐ AXDL1416	☐ AXDL1616	☐ AXDL1816	☐ AXDL2016	
		18"	☐ AXDL1218	☐ AXDL1418	☐ AXDL1618	☐ AXDL1818	☐ AXDL2018	☐ AXDL2218
	Tall	20"	☐ AXDL1220	☐ AXDL1420	☐ AXDL1620	☐ AXDL1820	☐ AXDL2020	
		22"		☐ AXDL1422	☐ AXDL1622	☐ AXDL1822	☐ AXDL2022	☐ AXDL2222



Axiom Deep Lateral Back - Mid



Axiom Deep Lateral Back - Tall


Axiom Deep Lateral Back Measurement Reference									
		Back Height	A	B ¹					Shell Contour Depth
Height	Mid	14"	16"	8.25"	8.25"	8.25"	8.25"	8.25"	6"
	Tall	16"	18"	9.25"	9.25"	9.25"	9.25"	9.25"	
		18"	20"	10.25"	10.25"	10.25"	10.25"	10.25"	
		20"	22"	9.5"	9.75"	10"	10.25"	10.25"	
		22"	24"	10.75"	11"	11.25"	11.5"	11.75"	

¹ Measurement based on back width

Axiom Back Accessories

Headrest and Hardware

Choose headrest and hardware combination in chart below. Not available with Low Backs.



			Hardware					
			Standard Straight Back Mount Hardware with Ball		Multi Axis Removable Hardware with Horizontal Bar		Headrest Hardware with Ball	
			HEDHDW101 / HEDHDW117		HEDHDW201 / HEDHDW217		HEDHDW401 / HEDHDW417	
	Description	Item Code	Kio Code	Retail	Kio Code	Retail	Kio Code	Retail
Headrest Pad	Comfort Plus Headrest 8"	HEDPAD408	☐ HDR001	£233	☐ HDR009	£263	☐ HDR017	£326
	Comfort Plus Headrest 10"	HEDPAD410	☐ HDR002	£233	☐ HDR010	£263	☐ HDR018	£326
	Comfort Plus Headrest 14"	HEDPAD414	☐ HDR003	£233	☐ HDR011	£263	☐ HDR019	£326
	Small All Positioning Headrest Pad and Ring Only	HEDPAD201	☐ HDR004	£233	☐ HDR012	£263	☐ HDR020	£326
	Medium All Positioning Headrest Pad and Ring Only	HEDPAD202	☐ HDR005	£233	☐ HDR013	£263	☐ HDR021	£326
	Large All Positioning Headrest Pad and Ring Only	HEDPAD203	☐ HDR006	£233	☐ HDR014	£263	☐ HDR022	£326
	Tri-Comfort Plus Headrest 14"	HEDPAD502	☐ HDR007	£242	☐ HDR015	£272	☐ HDR023	£335
	Tri-Comfort Plus Headrest 19"	HEDPAD503	☐ HDR008	£242	☐ HDR016	£272	☐ HDR024	£335

Headrest Hardware Options

	Description	Item Code	Retail
☐	Headrest Flip Down Mount <i>Available only with HEDHDW401/HEDHDW417</i>	HDR100	£147
☐	Universal Mounting Plate <i>For use with other manufacturer headrests</i>	UMP	£37

Axiom Back Accessories continued

Lateral Support						
Pair	Left	Right	Description	Item Code	Pair	Each
Laterals for Posterior Back						
☐	☐	☐	3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back	LF3x4P	£256	£128
☐	☐	☐	4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back	LF4x5P	£256	£128
☐	☐	☐	3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back	LS3x4P	£546	£273
☐	☐	☐	4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back	LS4x5P	£546	£273
Laterals for Posterior Lateral Back						
☐	☐	☐	3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back	LF3x4L	£256	£128
☐	☐	☐	4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back	LF4x5L	£256	£128
☐	☐	☐	3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back	LS3x4L	£546	£273
☐	☐	☐	4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back	LS4x5L	£546	£273

Axiom Lateral Support Modifications

Extended lead times will apply. Modified Lateral Supports are not eligible for refund.

	Description	Item Code	Retail
☐	Modification for Chest Strap attached to Lateral Supports	AXMLATMT	N/C

Privacy Flap

	Description	Retail
☐	Privacy Flap	£38

Shoulder Guides

	Description	Item Code	Retail
☐	Shoulder Strap Guides and Mount <i>Selection of a modification from the Axiom Back Modifications required in order to mount Shoulder Guides. Not available with Low Backs.</i>	AXBSG	£37

Therafin Posture Support

Selection of a modification from the Axiom Back Modifications section required in order to mount posture supports. Not available with Low Backs.

	Description	Item Code	Retail
☐	Small Static Therafit 2-Piece Chest Strap with H & L Adj and 1" Safety Strap	60153T	£121
☐	Medium Static Therafit 2-Piece Chest Strap with H & L Adj and 1" Safety Strap	60154T	£121
☐	Large Static Therafit 2-Piece Chest Strap with H & L Adj and 1" Safety Strap	60155T	£121
☐	Small Dynamic Therafit Chest Strap with Dual D-Ring Adjustment	60110T	£95
☐	Medium Dynamic Therafit Chest Strap with Dual D-Ring Adjustment	60111T	£95
☐	Large Dynamic Therafit Chest Strap with Dual D-Ring Adjustment	60112T	£95
☐	Small Dynamic Therafit Comfortfit Full Vest	30589T	£140
☐	Medium Dynamic Therafit Comfortfit Full Vest	30590T	£140
☐	Large Dynamic Therafit Comfortfit Full Vest	30591T	£140
☐	XLarge Dynamic Therafit Comfortfit Full Vest	30592T	£140
☐	XXLarge Dynamic Therafit Comfortfit Full Vest	30593T	£140
☐	Small Dynamic Therafit Comfortfit Trim Vest	30618T	£142
☐	Medium Dynamic Therafit Comfortfit Trim Vest	30619T	£142
☐	Large Dynamic Therafit Comfortfit Trim Vest	30620T	£142
☐	XLarge Dynamic Therafit Comfortfit Trim Vest	30621T	£142
☐	XXLarge Dynamic Therafit Comfortfit Trim Vest	30622T	£142

Axiom Back Modifications

Extended lead times will apply for all modifications. Makes back non-refundable.

Axiom Back Modifications

Choose modification to add holes to the back shell in order to mount posture supports. Not available with Low Backs.

	Description	Item Code	Retail
☐	Modification for 4 Point Harness	AXM4PT	N/C
☐	Modification for 2 Point Chest Strap	AXM2PT	N/C
☐	Modification for Shoulder Guides	AXMSG	N/C

Painted Back Shell Modification

If no selection is made, standard Shadow Black will be used with no additional lead time

	Gloss Finish	Item Code	Retail
☐	Black	AXMBLK	N/C
☐	Blue Skies	AXMSKY	N/C
☐	Bubblegum Pink	AXMPNK	N/C
☐	Candy Blue	AXMCBL	N/C
☐	Candy Purple	AXMPUR	N/C
☐	Candy Red	AXMRED	N/C
☐	Copper	AXMCOP	N/C
☐	Granny Smith	AXMGSM	N/C
☐	Metallic Orange	AXMMOR	N/C
☐	Silver	AXMSIL	N/C
☐	Teal	AXMTEA	N/C
☐	Yellow	AXMYLW	N/C

	Matte Finish	Item Code	Retail
☐	Charcoal Gray	AXMCGM	N/C
☐	Emerald Green	AXMEGM	N/C
☐	Iridescent Burgundy	AXMIBM	N/C
☐	Iridescent Navy	AXMINM	N/C
☐	Military Green	AXMAGM	N/C
☐	Retro Red	AXMRRM	N/C
☐	Royal Blue	AXMRBM	N/C
☐	Smooth White	AXMSWM	N/C