

Backrest

Back Post / Height Range / Back Height						
Back Post	Height Range	Back Height				
Height Adjustable (Standard or 1" Inward Offset)	Extra Short [CLBOP024] 1,2	9"	10"	11"	12"	
	Short [CLBOP025] 2	11"	12"	13"	14" 1	
	Medium [CLBOP026] 2	13"	14"	15"	16"	
	Tall [CLBOP027] 2	15"	16"	17"	18"	
Fixed Height (Straight or 1" Inward Offset)	Fixed [CLBOP011] 3	9" 1	10" 1	11" 1	12"	13"
		14"	15"	16"	17"	18"
Fixed Height 2" Offset	Fixed Offset 2" 3 [CLBOP022]	11" 1	12" 1	13" 1	14" 1	
		15"	16"	17"		

- 1 Not available with Fold Down Push Handle
- 2 Each Offset Back Post is offset 1/2" creating 1" overall offset. Only Available with Standard or 1" Inward Offset.
- 3 Uses 1" tubes compatible with after-market backs. Not available with back upholstery. 2" Inward or Outward Offset Back Posts are offset 1" creating 2" overall offset. 1" Offset Back Post is offset 1/2" creating 1" overall offset

Push Handles

Omit Push Handle [CLBCK008]	N/C
Standard Push Handle [CLBCK004]	N/C
Fold Down Push Handle [CLBCK006]	£164
<i>Not available on Height Adjustable Back Posts with 9"-11" Back Heights or Offset Back Posts with 11" - 14" Back Heights.</i>	
Bolt-On Push Handle [CLBOP005]	£182
<i>Handles will be 1" - 2" above back height</i>	

Stroller Handles

Stroller Handles are removable and height adjustable.

Ergo Stroller Handle [CLBCK009]	£182
Single Stroller Handle [CLBOP007]	£164
Double Stroller Handle [CLBOP006]	£218
<i>Not available with 8"-10" Seat Widths.</i>	

Backrest Release Cable Color

Select desired color. If no selection is made, black will be used.

Black Backrest Cable [CLSBC008]	Red Backrest Cable [CLSBC005]
White Backrest Cable [CLSBC007]	Blue Backrest Cable [CLSBC006]
Yellow Backrest Cable [CLSBC004]	Orange Backrest Cable [CLSBC003]
Green Backrest Cable [CLSBC002]	Pink Backrest Cable [CLSBC001] N/C

Back Upholstery Options

Available only with Height Adjustable Back Posts. If selecting both upholstery and Axiom Back, the upholstery will be installed.

Padded Velcro Adjustable with Standard Cover [CLBUP001]	N/C
<i>Not available with 17" or 18" Back Heights in combination with 8" or 9" Seat Widths.</i>	
Padded Velcro Adjustable with Air Mesh Breathable Cover [CLBUP006]	£73
<i>Not available on chair widths less than 12"</i>	
Tension Adjustable with Standard Cover [CLBUP002]	N/C
<i>Not available with 8" and 9" Seat Widths. Not available on 10" Seat Width and 1" Inward Offset.</i>	
Tension Adjustable with Air Mesh Breathable Cover [CLBUP007]	£182
<i>Not available on chair widths less than 12".</i>	
Omit [CLBUP003]	N/C

Seat Back Angle

Standard is 90° (perpendicular to floor). Minus (-°) closes back angle. Plus (+°) opens back angle.

Standard [CLSBA001]	2° [CLSBA009]	4° [CLSBA010]
6° [CLSBA011]	8° [CLSBA012]	10° [CLSBA013]
-2° [CLSBA014]	-4° [CLSBA015]	-6° [CLSBA016]
-8° [CLSBA017]	-10° [CLSBA018]	N/C

Front Seat Height and Casters

Caster Size/ Front Seat Height

Select Front Seat Height in row corresponding to desired Caster Size.
Front Seat Height is measured from floor to top of seat tube at front of seat sling. Adjustable up or down by 1/2". See reference E on page 7.

Caster Size	Front Seat Height														
3"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"
4"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"
5"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"
6"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"		

Casters

		Caster Size				
		Price	3"	4"	5"	6"
Caster Type	.75" Rollerblade	N/C	[CLCAS002]	[CLCAS004]		
	.75" Lighted Rollerblade	£81	[CLCAS001]	[CLCAS005]	[CLCAS033]	
	1" Poly	N/C		[CLCAS006]	[CLCAS010]	[CLCAS014]
	1" Poly Aluminum	N/C		[CLCAS007]	[CLCAS011]	
	1" Pneumatic	£103				[CLCAS017]
	1.5" Poly	N/C		[CLCAS032]	[CLCAS025]	[CLCAS026]
	1.5" Soft Roll Aluminum	£97	[CLCAS003]	[CLCAS009]	[CLCAS013]	[CLCAS016]

Caster Options

Standard Fork [CLFRK001]	N/C
Single Sided Fork [CLFRK003]	£121
Frog Legs Suspension Fork [CLCST013]	£344
<i>Not included in Colored Anodizing Package.</i>	

User Weight Under 80 lbs [CLCST004]	N/C
User Weight 80-120 lbs [CLCST008]	N/C
User Weight 120-160 lbs [CLCST009]	N/C

Footrest

Footplates

Aluminum Angle Adjustable [CLFPL006]	£151
<i>Not available with Front Seat Heights less than 14.5".</i>	
Aluminum Angle Adjustable with Riser [CLFPL007]	£151
<i>Not available with Front Seat Heights less than 14.5".</i>	
High Mount Aluminum Angle Adjustable [CLFPL008]	£151
Angle Adjustable Flip Under [CLFPL012]	N/C
<i>Not available with footrest widths less than 9" or with Front Seat Heights less than 15". See reference F on page 8.</i>	
High Mount Angle Adjustable Flip Under [CLFPL011]	£151
<i>Not available with 16" Seat Width in combination with Offset Frame.</i>	

Seat to Footrest Length

Measured from front edge of seat sling to top, rear of footrest. Footrest length at least 2.5" shorter than Front Seat Height recommended. See reference G on page 7.

5.5"	6"	6.5"	7"	7.5"	8"
8.5"	9"	9.5"	10"	10.5"	11"
11.5"	12"	12.5"	13"	13.5"	14"
14.5"	15"	15.5"	16"	16.5"	17"
17.5"	18"				

Rear Seat Height and Wheels

Wheel Size / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size. For proper setup, selection is required when omitting rear wheels. Measured from floor to top of seat tube at front of back post. Rear Seat Height is custom to needs. See reference H on page 7.

Wheel Size	Rear Seat Height															
18"	11.5" ¹	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"
20"	12.5" ¹	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
22"	13.5" ¹	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"		
24"	14.5" ¹	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"				

¹ Available only with Low Profile tires

One Arm Drive

One Arm Drive - Quick Release [CLDRV007]	£846
<i>Available only with 20", 22" and 24" spoke wheels and Quick Release axles. Spoke wheel style is exclusive to OAD. Standard handrim will be Aluminum anodized. Rear Wheel Spacing limitations exist based on camber and armrest selections. Not available with 10" wide chairs.</i>	

Left Side Drive [CLDRV002]	N/C
Right Side Drive [CLDRV003]	N/C
Aluminum Drive Handrim [CLDRV004]	N/C
Poly Coated Drive Handrim [CLDRV005]	£182

Rear Wheels

No selection needed with One Arm Drive.

		Rear Wheel Size				
		Price	18"	20"	22"	24"
Rear Wheel Type	MAXX Spoke	N/C	[CLWH1118]	[CLWH1120]	[CLWH1122]	[CLWH1124]
	Mag	£121		[CLWH1220]	[CLWH1222]	[CLWH1224]
	Colored Mag ¹	£121		[CLWH1320]	[CLWH1322]	
	Spinergy Spox	£665		[CLWHE420]	[CLWHE422]	[CLWHE424]
	Spinergy LX	£937			[CLWHE522]	[CLWHE524]

¹ Color selection on page 6

Spinergy Options

Select spoke color and handrim tab for Spinergy wheels

Spoke Color

Black Spokes [CLSSC001]	Red Spokes [CLSSC002]	
White Spokes [CLSSC003]	Blue Spokes [CLSSC005]	N/C

Rear Tire Type

Full Profile

Street Pneumatic 1 3/8" [CLTIR001]	£110
Full Poly 1 3/8" [CLTIR002]	N/C
Pneumatic with Airless Insert 1 3/8" [CLTIR004]	£110

Low Profile

Low Poly 1" [CLTIR003]	N/C
Sentinel Solid Tire 1" [CLTIR005]	£110
<i>Not available on 18" wheel size or 20" and 22" Mag wheels</i>	
High Pressure 1" [CLTIR006]	N/C
<i>Not available on Mag wheel or 18" wheel size</i>	
Schwalbe Marathon Plus Evolution 1" [CLTIR010]	£110
<i>Available only on 24" MAXX Spoke or Spinergy wheels</i>	

Handrims

No selection needed with One Arm Drive.

Standard Handrims

Aluminum Anodized [CLHRM001]	N/C
Plastic Coated [CLHRM002]	£121
Projection (8 Vertical Projections) [CLHRM003]	£272
<i>Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels</i>	
Omit Handrim [CLHRM012]	N/C

Ergonomic Handrims

Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.

Ergo Lite Regular Grip [CLHRM005]	£272
Ergo Standard Regular Grip [CLHRM008]	£272
Ergo Lite Extra Grip [CLHRM006]	£272
Ergo Standard Extra Grip [CLHRM009]	£272
Ergo Lite Omit Grip [CLHRM004]	£272
Ergo Standard Omit Grip [CLHRM007]	£272

Handrim Options

3M Micro Grip [CLACC006]	£91
<i>Not available on Plastic Coated or Projection Handrims.</i>	

Handrim Mount

If no selection is made, Standard Mount will be used

Close Mount [CLTAB001]	N/C
------------------------	-----

Axles

No selection needed with Drum Brake.

Quick Release [CLAXO003]	N/C
Quad Quick Release [CLAXO001]	£182
<i>Not available with Mag or Spinergy wheels.</i>	

Center of Gravity presets

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post. Center of Gravity setup may be limited with certain configurations. 2" minimum Center of Gravity with 5th Wheel and One Arm Drive and 0" minimum with Transit Option. See reference I on page 7.

Amputee Positions				
-1.5" [CLCGA150]	-1" [CLCGA100]	-0.5" [CLCGA050]	0" [CLCGS000]	.25" [CLCGS025]
Standard Positions				
.5" [CLCGS050]	.75" [CLCGS075]	1" [CLCGS100]	1.25" [CLCGS125]	1.5" [CLCGS150]
1.75" [CLCGS175]	2" [CLCGS200]	2.25" [CLCGS225]	2.5" [CLCGS250]	2.75" [CLCGS275]
3" [CLCGS300]	3.25" [CLCGS325]	3.5" [CLCGS350]	3.75" [CLCGS375]	4" [CLCGS400]
4.25" [CLCGS425]				

Rear Wheel Spacing

Adjustable up to .5" outward. Minimum spacing with offset frame is 1.25" with 0° camber, 1" with 2° or 4° camber and .75" with 6° or 8° camber. .75" minimum with 0° camber on standard frame. 0° Camber and Mag wheel combination requires 1" minimum Rear Wheel Spacing. No selection needed with Wheelchair Drive Ready option. .5 and .75" not available with One Arm Drive. See Reference J on page 7.

0.5" [CLWS050]	0.75" [CLWS075]	1" [CLWS100]	
1.25" [CLWS125]	1.5" [CLWS150]	1.75" [CLWS175]	N/C

Rear Wheel Camber

The following apply if One Arm Drive is selected: 6° and 8° not available. 4° not available with T-Arms. 0° required with the combination of a 5th Wheel and Seat Widths less than 16"

0° [CLWC00]	2° [CLWC02]	4° [CLWC04]	6° [CLWC06]	
8° [CLWC08]				N/C

Wheel Locks

Wheel Lock

Push to Lock [CLWLK002]	N/C
Push to Lock Flush Mount [CLWLK001]	N/C
<i>Not available with Seat Depths less than 11" or Offset Frame.</i>	
Pull to Lock [CLWLK003]	£67
Low Profile Scissor Lock [CLWLK012]	£55
<i>Not available with Low Poly, Full Poly, Sentinel tires or Offset Hanger.</i>	
Short Thro Scissor Lock [CLWLK006]	£121
Drum Brake [CLWLK015]	£532
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Fold Down Push Handles or Ergo and Single Stroller Handles without integrated push handles. Axle is Quick Release.</i>	
Dual - Drum Brake with Push to Lock [CLWLK016]	£653
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Fold Down Push Handles or Ergo and Single Stroller Handles without integrated push handles. Axle is Quick Release.</i>	
Dual - Drum Brake with Pull to Lock [CLWLK017]	£653
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Fold Down Push Handles or Ergo and Single Stroller Handles without integrated push handles. Axle is Quick Release.</i>	
Omit [CLWLK008]	N/C

Wheel Lock Options

Extension Handles [CLWLK004]	£49
<i>Available only with Push to Lock or Pull to Lock Wheel Locks.</i>	
Grade Aids [CLWLK009]	£145
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic w/ Airless Insert tires.</i>	

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") ¹	Desk [CLARM039]	£266
	Full [CLARM040]	
Height Adjustable T-Arm (8.5"-13.5") ¹	Desk [CLARM003]	£266
	Full [CLARM004]	
Height Adjustable T-Arm Low (6.5"-10.5") ²	Desk [CLARM022]	£266
	Full [CLARM023]	
Angle Adjustable Locking Extendable Flip Up ³	Desk [CLARM010]	£266
	Full [CLARM011]	

Armrest Type		Price
Height Adjustable Pediatric Desk (6.5" - 10") ⁴	[CLARM014]	£242
Tubular Flip Up ⁵	[CLARM015]	N/C
Swing Away (7"-9") ⁶	[CLARM001]	£194
Tall Swing Away (9.25" - 11.25") ⁷	[CLARM021]	£194
Omit Armrests	[CLARM007]	N/C

- ¹ Not available with Seat Depths less than 13". 16" minimum Seat Depth with Offset Frame. 1.25" minimum Rear Wheel Spacing.
- ² Not available with Seat Depths less than 12". 14" minimum Seat Depth with Offset Frame. 1.25" minimum Rear Wheel Spacing. If requested configuration is unavailable, Standard Height T-Arm will be used.
- ³ Requires Fixed Height Backrest. Not available with backrest heights less than 12"
- ⁴ Seat Depth plus Sling Position must be at least 10". Seat Depth plus Sling Position must be at least 12" with Offset Frame 1.25" minimum Rear Wheel Spacing.
- ⁵ Requires Fixed Height Backrest. Not available with Arm Pads, Back Heights less than 13" or Back Heights less than 14" with Offset.
- ⁶ Not available with Arm Pads or Outward Offset Back Posts. If selection is not available with rear wheel and rear seat height configuration, then Tall Swing Away armrest will be used.
- ⁷ Not available with Arm Pads or Outward Offset Back Posts

Arm Pads

Not available with Tubular Flip Up or Swing Away Armrests.

Standard Arm Pad [CLARM026]	N/C
Waterfall Arm Pad [CLARM020]	£30
<i>Not available with Angle Adjustable Locking Extendable Flip Up Armrest.</i>	
Foam Grip [CLARM025]	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest.</i>	

Side Guards

Direct Mount Side Guards

Not available with T-Arms or Height Adjustable Pediatric Desk Armrests.

Composite Side Guards - Adult - Direct Mount Receiver [CLSDG002]	£140
<i>Seat Depth plus Sling Position must be at least 12".</i>	
Composite Side Guards - Pediatric - Direct Mount Receiver [CLSDG001]	£140
<i>Seat Depth plus Sling Position must be at least 10".</i>	

Clamp-On Mount Side Guards

10" minimum Seat Depth. Not available with T-Arms or Height Adjustable Pediatric Desk Armrests.

Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver - Black [CLSDG017]	£296
Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver [CLSDG018]	£357
Carbon Fiber Adjustable Fender - Clamp-On Receiver [CLSDG016]	£393

Accessories

5th Wheel

Rotating removable rear caster.

Standard 5th Wheel [CLATP002]	£206
Dynamic 5th Wheel [CLATP004]	£206
Soft Elastomer [CLELA001] Firm Elastomer [CLELA002]	N/C

Both elastomers are included; the selected elastomer will be installed.

Anti-Tip

Rear Anti-Tips [CLATP001]	£212
Omit [CLATP000]	N/C

Calf Strap

Not available with 8" or 9" wide chairs.

Velcro Adjustable [CLACC002]	N/C
Padded Velcro Adjustable [CLACC003]	£67
Bodypoint Padded AeroMesh [CLACC010]	£103

Pelvic Positioning Belts

1" Auto Buckle Pelvic Positioning Belt [CLPOB004]	£55
1" Padded Auto Buckle Pelvic Positioning Belt [CLPOB005]	£151
1 1/2" Auto Buckle Pelvic Positioning Belt [CLPOB001]	£55
1 1/2" Padded Auto Buckle Pelvic Positioning Belt [CLPOB002]	£151
Small Bodypoint Evoflex - Pelvic Stabilizer (16-24") [CLPOB007]	£188
Medium Bodypoint Evoflex - Pelvic Stabilizer (25-37") [CLPOB008]	£188
Large Bodypoint Evoflex - Pelvic Stabilizer (25-37") [CLPOB009]	£188

Belt Mounting Kit


Ki Mobility Belt Mounting Clamp [CLACC023]	£37
--	-----

Other Accessories

Tool Kit [CLACC005]	£20
Removable Underseat Pouch [CLACC015]	£30
Pediatric Ki Mobility Backpack [CLACC026]	£100
<i>Best fit on chairs more than 8" wide.</i>	
Ki Mobility Backpack [CLACC009]	£100
<i>Best fit on chairs more than 10" wide.</i>	
Luggage Carrier [CLACC025]	£73
Cane and Crutch Holder [CLACC007]	£91
Canopy with Hardware (K0108) [CLACC042]	£272
<i>See image to the right.</i>	

Spoke Guards

Clear Spoke Guards [CLSPG001]	£151
<i>Not available with Mag wheel.</i>	
Black Spoke Guards [CLSPG002]	£151
<i>Not available with 18" or 20" wheel sizes or Mag wheel.</i>	
Blue Monster Spoke Guards [CLSPG003]	£182
<i>Available only with 20" and 22" wheel sizes. Not available with Mag wheel.</i>	
Robo Fun Spoke Guards [CLSPG004]	£182
<i>Available only with 20" and 22" wheel sizes. Not available with Mag wheel.</i>	
Tutu Cute Spoke Guards [CLSPG005]	£182
<i>Available only with 20" and 22" wheel sizes. Not available with Mag wheel.</i>	
Tuxedo Time Spoke Guards [CLSPG008]	£182
<i>Available only with 20" and 22" wheel sizes. Not available with Mag wheel.</i>	



Color

Color

Frame/Accent: Select desired frame and accent colors. Accent color is used on backrest rigidizer bar and front caster wings; if no accent color is selected, frame color will be used.

Hydrographics: Available on frame only; please choose accent color above. Does not include front caster wing

Options	Price	Colors																			
		Gloss												Matte							
		Black	Blue Skies	Bubblegum Pink	Candy Blue	Candy Purple	Candy Red	Copper	Granny Smith	Metallic Orange	Silver	Teal	Yellow	Charcoal Gray	Emerald Green	Iridescent Burgundy	Iridescent Navy	Military Green	Retro Red	Royal Blue	Shadow Black
Frame	N/C																				
Accent	N/C																				
Mag Color																					

Color Anodizing Options

Select desired color. If no selection is made, black will be used.

	Color	Black ¹	Blue	Bright Green	Gold	Orange	Pewter	Pink	Purple	Red
Retail		N/C	£35	£35	£35	£35	£35	£35	£35	£35
Option	Back Plates ²									
	Camber Tube ³									
	Towers									
	Side Guard Receiver									
	Caster Fork ³									
	CG Clamp									
	5th Wheel Split Clamp									

¹ Black is no charge

² Not available with Transit Option

³ Available only with Standard or Single Sided

References

A. Seat Width See Diagram 1

Measured from outside of frame tube on one side to outside of frame tube on other side.

B. Seat Depth See Diagram 3

Measured from front of back posts to front edge of seat sling. Seat sling starts at beginning of bend at front of frame.

C. Front Frame Bend See Diagram 3

Angle between front seat tube and ground.

D. Sling Position See Diagram 2

With 0 Sling Position, sling is approximately 4" from the front tube. The 1" and 2" Sling Positions extend the frame in front of the sling by that amount. -1 pushes the seat 1" farther forward. See Charts 1 and 2 for more detail.

E. Front Seat Height See Diagram 1

Measured from floor to top of seat tube at front of seat sling.

F. Footrest Width See Diagram 1

Measured from inside of front tube to inside of front tube on other side. Footrest width is listed below for each seat width. See Chart 1 for more detail.

G. Seat to Footrest Length See Diagram 3

Measured from front edge of seat sling to top of rear of footrest. Footrest length of at least 2.5" shorter than front seat height recommended.

H. Rear Seat Height See Diagram 3

Measured from floor to top of seat tube at front of back post. Rear seat height is custom to needs.

I. Center of Gravity Preset See Diagram 3

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post.

J. Rear Wheel Spacing See Diagram 1

Measured from outside of seat back to inside of rear tire. Adjustable up to .5" outward.

Diagram 1




Diagram 2




Diagram 3




Chart 1

Seat Width	F. Footrest Width	
	Straight Frame (Diagram 1)	Offset Frame
8"	5.75"	9"
9"	6.75"	10"
10"	7.75"	11"
11"	8.75"	12"
12"	9.75"	13"
13"	10.75"	14"
14"	11.75"	15"
15"	12.75"	16"
16"	13.75"	17"

Chart 2

Sling Position / Adjustable Depth / Seat Depth Compatibility										
Seat Depth	Sling Position	-1", 0"			1"			2"		
		1"	2"	3"	1"	2"	3"	1"	2"	3"
14"	14"									
	15"						X		X	X
	16"			X		X	X	X	X	X

*X" indicates Sling Position, Adjustable Depth and Seat Depth combination is not available.
-1" Sling Position, 0" Adjustable Depth and 8" - 13" Seat Depths have no restrictions.

Axiom G

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover														
Retail		£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
	22"																

Axiom S

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover														
Retail		£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
	22"																

Axiom P

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover														
Retail		£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
	22"																

Axiom SP - Visco

		Width																										
		14"			15"			16"			17"			18"			19"			20"			21"			22"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover				
Retail		£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70			
Height	14"																											
	15"																											
	16"																											
	17"																											
	18"																											
	19"																											
	20"																											
	21"																											
	22"																											

Axiom SP - Fluid

		Width																										
		14"			15"			16"			17"			18"			19"			20"			21"			22"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover				
Retail		£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70			
Height	14"																											
	15"																											
	16"																											
	17"																											
	18"																											
	19"																											
	20"																											
	21"																											
	22"																											


Solid Seat Insert

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
		Cushion								
Retail		£55	£55	£55	£55	£55	£55	£55	£55	£55
Height	14"									
	15"									
	16"									
	17"									
	18"									
	19"									
	20"									
	21"									
	22"									

Axiom Back Width	Fits Back Frame Widths	
	Inner Slots	Outer Slots
12" (30.5cm)	11" (28.0cm)	12" -14" (30.5cm -35.5cm)
14" (35.5cm)	12" -13" (30.5cm -33.0cm)	14" -16" (35.5cm -40.5cm)
16" (40.5cm)	14" -15" (35.5cm -38.0cm)	16" -18" (40.5cm -45.5cm)
18" (45.5cm)	16" - 17" (40.5cm - 43.0cm)	18" - 20" (45.5cm - 50.5cm)
20" (51.0cm)	18" -19" (46.0cm -48.5cm)	20" -22" (50.5cm -55.5cm)
22" (56.0cm)	20" -21" (51.0cm -53.5cm)	22" -24" (56.0cm -61.0cm)

Hardware Description	Code
Improved DEX Hardware, removable adjustable hardware (300 lb max user weight)	DEX
New BTX Hardware, removable adjustable hardware not recommended for, or installed on Ultralight Rigid Frame Wheelchairs (400 lb max user weight)	BTX
New Fixed Hardware, non-removeable hardware up to (300 lb max user weight)	FXT



Axiom Back

Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports. Weight capacity: 300 lbs.

Axiom Posterior Back (PB)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPB1208]	[AXPB1408]	[AXPB1608]	[AXPB1808]	[AXPB2008]	
		10"		[AXPB1410]	[AXPB1610]	[AXPB1810]	[AXPB2010]	[AXPB2210]
		12"	[AXPB1212]	[AXPB1412]	[AXPB1612]	[AXPB1812]	[AXPB2012]	
	Mid	14"		[AXPB1414]	[AXPB1614]	[AXPB1814]	[AXPB2014]	[AXPB2214]
		16"	[AXPB1216]	[AXPB1416]	[AXPB1616]	[AXPB1816]	[AXPB2016]	
		18"		[AXPB1418]	[AXPB1618]	[AXPB1818]	[AXPB2018]	[AXPB2218]

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support for AXPB Backs

Lateral Trunk Supports for AXPB Backs							
Pair	Left	Right	Description	Item Code	Pair	Each	
			3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back	PBLF3x4P	£295	£148	
			4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back	PBLF4x5P	£295	£148	
			3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back	PBLS3x4P	£628	£314	
			4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back	PBLS4x5P	£628	£314	
			3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware	PBLAS3x4P	£628	£314	
			4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware	PBLAS4x5P	£628	£314	

Axiom Posterior Lateral Back (PL)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPL1208]	[AXPL1408]	[AXPL1608]	[AXPL1808]	[AXPL2008]	
		10"		[AXPL1410]	[AXPL1610]	[AXPL1810]	[AXPL2010]	[AXPL2210]
		12"	[AXPL1212]	[AXPL1412]	[AXPL1612]	[AXPL1812]	[AXPL2012]	
	Mid	14"		[AXPL1414]	[AXPL1614]	[AXPL1814]	[AXPL2014]	[AXPL2214]
		16"	[AXPL1216]	[AXPL1416]	[AXPL1616]	[AXPL1816]	[AXPL2016]	
		18"		[AXPL1418]	[AXPL1618]	[AXPL1818]	[AXPL2018]	[AXPL2218]
Tall	20"	[AXPL1220]	[AXPL1420]	[AXPL1620]	[AXPL1820]	[AXPL2020]		
	22"		[AXPL1422]	[AXPL1622]	[AXPL1822]	[AXPL2022]	[AXPL2222]	

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support for AXPL

Lateral Trunk Supports for AXPL & AXLB Backs						
Pair	Left	Right	Description	Item Code	Pair	Each
			3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back	PLLF3x4P	£295	£295
			4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back	PLLF4x5P	£295	£295
			3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back	PLLS3x4P	£628	£314
			4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back	PLLS4x5P	£628	£314
			3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware	PLLAS3x4P	£628	£314
			4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware	PLLAS4x5P	£628	£314

Axiom Deep Lateral Back (DL)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
HCPC		E2620	E2620	E2620	E2620	E2620	E2621	
Height	Mid	14"	[AXDL1214]	[AXDL1414]	[AXDL1614]	[AXDL1814]	[AXDL2014]	[AXDL2214]
		16"		[AXDL1416]	[AXDL1616]	[AXDL1816]	[AXDL2016]	
		18"	[AXDL1218]	[AXDL1418]	[AXDL1618]	[AXDL1818]	[AXDL2018]	[AXDL2218]
	Tall	20"	[AXDL1220]	[AXDL1420]	[AXDL1620]	[AXDL1820]	[AXDL2020]	
		22"		[AXDL1422]	[AXDL1622]	[AXDL1822]	[AXDL2022]	[AXDL2222]

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Headrest

Description	Std. Straight Removable Adjustable Hardware		Multi-Axis Removable Adjustable Hardware		Posalinc Removable Adjustable Hardware		HCPCS
	Item Code	Retail	Item Code	Retail	Item Code	Retail	
Comfort Soft Headrest 8"	[CS08SS]	£268	[CS08MA]	£268	[CS08PL]	£268	E0955 / E1028
Comfort Soft Headrest 10"	[CS10SS]	£268	[CS10MA]	£268	[CS10PL]	£268	E0955 / E1028
Comfort Soft Headrest 14"	[CS14SS]	£268	[CS14MA]	£268	[CS14PL]	£268	E0955 / E1028
Contour Positioning Pad, 8" X 6.25"	[CP08SS]	£303	[CP08MA]	£303	[CP08PL]	£303	E0955 / E1028
Contour Positioning Pad, 9" X 6.5"	[CP09SS]	£303	[CP09MA]	£303	[CP09PL]	£303	E0955 / E1028
Contour Positioning Pad, 10" X 7.5"	[CP10SS]	£303	[CP10MA]	£303	[CP10PL]	£303	E0955 / E1028
Tri-soft Headrest 10"	[TS10SS]	£386	[TS10MA]	£386	[TS10PL]	£386	E0955 / E1028
Tri-soft Headrest 14"	[TS14SS]	£386	[TS14MA]	£386	[TS14PL]	£386	E0955 / E1028

Description	Item Code	Retail	HCPCS
Universal Headrest Mounting Plate	[UMP]	£43	K0108

Not required for Axiom Headrests

Anterior Thoracic Supports

Selection of a modification from the Axiom Back Modifications section required in order to mount posture supports. Not available with Low Backs.

Anterior Thoracic Supports

Description	Item Code	HCPCS	Retail
Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small	[14941]	E0960	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium	[14942]	E0960	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large	[14943]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 20" Width, Small	[14924]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	[14927]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	[14928]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small	[14851]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Small	[14852]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Medium	[14853]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Large	[14854]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small	[14856]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small	[14857]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14858]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large	[14859]	E0960	£100
AirLogic Anterior Trunk Support: Standard, Structure, Small	[14282]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, Medium	[14283]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, Large	[14284]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, X-Large	[14285]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	[14257]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14258]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Large	[14259]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, X-Large	[14260]	E0960	£300

Accessories

Adjustable Shoulder Guides (K0108) [AXBSG] <i>Includes Anterior Thoracic Mounting Modification. Not available with Low Backs.</i>	£43
Extra Lumbar Pad (K0108) [LUMPAD] <i>Available only with Axiom Lateral Back and Axiom Deep Lateral back</i>	N/C
Privacy Flap (K0108) [PRFLAP]	£44

Axiom Back Modifications

Axiom Back Modifications

Extended lead times will apply for all modifications. Makes back non-refundable.

Anterior Thoracic Mounting Modification [AXSMM] <i>Provides holes for Anterior Thoracic Supports and Shoulder Guides</i>	N/C
---	-----

Gloss Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

Black Back Shell [AXMBLK]	N/C
Blue Skies Back Shell [AXMSKY]	N/C
Bubblegum Pink Back Shell [AXMPNK]	N/C
Candy Blue Back Shell [AXMCBL]	N/C
Candy Purple Back Shell [AXMPUR]	N/C
Candy Red Back Shell [AXMRED]	N/C
Copper Back Shell [AXMCOP]	N/C
Granny Smith Back Shell [AXMGSM]	N/C
Metallic Orange Back Shell [AXMMOR]	N/C
Yellow Back Shell [AXMYLW]	N/C
Silver Back Shell [AXMSIL]	N/C
Teal Back Shell [AXMTEA]	N/C

Matte Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

Military Green Back Shell [AXMAGM]	N/C
Charcoal Gray Back Shell [AXMCGM]	N/C
Emerald Green Back Shell [AXMEGM]	N/C
Iridescent Burgundy Back Shell [AXMIBM]	N/C
Iridescent Navy Back Shell [AXMINM]	N/C
Retro Red Back Shell [AXMRRM]	N/C
Royal Blue Back Shell [AXMRBM]	N/C
Smooth White Back Shell [AXMSWM]	N/C