

Loop Family



*Growable tiltable wheelchair
and seat shell base frame*

The new world of positioning




*Growable seat shell base
frame*



*Tiltable wheelchair for
individual positioning*

The seat shell base

German aid index N° 26.99.01.3034



back height
43 or 58 cm



seat depth
from 32 cm
(+6 cm)



weight capacity
120 kg
(+4 cm)



wheel sizes
12", 20" bis 24"



caster sizes
5", 5,5", 6", 7"



tilting angle
-5° till +35°
or
+2,5° till + 40°



Overview

- Seat shell base frame that grows with the patient
- Tilttable from -5° to +35° or from +2,5° to +40°
- Parameters such as seat height, seat position/centre of gravity, distance of foot plate/s etc. can be adjusted precisely to the user
- alternatively available with 12" wheels
- suitable as a seat for transport in a vehicle according to ISO 7176-19




Loop^{SORG} RS


The new development



seat width
30-48 cm
(+ 4 cm)




back height
30-60 cm
(+5 cm)



seat depth
30-48 cm
(max. +6 cm)



weight capacity
120 kg



wheel sizes
20"-24"



tilting angle
-3° to max. +34°



Overview

- Tiltable active wheelchairs
- The optionally available ERGO seat and back system offers an individual and best possible body guidance
- The 3-part ERGO back is growing in the width and offers space for a seat designed individually by a specialist shop
- suitable as a seat for transport in a vehicle according to ISO 7176-19

Unique people, unique solutions.

One name - two models

Family overview:


	Loop ^{SORG}	Loop ^{SORG} RS
seat depth	from 32 cm	30 - 48 cm
back height	43 or 58 cm	30 - 60 cm
max. weight capacity	120 kg	120 kg
wheel sizes	12", 20" - 24"	20" - 24"
caster sizes	5", 5,5", 6", 7"	5", 5,5", 6", 7"
tilting angle	✓	✓
growable in seat width, seat height and back height	✓	✓
suitable as a seat for transport in a vehicle according to ISO 7176-19	✓	✓
specialy adapted seat shell from the specialized supplier necessary	✓	

Equipment examples

All equipment options are on the order form of the corresponding wheelchair or sheat shell base




12" wheel



detachable leg supports



seat shell adapter with wedge
(only for Loop^{SORG})



drum brake



back can be folded down to the front



Outdoor Front End



ERGO-back



Push bail with synthetic hinge



ERGO-side guards with upholstery

And much more ...