

Medical Devices Register of the Health Ministry. TUV Certification No. 15-PEL-0010022-TIC







48 V



max 1500 W



12,3 Kg⁽¹⁾



up to 80 Km⁽²⁾



High Torque*

TECHNICAL FEATURES	EN
Max Power 1500W	48 V Lithium-ion battery single or double
Electronics by Triride	Battery switch (1-2)
Assisted Load & Reverse	Powerful disc brake
Intelligent Cruise Control ICC	City Tyre
Intelligent Braking System IBS	Stainless steel frame
Weight: 12,3 Kg	Carbon fender
16"x3" driving wheel	5 Speed levels
LCD Display + USB Port	Customized frame color

¹ Weight refers to device excluding battery and arms for wheelchair connection.
² Autonomy refers to double 48 V battery on a flat route at medium speed of 6 km/h (level 2) with user weight 70 kg. The autonomy depends on user weight, road conditions, battery model, speed/driving mode and on the general state of the wheelchair and tires.

Triride srl reserves the right to modify the models with updates and innovations and therefore the images may differ from the final product sold.